

epServices
for premium performance



BioFlo[®]/CelliGen[®] 510

BioFlo/CelliGen 510 Performance Plans

Bioprocess systems include several sophisticated technologies. Peak performance requires smooth interplay of the subsystems and fully functioning consumable parts. Regular maintenance by qualified engineers helps to ensure reliable operation, while deferring it can cause unreliable results, expensive repairs, and prolonged downtime.

The BioFlo/CelliGen 510 is a sterilize-in-place pilot/production-scale fermentor/bioreactor that can be placed

on a laboratory bench or on a mobile cart. The system can be controlled by the Eppendorf RPC controller or an Allen Bradley Programmable Logic Controller (PLC).

Eppendorf Bioprocess Performance Plans help you to maintain the optimal performance of your equipment over years of frequent use.

Bioprocess maintenance —
see the advantage!

eppendorf

Performance tested on

Model: _____

Serial no.: _____

Serviced by: _____

Service no.: _____

Next service: _____

Date: _____

epServices
for premium performance

BioFlo®/CelliGen® 510 Performance Plans

Service Operation	ESSENTIAL CHECK	ADVANCED MAINTENANCE	PREMIUM SERVICE	Installation Qualification (IQ)	Operational Qualification (OQ)	IQ/OQ
Catalog Number	0082140293	0082140294	0082140295	0082140297	0082140298	0082140299
List USD	\$9,216.00	\$11,650.00	\$13,750.00	\$1,969.00	\$3,275.00	\$4,900.00
External Inspection and Maintenance						
Check installation environment	■	■	■	■		■
Verify utility connections	■	■	■	■		■
Inspect addition pump heads and motors	■	■	■	■		■
Inspect motor cables and sensor cables	■	■	■	■		■
Check vessel for damage and leaks	■	■	■	■		■
Verify loop PI settings	■	■	■	■		■
Check valve operation	■	■	■	■		■
Check hoses	■	■	■	■		■
Internal Equipment						
Check motor and bearing housing for speed and operation	■	■	■		■	■
Rebuild bearing housing		■	■			
Replace o-rings, gaskets, and ferrules	■	■	■			
Rebuild solenoid, pneumatic, hand-operated, addition, sample, and harvest valves	■	■	■			
Inspect angle seat valves, gas sequencing valves, and circulating pump	■	■	■			
Replace rupture disk	■	■	■			
Test TMFC for operation	■	■	■		■	■
System Checks						
Calibrate pH and DO loops with simulator	■	■	■	■	■	■
Verify motor speed	■	■	■	■	■	■
Verify temperature	■	■	■			
Perform pressure hold	■	■	■	■	■	■
Perform test run	■	■	■	■	■	■
Check for latest software revision	■	■	■	■	■	■
Verify valve sequence	■	■	■		■	■
Documentation						
Check list provided	■	■	■			
Dated service sticker to confirm Eppendorf service	■	■	■	■	■	■
IQ report and signed documentation				■		■
OQ report and signed documentation					■	■
Supporting Information						
Travel	not included	not included	included	included	not included	included
Time/Labor	included	included	included	included	included	included
Additional parts	not included	not included	included	not included	not included	not included

Advantage OS

Description	Catalog No.	List USD
BioFlo/CelliGen 510	ADSERV-510	\$12,400.00

Advantage OS

An easy to use value packaged service option to help in advance and protect the system. Bioprocess Advantage OS provides labor and travel for up to two visits; an Advanced Maintenance and a single repair visit within the year. The plan also includes the Advanced Maintenance parts kit and a 20 % discount on parts utilized on a repair visit.

Eppendorf North America, Inc.
 Phone: 800-645-3050
 Email: info@eppendorf.com

www.eppendorf.com/epServices