

Pipette Calibration Service

Calibration Services – Ordering information

| Service | Calibration measurements »As Is« | Calibration measurements »As Left« | Outside decontamination ² | Checking, cleaning, lubricating ⁴ before calibration | Single channel pipette/dispenser | 4-channel pipette | 6-channel pipette | 8-channel pipette | 12-channel pipette | 16-channel pipette | 24-channel pipette |
|--|----------------------------------|------------------------------------|--------------------------------------|---|----------------------------------|-------------------|-------------------|-------------------|--------------------|--------------------|--------------------|
| ISO 8655 including PM⁵ | | 10 per test volume ¹ | included | included | 0082000142 | 0082000342 | 0082000343 | 0082000165 | 0082000166 | 0082000185 | 0082000186 |
| 3 x 4 including PM⁵ | | 4 per test volume ¹ | included | included | 0082000141 | 0082000340 | 0082000341 | 0082000163 | 0082000164 | 0082000183 | 0082000184 |

Calibration Service Packages – Ordering information

| Service | Calibration measurements »As Is« | Calibration measurements »As Left« | Outside decontamination ² checking, cleaning, lubricating ⁴ | ISO 17025 Accreditation calibration certificate ³ | Single channel pipette/dispenser | 4-channel pipette | 6-channel pipette | 8-channel pipette | 12-channel pipette | 16-channel pipette | 24-channel pipette |
|---|----------------------------------|------------------------------------|---|--|----------------------------------|-------------------|-------------------|-------------------|--------------------|--------------------|--------------------|
| ISO 8655 Verified | 4 per test volume ¹ | 10 per test volume ¹ | included | | 0082000234 | 0082000324 | 0082000325 | 0082000235 | 0082000236 | 0082000237 | 0082000238 |
| ISO 17025 Without PM⁵ | 10 per test volume ¹ | | | included | 0082000304 | 0082000305 | 0082000306 | 0082000300 | 0082000301 | 0082000302 | 0082000303 |
| ISO 17025 Advanced | 4 per test volume ¹ | 10 per test volume ¹ | included | included | 0082000227 | 0082000322 | 0082000323 | 0082000228 | 0082000229 | 0082000230 | 0082000231 |
| ISO 17025 Professional | 10 per test volume ¹ | 10 per test volume ¹ | included | included | 0082000220 | 0082000320 | 0082000321 | 0082000221 | 0082000222 | 0082000223 | 0082000224 |

Additional Services – Ordering information

Advanced Documentation, Excel file for calibration information on all pipettes/dispensers per order 0082000212

On-site verification services in your own laboratory Inquire

Additional Services Charges

Labor Charge of repair services Assessment fee will apply if repairing options are not accepted Inquire

¹ 10%, 50% and 100% of nominal volumes for variable volume instruments. ² Pipettes and dispensers and package equipment need to be decontaminated before sending in. ³ ISO 17025 accreditation calibration certificate is included in the second calibration after the maintenance where applicable. ⁴ Internal parts. ⁵ Preventative Maintenance. ⁶ Only available for new pipettes accredited calibration. [^] All prices in AUD and exclude GST.



Qualification and
Certification

ISO 8655

ISO 8655 defines the standard criteria by which a pipette should be calibrated and which maximum permissible error limits should be used at each volume. ISO 8655 criteria includes measurements at 3 volumes, 10 measurements at each volume, and controls for room environment.

ISO 17025

The ISO 17025 accreditation certifies the high level of competence and compliance with strict quality rules. ISO 17025 criteria include a traceability chain for measurement equipment, measurement uncertainty values, and approved quality standards.

Download a decontamination form at www.eppendorf.com/epServices.

Your local distributor: www.eppendorf.com/contact

Eppendorf South Pacific Pty. Ltd. · Level 1, 97 Waterloo Road · Macquarie Park, NSW 2113 · Australia info@eppendorf.com.au · www.eppendorf.com

www.eppendorf.com/epservices