

VisioNize® sense

Eppendorf Product Catalog p. 20	Order no.	Price
VisioNize® sense gateway router, routes the collected parameter data to VisioNize® Lab Suite	ESS1000001	609,00 €
VisioNize® sense environmental monitoring sensors, a wireless, rechargeable multi-sensor unit (ambient temperature, humidity, pressure & light)	ESS1000002	365,00 €
VisioNize® sense external temperature sensor, PT-100 extension that can measure temperature of any device (including ULT Freezers) within the range of +60 to -90 °C.	ESS1000003	609,00 €
VisioNize® sense monitoring kit , all-in-one package including temperature sensor, evironmental monitoring sensors and gateway router.	ESS1000004	1.583,00 €

Connected Electronic Pipettes

Eppendorf Product Catalog p. 31	Order no.	Price
VisioNize® pipette manager, an external touch server enabling the communication with connected electronic pipettes.	1004 000 001	945,00 €
Eppendorf Xplorer® connect , WiFi module to connect pipette to the VisioNize® pipette manager	4861 000 970	78,80 €

Move It® mechanical multi-channel pipettes, adjustable cone distance

Eppendorf Product Catalog p. 35	Order no.	Price
Eppendorf Research® plus Move It®, 4-channel, 9 mm cone distance, variable, incl. epT.I.P.S.® Box 2.0		
30 – 300 µL, orange	3125 000 150	1.379,06 €
Eppendorf Research® plus Move It®, 4-channel, 9 mm cone distance, variable		
120 – 1,200 µL, dark green	3125 000 184	1.379,06 €
Eppendorf Research® plus Move It®, 6-channel, 9 mm cone distance, variable, incl. epT.I.P.S.® Box 2.0		
30 – 300 µL, orange	3125 000 168	1.435,24 €
Eppendorf Research® plus Move It®, 6-channel, 9 mm cone distance, variable		
120 – 1,200 µL, dark green	3125 000 192	1.435,24 €
Eppendorf Research® plus Move It®, 8-channel, 4.5 mm cone distance, variable, incl. epT.I.P.S.® Box 2.0		
New1 – 20 µL, light pink	3125 000 117	1.481,00 €
New5 – 100 µL, light yellow	3125 000 133	1.481,00 €
Eppendorf Research® plus Move It®, 8-channel, 9 mm cone distance, variable, incl. epT.I.P.S.® Box 2.0		
30 – 300 µL, orange	3125 000 176	1.494,60 €
Eppendorf Research® plus Move It®, 8-channel, 9 mm cone distance, variable		
120 – 1,200 µL, dark green	3125 000 206	1.494,60 €
Eppendorf Research® plus Move It®, 12-channel, 4.5 mm cone distance, variable, incl. epT.I.P.S.® Box 2.0		
New1 – 20 µL, light pink	3125 000 125	1.594,00 €
New5 – 100 µL, light yellow	3125 000 141	1.594,00 €

Move It® electronic multi-channel pipettes, adjustable cone distance

Eppendorf Product Catalog p. 35	Order no.	Price
Eppendorf Xplorer® plus Move It®, 4-channel, 9 mm cone distance, variable		
15 – 300 µL, orange	4861 000 816	1.835,00 €
50 – 1,200 µL, green	4861 000 833	1.835,00 €
Eppendorf Xplorer® plus Move It®, 6-channel, 9 mm cone distance, variable		
15 – 300 µL, orange	4861 000 817	1.892,00 €
50 – 1,200 µL, green	4861 000 834	1.892,00 €
Eppendorf Xplorer® plus Move It®, 8-channel, 4.5 mm cone distance, variable		
1 – 20 µL, light pink	4861 000 781	1.949,00 €
5 – 100 µL, light yellow	4861 000 794	1.949,00 €
Eppendorf Xplorer® plus Move It®, 8-channel, 9 mm cone distance, variable		
15 – 300 µL, orange	4861 000 818	1.949,00 €
50 – 1,200 µL, green	4861 000 835	1.949,00 €
Eppendorf Xplorer® plus Move It®, 12-channel, 4.5 mm cone distance, variable		
1 – 20 µL, light pink	4861	

Eppendorf Research® plus, 1-channel, fixed

Eppendorf Product Catalog p. 36-38	Order no.	Price
Eppendorf Research® plus , 1-channel, fixed		
10 µL, medium gray	3124 000 016	169,00 €
10 µL, yellow	3124 000 024	169,00 €
20 µL, light gray	3124 000 032	169,00 €
20 µL, yellow	3124 000 040	169,00 €
25 µL, yellow	3124 000 059	169,00 €
50 µL, yellow	3124 000 067	169,00 €
100 µL, yellow	3124 000 075	169,00 €
200 µL, yellow	3124 000 083	169,00 €
200 µL, blue	3124 000 091	169,00 €
250 µL, blue	3124 000 105	169,00 €
500 µL, blue	3124 000 113	169,00 €
1,000 µL, blue	3124 000 121	169,00 €

Eppendorf Research® plus, 1-channel, variable

Eppendorf Product Catalog p. 36-38	Order no.	Price
Eppendorf Research® plus , 1-channel, variable, incl. epT.I.P.S.® Box		
0,1 – 2,5 µL, dark gray	3123 000 012	270,00 €
0,5 – 10 µL, medium gray	3123 000 020	270,00 €
2 – 20 µL, light gray	3123 000 098	270,00 €
2 – 20 µL, yellow	3123 000 039	270,00 €
10 – 100 µL, yellow	3123 000 047	270,00 €
20 – 200 µL, yellow	3123 000 055	270,00 €
30 – 300 µL, orange	3123 000 101	270,00 €
100 – 1,000 µL, blue	3123 000 063	270,00 €
Eppendorf Research® plus , 1-channel, variable, incl. epT.I.P.S.® sample bag		
0,25 – 2,5 mL, red	3123 000 144	270,00 €
0,5 – 5 mL, violet	3123 000 07	

Eppendorf Reference® 2, 1-channel, fixed

Eppendorf Product Catalog p. 40-41	Order no.	Price
Eppendorf Reference® 2 , 1-channel, fixed		
1 µL, dark gray	4925 000 014	200,00 €
2 µL, dark gray	4925 000 022	200,00 €
5 µL, medium gray	4925 000 030	200,00 €
10 µL, medium gray	4925 000 049	200,00 €
10 µL, yellow	4925 000 057	200,00 €
20 µL, light gray	4925 000 065	200,00 €
20 µL, yellow	4925 000 073	200,00 €
25 µL, yellow	4925 000 081	200,00 €
50 µL, yellow	4925 000 090	200,00 €
100 µL, yellow	4925 000 103	200,00 €
200 µL, yellow	4925 000 111	200,00 €
200 µL, blue	4925 000 120	200,00 €
250 µL, blue	4925 000 138	200,00 €
500 µL, blue	4925 000 146	200,00 €
1,000 µL, blue	4925 000 154	200,00 €
2 mL, red	4925 000 162	200,00 €
2,5 mL, red	4925 000 170	200,00 €

Eppendorf Reference® 2, 1-channel, variable

Eppendorf Product Catalog p. 40-41	Order no.	Price
Eppendorf Reference® 2 , 1-channel, variable, incl. epT.I.P.S.® Box		
0,1 – 2,5 µL, dark gray	4924 000 010	352,00 €
0,5 – 10 µL, medium gray	4924 000 029	352,00 €
2 – 20 µL, light gray	4924 000 037	352,00 €
2 – 20 µL, yellow	4924 000 045	352,00 €
10 – 100 µL, yellow	4924 000 053	352,00 €
20 – 200 µL, yellow	4924 000 061	352,00 €

Eppendorf Xplorer®, 12-channel, variable

Eppendorf Product Catalog p. 42-44	Order no.	Price
Eppendorf Xplorer®, 12-channel, variable		
0,5 – 10 µL, medium gray	4861 000 112	1.076,00 €
5 – 100 µL, yellow	4861 000 139	1.076,00 €
15 – 300 µL, orange	4861 000 155	1.076,00 €
50 – 1,200 µL, green	4861 000 171	1.076,00 €

Eppendorf Xplorer® plus, 1-channel, variable

Eppendorf Product Catalog p. 42-44	Order no.	Price
Eppendorf Xplorer® plus, 1-channel, variable		
0,5 – 10 µL, medium gray	4861 000 708	752,00 €
1 – 20 µL, light gray	4861 000 710	752,00 €
5 – 100 µL, yellow	4861 000 716	752,00 €
10 – 200 µL, yellow	4861 000 720	752,00 €
15 – 300 µL, orange	4861 000 724	752,00 €
50 – 1,000 µL, blue	4861 000 732	752,00 €
0,1 – 2,5 mL, red	4861 000 736	752,00 €
0,2 – 5 mL, violet	4861 000 740	752,00 €
0,5 – 10 mL, turquoise	4861 000 759	752,00 €

Eppendorf Xplorer® plus, 8-channel, variable

Eppendorf Product Catalog p. 42-44	Order no.	Price
Eppendorf Xplorer® plus, 8-channel, variable		
0,5 – 10 µL, medium gray	4861 000 767	1.107,00 €
5 – 100 µL, yellow	4861 000 783	1.107,00 €
15 – 300 µL, orange	4861 000 805	1.107,00 €
50 – 1,200 µL, green	4861 000 821	1.107,00 €

Eppendorf Xplorer® plus, 12-channel, variable

epT.I.P.S.® Standard

Eppendorf Product Catalog p. 54-61	Order no.	Price
epT.I.P.S.® Standard		
Eppendorf Quality™, 0,1 – 10 µL, 34 mm, dark gray, colorless tips, 1.000 tips (2 bags × 500 tips)	0030 000 811	63,20 €
Eppendorf Quality™, 0,1 – 20 µL, 40 mm, medium gray, colorless tips, 1.000 tips (2 bags × 500 tips)	0030 000 838	63,20 €
Eppendorf Quality™, 0,5 – 20 µL L, 46 mm, light gray, colorless tips, 1.000 tips (2 bags × 500 tips)	0030 000 854	60,10 €
Eppendorf Quality™, 2 – 200 µL, 53 mm, yellow, colorless tips, 1.000 tips (2 bags × 500 tips)	0030 000 889	38,90 €
Eppendorf Quality™, 2 – 200 µL, 53 mm, yellow, yellow tips, 1.000 tips (2 bags × 500 tips)	0030 000 870	38,90 €
Eppendorf Quality™, 20 – 300 µL, 55 mm, orange, colorless tips, 1.000 tips (2 bags × 500 tips)	0030 000 900	38,90 €
Eppendorf Quality™, 20 – 300 µL, 55 mm, orange, yellow tips, 1.000 tips (2 bags × 500 tips)	0030 000 897	38,90 €
Eppendorf Quality™, 50 – 1,000 µL, 71 mm, blue, colorless tips, 1.000 tips (2 bags × 500 tips)	0030 000 927	38,90 €
Eppendorf Quality™, 50 – 1,000 µL, 71 mm, blue, blue tips, 1.000 tips (2 bags × 500 tips)	0030 000 919	38,90 €
Eppendorf Quality™, 50 – 1,250 µL, 76 mm, green, colorless tips, 1.000 tips (4 bags × 250 tips)	0030 000 935	66,30 €
Eppendorf Quality™, 50 – 1,250 µL L, 103 mm, dark green, colorless tips, 1.000 tips (4 bags × 250 tips)	0030 000 730	67,20 €
Eppendorf Quality™, 0,25 – 2,5 mL, 115 mm, red, colorless tips, 500 tips (5 bags × 100 tips)	0030 000 951	45,60 €
Eppendorf Quality™, 0,1 – 5 mL, 120 mm, violet, colorless tips, 500 tips (5 bags × 100 tips)	0030 000 978	53,20 €
Eppendorf Quality™, 0,2 – 5 mL L, 175 mm, violet, colorless tips, 300 tips (3 bags × 100 tips)	0030 000 650	36,30 €
Eppendorf Quality™, 0,5 – 10 mL, 165 mm, turquoise, colorless tips, 200 tips		

epT.I.P.S.® LoRetention® Sets

Eppendorf Product Catalog p. 54-61	Order no.	Price
epT.I.P.S.® LoRetention® Sets		
Eppendorf Quality™, 0,1 – 10 µL S, 34 mm, dark gray, colorless tips, 480 tips (5 trays × 96 tips), 1 autoclavable reusable box	0030 072 367	POA
Eppendorf Quality™, 0,5 – 20 µL L, 46 mm, light gray, colorless tips, 480 tips (5 trays × 96 tips), 1 autoclavable reusable box	0030 072 375	POA
Eppendorf Quality™, 2 – 200 µL, 53 mm, yellow, colorless tips, 480 tips (5 trays × 96 tips), 1 autoclavable reusable box	0030 072 383	POA
Eppendorf Quality™, 50 – 1,000 µL, 71 mm, blue, colorless tips, 480 tips (5 trays × 96 tips), 1 autoclavable reusable box	0030 072 391	POA

epT.I.P.S.® Singles

Eppendorf Product Catalog p. 54-61	Order no.	Price
epT.I.P.S.® Singles		
Biopur®, 0,1 – 20 µL, 40 mm, medium gray, colorless tips, 100 tips, individually packed	0030 010 027	30,60 €
Biopur®, 2 – 200 µL, 53 mm, yellow, colorless tips, 100 tips, individually packed	0030 010 043	30,60 €
Biopur®, 50 – 1,000 µL, 71 mm, blue, colorless tips, 100 tips, individually packed	0030 010 060	30,60 €

epT.I.P.S.® Racks

Eppendorf Product Catalog p. 54-61	Order no.	Price
epT.I.P.S.® Racks		
Biopur®, 0,1 – 20 µL, 40 mm, medium gray, colorless tips, 480 tips (5 racks × 96 tips)	0030 075 226	58,90 €
Biopur®, 2 – 200 µL, 53 mm, yellow, colorless tips, 480 tips (5 racks × 96 tips)	0030 075 234	54,80 €
Biopur®, 20 – 300 µL, 55 mm, orange, colorless tips, 480 tips (5 racks × 96 tips)	0030 075 242	54,80 €
Biopur®, 50 – 1,000 µL, 71 mm, blue, colorless tips, 480 tips (5 racks × 96 tips)	0030 075 250	54,80 €

ep Dualfilter T.I.P.S.® SealMax, PCR clean/sterile (sterile and pyrogen-free)

Eppendorf Product Catalog p. 54-61	Order no.	Price
ep Dualfilter T.I.P.S.® SealMax		
PCR clean and sterile, 0,1 – 10 µL S, 34 mm, dark gray, colorless tips, 960 tips (10 racks × 96 tips)	0030 078 691	183,00 €
PCR clean and sterile, 0,5 – 20 µL L, 46 mm, light gray, colorless tips, 960 tips (10 racks × 96 tips)	0030 078 705	175,00 €
PCR clean and sterile, 2 – 100 µL, 53 mm, yellow, colorless tips, 960 tips (10 racks × 96 tips)	0030 078 713	175,00 €
PCR clean and sterile, 2 – 200 µL, 55 mm, yellow, colorless tips, 960 tips (10 racks × 96 tips)	0030 078 721	182,00 €
PCR clean and sterile, 20 – 300 µL, 55 mm, orange, colorless tips, 960 tips (10 racks × 96 tips)	0030 078 730	175,00 €
PCR clean and sterile, 50 – 1,000 µL, 76 mm, blue, colorless tips, 960 tips (10 racks × 96 tips)	0030 078 748	190,00 €

GELoader®

Eppendorf Product Catalog p. 54-61	Order no.	Price
GELoader® , Eppendorf Quality™, 0,5 – 20 µL, 62 mm, light gray, 192 tips (2 racks × 96 tips)	0030 001 222	69,10 €

Microloader™

Eppendorf Product Catalog p. 54-61	Order no.	Price
Microloader™ , tip for filling Femtotips® and other glass microcapillaries (for research use only), 0,5 – 20 µL, 100 mm, light gray, 192 pcs. (2 racks × 96 pcs.)	5242 956 003	94,40 €

epT.I.P.S.® 384

Eppendorf Product Catalog p. 54-61	Order no.	Price
epT.I.P.S.® 384 Reloads		
PCR clean, 0,1 – 20 µL, 42 mm, light pink, colorless tips, 3.840 tips (10 trays × 384 tips)	0030 076 001	525,00 €
PCR clean, 5 – 100 µL,		

Tips

Combitips® advanced

Eppendorf Product Catalog p. 66-67	Order no.	Price
Combitips® advanced		
Eppendorf Quality™, 5,0 mL, blue, colorless tips, 100 pcs. (4 bags × 25 pcs.)	0030 089 456	96,10 €
Biopur®, 5,0 mL, blue, colorless tips, 100 pcs., individually blister-wrapped	0030 089 669	176,00 €
Forensic DNA Grade, 5,0 mL, blue, colorless tips, 100 pcs., individually blister-wrapped	0030 089 871	135,00 €
PCR clean, 5,0 mL, blue, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)	0030 089 812	121,00 €
Eppendorf Quality™, 10 mL, orange, colorless tips, 100 pcs. (4 bags × 25 pcs.)	0030 089 464	96,10 €
Biopur®, 10 mL, orange, colorless tips, 100 pcs., individually blister-wrapped	0030 089 677	176,00 €
PCR clean, 10 mL, orange, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)	0030 089 820	121,00 €
Eppendorf Quality™, 25 mL, red, colorless tips, 100 pcs. (4 boxes × 25 pcs.), each box containing one adapter	0030 089 472	113,00 €
Biopur®, 25 mL, red, colorless tips, 100 pcs. (4 bags × 25 pcs.), individually blister-wrapped, each box containing one adapter	0030 089 685	207,00 €
PCR clean, 25 mL, red, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.), each box containing one adapter	0030 089 839	143,00 €
Eppendorf Quality™, 50 mL, light gray, colorless tips, 100 pcs. (4 boxes × 25 pcs.), each box containing one adapter	0030 089 480	125,00 €
Biopur®, 50 mL, light gray, colorless tips, 100 pcs. (4 bags × 25 pcs.), individually blister-wrapped, each box containing one adapter	0030 089 693	223,00 €
PCR clean, 50 mL, light gray, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.), each box containing one adapter	0030 089 847	154,00 €
Adapter advanced, 25 mL, autoclavable		
Eppendorf Quality™, red, 1 piece	0030 089 715	10,60 €
Biopur®, red, 7 pcs., individually wrapped	0030 089 731	77,20 €
Adapter advanced, 50 mL, autoclavable		
Eppendorf Quality™, light gray, 1 piece	0030 089 723	10,60 €
Biopur®, light gray, 7 pcs., individually wrapped	0030 089 740	77,20 €</

Tips

Varispenser® 2/Varispenser® 2x

Eppendorf Product Catalog p. 70-71	Order no.	Price
Varispenser® 2 , 1-channel, bottle top dispenser with adapters GL 25, GL 28/S 28, GL 32, GL 38, S 40 and telescopic intake tube (length 125 – 240 mm)		
0.2 – 2 mL	4966 000 010	286,00 €
0.5 – 5 mL	4966 000 029	286,00 €
1 – 10 mL	4966 000 037	286,00 €
Varispenser® 2 , 1-channel, bottle top dispenser with adapters GL 32, GL 38, S 40 and telescopic intake tube (length 170 – 330 mm)		
2.5 – 25 mL	4966 000 045	380,00 €
5 – 50 mL	4966 000 053	418,00 €
10 – 100 mL	4966 000 061	550,00 €
Varispenser® 2x , 1-channel, bottle top dispenser with recirculation valve with adapters GL 25, GL 28/S 28, GL 32, GL 38, S 40 and telescopic intake tube (length 125 – 240 mm)		
0.2 – 2 mL	4967 000 014	338,00 €
0.5 – 5 mL	4967 000 022	338,00 €
1 – 10 mL	4967 000 030	338,00 €
Varispenser® 2x , 1-channel, bottle top dispenser with recirculation valve with adapters GL 32, GL 38, S 40 and telescopic intake tube (length 170 – 330 mm)		
2.5 – 25 mL	4967 000 049	451,00 €
5 – 50 mL	4967 000 057	488,00 €
10 – 100 mL	4967 000 065	612,00 €

Easypet® 3

Eppendorf Product Catalog p. 72	Order no.	Price
Easypet® 3 , 1-channel, incl. mains/power supply device, wall mounting device, shelf stand, 2 membrane filters 0.45 µm, 0,1 – 100 mL	4430 000 018	494,00 €
Membrane filter , for Easypet®3 , 0.2 µm, sterile, 5 pcs.	4430 606 005	103,00 €
Membrane filter , for Easypet® 3, 0.45 µm, sterile, 5 pcs.	4421 601 009	103,00 €
Lithium polymer		

Automated Pipetting

epMotion® 5075I

Eppendorf Product Catalog p. 84

epMotion® 5075I, basic device incl. epBlue™ software, mouse, waste box, 100 – 240 V ±10 %/50 – 60 Hz ±5 %, 0,2 µL – 1 mL

Order no.	Price
5075 000 301	POA

epMotion® 5075t

Eppendorf Product Catalog p. 85

epMotion® 5075t, basic device incl. Eppendorf ThermoMixer®, epBlue™ software, mouse, waste box, 100 – 240 V ±10 %/50 – 60 Hz ±5 %, 0,2 µL – 1 mL

Order no.	Price
5075 000 302	POA

epMotion® 5075v/5075vt

Eppendorf Product Catalog p. 85

epMotion® 5075v, basic device incl. vacuum system, gripper, vac frame 2, vac frame holder, epBlue™ software, mouse, waste box, 100 – 240 V ±10 %/50 – 60 Hz ±5 %, 0,2 µL – 1 mL

epMotion® 5075vt, basic device incl. vacuum system, gripper, vac frame 2, vac frame holder, Eppendorf ThermoMixer®, epBlue™ software, mouse, waste box, 100 – 240 V ±10 %/50 – 60 Hz ±5 %, 0,2 µL – 1 mL

Order no.	Price
5075 000 303	POA
5075 000 304	POA

epMotion® 5075 NGS solution

Eppendorf Product Catalog p. 86

epMotion® 5075t NGS solution, package with completely contained housing, MultiCon PC, Enhanced feature set 1, C2 thermal module, dispensing tools, plus NGS specific accessories, plus NGS specific consumables to start automated library preparation, 100 – 240 V ±10 %/50 – 60 Hz ±5 % (EU), 0,2 µL – 1 mL

epMotion® 5075tc NGS solution, package with CleanCap, MultiCon PC, Enhanced feature set 1, C2 thermal module, dispensing tools, plus NGS specific accessories, waste bag holder, plus NGS specific consumables to start automated library preparation,

Eppendorf ThermoMixer® C

Eppendorf Product Catalog p. 98

	Order no.	Price
Eppendorf ThermoMixer® C, basic device without thermoblock, 220 – 240 V/50 – 60 Hz (EU)	5382 000 015	2.804,00 €

Eppendorf ThermoStat™ C

Eppendorf Product Catalog p. 99

	Order no.	Price
Eppendorf ThermoStat™ C, basic device without thermoblock, 220 – 240 V/50 – 60 Hz (EU)	5383 000 019	1.917,00 €
ThermoTop®, with condens.protect technology	5308 000 003	491,00 €
SmartExtender, 1.5 mL	5322 000 008	764,00 €

Eppendorf ThermoMixer® F

Eppendorf Product Catalog p. 100

	Order no.	Price
Eppendorf ThermoMixer® FP, with thermoblock for Microplates and Deepwell Plates, including lid, 220 – 240 V/50 – 60 Hz (EU)	5385 000 016	2.863,00 €
Eppendorf ThermoMixer® F0.5, with thermoblock for 24 reaction vessels 0.5 mL, 200 – 240 V/50 – 60 Hz (EU)	5386 000 010	2.749,00 €
Eppendorf ThermoMixer® F1.5, with thermoblock for 24 reaction vessels 1.5 mL, 220 – 240 V/50 – 60 Hz (EU)	5384 000 012	2.682,00 €
Eppendorf ThermoMixer® F2.0, with thermoblock for 24 reaction vessels 2.0 mL, 220 – 240 V/50 – 60 Hz (EU)	5387 000 013	2.682,00 €

MixMate®

Eppendorf Product Catalog p. 101

	Order no.	Price
MixMate®, basic device without tube holders, 230 V/50 – 60 Hz (EU)	5353 000 510	1.857,00 €
Tube Holder 0.5 mL, for 24 × 0.5 mL reaction vessels	5353 040 121	98,30 €
Tube Holder 1.5/2.0 mL, for 24 × 1		

CellXpert® C170

Eppendorf Product Catalog p. 264	Order no.	Price
CellXpert® C170, non-segmented inner door, handle right side, 220 – 240 V/50 – 60 Hz (EU)	6734 000 011	8.232,17 €
CellXpert® C170, inner door with 4 door segments, handle right side, 220 – 240 V/50 – 60 Hz (EU)	6734 000 041	10.002,46 €

New Brunswick™ S41i

Eppendorf Product Catalog p. 265	Order no.	Price
New Brunswick™ S41i, 170 L, CO ₂ incubator shaker with inner shelf and touch screen control, 1 shelf, stackable up to 2-fold, 1 (2. optional) shelves, 230 V/50 Hz, orbit 2.5 cm (1 in)	S41I230011	15.767,50 €
Interchangeable platform for New Brunswick™ S41i, 36 × 61 cm (14 × 24 in)		
Universal platform	M1334-9920	684,00 €
125 mL Erlenmeyer flask dedicated platform, comes complete with 32 clamps installed	M1334-9921	785,00 €
250 mL Erlenmeyer flask dedicated platform, comes complete with 24 clamps installed	M1334-9922	736,00 €
500 mL Erlenmeyer flask dedicated platform, comes complete with 15 clamps installed	M1334-9923	724,00 €
1 L Erlenmeyer flask dedicated platform, comes complete with 11 clamps installed	M1334-9924	715,00 €
2 L Erlenmeyer flask dedicated platform, comes complete with 6 clamps installed	M1334-9925	715,00 €
2.8 L Erlenmeyer flask dedicated platform, comes complete with 3 clamps installed	M1334-9926	715,00 €
4 L Erlenmeyer flask dedicated platform, comes complete with 3 clamps installed	M1334-9927	715,00 €
Stacking kit, allows up to 2 each New Brunswick™ S41i shakers to be stacked one on top of the other	P0628-6502	1.255,00 €
Additional shelf for New Brunswick™ S41i, perforated, can easily be added to the New Brunswick™ S41i for static cell culture of adherent cells	P0628-6181	258,00 €
Sticky pad adapter kit, 2 anodized aluminum sheets, each 15.2 × 45.7 cm (6 × 18 in)	M1250-9504	191,00 €
Sticky pad, 20 × 20 cm (8 × 8 in), Compatible platform or adapter kit necessary	M1250-9700	79,50 €

Galaxy® 48 R

Eppendorf Product Catalog p

Clamps and Accessories for Universal Platforms

Eppendorf Product Catalog p. 283	Order no.	Price
Test tube rack , small, 30 × 8 × 12.9 cm		
diameter 8 – 11 mm, 48 tubes	M1289-0001	321,00 €
diameter 12 – 15 mm, 34 tubes	M1289-0002	321,00 €
diameter 15 – 18 mm, 24 tubes	M1289-0003	321,00 €
diameter 18 – 21 mm, 18 tubes	M1289-0004	321,00 €
diameter 22 – 26 mm, 13 tubes	M1289-0005	321,00 €
diameter 26 – 30 mm, 12 tubes	M1289-0006	321,00 €
Test tube rack , medium, 36 × 8 × 12.9 cm		
diameter 8 – 11 mm, 60 tubes	M1289-0010	321,00 €
diameter 12 – 15 mm, 44 tubes	M1289-0020	321,00 €
diameter 15 – 18 mm, 31 tubes	M1289-0030	321,00 €
diameter 18 – 21 mm, 23 tubes	M1289-0040	321,00 €
diameter 22 – 26 mm, 16 tubes	M1289-0050	321,00 €
diameter 26 – 30 mm, 16 tubes	M1289-0060	321,00 €
Test tube rack , large, 45 × 8 × 12.9 cm		
diameter 8 – 11 mm, 80 tubes	M1289-0100	321,00 €
diameter 12 – 15 mm, 60 tubes	M1289-0200	321,00 €
diameter 15 – 18 mm, 42 tubes	M1289-0300	321,00 €
diameter 18 – 21 mm, 30 tubes	M1289-0400	321,00 €
diameter 22 – 26 mm, 22 tubes	M1289-0500	321,00 €
diameter 26 – 30 mm, 20 tubes	M1289-0600	321,00 €
Clamp , for 1 L media bottles	ACSB-1000S	102,00 €
Clamp , for 500 mL media bottles and TPP Bioreactors 450 mL or 600 mL	ACSB-500S	102,00 €
Angled test tube rack holder , for user-supplied test tube racks and microplates	TTR-210	335,00 €
Angled Test Tube Rack Spacer Bar , for use with TTR-210		
to accommodate test tubes racks that are less than 13 mm (5 in) wide		

Application and Technical support



North America:

+1 800 645 3050, menu option 2
techserv@eppendorf.com

Asia Pacific Africa:

+603 8023 6869
support_asiapacific@eppendorf.com



Your local distributor:

www.eppendorf.com/contact



Eppendorf Online Shop

www.eppendorf.com/eshop

Eppendorf® and the Eppendorf Brand Design are registered trademarks of Eppendorf SE, Germany. U.S. Design Patents are listed on www.eppendorf.com/ip. All rights reserved, including graphics and images. EN1/INT/Price-List_2022/032022/CCO_MEP/
Konradin Heckel. Copyright© 2022 by Eppendorf SE.

Eppendorf Bioprocessing-Catalog. More info at www.eppendorf.com/catalog



Your local distributor: www.eppendorf.com/contact

Eppendorf SE · Barkhausenweg 1 22339 · Hamburg · Germany

Eppendorf North America, Inc. · USA · Phone: +1 800 645 3050 · Email: info@eppendorf.com