

Test Report

– Eppendorf Forensic DNA Grade –

according to ISO 18385

Blank example certificate

• **Human DNA**

Method: *real-time* PCR acc. to SOP 111141/V1
Limit: < 0.5 pg/ μ L

• **PCR inhibition**

Method: *real-time* PCR acc. to SOP 100559/V1
Limit: Less than 10 targets amplifiable

• **DNase**

Method: DNA Digestion acc. to SOP 100570/V1
Limit: Not detectable
(LOD: 1.0×10^{-6} Kunitz units)

• **RNase**

Method: RNA Digestion acc. to SOP 100571/V1
Limit: Not detectable
(LOD: 1.0×10^{-9} Kunitz units)

Customer: Eppendorf SE, Barkhausenweg 1, D-22339 Hamburg

Date of sample receipt / examination:

Testing of the above described lot showed conformity within the limits of detection.

.....
Date of release

.....
Registration number

.....
Dipl.-Biol. Thomas Weigel (Head)

The test results refer exclusively to the items tested. The certificate must not be copied partially without the approval of LADR GmbH.