

Test Report

– Eppendorf Sterile –

Blank example certificate

• **Sterility**

Method:

Sterilization with SAL 10^{-6} by Irradiation according DIN EN ISO 11137:2015* or ethylene oxide*

*please see attached label for the standard method of sterilization for this product

Test method:

USP, Ph. Eur. 8.0 / 2.6.12

• **Endotoxin**

Test method:

Kinetic-chromogenic LAL test (Ph. Eur. 8.0 / 2.6.14)

Limit:

< 0.001 EU/mL

Customer: Eppendorf SE, Barkhausenweg 1, D-22339 Hamburg

Date of sample receipt / examination:

Testing of the above described lot showed conformity within the limits of detection.

.....
Date of release

.....
Registration number

.....
Dipl.-Biol. Thomas Weigel (Head)

The test results refer exclusively to the items tested. The certificate must not be copied partially without the approval of LADR GmbH.