


Test Report

– Eppendorf Sterile –

Blank example certificate

• **Sterility**

Method: Sterilization with SAL 10^{-6} by  Irradiation according DIN EN ISO 11137
Test method: USP, Ph. Eur. 10.0 / 2.6.12

• **Endotoxin**

Test method: Kinetic-chromogenic LAL test (Ph. Eur. 10.0 / 2.6.14)
Limit: < 0.001 EU/mL

Customer: Eppendorf SE, Barkhausenweg 1, D-22339 Hamburg

Date of sample receipt / examination:

Testing of the rinse solution of the above described lot showed conformity within the limits of detection.

.....
Date of release

.....
Registration number

.....
Dr. rer. nat. Sybille Kuhnert (Head)

The test results refer exclusively to the items tested. The test report must not be copied partially without the approval of LADR GmbH.