



# Price List 2023

## Digital Lab Management

### VisioNize® sense

Eppendorf Product Catalog p. 20

|   | <b>Order no.</b> | <b>Price</b> |
|---|------------------|--------------|
| <b>VisioNize® sense gateway router</b> , routes the collected parameter data to VisioNize® Lab Suite  | ESS1000001       | 649,00 €     |
| <b>VisioNize® sense environmental monitoring sensors</b> , a wireless, rechargeable multi-sensor unit (ambient temperature, humidity, pressure & light)                       | ESS1000002       | 389,00 €     |
| <b>VisioNize® sense external temperature sensor</b> , PT-100 extension that can measure temperature of any device (including ULT Freezers) within the range of +60 to -90 °C. | ESS1000003       | 649,00 €     |
| <b>VisioNize® sense monitoring kit</b> , all-in-one package including temperature sensor, environmental monitoring sensors and gateway router.                                | ESS1000004       | 1.687,00 €   |

## Connected Electronic Pipettes

Eppendorf Product Catalog p. 31

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Pipette Manager</b> , an external touch server enabling the communication with connected electronic pipettes. | 1004 000 001 | 1.009,00 € |
| <b>Eppendorf Xplorer® connect</b> , WiFi module to connect pipette to the Pipette Manager                        | 4861 000 970 | 84,20 €    |

## Move It® mechanical multi-channel pipettes, adjustable cone distance

Eppendorf Product Catalog p. 35

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Eppendorf Research® plus Move It®</b> , 4-channel, 9 mm cone distance, variable, incl. epT.I.P.S.® Box 2.0    |              |            |
| 30 – 300 µL, orange  | 3125 000 150 | 1.473,00 € |
| <b>Eppendorf Research® plus Move It®</b> , 4-channel, 9 mm cone distance, variable                               |              |            |
| 120 – 1,200 µL, dark green   | 3125 000 184 | 1.473,00 € |
| <b>Eppendorf Research® plus Move It®</b> , 6-channel, 9 mm cone distance, variable, incl. epT.I.P.S.® Box 2.0    |              |            |
| 30 – 300 µL, orange  | 3125 000 168 | 1.533,00 € |
| <b>Eppendorf Research® plus Move It®</b> , 6-channel, 9 mm cone distance, variable                               |              |            |
| 120 – 1,200 µL, dark green   | 3125 000 192 | 1.533,00 € |
| <b>Eppendorf Research® plus Move It®</b> , 8-channel, 4.5 mm cone distance, variable, incl. epT.I.P.S.® Box 2.0  |              |            |
| 1 – 20 µL, light pink  | 3125 000 117 | 1.582,00 € |
| 5 – 100 µL, light yellow   | 3125 000 133 | 1.582,00 € |
| <b>Eppendorf Research® plus Move It®</b> , 8-channel, 9 mm cone distance, variable, incl. epT.I.P.S.® Box 2.0    |              |            |
| 30 – 300 µL, orange  | 3125 000 176 | 1.596,00 € |
| <b>Eppendorf Research® plus Move It®</b> , 8-channel, 9 mm cone distance, variable                               |              |            |
| 120 – 1,200 µL, dark green   | 3125 000 206 | 1.596,00 € |
| <b>Eppendorf Research® plus Move It®</b> , 12-channel, 4.5 mm cone distance, variable, incl. epT.I.P.S.® Box 2.0 |              |            |
| 1 – 20 µL, light pink  | 3125 000 125 | 1.702,00 € |
| 5 – 100 µL, light yellow   | 3125 000 141 | 1.702,00 € |

## Move It® electronic multi-channel pipettes, adjustable cone distance

Eppendorf Product Catalog p. 35

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Eppendorf Xplorer® plus Move It®</b> , 4-channel, 9 mm cone distance, variable    |              |            |
| 15 – 300 µL, orange  | 4861 000 816 | 1.960,00 € |
| 50 – 1,200 µL, green   | 4861 000 833 | 1.960,00 € |
| <b>Eppendorf Xplorer® plus Move It®</b> , 6-channel, 9 mm cone distance, variable    |              |            |
| 15 – 300 µL, orange  | 4861 000 817 | 2.021,00 € |
| 50 – 1,200 µL, green   | 4861 000 834 | 2.021,00 € |
| <b>Eppendorf Xplorer® plus Move It®</b> , 8-channel, 4.5 mm cone distance, variable  |              |            |
| 1 – 20 µL, light pink  | 4861 000 781 | 2.082,00 € |
| 5 – 100 µL, light yellow   | 4861 000 794 | 2.082,00 € |
| <b>Eppendorf Xplorer® plus Move It®</b> , 8-channel, 9 mm cone distance, variable    |              |            |
| 15 – 300 µL, orange  | 4861 000 818 | 2.082,00 € |
| 50 – 1,200 µL, green   | 4861 000 835 | 2.082,00 € |
| <b>Eppendorf Xplorer® plus Move It®</b> , 12-channel, 4.5 mm cone distance, variable |              |            |
| 1 – 20 µL, light pink  | 4861 000 782 | 2.204,00 € |
| 5 – 100 µL, light yellow   | 4861 000 795 | 2.204,00 € |

## Pipetting

### Eppendorf Research® plus, 1-channel, fixed

Eppendorf Product Catalog p. 36-38

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Eppendorf Research® plus, 1-channel, fixed</b> |              |          |
| 10 µL, medium gray                                | 3124 000 016 | 180,00 € |
| 10 µL, yellow                                     | 3124 000 024 | 180,00 € |
| 20 µL, light gray                                 | 3124 000 032 | 180,00 € |
| 20 µL, yellow                                     | 3124 000 040 | 180,00 € |
| 25 µL, yellow                                     | 3124 000 059 | 180,00 € |
| 50 µL, yellow                                     | 3124 000 067 | 180,00 € |
| 100 µL, yellow                                    | 3124 000 075 | 180,00 € |
| 200 µL, yellow                                    | 3124 000 083 | 180,00 € |
| 200 µL, blue                                      | 3124 000 091 | 180,00 € |
| 250 µL, blue                                      | 3124 000 105 | 180,00 € |
| 500 µL, blue                                      | 3124 000 113 | 180,00 € |
| 1,000 µL, blue                                    | 3124 000 121 | 180,00 € |

### Eppendorf Research® plus, 1-channel, variable

Eppendorf Product Catalog p. 36-38

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Eppendorf Research® plus, 1-channel, variable, incl. epT.I.P.S.® Box 2.0</b>    |              |          |
| 0,1 – 2,5 µL, dark gray  | 3123 000 012 | 288,00 € |
| 0,5 – 10 µL, medium gray   | 3123 000 020 | 288,00 € |
| 2 – 20 µL, light gray  | 3123 000 098 | 288,00 € |
| 2 – 20 µL, yellow  | 3123 000 039 | 288,00 € |
| 10 – 100 µL, yellow  | 3123 000 047 | 288,00 € |
| 20 – 200 µL, yellow  | 3123 000 055 | 288,00 € |
| 30 – 300 µL, orange  | 3123 000 101 | 288,00 € |
| 100 – 1,000 µL, blue   | 3123 000 063 | 288,00 € |
| <b>Eppendorf Research® plus, 1-channel, variable, incl. epT.I.P.S.® sample bag</b> |              |          |
| 0,25 – 2,5 mL, red   | 3123 000 144 | 288,00 € |
| 0,5 – 5 mL, violet   | 3123 000 071 | 288,00 € |
| 1 – 10 mL, turquoise   | 3123 000 080 | 288,00 € |

### Eppendorf Research® plus, 8-channel, variable

Eppendorf Product Catalog p. 36-38

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Eppendorf Research® plus, 8-channel, variable, incl. epT.I.P.S.® Box 2.0</b> |              |          |
| 0,5 – 10 µL, medium gray  | 3125 000 010 | 836,00 € |
| 10 – 100 µL, yellow   | 3125 000 036 | 836,00 € |
| 30 – 300 µL, orange   | 3125 000 052 | 836,00 € |
| <b>Eppendorf Research® plus, 8-channel, variable</b>                            |              |          |
| 120 – 1,200 µL, dark green  | 3125 000 214 | 836,00 € |

### Eppendorf Research® plus, 12-channel, variable

Eppendorf Product Catalog p. 36-38

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Eppendorf Research® plus, 12-channel, variable, incl. epT.I.P.S.® Box 2.0</b> |              |          |
| 0,5 – 10 µL, medium gray   | 3125 000 028 | 915,00 € |
| 10 – 100 µL, yellow  | 3125 000 044 | 915,00 € |
| 30 – 300 µL, orange  | 3125 000 060 | 915,00 € |
| <b>Eppendorf Research® plus, 12-channel, variable</b>                            |              |          |
| 120 – 1,200 µL, dark green   | 3125 000 222 | 915,00 € |

### Eppendorf Research® plus, 16-channel, variable

Eppendorf Product Catalog p. 36-38

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Eppendorf Research® plus</b>   |              |          |
| 16-channel, variable, incl. epT.I.P.S.® Box 2.0, 1 – 20 µL, light pink    | 3125 000 079 | 979,00 € |
| 16-channel, variable, incl. epT.I.P.S.® Box 2.0, 5 – 100 µL, light yellow | 3125 000 095 | 979,00 € |

## Eppendorf Research® plus, 24-channel, variable

| Eppendorf Product Catalog p. 36-38  | Order no.    | Price      |
|---|--------------|------------|
| <b>Eppendorf Research® plus</b>   |              |            |
| 24-channel, variable, incl. epT.I.P.S.® Box 2.0, 1 – 20 µL, light pink    | 3125 000 087 | 1.040,00 € |
| 24-channel, variable, incl. epT.I.P.S.® Box 2.0, 5 – 100 µL, light yellow | 3125 000 109 | 1.040,00 € |

## Eppendorf Research® plus Multi-packs

| Eppendorf Product Catalog p. 36-38   | Order no.    | Price      |
|--|--------------|------------|
| <b>Eppendorf Research® plus, 3-pack</b> , 1-channel, variable, incl. 3x epT.I.P.S.® Box 2.0 and ballpoint pen                            |              |            |
| <b>Option 1:</b> 0.5 – 10 µL, 10 – 100 µL, 100 – 1,000 µL  | 3123 000 900 | 708,00 €   |
| <b>Option 2:</b> 2 – 20 µL yellow, 20 – 200 µL, 100 – 1,000 µL   | 3123 000 918 | 708,00 €   |
| <b>Eppendorf Research® plus, 3-pack</b> , 1-channel, variable, incl. 1x epT.I.P.S.® Box 2.0, 2x epT.I.P.S. sample bags and ballpoint pen |              |            |
| <b>Option 3:</b> 100 – 1,000 µL, 0.5 – 5 mL, 1 – 10 mL   | 3123 000 926 | 708,00 €   |
| <b>Eppendorf Research® plus, 6-pack with Pipette Carousel 2</b> , 1-channel, incl. 6x epT.I.P.S.® Box 2.0 and ballpoint pen              |              |            |
| 0.1 – 2.5 µL, 0.5 – 10 µL, 2 – 20 µL (yellow), 10 – 100 µL, 20 – 200 µL, 100 – 1,000 µL  | 3123 000 942 | 1.502,00 € |

## Accessories for Eppendorf Research® plus

| Eppendorf Product Catalog p. 36-38   | Order no.    | Price   |
|--|--------------|---------|
| <b>Reagent reservoir Tip-Tub</b> , autoclavable reservoir for the aspiration of liquids with multi-channel pipettes, 1 set = 10 reservoirs and 10 lids                       | 0030 058 607 | 27,20 € |
| <b>Grease for pipettes</b> , incl. lint-free applicators, to relubricate the piston or cylinder in pipette lower parts   | 0013 022 153 | 22,30 € |
| <b>Protection filter</b> , 2,5 mL, for Eppendorf Research® plus, Reference® 2 and Xplorer® (plus), 10 filters with 1 filter sleeve, for 2.5 mL pipettes, color code: red     | 4920 622 009 | 15,40 € |
| <b>Protection filter</b> , 5 mL, for Eppendorf Research® plus, Reference® 2 and Xplorer® (plus), 10 filters with 1 filter sleeve, for 5 mL pipettes, color code: violet      | 4920 623 005 | 15,40 € |
| <b>Protection filter</b> , 10 mL, for Eppendorf Research® plus, Reference® 2 and Xplorer® (plus), 10 filters with 1 filter sleeve, for 10 mL pipettes, color code: turquoise | 4920 624 001 | 15,40 € |
| <b>Safety plug tools</b> , for Eppendorf Research® plus, to change the factory adjustment, 5 pcs.  | 3120 639 004 | 7,80 €  |
| <b>Adjustment tool</b> , for Eppendorf Research® plus and Reference® 2, to change factory or temporary adjustment  | 3120 633 006 | 7,00 €  |
| <b>Adjustment seal ADJ</b> , for Eppendorf Research® plus, red, 5 pcs.   | 3120 636 005 | 2,10 €  |
| <b>Pipette wrench 2 mL – 10 mL</b> , to open multi-channel pipette lower parts up to 10 mL   | 3120 634 002 | 7,00 €  |
| <b>Multi-channel tool</b> , for cutting the O-rings on tip cones of multi-channel pipettes, for 5 – 100 µL or 15 – 300 µL pipettes   | 3122 610 003 | 7,00 €  |
| <b>Multi-channel tool</b> , for cutting the O-rings on tip cones of multi-channel pipettes, for 50 – 1,200 µL pipettes   | 4860 736 000 | 18,40 € |
| <b>Unlocking tool</b> , to open multi-channel lower parts of 1,200 µL pipettes   | 4861 605 006 | 10,20 € |
| <b>Locking ring for single-channel pipettes</b> , to prevent spring action in single-channel pipettes, for pipettes from 2.5 µL – 1,000 µL                                   | 3120 635 009 | 3,50 €  |
| <b>Locking clip for 8-channel pipettes</b> , to prevent spring action in multi-channel pipettes, 2 pieces  | 3122 612 006 | 6,50 €  |
| <b>Locking clip for 12-channel pipettes</b> , to prevent spring action in multi-channel pipettes, 3 pieces   | 3122 613 002 | 8,00 €  |
| <b>O-ring</b> , red, for multi-channel pipette lower parts, 5 – 100 µL, 15 – 300 µL (1 set = 24 pcs)   | 3122 611 000 | 18,40 € |
| <b>O-ring</b> , red, for multi-channel pipette lower parts, 50 – 1,200 µL (1 set = 24 pcs)   | 4860 716 009 | 21,30 € |
| <b>Lower part made of PTFE (higher chemical resistance)</b>  |              |         |
| for Eppendorf Research® plus, single-channel, 5 mL, color code: violet   | 3120 642 005 | 93,90 € |
| for Eppendorf Research® plus, single-channel, 1,000 µL, color code: blue   | 3120 641 009 | 44,40 € |
| for Eppendorf Research® plus, single-channel, 10 mL, color code: turquoise   | 3120 646 000 | 93,90 € |

## Pipetting

### Eppendorf Reference® 2, 1-channel, fixed

Eppendorf Product Catalog p. 40-41

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Eppendorf Reference® 2, 1-channel, fixed</b> |              |          |
| 1 µL, dark gray                                 | 4925 000 014 | 214,00 € |
| 2 µL, dark gray                                 | 4925 000 022 | 214,00 € |
| 5 µL, medium gray                               | 4925 000 030 | 214,00 € |
| 10 µL, medium gray                              | 4925 000 049 | 214,00 € |
| 10 µL, yellow                                   | 4925 000 057 | 214,00 € |
| 20 µL, light gray                               | 4925 000 065 | 214,00 € |
| 20 µL, yellow                                   | 4925 000 073 | 214,00 € |
| 25 µL, yellow                                   | 4925 000 081 | 214,00 € |
| 50 µL, yellow                                   | 4925 000 090 | 214,00 € |
| 100 µL, yellow                                  | 4925 000 103 | 214,00 € |
| 200 µL, yellow                                  | 4925 000 111 | 214,00 € |
| 200 µL, blue                                    | 4925 000 120 | 214,00 € |
| 250 µL, blue                                    | 4925 000 138 | 214,00 € |
| 500 µL, blue                                    | 4925 000 146 | 214,00 € |
| 1,000 µL, blue                                  | 4925 000 154 | 214,00 € |
| 2 mL, red                                       | 4925 000 162 | 214,00 € |
| 2,5 mL, red                                     | 4925 000 170 | 214,00 € |

### Eppendorf Reference® 2, 1-channel, variable

Eppendorf Product Catalog p. 40-41

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Eppendorf Reference® 2, 1-channel, variable, incl. epT.I.P.S.® Box 2.0</b>    |              |          |
| 0,1 – 2,5 µL, dark gray  | 4924 000 010 | 376,00 € |
| 0,5 – 10 µL, medium gray   | 4924 000 029 | 376,00 € |
| 2 – 20 µL, light gray  | 4924 000 037 | 376,00 € |
| 2 – 20 µL, yellow  | 4924 000 045 | 376,00 € |
| 10 – 100 µL, yellow  | 4924 000 053 | 376,00 € |
| 20 – 200 µL, yellow  | 4924 000 061 | 376,00 € |
| 30 – 300 µL, orange  | 4924 000 070 | 376,00 € |
| 100 – 1,000 µL, blue   | 4924 000 088 | 376,00 € |
| <b>Eppendorf Reference® 2, 1-channel, variable, incl. epT.I.P.S.® sample bag</b> |              |          |
| 0,25 – 2,5 mL, red   | 4924 000 096 | 376,00 € |
| 0,5 – 5 mL, violet   | 4924 000 100 | 376,00 € |
| 1 – 10 mL, turquoise   | 4924 000 118 | 376,00 € |

### Eppendorf Reference® 2, multi-channel, variable

Eppendorf Product Catalog p. 40-41

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Eppendorf Reference® 2, 8-channel, variable, incl. epT.I.P.S.® Box 2.0</b>  |              |            |
| 0,5 – 10 µL, medium gray   | 4926 000 018 | 1.037,00 € |
| 10 – 100 µL, yellow  | 4926 000 034 | 1.037,00 € |
| 30 – 300 µL, orange  | 4926 000 050 | 1.037,00 € |
| <b>Eppendorf Reference® 2, 12-channel, variable, incl. epT.I.P.S.® Box 2.0</b> |              |            |
| 0,5 – 10 µL, medium gray   | 4926 000 026 | 1.135,00 € |
| 10 – 100 µL, yellow  | 4926 000 042 | 1.135,00 € |
| 30 – 300 µL, orange  | 4926 000 069 | 1.135,00 € |

## Eppendorf Reference® 2 Multi-packs

Eppendorf Product Catalog p. 40-41

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Eppendorf Reference® 2, 3-pack</b> , 1-channel, variable, incl. 3x epT.I.P.S.® Box 2.0 and ballpoint pen                            |              |            |
| <b>Option 1:</b> 0.5 – 10 µL, 10 – 100 µL, 100 – 1,000 µL  | 4924 000 908 | 882,00 €   |
| <b>Option 2:</b> 2 – 20 µL yellow, 20 – 200 µL, 100 – 1,000 µL   | 4924 000 916 | 882,00 €   |
| <b>Eppendorf Reference® 2, 3-pack</b> , 1-channel, variable, incl. 1x epT.I.P.S.® Box 2.0, 2x epT.I.P.S. sample bags and ballpoint pen |              |            |
| <b>Option 3:</b> 100 – 1,000 µL, 0.5 – 5 mL, 1 – 10 mL   | 4924 000 924 | 882,00 €   |
| <b>Eppendorf Reference® 2, 6-pack with Pipette Carousel 2</b> , 1-channel, incl. 6x epT.I.P.S.® Box 2.0 and ballpoint pen              |              |            |
| 0.1 – 2.5 µL, 0.5 – 10 µL, 2 – 20 µL (yellow), 10 – 100 µL, 20 – 200 µL, 100 – 1,000 µL  | 4924 000 940 | 1.926,00 € |

## Accessories for Eppendorf Reference® 2

Eppendorf Product Catalog p. 40-41

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Reagent reservoir Tip-Tub</b> , autoclavable reservoir for the aspiration of liquids with multi-channel pipettes, 1 set = 10 reservoirs and 10 lids                       | 0030 058 607 | 27,20 €  |
| <b>Grease for pipettes</b> , incl. lint-free applicators, to relubricate the piston or cylinder in pipette lower parts   | 0013 022 153 | 22,30 €  |
| <b>Protection filter</b> , 2,5 mL, for Eppendorf Research® plus, Reference® 2 and Xplorer® (plus), 10 filters with 1 filter sleeve, for 2.5 mL pipettes, color code: red     | 4920 622 009 | 15,40 €  |
| <b>Protection filter</b> , 5 mL, for Eppendorf Research® plus, Reference® 2 and Xplorer® (plus), 10 filters with 1 filter sleeve, for 5 mL pipettes, color code: violet      | 4920 623 005 | 15,40 €  |
| <b>Protection filter</b> , 10 mL, for Eppendorf Research® plus, Reference® 2 and Xplorer® (plus), 10 filters with 1 filter sleeve, for 10 mL pipettes, color code: turquoise | 4920 624 001 | 15,40 €  |
| <b>Adjustment seal ADJ, red</b> , for Eppendorf Reference® 2, incl. pin for removing the inserted seal, 1 set = 5 pieces   | 4920 626 004 | 9,50 €   |
| <b>Adjustment tool</b> , for Eppendorf Research® plus and Reference® 2, to change factory or temporary adjustment  | 3120 633 006 | 7,00 €   |
| <b>Pipette wrench 2 mL – 10 mL</b> , to open multi-channel pipette lower parts up to 10 mL   | 3120 634 002 | 7,00 €   |
| <b>Multi-channel tool</b> , for cutting the O-rings on tip cones of multi-channel pipettes, for 5 – 100 µL or 15 – 300 µL pipettes   | 3122 610 003 | 7,00 €   |
| <b>Lower part made of PTFE (higher chemical resistance)</b>  |              |          |
| for Eppendorf Reference® 2, single-channel, 1,000 µL, color code: blue   | 4920 627 000 | 44,40 €  |
| for Eppendorf Reference® 2, fixed, single-channel, 2 mL, color code: red   | 4921 610 004 | 108,00 € |
| for Eppendorf Reference® 2, single-channel, 2.5 mL, color code: red  | 4920 629 003 | 108,00 € |
| for Eppendorf Reference® 2, single-channel, 5 mL, color code: violet   | 4920 633 000 | 108,00 € |
| for Eppendorf Reference® 2, single-channel, 10 mL, color code: turquoise   | 4920 636 000 | 122,00 € |
| <b>Locking ring for single-channel pipettes</b> , to prevent spring action in single-channel pipettes, for pipettes from 2.5 µL – 1,000 µL                                   | 3120 635 009 | 3,50 €   |
| <b>O-ring</b> , red, for multi-channel pipette lower parts, 5 – 100 µL, 15 – 300 µL (1 set = 24 pcs)   | 3122 611 000 | 18,40 €  |

## Eppendorf Xplorer®, 1-channel, variable

Eppendorf Product Catalog p. 42-44

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Eppendorf Xplorer®</b> , 1-channel, variable |              |          |
| 0,5 – 10 µL, medium gray                        | 4861 000 015 | 664,00 € |
| 1 – 20 µL, light gray                           | 4861 000 017 | 664,00 € |
| 5 – 100 µL, yellow                              | 4861 000 023 | 664,00 € |
| 10 – 200 µL, yellow                             | 4861 000 027 | 664,00 € |
| 15 – 300 µL, orange                             | 4861 000 031 | 664,00 € |
| 50 – 1,000 µL, blue                             | 4861 000 040 | 664,00 € |
| 0,1 – 2,5 mL, red                               | 4861 000 044 | 664,00 € |
| 0,2 – 5 mL, violet                              | 4861 000 058 | 664,00 € |
| 0,5 – 10 mL, turquoise                          | 4861 000 066 | 664,00 € |

## Eppendorf Xplorer®, 8-channel, variable

Eppendorf Product Catalog p. 42-44

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Eppendorf Xplorer®</b> , 8-channel, variable |              |          |
| 0,5 – 10 µL, medium gray                        | 4861 000 104 | 980,00 € |
| 5 – 100 µL, yellow                              | 4861 000 120 | 980,00 € |
| 15 – 300 µL, orange                             | 4861 000 147 | 980,00 € |
| 50 – 1,200 µL, green                            | 4861 000 163 | 980,00 € |

## Pipetting

### Eppendorf Xplorer®, 12-channel, variable

Eppendorf Product Catalog p. 42-44

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Eppendorf Xplorer®</b> , 12-channel, variable |              |            |
| 0,5 – 10 µL, medium gray                         | 4861 000 112 | 1.149,00 € |
| 5 – 100 µL, yellow                               | 4861 000 139 | 1.149,00 € |
| 15 – 300 µL, orange                              | 4861 000 155 | 1.149,00 € |
| 50 – 1,200 µL, green                             | 4861 000 171 | 1.149,00 € |

### Eppendorf Xplorer® plus, 1-channel, variable

Eppendorf Product Catalog p. 42-44

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Eppendorf Xplorer® plus</b> , 1-channel, variable |              |          |
| 0,5 – 10 µL, medium gray                             | 4861 000 708 | 803,00 € |
| 1 – 20 µL, light gray                                | 4861 000 710 | 803,00 € |
| 5 – 100 µL, yellow                                   | 4861 000 716 | 803,00 € |
| 10 – 200 µL, yellow                                  | 4861 000 720 | 803,00 € |
| 15 – 300 µL, orange                                  | 4861 000 724 | 803,00 € |
| 50 – 1,000 µL, blue                                  | 4861 000 732 | 803,00 € |
| 0,1 – 2,5 mL, red                                    | 4861 000 736 | 803,00 € |
| 0,2 – 5 mL, violet                                   | 4861 000 740 | 803,00 € |
| 0,5 – 10 mL, turquoise                               | 4861 000 759 | 803,00 € |

### Eppendorf Xplorer® plus, 8-channel, variable

Eppendorf Product Catalog p. 42-44

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Eppendorf Xplorer® plus</b> , 8-channel, variable |              |            |
| 0,5 – 10 µL, medium gray                             | 4861 000 767 | 1.182,00 € |
| 5 – 100 µL, yellow                                   | 4861 000 783 | 1.182,00 € |
| 15 – 300 µL, orange                                  | 4861 000 805 | 1.182,00 € |
| 50 – 1,200 µL, green                                 | 4861 000 821 | 1.182,00 € |

### Eppendorf Xplorer® plus, 12-channel, variable

Eppendorf Product Catalog p. 42-44

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Eppendorf Xplorer® plus</b> , 12-channel, variable |              |            |
| 0,5 – 10 µL, medium gray                              | 4861 000 775 | 1.388,00 € |
| 5 – 100 µL, yellow                                    | 4861 000 791 | 1.388,00 € |
| 15 – 300 µL, orange                                   | 4861 000 813 | 1.388,00 € |
| 50 – 1,200 µL, green                                  | 4861 000 830 | 1.388,00 € |

### Eppendorf Xplorer® plus, 16-channel, variable, for 384-well plates

Eppendorf Product Catalog p. 42-44

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Eppendorf Xplorer® plus</b> , 16-channel, variable |              |            |
| 1 – 20 µL, light pink                                 | 4861 000 778 | 1.550,00 € |
| 5 – 100 µL, light yellow                              | 4861 000 792 | 1.550,00 € |



## Eppendorf Xplorer® plus, 24-channel, variable, for 384-well plates

Eppendorf Product Catalog p. 42-44

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Eppendorf Xplorer® plus, 24-channel, variable</b> |              |            |
| 1 – 20 µL, light pink                                | 4861 000 779 | 1.714,00 € |
| 5 – 100 µL, light yellow                             | 4861 000 793 | 1.714,00 € |

## Accessories for Eppendorf Xplorer®/Eppendorf Xplorer® plus

Eppendorf Product Catalog p. 42-44

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Reagent reservoir Tip-Tub</b> , autoclavable reservoir for the aspiration of liquids with multi-channel pipettes, 1 set = 10 reservoirs and 10 lids | 0030 058 607 | 27,20 € |
| <b>Lithium polymer rechargeable battery</b> , for Eppendorf Xplorer® pipette family  | 4861 601 000 | 20,10 € |
| <b>Lithium polymer rechargeable battery (small)</b> , for Eppendorf Xplorer® connect   | 4861 618 000 | 13,10 € |
| <b>Grease for pipettes</b> , incl. lint-free applicators, to relubricate the piston or cylinder in pipette lower parts                                 | 0013 022 153 | 22,30 € |
| <b>Power supply with power plug adapters</b> , for Eppendorf pipettes, multi-dispenser and charger stand, 100 – 240 VAC / 5 VDC                        | 4986 603 005 | 21,30 € |
| <b>Pipette wrench 2 mL – 10 mL</b> , to open multi-channel pipette lower parts up to 10 mL   | 3120 634 002 | 7,00 €  |
| <b>Multi-channel tool</b> , for cutting the O-rings on tip cones of multi-channel pipettes, for 5 – 100 µL or 15 – 300 µL pipettes                     | 3122 610 003 | 7,00 €  |
| <b>Multi-channel tool</b> , for cutting the O-rings on tip cones of multi-channel pipettes, for 50 – 1,200 µL pipettes                                 | 4860 736 000 | 18,40 € |
| <b>Unlocking tool</b> , to open multi-channel lower parts of 1,200 µL pipettes   | 4861 605 006 | 10,20 € |
| <b>Locking ring for single-channel pipettes</b> , to prevent spring action in single-channel pipettes, for pipettes from 2.5 µL – 1,000 µL             | 3120 635 009 | 3,50 €  |
| <b>Locking clip for 8-channel pipettes</b> , to prevent spring action in multi-channel pipettes, 2 pieces  | 3122 612 006 | 6,50 €  |
| <b>Locking clip for 12-channel pipettes</b> , to prevent spring action in multi-channel pipettes, 3 pieces   | 3122 613 002 | 8,00 €  |
| <b>O-ring, red</b> , for multi-channel pipette lower parts, 5 – 100 µL, 15 – 300 µL (1 set = 24 pcs)   | 3122 611 000 | 18,40 € |
| <b>O-ring, red</b> , for multi-channel pipette lower parts, 50 – 1,200 µL (1 set = 24 pcs)   | 4860 716 009 | 21,30 € |

## Pipette carousels, charger stands and wall mounting devices

Eppendorf Product Catalog p. 46-47

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Pipette Carousel 2</b> , with holders for 6 Eppendorf Research®, Eppendorf Research® plus, Eppendorf Reference®, Eppendorf Reference® 2 or Biomaster®, additional pipette holders are optionally available                     | 3116 000 015 | 156,00 € |
| <b>Charger Carousel 2</b> , with charger shells for 6 Eppendorf Xplorer®/Xplorer® plus, mains/power adapter included, additional charger shells and pipette holders are optionally available                                      | 3116 000 023 | 493,00 € |
| <b>Charger Stand 2</b> , for one Eppendorf Xplorer®/Xplorer® plus   | 3116 000 031 | 261,00 € |
| <b>Charger Stand 2</b> , for one Multipipette® E3/E3x or Multipipette® stream/Xstream, operated with mains/power adapter supplied with Multipipette® E3/E3x or Multipipette® stream/Xstream                                       | 3116 000 040 | 261,00 € |
| <b>Pipette Stand 2</b> , for one Multipipette® M4, without charging functionality, additional pipette holders are optionally available  | 3116 000 058 | 121,00 € |
| <b>Pipette Holder 2</b> , for one Multipipette® E3/E3x or Multipipette® stream/Xstream, for Pipette Carousel 2 or wall mounting, sticky tape included, without charging functionality   | 3116 000 139 | 41,70 €  |
| <b>Pipette Holder 2</b> , for one Eppendorf Research®, Eppendorf Research® plus, Eppendorf Reference®, Eppendorf Reference® 2 or Biomaster®, for Pipette Carousel 2 and Charger Carousel 2 or wall mounting, sticky tape included | 3116 000 112 | 32,20 €  |
| <b>Pipette Holder 2</b> , for one Eppendorf Xplorer®/Xplorer® plus, for Pipette Carousel 2 or wall mounting, sticky tape included, without charging functionality   | 3116 000 120 | 41,70 €  |
| <b>Pipette Holder 2</b> , for one Multipipette® M4, for Pipette Carousel 2 and Charger Carousel 2 or wall mounting, sticky tape included  | 3116 000 147 | 41,70 €  |
| <b>Charger Shell 2</b> , for one Eppendorf Xplorer®/Xplorer® plus, for Charger Carousel 2 (with charging functionality)   | 3116 602 007 | 54,70 €  |
| <b>Charger Shell 2</b> , for one Multipipette® E3/E3x or Multipipette® stream/Xstream, for Charger Carousel 2 (with charging functionality)   | 3116 603 003 | 54,70 €  |

## Tips

### epT.I.P.S.<sup>®</sup> Standard

Eppendorf Product Catalog p. 56-63

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>epT.I.P.S.<sup>®</sup> Standard</b>   |              |         |
| Eppendorf Quality, 0,1 – 10 µL, 34 mm, dark gray, colorless tips, 1.000 tips (2 bags × 500 tips)       | 0030 000 811 | 67,20 € |
| Eppendorf Quality, 0,1 – 20 µL, 40 mm, medium gray, colorless tips, 1.000 tips (2 bags × 500 tips)     | 0030 000 838 | 67,20 € |
| Eppendorf Quality, 0,5 – 20 µL L, 46 mm, light gray, colorless tips, 1.000 tips (2 bags × 500 tips)    | 0030 000 854 | 63,90 € |
| Eppendorf Quality, 2 – 200 µL, 53 mm, yellow, colorless tips, 1.000 tips (2 bags × 500 tips)           | 0030 000 889 | 41,40 € |
| Eppendorf Quality, 2 – 200 µL, 53 mm, yellow, yellow tips, 1.000 tips (2 bags × 500 tips)              | 0030 000 870 | 41,40 € |
| Eppendorf Quality, 20 – 300 µL, 55 mm, orange, colorless tips, 1.000 tips (2 bags × 500 tips)          | 0030 000 900 | 41,40 € |
| Eppendorf Quality, 20 – 300 µL, 55 mm, orange, yellow tips, 1.000 tips (2 bags × 500 tips)             | 0030 000 897 | 41,40 € |
| Eppendorf Quality, 50 – 1,000 µL, 71 mm, blue, colorless tips, 1.000 tips (2 bags × 500 tips)          | 0030 000 927 | 41,40 € |
| Eppendorf Quality, 50 – 1,000 µL, 71 mm, blue, blue tips, 1.000 tips (2 bags × 500 tips)               | 0030 000 919 | 41,40 € |
| Eppendorf Quality, 50 – 1,250 µL, 76 mm, green, colorless tips, 1.000 tips (4 bags × 250 tips)         | 0030 000 935 | 70,50 € |
| Eppendorf Quality, 50 – 1,250 µL L, 103 mm, dark green, colorless tips, 1.000 tips (4 bags × 250 tips) | 0030 000 730 | 71,40 € |
| Eppendorf Quality, 0,25 – 2,5 mL, 115 mm, red, colorless tips, 500 tips (5 bags × 100 tips)            | 0030 000 951 | 48,50 € |
| Eppendorf Quality, 0,1 – 5 mL, 120 mm, violet, colorless tips, 500 tips (5 bags × 100 tips)            | 0030 000 978 | 56,60 € |
| Eppendorf Quality, 0,2 – 5 mL L, 175 mm, violet, colorless tips, 300 tips (3 bags × 100 tips)          | 0030 000 650 | 38,60 € |
| Eppendorf Quality, 0,5 – 10 mL, 165 mm, turquoise, colorless tips, 200 tips (2 bags × 100 tips)        | 0030 000 765 | 58,20 € |
| Eppendorf Quality, 0,5 – 10 mL L, 243 mm, turquoise, colorless tips, 200 tips (2 bags × 100 tips)      | 0030 000 781 | 82,90 € |

### epT.I.P.S.<sup>®</sup> Reloads Eppendorf Quality

Eppendorf Product Catalog p. 56-63

|   | Order no.    | Price   |
|---|--------------|---------|
| <b>epT.I.P.S.<sup>®</sup> Reloads</b>   |              |         |
| Eppendorf Quality, 0,1 – 10 µL, 34 mm, dark gray, colorless tips, 960 tips (10 trays × 96 tips)       | 0030 073 371 | 80,40 € |
| Eppendorf Quality, 0,1 – 20 µL, 40 mm, medium gray, colorless tips, 960 tips (10 trays × 96 tips)     | 0030 073 398 | 80,40 € |
| Eppendorf Quality, 0,5 – 20 µL L, 46 mm, light gray, colorless tips, 960 tips (10 trays × 96 tips)    | 0030 073 410 | 72,60 € |
| Eppendorf Quality, 2 – 200 µL, 53 mm, yellow, colorless tips, 960 tips (10 trays × 96 tips)           | 0030 073 436 | 68,80 € |
| Eppendorf Quality, 20 – 300 µL, 55 mm, orange, colorless tips, 960 tips (10 trays × 96 tips)          | 0030 073 452 | 68,80 € |
| Eppendorf Quality, 50 – 1,000 µL, 71 mm, blue, colorless tips, 960 tips (10 trays × 96 tips)          | 0030 073 479 | 68,80 € |
| Eppendorf Quality, 50 – 1,250 µL, 76 mm, green, colorless tips, 960 tips (10 trays × 96 tips)         | 0030 073 495 | 83,30 € |
| Eppendorf Quality, 50 – 1,250 µL L, 103 mm, dark green, colorless tips, 960 tips (10 trays × 96 tips) | 0030 073 630 | 84,60 € |
| Eppendorf Quality, 0,25 – 2,5 mL, red, colorless tips, 480 tips (10 trays × 48 tips)                  | 0030 073 517 | 57,30 € |

### epT.I.P.S.<sup>®</sup> Reloads PCR clean

Eppendorf Product Catalog p. 56-63

|   | Order no.    | Price   |
|---|--------------|---------|
| <b>epT.I.P.S.<sup>®</sup> Reloads</b>   |              |         |
| PCR clean, 0,1 – 10 µL, 34 mm, dark gray, colorless tips, 960 tips (10 trays × 96 tips)       | 0030 073 754 | 84,80 € |
| PCR clean, 0,1 – 20 µL, 40 mm, medium gray, colorless tips, 960 tips (10 trays × 96 tips)     | 0030 073 770 | 87,00 € |
| PCR clean, 0,5 – 20 µL L, 46 mm, light gray, colorless tips, 960 tips (10 trays × 96 tips)    | 0030 073 797 | 80,70 € |
| PCR clean, 2 – 200 µL, 53 mm, yellow, colorless tips, 960 tips (10 trays × 96 tips)           | 0030 073 819 | 76,20 € |
| PCR clean, 20 – 300 µL, 55 mm, orange, colorless tips, 960 tips (10 trays × 96 tips)          | 0030 073 835 | 76,20 € |
| PCR clean, 50 – 1,000 µL, 71 mm, blue, colorless tips, 960 tips (10 trays × 96 tips)          | 0030 073 851 | 76,20 € |
| PCR clean, 50 – 1,250 µL, 76 mm, green, colorless tips, 960 tips (10 trays × 96 tips)         | 0030 073 878 | 93,10 € |
| PCR clean, 50 – 1,250 µL L, 103 mm, dark green, colorless tips, 960 tips (10 trays × 96 tips) | 0030 073 649 | 94,40 € |
| PCR clean, 0,25 – 2,5 mL, 115 mm, red, colorless tips, 480 tips (10 trays × 48 tips)          | 0030 073 894 | 61,90 € |

## epT.I.P.S.<sup>®</sup> Box, Eppendorf Quality™

Eppendorf Product Catalog p. 56-63

Order no.

Price

### epT.I.P.S.<sup>®</sup> Box 2.0

|  |              |         |
|--|--------------|---------|
| Eppendorf Quality, 0,1 – 10 µL, dark gray, colorless tips, 1 reusable box × 96 tips. Replaces order no. 0030073010.              | 0030 076 125 | 40,10 € |
| Eppendorf Quality, 0,1 – 20 µL, medium gray, colorless tips, 1 reusable box × 96 tips. Replaces order no. 0030073037.            | 0030 076 133 | 40,10 € |
| Eppendorf Quality, 0,5 – 20 µL L, 46 mm, light gray, colorless tips, 1 reusable box × 96 tips. Replaces order no. 0030073053.    | 0030 076 141 | 35,50 € |
| Eppendorf Quality, 2 – 200 µL, yellow, colorless tips, 1 reusable box × 96 tips. Replaces order no. 0030073070.                  | 0030 076 150 | 35,50 € |
| Eppendorf Quality, 20 – 300 µL, orange, colorless tips, 1 reusable box × 96 tips. Replaces order no. 0030073096.                 | 0030 076 168 | 35,50 € |
| Eppendorf Quality, 50 – 1,000 µL, blue, colorless tips, 1 reusable box × 96 tips. Replaces order no. 0030073118.                 | 0030 076 176 | 35,50 € |
| Eppendorf Quality, 50 – 1,250 µL, green, colorless tips, 1 reusable box × 96 tips. Replaces order no. 0030073134.                | 0030 076 184 | 40,10 € |
| Eppendorf Quality, 50 – 1,250 µL L, 103 mm, dark green, colorless tips, 1 reusable box × 96 tips. Replaces order no. 0030073657. | 0030 076 192 | 40,70 € |
| Eppendorf Quality, 0,25 – 2,5 mL, red, colorless tips, 1 reusable box × 48 tips. Replaces order no. 0030073150.                  | 0030 076 206 | 47,20 € |
| Eppendorf Quality, 0,1 – 5 mL, 120 mm, violet, colorless tips, 1 reusable box × 24 tips. Replaces order no. 0030073177.          | 0030 076 214 | 42,10 € |

## epT.I.P.S.<sup>®</sup> Set

Eppendorf Product Catalog p. 56-63

Order no.

Price

### epT.I.P.S.<sup>®</sup> Set

|   |              |         |
|---|--------------|---------|
| Eppendorf Quality, 0,1 – 10 µL, dark gray, colorless tips, 480 tips (5 trays × 96 tips), 1 reusable box. Replaces order no. 0030073215.           | 0030 076 290 | 65,50 € |
| Eppendorf Quality, 0,1 – 20 µL, medium gray, colorless tips, 480 tips (5 trays × 96 tips), 1 reusable box. Replaces order no. 0030073231.         | 0030 076 303 | 65,50 € |
| Eppendorf Quality, 0,5 – 20 µL L, 46 mm, light gray, colorless tips, 480 tips (5 trays × 96 tips), 1 reusable box. Replaces order no. 0030073258. | 0030 076 311 | 60,20 € |
| Eppendorf Quality, 2 – 200 µL, yellow, colorless tips, 480 tips (5 trays × 96 tips), 1 reusable box. Replaces order no. 0030073274.               | 0030 076 320 | 58,00 € |
| Eppendorf Quality, 20 – 300 µL, orange, colorless tips, 480 tips (5 trays × 96 tips), 1 reusable box. Replaces order no. 0030073290.              | 0030 076 338 | 58,00 € |
| Eppendorf Quality, 50 – 1,000 µL, blue, colorless tips, 480 tips (5 trays × 96 tips), 1 reusable box. Replaces order no. 0030073312.              | 0030 076 346 | 61,70 € |
| Eppendorf Quality, 50 – 1,250 µL, green, colorless tips, 480 tips (5 trays × 96 tips), 1 reusable box. Replaces order no. 0030073339.             | 0030 076 354 | 80,60 € |
| Eppendorf Quality, 0,25 – 2,5 mL, red, colorless tips, 240 tips (5 trays × 48 tips), 1 reusable box. Replaces order no. 0030073355.               | 0030 076 362 | 56,60 € |

## Tips

### epT.I.P.S.<sup>®</sup> Singles

Eppendorf Product Catalog p. 56-63

Order no.

Price

#### epT.I.P.S.<sup>®</sup> Singles

|  |              |         |
|--|--------------|---------|
| Biopur <sup>®</sup> , 0,1 – 20 µL, 40 mm, medium gray, colorless tips, 100 tips, individually packed | 0030 010 027 | 32,50 € |
| Biopur <sup>®</sup> , 2 – 200 µL, 53 mm, yellow, colorless tips, 100 tips, individually packed       | 0030 010 043 | 32,50 € |
| Biopur <sup>®</sup> , 50 – 1,000 µL, 71 mm, blue, colorless tips, 100 tips, individually packed      | 0030 010 060 | 32,50 € |

### epT.I.P.S.<sup>®</sup> Racks

Eppendorf Product Catalog p. 56-63

Order no.

Price

#### epT.I.P.S.<sup>®</sup> Racks

|   |              |         |
|---|--------------|---------|
| Biopur <sup>®</sup> , 0,1 – 20 µL, 40 mm, medium gray, colorless tips, 480 tips (5 racks × 96 tips)     | 0030 075 226 | 62,60 € |
| Biopur <sup>®</sup> , 2 – 200 µL, 53 mm, yellow, colorless tips, 480 tips (5 racks × 96 tips)           | 0030 075 234 | 58,30 € |
| Biopur <sup>®</sup> , 20 – 300 µL, 55 mm, orange, colorless tips, 480 tips (5 racks × 96 tips)          | 0030 075 242 | 58,30 € |
| Biopur <sup>®</sup> , 50 – 1,000 µL, 71 mm, blue, colorless tips, 480 tips (5 racks × 96 tips)          | 0030 075 250 | 58,30 € |
| Biopur <sup>®</sup> , 50 – 1,250 µL, 76 mm, green, colorless tips, 480 tips (5 racks × 96 tips)         | 0030 075 269 | 58,30 € |
| Biopur <sup>®</sup> , 50 – 1,250 µL L, 103 mm, dark green, colorless tips, 480 tips (5 racks × 96 tips) | 0030 075 277 | 58,30 € |
| Biopur <sup>®</sup> , 0,25 – 2,5 mL, 115 mm, red, colorless tips, 240 tips (5 racks × 48 tips)          | 0030 075 285 | 52,00 € |
| Biopur <sup>®</sup> , 0,1 – 5 mL, 120 mm, violet, colorless tips, 120 tips (5 racks × 24 tips)          | 0030 075 293 | 58,00 € |
| Biopur <sup>®</sup> , 0,2 – 5 mL L, 175 mm, violet, colorless tips, 120 tips (5 racks × 24 tips)        | 0030 075 307 | 67,30 € |
| Biopur <sup>®</sup> , 0,5 – 10 mL, 165 mm, turquoise, colorless tips, 120 tips (5 racks × 24 tips)      | 0030 075 315 | 69,40 € |

### epT.I.P.S.<sup>®</sup> Racks

Eppendorf Product Catalog p. 56-63

Order no.

Price

#### epT.I.P.S.<sup>®</sup> Racks

|   |              |         |
|---|--------------|---------|
| Eppendorf Quality, 0,1 – 5 mL, 120 mm, violet, 120 tips (5 racks × 24 tips)     | 0030 071 638 | 43,30 € |
| PCR clean, 0,1 – 5 mL, 120 mm, violet, 120 tips (5 racks × 24 tips)             | 0030 071 662 | 51,20 € |
| Eppendorf Quality, 0,2 – 5 mL L, 175 mm, violet, 120 tips (5 racks × 24 tips)   | 0030 071 646 | 52,90 € |
| Eppendorf Quality, 0,5 – 10 mL, 165 mm, turquoise, 120 tips (5 racks × 24 tips) | 0030 071 654 | 64,80 € |

### epT.I.P.S.<sup>®</sup> Reloads

New

Eppendorf Product Catalog p. 56-63

Order no.

Price

#### epT.I.P.S.<sup>®</sup> Reloads

|   |              |         |
|---|--------------|---------|
| Biopur <sup>®</sup> , 0,1 – 20 µL, 40 mm, medium gray, colorless tips, 480 tips (5 trays × 96 tips)     | 0030 075 420 | Inquire |
| Biopur <sup>®</sup> , 2 – 200 µL, 53 mm, yellow, colorless tips, 480 tips (5 trays × 96 tips)           | 0030 075 439 | Inquire |
| Biopur <sup>®</sup> , 20 – 300 µL, 55 mm, orange, colorless tips, 480 tips (5 trays × 96 tips)          | 0030 075 447 | Inquire |
| Biopur <sup>®</sup> , 50 – 1,000 µL, 71 mm, blue, colorless tips, 480 tips (5 trays × 96 tips)          | 0030 075 455 | Inquire |
| Biopur <sup>®</sup> , 50 – 1,250 µL, 76 mm, green, colorless tips, 480 tips (5 trays × 96 tips)         | 0030 075 463 | Inquire |
| Biopur <sup>®</sup> , 50 – 1,250 µL L, 103 mm, dark green, colorless tips, 480 tips (5 trays × 96 tips) | 0030 075 471 | Inquire |

## ep Dualfilter T.I.P.S.<sup>®</sup> PCR clean/sterile (sterile and pyrogen-free) and Forensic DNA Grade

Eppendorf Product Catalog p. 54-61

Order no.

Price

### ep Dualfilter T.I.P.S.<sup>®</sup>

|  |              |          |
|--|--------------|----------|
| PCR clean and sterile, 0,1 – 10 µL S, 34 mm, dark gray, colorless tips, 960 tips (10 racks × 96 tips)    | 0030 078 500 | 188,00 € |
| PCR clean and sterile, 0,1 – 10 µL M, 40 mm, medium gray, colorless tips, 960 tips (10 racks × 96 tips)  | 0030 078 519 | 188,00 € |
| PCR clean and sterile, 0,5 – 20 µL L, 46 mm, light gray, colorless tips, 960 tips (10 racks × 96 tips)   | 0030 078 527 | 181,00 € |
| PCR clean and sterile, 2 – 20 µL, 53 mm, yellow, colorless tips, 960 tips (10 racks × 96 tips)           | 0030 078 535 | 181,00 € |
| PCR clean and sterile, 2 – 100 µL, 53 mm, yellow, colorless tips, 960 tips (10 racks × 96 tips)          | 0030 078 543 | 181,00 € |
| PCR clean and sterile, 2 – 200 µL, 55 mm, yellow, colorless tips, 960 tips (10 racks × 96 tips)          | 0030 078 551 | 181,00 € |
| PCR clean and sterile, 20 – 300 µL, 55 mm, orange, colorless tips, 960 tips (10 racks × 96 tips)         | 0030 078 560 | 181,00 € |
| PCR clean and sterile, 50 – 1,000 µL, 76 mm, blue, colorless tips, 960 tips (10 racks × 96 tips)         | 0030 078 578 | 195,00 € |
| PCR clean and sterile, 50 – 1,250 µL L, 103 mm, dark green, colorless tips, 480 tips (5 racks × 96 tips) | 0030 078 594 | 100,00 € |
| PCR clean and sterile, 0,25 – 2,5 mL, 115 mm, red, colorless tips, 240 tips (5 racks × 48 tips)          | 0030 078 586 | 49,50 €  |
| PCR clean and sterile, 0,1 – 5 mL, 120 mm, violet, colorless tips, 120 tips (5 racks × 24 tips)          | 0030 078 616 | 98,70 €  |
| PCR clean and sterile, 0,2 – 5 mL L, 175 mm, violet, colorless tips, 120 tips (5 racks × 24 tips)        | 0030 078 624 | 97,70 €  |
| Forensic DNA Grade, 0,1 – 10 µL, medium gray, colorless tips, 960 tips (10 racks × 96 tips)              | 0030 078 810 | 213,00 € |
| Forensic DNA Grade, 2 – 20 µL, yellow, colorless tips, 960 tips (10 racks × 96 tips)                     | 0030 078 829 | 205,00 € |
| Forensic DNA Grade, 2 – 200 µL, yellow, colorless tips, 960 tips (10 racks × 96 tips)                    | 0030 078 837 | 205,00 € |
| Forensic DNA Grade, 50 – 1,000 µL, blue, colorless tips, 960 tips (10 racks × 96 tips)                   | 0030 078 845 | 215,00 € |
| PCR clean and sterile, 0,5 – 10 mL L, 243 mm, turquoise, colorless tips, 100 tips, individually packed   | 0030 077 598 | 81,50 €  |

### ep Dualfilter T.I.P.S.<sup>®</sup> Reloads

|            |   |              |         |
|------------|---|--------------|---------|
| <b>New</b> | PCR clean and sterile, 0,1 – 10 µL S, 34 mm, dark gray, colorless tips, 960 tips (10 racks × 96 tips)     | 0030 081 021 | Inquire |
| <b>New</b> | PCR clean and sterile, 0,1 – 10 µL M, 40 mm, medium gray, colorless tips, 960 tips (10 trays × 96 tips)   | 0030 081 030 | Inquire |
| <b>New</b> | PCR clean and sterile, 0,5 – 20 µL L, 46 mm, light gray, colorless tips, 960 tips (10 trays × 96 tips)    | 0030 081 048 | Inquire |
| <b>New</b> | PCR clean and sterile, 2 – 20 µL, 53 mm, yellow, colorless tips, 960 tips (10 trays × 96 tips)            | 0030 081 056 | Inquire |
| <b>New</b> | PCR clean and sterile, 2 – 100 µL, 53 mm, yellow, colorless tips, 960 tips (10 trays × 96 tips)           | 0030 081 064 | Inquire |
| <b>New</b> | PCR clean and sterile, 2 – 200 µL, 55 mm, yellow, colorless tips, 960 tips (10 trays × 96 tips)           | 0030 081 072 | Inquire |
| <b>New</b> | PCR clean and sterile, 20 – 300 µL, 55 mm, orange, colorless tips, 960 tips (10 trays × 96 tips)          | 0030 081 080 | Inquire |
| <b>New</b> | PCR clean and sterile, 50 – 1,000 µL, 76 mm, blue, colorless tips, 960 tips (10 trays × 96 tips)          | 0030 081 099 | Inquire |
| <b>New</b> | PCR clean and sterile, 50 – 1,250 µL L, 103 mm, dark green, colorless tips, 960 tips (10 trays × 96 tips) | 0030 081 102 | Inquire |

## Tips

### ep Dualfilter T.I.P.S.<sup>®</sup> SealMax, PCR clean/sterile (sterile and pyrogen-free)

Eppendorf Product Catalog p. 54-61

Order no.

Price

#### ep Dualfilter T.I.P.S.<sup>®</sup> SealMax

|  |              |          |
|--|--------------|----------|
| PCR clean and sterile, 0,1 – 10 µL S, 34 mm, dark gray, colorless tips, 960 tips (10 racks × 96 tips)  | 0030 078 691 | 195,00 € |
| PCR clean and sterile, 0,5 – 20 µL L, 46 mm, light gray, colorless tips, 960 tips (10 racks × 96 tips) | 0030 078 705 | 186,00 € |
| PCR clean and sterile, 2 – 100 µL, 53 mm, yellow, colorless tips, 960 tips (10 racks × 96 tips)        | 0030 078 713 | 186,00 € |
| PCR clean and sterile, 2 – 200 µL, 55 mm, yellow, colorless tips, 960 tips (10 racks × 96 tips)        | 0030 078 721 | 193,00 € |
| PCR clean and sterile, 20 – 300 µL, 55 mm, orange, colorless tips, 960 tips (10 racks × 96 tips)       | 0030 078 730 | 186,00 € |
| PCR clean and sterile, 50 – 1,000 µL, 76 mm, blue, colorless tips, 960 tips (10 racks × 96 tips)       | 0030 078 748 | 202,00 € |

### ep Dualfilter T.I.P.S.<sup>®</sup> SealMax, PCR clean/sterile (sterile and pyrogen-free)

New

Eppendorf Product Catalog p. 54-61

Order no.

Price

#### ep Dualfilter T.I.P.S.<sup>®</sup> SealMax<sup>®</sup> Reloads

|   |              |         |
|---|--------------|---------|
| Biopur <sup>®</sup> , 0,1 – 10 µL S, 34 mm, dark gray, colorless tips, 960 tips (10 trays × 96 tips)  | 0030 081 226 | Inquire |
| Biopur <sup>®</sup> , 0,5 – 20 µL L, 46 mm, light gray, colorless tips, 960 tips (10 trays × 96 tips) | 0030 081 234 | Inquire |
| Biopur <sup>®</sup> , 2 – 100 µL, 53 mm, yellow, colorless tips, 960 tips (10 trays × 96 tips)        | 0030 081 242 | Inquire |
| Biopur <sup>®</sup> , 2 – 200 µL, 55 mm, yellow, colorless tips, 960 tips (10 trays × 96 tips)        | 0030 081 250 | Inquire |
| Biopur <sup>®</sup> , 20 – 300 µL, 55 mm, orange, colorless tips, 960 tips (10 trays × 96 tips)       | 0030 081 269 | Inquire |
| Biopur <sup>®</sup> , 50 – 1,000 µL, 76 mm, blue, colorless tips, 960 tips (10 trays × 96 tips)       | 0030 081 277 | Inquire |

### GEloader<sup>®</sup>

Eppendorf Product Catalog p. 54-61

Order no.

Price

**GEloader<sup>®</sup>**, Eppendorf Quality, 0,5 – 20 µL, 62 mm, light gray, 192 tips (2 racks × 96 tips)

0030 001 222

73,50 €

### Microloader

Eppendorf Product Catalog p. 54-61

Order no.

Price

**Microloader**, tip for filling Femtotips<sup>®</sup> and other glass microcapillaries (for research use only), 0,5 – 20 µL, 100 mm, light gray, 192 pcs. (2 racks × 96 pcs.)

5242 956 003

102,00 €

### epT.I.P.S.<sup>®</sup> 384

Eppendorf Product Catalog p. 54-61

Order no.

Price

#### epT.I.P.S.<sup>®</sup> 384 Reloads

|  |              |          |
|--|--------------|----------|
| PCR clean, 0,1 – 20 µL, 42 mm, light pink, colorless tips, 3.840 tips (10 trays × 384 tips)          | 0030 076 001 | 558,00 € |
| PCR clean, 5 – 100 µL, 53 mm, light yellow, colorless tips, 3.840 tips (10 trays × 384 tips)         | 0030 076 010 | 558,00 € |
| Eppendorf Quality, 0,1 – 20 µL, 42 mm, light pink, colorless tips, 3.840 tips (10 trays × 384 tips)  | 0030 076 044 | 503,00 € |
| Eppendorf Quality, 5 – 100 µL, 53 mm, light yellow, colorless tips, 3.840 tips (10 trays × 384 tips) | 0030 076 052 | 503,00 € |

#### epT.I.P.S.<sup>®</sup> 384 Set

|   |              |          |
|---|--------------|----------|
| Eppendorf Quality, 0,1 – 20 µL, 42 mm, light pink, colorless tips, 1.920 tips (5 trays × 384 tips), 1 reusable box. Replaces order no. 0030076028.  | 0030 076 273 | 285,00 € |
| Eppendorf Quality, 5 – 100 µL, 53 mm, light yellow, colorless tips, 1.920 tips (5 trays × 384 tips), 1 reusable box. Replaces order no. 0030076036. | 0030 076 281 | 285,00 € |

### Dualfilter T.I.P.S.<sup>®</sup> 384

Eppendorf Product Catalog p. 54-61

Order no.

Price

#### ep Dualfilter T.I.P.S.<sup>®</sup> 384

|  |              |          |
|--|--------------|----------|
| PCR clean and sterile, 0,1 – 20 µL, 42 mm, light pink, colorless tips, 3.840 tips (10 racks × 384 tips)  | 0030 078 853 | 848,00 € |
| PCR clean and sterile, 5 – 100 µL, 53 mm, light yellow, colorless tips, 3.840 tips (10 racks × 384 tips) | 0030 078 861 | 848,00 € |

## ep Dualfilter T.I.P.S.® PCR clean/sterile (sterile and pyrogen-free) and Forensic DNA Grade

| Eppendorf Product Catalog p. 54-61  | Order no.    | Price   |
|---|--------------|---------|
| <b>ep Dualfilter T.I.P.S.® 384 Reloads</b>  |              |         |
| PCR clean and sterile, 0,1 – 20 µL, 42 mm, light pink, 3.840 tips (10 rack × 384 tips)  | 0030 081 005 | Inquire |
| PCR clean and sterile, 5 – 100 µL, 53 mm, light yellow, 3.840 tips (10 rack × 384 tips) | 0030 081 013 | Inquire |

## Multipette® M4

| Eppendorf Product Catalog p. 66  | Order no.    | Price    |
|--|--------------|----------|
| <b>Multipette® M4</b> , 1-channel, incl. 1x Combitips® advanced 2.5 mL (green) tip and Pipette Holder 2, 1 µL – 10 mL  | 4982 000 012 | 410,00 € |
| <b>Multipette® M4 Starter Kit</b> , 1-channel, incl. Combitips® advanced Rack, Combitips® advanced assortment pack (1 tip of each size) and Pipette Holder 2, 1 µL – 10 mL | 4982 000 314 | 422,00 € |
| <b>Pipette Holder 2</b> , for one Multipette® M4, for Pipette Carousel 2 and Charger Carousel 2 or wall mounting, sticky tape included                                     | 3116 000 147 | 41,70 €  |
| <b>Battery CR2032</b> , for Multipette® M4, 1 piece  | 4980 215 003 | 7,40 €   |

## Multipette® E3/Multipette® E3x

| Eppendorf Product Catalog p. 67   | Order no.    | Price    |
|---|--------------|----------|
| <b>Multipette® E3</b> , 1-channel, incl. charging cable and Combitips® advanced assortment pack (1 tip of each size), 1 µL – 50 mL  | 4987 000 010 | 711,00 € |
| <b>Multipette® E3 bundle with charger stand</b> , 1-channel, incl. charging cable, Charger Stand 2 and Combitips® advanced assortment pack (1 tip of each size), 1 µL – 50 mL     | 4987 000 371 | 803,00 € |
| <b>Multipette® E3x</b> , 1-channel, incl. charging cable and Combitips® advanced assortment pack (1 tip of each size), 1 µL – 50 mL   | 4987 000 029 | 908,00 € |
| <b>Multipette® E3x bundle with charger stand</b> , 1-channel, incl. charging cable, Charger Stand 2 and Combitips® advanced assortment pack (1 tip of each size), 1 µL – 50 mL    | 4987 000 380 | 978,00 € |
| <b>Pipette Holder 2</b> , for one Multipette® E3/E3x or Multipette® stream/Xstream, for Pipette Carousel 2 or wall mounting, sticky tape included, without charging functionality | 3116 000 139 | 41,70 €  |
| <b>Lithium ion rechargeable battery</b> , for Multipette® E3/E3x and Multipette (X)stream   | 4986 602 009 | 18,40 €  |
| <b>Power supply with power plug adapters</b> , for Eppendorf pipettes, multi-dispenser and charger stand, 100 – 240 VAC / 5 VDC   | 4986 603 005 | 21,30 €  |

## Combitips® advanced

| Eppendorf Product Catalog p. 66-67   | Order no.    | Price    |
|--|--------------|----------|
| <b>Combitips® advanced</b>   |              |          |
| Eppendorf Quality, 0,1 mL, white, colorless tips, 100 pcs. (4 bags × 25 pcs.)              | 0030 089 405 | 103,00 € |
| Biopur®, 0,1 mL, white, colorless tips, 100 pcs., individually blister-wrapped             | 0030 089 618 | 188,00 € |
| PCR clean, 0,1 mL, white, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)           | 0030 089 766 | 129,00 € |
| Eppendorf Quality, 0,2 mL, light blue, colorless tips, 100 pcs. (4 bags × 25 pcs.)         | 0030 089 413 | 103,00 € |
| Biopur®, 0,2 mL, light blue, colorless tips, 100 pcs., individually blister-wrapped        | 0030 089 626 | 188,00 € |
| PCR clean, 0,2 mL, light blue, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)      | 0030 089 774 | 129,00 € |
| Eppendorf Quality, 0,5 mL, violet, colorless tips, 100 pcs. (4 bags × 25 pcs.)             | 0030 089 421 | 103,00 € |
| Biopur®, 0,5 mL, violet, colorless tips, 100 pcs., individually blister-wrapped            | 0030 089 634 | 188,00 € |
| PCR clean, 0,5 mL, violet, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)          | 0030 089 782 | 129,00 € |
| Eppendorf Quality, 1,0 mL, yellow, colorless tips, 100 pcs. (4 bags × 25 pcs.)             | 0030 089 430 | 103,00 € |
| Biopur®, 1,0 mL, yellow, colorless tips, 100 pcs., individually blister-wrapped            | 0030 089 642 | 188,00 € |
| Forensic DNA Grade, 1,0 mL, yellow, colorless tips, 100 pcs., individually blister-wrapped | 0030 089 855 | 144,00 € |
| PCR clean, 1,0 mL, yellow, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)          | 0030 089 790 | 129,00 € |
| Eppendorf Quality, 2,5 mL, green, colorless tips, 100 pcs. (4 bags × 25 pcs.)              | 0030 089 448 | 103,00 € |
| Biopur®, 2,5 mL, green, colorless tips, 100 pcs., individually blister-wrapped             | 0030 089 650 | 188,00 € |
| Forensic DNA Grade, 2,5 mL, green, colorless tips, 100 pcs., individually blister-wrapped  | 0030 089 863 | 144,00 € |
| PCR clean, 2,5 mL, green, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)           | 0030 089 804 | 129,00 € |

## Tips

### Combitips® advanced

| Eppendorf Product Catalog p. 66-67   | Order no.    | Price    |
|--|--------------|----------|
| <b>Combitips® advanced</b>   |              |          |
| Eppendorf Quality, 5,0 mL, blue, colorless tips, 100 pcs. (4 bags × 25 pcs.)   | 0030 089 456 | 103,00 € |
| Biopur®, 5,0 mL, blue, colorless tips, 100 pcs., individually blister-wrapped  | 0030 089 669 | 188,00 € |
| Forensic DNA Grade, 5,0 mL, blue, colorless tips, 100 pcs., individually blister-wrapped   | 0030 089 871 | 144,00 € |
| PCR clean, 5,0 mL, blue, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)  | 0030 089 812 | 129,00 € |
| Eppendorf Quality, 10 mL, orange, colorless tips, 100 pcs. (4 bags × 25 pcs.)  | 0030 089 464 | 103,00 € |
| Biopur®, 10 mL, orange, colorless tips, 100 pcs., individually blister-wrapped   | 0030 089 677 | 188,00 € |
| PCR clean, 10 mL, orange, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.)   | 0030 089 820 | 129,00 € |
| Eppendorf Quality, 25 mL, red, colorless tips, 100 pcs. (4 boxes × 25 pcs.), each box containing one adapter                           | 0030 089 472 | 121,00 € |
| Biopur®, 25 mL, red, colorless tips, 100 pcs. (4 bags × 25 pcs.), individually blister-wrapped, each box containing one adapter        | 0030 089 685 | 221,00 € |
| PCR clean, 25 mL, red, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.), each box containing one adapter                         | 0030 089 839 | 153,00 € |
| Eppendorf Quality, 50 mL, light gray, colorless tips, 100 pcs. (4 boxes × 25 pcs.), each box containing one adapter                    | 0030 089 480 | 134,00 € |
| Biopur®, 50 mL, light gray, colorless tips, 100 pcs. (4 bags × 25 pcs.), individually blister-wrapped, each box containing one adapter | 0030 089 693 | 238,00 € |
| PCR clean, 50 mL, light gray, colorless tips, 100 pcs. (4 reclosable bags × 25 pcs.), each box containing one adapter                  | 0030 089 847 | 164,00 € |
| <b>Adapter advanced, 25 mL, autoclavable</b>   |              |          |
| Eppendorf Quality, red, 1 piece  | 0030 089 715 | 11,30 €  |
| Biopur®, red, 7 pcs., individually wrapped   | 0030 089 731 | 82,50 €  |
| <b>Adapter advanced, 50 mL, autoclavable</b>   |              |          |
| Eppendorf Quality, light gray, 1 piece   | 0030 089 723 | 11,30 €  |
| Biopur®, light gray, 7 pcs., individually wrapped  | 0030 089 740 | 82,50 €  |
| <b>Combitips® advanced Rack, for 8 Combitips® advanced (0.1 – 10 mL), Eppendorf Quality, 1 piece</b>                                   | 0030 089 758 | 40,10 €  |
| <b>Combitips® advanced Assortment Pack, contains one Combitips® advanced sample of each size and Adapter advanced 25 and 50 mL</b>     |              |          |
| Eppendorf Quality  | 0030 089 936 | 15,10 €  |

### ViscoTip®

| Eppendorf Product Catalog p. 66-67   | Order no.    | Price    |
|--|--------------|----------|
| <b>ViscoTip®</b> , Eppendorf Quality, 10 mL, orange, colorless tips, 100 pcs. (4 bags × 25 pcs.) | 0030 089 502 | 113,00 € |

### Varipette® 4720 with Varitips® P and Varitips® S

| Eppendorf Product Catalog p. 70   | Order no.    | Price    |
|---|--------------|----------|
| <b>Varipette® 4720</b> , 1-channel, 1 – 10 mL (continuous volume selection)   | 4720 000 011 | 551,00 € |
| <b>Eppendorf Varitips® P</b> , to remove liquid from smaller vessels, 100 pcs.  | 0030 048 130 | 122,00 € |
| <b>Eppendorf Varitips® S Starter Kit</b> , to remove liquid from narrow-necked vessels and measuring flasks, consisting of 100 Maxitips, 10 dispensing parts, 10 valves | 0030 050 525 | 92,20 €  |
| <b>Eppendorf Varitips® S</b> , graduated, 200 pcs.  | 0030 050 568 | 87,30 €  |
| <b>Eppendorf Varitips® S valve</b> , 100 pcs.   | 0030 050 541 | 71,10 €  |
| <b>Eppendorf Varitips® S dispensing part</b> , 30 pcs.  | 0030 050 533 | 87,30 €  |
| <b>Wall mount for Varipette® 4720 &amp; Multipette® plus</b> , with adhesive surface  | 4981 203 009 | 22,80 €  |



## Eppendorf Top Buret®

Eppendorf Product Catalog p. 69

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Eppendorf Top Buret M</b> , manual bottle-top buret with discharge valve, telescopic aspiration tube and variable discharge tube, 2.5 mL per turn, digital display ranging from 0.01 – 999.9 mL, for 45 mm bottle thread, incl. 3 PP adapters for 40/38/32 mm threads, 2 micro batteries 1.5 V, certificate of conformity | 4965 000 017 | 928,00 € |
| <b>Eppendorf Top Buret H</b> , manual bottle-top buret with discharge valve, telescopic aspiration tube and variable discharge tube, 5 mL per turn, digital display ranging from 0.01 – 999.9 mL, for 45 mm bottle thread, incl. 3 PP adapters for 40/38/32 mm threads, 2 micro batteries 1.5 V, certificate of conformity   | 4965 000 025 | 928,00 € |

## Accessories for Eppendorf Top Buret®

Eppendorf Product Catalog p. 69

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Drying tube</b> , for Eppendorf Top Buret, without drying agent   | 4960 851 000 | 42,80 € |
| <b>Telescopic aspirating tube</b> , for Eppendorf Top Buret, for 2.5 mL, 5 mL, 10 mL, 25 mL, 50 mL, 100 mL | 4960 805 009 | 28,10 € |
| <b>Batteries (2 pcs.)</b> , for Eppendorf Top Buret, LR03/AAA, 1.5 V                                       | 4965 625 009 | 9,00 €  |
| <b>Bottle thread adapter</b> , PP, from GL 32 to GL 25   | 4960 800 040 | 9,00 €  |
| <b>Bottle thread adapter</b> , PP, from GL 32 to GL 27   | 4960 800 139 | 10,80 € |
| <b>Bottle thread adapter</b> , PP, from GL 32 to GL/S 28   | 4960 800 058 | 9,00 €  |
| <b>Bottle thread adapter</b> , PP, from GL 45 to GL 32   | 4960 800 120 | 13,80 € |
| <b>Bottle thread adapter</b> , PP, from GL 45 to GL 38   | 4960 800 155 | 4,40 €  |
| <b>Bottle thread adapter</b> , PP, from GL 45 to S 40  | 4960 800 147 | 4,10 €  |
| <b>Bottle thread adapter</b> , PP, from GL 32 to NS 19/26  | 4960 800 082 | 22,60 € |
| <b>Bottle thread adapter</b> , PP, from GL 32 to NS 24/29  | 4960 800 090 | 24,60 € |
| <b>Bottle thread adapter</b> , PP, from GL 32 to NS 29/32  | 4960 800 104 | 25,30 € |

## Tips

### Varispenser® 2/Varispenser® 2x

Eppendorf Product Catalog p. 72-73

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Varispenser® 2</b> , 1-channel, bottle-top dispenser with telescopic intake tube (length 125 – 240 mm)                          |              |          |
| 0,2 – 2 mL, incl. adapters GL 25, GL 28/S 28, GL 32, GL 38, S 40   | 4966 000 010 | 305,00 € |
| 0,5 – 5 mL, incl. adapters GL 25, GL 28/S 28, GL 32, GL 38, S 40   | 4966 000 029 | 305,00 € |
| 1 – 10 mL, incl. adapters GL 25, GL 28/S 28, GL 32, GL 38, S 40  | 4966 000 037 | 305,00 € |
| <b>Varispenser® 2</b> , 1-channel, bottle-top dispenser with telescopic intake tube (length 170 – 330 mm)                          |              |          |
| 2,5 – 25 mL, incl. adapters GL 32, GL 38, S 40   | 4966 000 045 | 406,00 € |
| 5 – 50 mL, incl. adapters GL 32, GL 38, S 40   | 4966 000 053 | 446,00 € |
| 10 – 100 mL, incl. adapters GL 32, GL 38, S 40   | 4966 000 061 | 587,00 € |
| <b>Varispenser® 2x</b> , 1-channel, bottle top dispenser with recirculation valve and telescopic intake tube (length 125 – 240 mm) |              |          |
| 0,2 – 2 mL, incl. adapters GL 25, GL 28/S 28, GL 32, GL 38, S 40   | 4967 000 014 | 361,00 € |
| 0,5 – 5 mL, incl. adapters GL 25, GL 28/S 28, GL 32, GL 38, S 40   | 4967 000 022 | 361,00 € |
| 1 – 10 mL, incl. adapters GL 25, GL 28/S 28, GL 32, GL 38, S 40  | 4967 000 030 | 361,00 € |
| <b>Varispenser® 2x</b> , 1-channel, bottle top dispenser with recirculation valve and telescopic intake tube (length 170 – 330 mm) |              |          |
| 2,5 – 25 mL, incl. adapters GL 32, GL 38, S 40   | 4967 000 049 | 482,00 € |
| 5 – 50 mL, incl. adapters GL 32, GL 38, S 40   | 4967 000 057 | 521,00 € |
| 10 – 100 mL, incl. adapters GL 32, GL 38, S 40   | 4967 000 065 | 654,00 € |

### Accessories for Varispenser® 2/Varispenser® 2x

Eppendorf Product Catalog p. 72-73

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Drying tube</b> , for Varispenser® 2/2x, without filling, including sealing washer (PTFE)              | 4966 509 002 | 46,40 €  |
| <b>Ventilation screw</b> , for microfilters with Luer connector, PP, with PTFE sealing washer             | 4966 511 007 | 41,90 €  |
| <b>Ventilation screw</b> , for Varispenser® 2/2x, PP  | 4966 601 006 | Inquire  |
| <b>Discharge tube</b> , for Varispenser® 2/2x, spiraled, for 2 mL, 5 mL, 10 mL                            | 4966 501 001 | 116,00 € |
| <b>Discharge tube</b> , for Varispenser® 2/2x, spiraled, for 25 mL, 50 mL, 100 mL                         | 4966 502 008 | 116,00 € |
| <b>Telescopic aspirating tube</b> , FEP; for Varispenser® 2/2x, 70 mm – 140 mm, for 2 mL, 5 mL, 10 mL     | 4966 503 004 | 33,90 €  |
| <b>Telescopic aspirating tube</b> , FEP; for Varispenser® 2/2x, 125 mm – 240 mm, for 2 mL, 5 mL, 10 mL    | 4966 504 000 | 33,90 €  |
| <b>Telescopic aspirating tube</b> , FEP; for Varispenser® 2/2x, 195 mm – 350 mm, for 2 mL, 5 mL, 10 mL    | 4966 505 007 | 33,90 €  |
| <b>Telescopic aspirating tube</b> , FEP; for Varispenser® 2/2x, 250 mm – 480 mm, for 2 mL, 5 mL, 10 mL    | 4966 506 003 | 33,90 €  |
| <b>Telescopic aspirating tube</b> , FEP; for Varispenser® 2/2x, 250 mm – 480 mm, for 25 mL, 50 mL, 100 mL | 4966 508 006 | 33,90 €  |
| <b>Telescopic aspirating tube</b> , FEP; for Varispenser® 2/2x, 170 mm – 330 mm, for 25 mL, 50 mL, 100 mL | 4966 507 000 | 33,90 €  |
| <b>Bottle thread adapter</b> , PP, from GL 32 to GL 25  | 4960 800 040 | 9,00 €   |
| <b>Bottle thread adapter</b> , PP, from GL 32 to GL 27  | 4960 800 139 | 10,80 €  |
| <b>Bottle thread adapter</b> , PP, from GL 32 to GL/S 28  | 4960 800 058 | 9,00 €   |
| <b>Bottle thread adapter</b> , PP, from GL 45 to GL 32  | 4960 800 120 | 13,80 €  |
| <b>Bottle thread adapter</b> , PP, from GL 45 to GL 38  | 4960 800 155 | 4,40 €   |
| <b>Bottle thread adapter</b> , PP, from GL 45 to S 40   | 4960 800 147 | 4,10 €   |
| <b>Bottle thread adapter</b> , ETFE, from GL 32 to GL 25  | 4966 614 000 | 33,90 €  |
| <b>Bottle thread adapter</b> , ETFE, from GL 32 to GL/S 28  | 4960 835 005 | 20,40 €  |
| <b>Bottle thread adapter</b> , ETFE, from GL 45 to GL 32  | 4966 615 007 | 33,90 €  |
| <b>Bottle thread adapter</b> , ETFE, from GL 45 to GL 38  | 4960 839 000 | 34,30 €  |
| <b>Bottle thread adapter</b> , PTFE, from GL 45 to S 40   | 4960 834 009 | 28,60 €  |
| <b>Thread adapter</b> , for jerrycan 5 L, ETFE, from 45 mm to 17/8" thread                                | 4960 832 006 | 58,70 €  |

### Easypet® 3

Eppendorf Product Catalog p. 74

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Easypet® 3</b> , 1-channel, incl. mains/power supply device, wall mounting device, shelf stand, 2 membrane filters<br>0.45 µm, 0,1 – 100 mL | 4430 000 018 | 528,00 € |
| <b>Membrane filter</b> , for Easypet® 3, 0.2 µm, sterile, 5 pcs.   | 4430 606 005 | 110,00 € |
| <b>Membrane filter</b> , for Easypet® 3, 0.45 µm, sterile, 5 pcs.  | 4421 601 009 | 110,00 € |
| <b>Lithium polymer rechargeable battery</b> , for Easypet® 3   | 4430 605 009 | 31,50 €  |

### Pipet Helper®

Eppendorf Product Catalog p. 75

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Pipet Helper®</b> , 1-channel, 0,1 – 100 mL                           | 4423 000 010 | 96,00 € |
| <b>Membrane filter</b> , for Pipet Helper®, 3 µm, non-sterile, 10 pieces | 4423 601 014 | 62,70 € |

## Eppendorf Serological Pipets

Eppendorf Product Catalog p. 75

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Eppendorf Serological Pipets</b> , sterile, free of detectable pyrogens, DNA, RNase and DNase. Non-cytotoxic                                      |              |          |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free. Non-cytotoxic, 1 mL, yellow, 800 pcs. (4 × 200 pcs.), individually blister-wrapped  | 0030 127 692 | 140,00 € |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free. Non-cytotoxic, 2 mL, green, 600 pcs. (4 × 150 pcs.), individually blister-wrapped   | 0030 127 706 | 118,00 € |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free. Non-cytotoxic, 5 mL, blue, 400 pcs. (4 × 100 pcs.), individually blister-wrapped    | 0030 127 714 | 100,00 € |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free. Non-cytotoxic, 10 mL, orange, 400 pcs. (4 × 100 pcs.), individually blister-wrapped | 0030 127 722 | 112,00 € |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free. Non-cytotoxic, 25 mL, red, 200 pcs. (4 × 50 pcs.), individually blister-wrapped     | 0030 127 730 | 106,00 € |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free. Non-cytotoxic, 50 mL, violet, 160 pcs. (4 × 40 pcs.), individually blister-wrapped  | 0030 127 749 | 247,00 € |

## epMotion® 5073t NGS Solution

Eppendorf Product Catalog p. 83

|  | Order no.    | Price       |
|--|--------------|-------------|
| <b>epMotion® 5073t NGS solution</b> , MultiCon PC, 3 tools, gripper, plus NGS specific accessories and consumables, waste bags and holder, 100 – 240 V ±10 %/50 – 60 Hz ±5 % (EU/AU/CN/GB/US), 0,2 µL – 1 mL | 5073 000 978 | 67.380,00 € |

## epMotion® 96 and 96xl

Eppendorf Product Catalog p. 84

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>epMotion® 96</b> , 96-channel, semi-automated electronic pipette for parallel 96 channel microplate processing, 100 – 240 V ±10 %/50 – 60 Hz ±5 % (EU), 0,5 – 300 µL  | 5069 000 012 | Inquire  |
| <b>epMotion® 96 with 2-position slider</b> , 96-channel, with 2-position slider, semi-automated electronic pipette for parallel 96 channel microplate processing, 100 – 240 V ±10 %/50 – 60 Hz ±5 % (EU), 0,5 – 300 µL | 5069 000 110 | Inquire  |
| <b>epMotion® 96xl</b> , 96-channel, semi-automated electronic pipette for parallel 96 channel microplate processing, 100 – 240 V ±10 %/50 – 60 Hz ±5 % (EU), 5 – 1,000 µL  | 5069 000 217 | Inquire  |
| <b>epMotion® 96xl with 2-position slider</b> , 96-channel, semi-automated electronic pipette for parallel 96 channel microplate processing, 100 – 240 V ±10 %/50 – 60 Hz ±5 % (EU), 5 – 1,000 µL                       | 5069 000 314 | Inquire  |
| <b>Upgrade set 2-position slider for epMotion® 96</b>  | 5069 074 008 | Inquire  |
| <b>epMotion® reservoir 195 mL</b> , low dead volume reservoir, single well with 96 cavities, PCR clean, 4 × 5 reservoirs, 195 mL, polypropylene (PP)   | 0030 126 556 | 360,00 € |

## Automated Pipetting

### epMotion® 5070/5070 PCR Solution

Eppendorf Product Catalog p. 85

**epMotion® 5070 MultiCon**, completely contained housing, system incl. MultiCon, epBlue software, keyboard, mouse, waste box, 100 – 240 V ±10 %/50 – 60 Hz ±5 % (EU), 0,2 µL – 1 mL

| Order no.    | Price       |
|--------------|-------------|
| 5070 000 282 | 23.650,00 € |

**epMotion® 5070 MultiCon PCR solution**, includes MultiCon PC, dispensing tool (TS 50), PCR specific accessories, 100 – 240 V ±10 %/50 – 60 Hz ±5 % (EU/AU/CN/GB/US), 0,2 µL – 1 mL

|              |             |
|--------------|-------------|
| 5070 000 948 | 19.060,00 € |
|--------------|-------------|

### epMotion® 5073I

Eppendorf Product Catalog p. 85

**epMotion® 5073I MultiCon**, completely contained housing system incl. MultiCon, epBlue software, keyboard, mouse, waste bags and holder, 100 – 240 V ±10 %/50 – 60 Hz ±5 % (EU), 0,2 µL – 1 mL

| Order no.    | Price       |
|--------------|-------------|
| 5073 000 590 | 43.210,00 € |

### epMotion® 5073t

Eppendorf Product Catalog p. 86

**epMotion® 5073t**, completely contained housing, system incl. MultiCon, Eppendorf ThermoMixer®, epBlue software, keyboard, mouse, waste bags and holder, 100 – 240 V ±10 %/50 – 60 Hz ±5 % (EU/AU/CN/GB/US), 0,2 µL – 1 mL

| Order no.    | Price       |
|--------------|-------------|
| 5073 000 345 | 45.640,00 € |

### epMotion® 5075I

Eppendorf Product Catalog p. 86

**epMotion® 5075I**, basic device incl. epBlue software, mouse, waste box, 100 – 240 V ±10 %/50 – 60 Hz ±5 % (EU/AU/CN/GB/US), 0,2 µL – 1 mL

| Order no.    | Price       |
|--------------|-------------|
| 5075 000 301 | 52.540,00 € |

### epMotion® 5075t

Eppendorf Product Catalog p. 87

**epMotion® 5075t**, basic device incl. Eppendorf ThermoMixer®, epBlue software, mouse, waste box, 100 – 240 V ±10 %/50 – 60 Hz ±5 %, 0,2 µL – 1 mL

| Order no.    | Price       |
|--------------|-------------|
| 5075 000 302 | 61.600,00 € |

### epMotion® 5075v/5075vt

Eppendorf Product Catalog p. 87

**epMotion® 5075v**, basic device incl. vacuum system, gripper, vac frame 2, vac frame holder, epBlue software, mouse, waste box, 100 – 240 V ±10 %/50 – 60 Hz ±5 %, 0,2 µL – 1 mL

| Order no.    | Price       |
|--------------|-------------|
| 5075 000 303 | 60.780,00 € |

**epMotion® 5075vt**, basic device incl. vacuum system, gripper, vac frame 2, vac frame holder, Eppendorf ThermoMixer®, epBlue software, mouse, waste box, 100 – 240 V ±10 %/50 – 60 Hz ±5 %, 0,2 µL – 1 mL

|              |             |
|--------------|-------------|
| 5075 000 304 | 72.450,00 € |
|--------------|-------------|

### epMotion® 5075 NGS solution

Eppendorf Product Catalog p. 88

**epMotion® 5075t NGS solution**, package with completely contained housing, MultiCon PC, Enhanced feature set 1, C2 thermal module, dispensing tools, plus NGS specific accessories, plus NGS specific consumables to start automated library preparation, 100 – 240 V ±10 %/50 – 60 Hz ±5 % (EU/AU/CN/GB/US), 0,2 µL – 1 mL

| Order no.    | Price        |
|--------------|--------------|
| 5075 000 962 | 103.500,00 € |

**epMotion® 5075tc NGS solution**, package with CleanCap, MultiCon PC, Enhanced feature set 1, C2 thermal module, dispensing tools, plus NGS specific accessories, waste bag holder, plus NGS specific consumables to start automated library preparation, 100 – 240 V ±10 %/50 – 60 Hz ±5 % (EU/AU/CN/GB/US), 0,2 µL – 1 mL

|              |              |
|--------------|--------------|
| 5075 000 963 | 111.000,00 € |
|--------------|--------------|

## epT.I.P.S.<sup>®</sup> Motion as Reload System

| Eppendorf Product Catalog p. 89                                       | Order no.    | Price    |
|---|--------------|----------|
| <b>epT.I.P.S.<sup>®</sup> Motion as Reload System, without filter</b> |              |          |
| Eppendorf Quality, 10 µL, 2.304 tips (24 trays × 96 tips)             | 0030 014 545 | 238,00 € |
| Eppendorf Quality, 50 µL, 2.304 tips (24 trays × 96 tips)             | 0030 014 421 | 238,00 € |
| Eppendorf Quality, 300 µL, 2.304 tips (24 trays × 96 tips)            | 0030 014 464 | 238,00 € |
| Eppendorf Quality, 1,000 µL, 2.304 tips (24 trays × 96 tips)          | 0030 014 502 | 238,00 € |
| <b>epT.I.P.S.<sup>®</sup> Motion as Reload System, with filter</b>    |              |          |
| PCR clean, 10 µL, 2.304 tips (24 trays × 96 tips)                     | 0030 014 553 | 456,00 € |
| PCR clean, 50 µL, 2.304 tips (24 trays × 96 tips)                     | 0030 014 430 | 456,00 € |
| PCR clean, 300 µL, 2.304 tips (24 trays × 96 tips)                    | 0030 014 472 | 456,00 € |
| PCR clean, 1,000 µL, 2.304 tips (24 trays × 96 tips)                  | 0030 014 510 | 456,00 € |
| PCR clean and sterile, 10 µL, 2.304 tips (24 trays × 96 tips)         | 0030 014 561 | 512,00 € |
| PCR clean and sterile, 50 µL, 2.304 tips (24 trays × 96 tips)         | 0030 014 529 | 512,00 € |
| PCR clean and sterile, 300 µL, 2.304 tips (24 trays × 96 tips)        | 0030 014 537 | 512,00 € |
| PCR clean and sterile, 1,000 µL, 2.304 tips (24 trays × 96 tips)      | 0030 014 570 | 512,00 € |
| <b>TipHolder, for epT.I.P.S.<sup>®</sup> Motion Reloads</b>           | 5075 751 399 | 457,00 € |
| <b>Clip for epT.I.P.S.<sup>®</sup> Motion Reloads</b>                 | 5075 751 070 | 151,00 € |
| <b>Box for epT.I.P.S.<sup>®</sup> Motion 1,000 µL Reload</b>          | 0030 014 677 | 52,20 €  |
| <b>Box for epT.I.P.S.<sup>®</sup> Motion 10/50/300 µL Reload</b>      | 0030 014 669 | 47,30 €  |

## epT.I.P.S.<sup>®</sup> Motion pipette tips

| Eppendorf Product Catalog p. 89                                       | Order no.    | Price    |
|---|--------------|----------|
| <b>epT.I.P.S.<sup>®</sup> Motion Pipette Tips, without filter</b>     |              |          |
| sterile, 10 µL, 960 tips (10 racks × 96 tips)                         | 0030 015 185 | 125,00 € |
| sterile, 50 µL, 960 tips (10 racks × 96 tips)                         | 0030 015 207 | 125,00 € |
| sterile, 300 µL, 960 tips (10 racks × 96 tips)                        | 0030 015 223 | 125,00 € |
| sterile, 1,000 µL, 960 tips (10 racks × 96 tips)                      | 0030 015 240 | 125,00 € |
| Eppendorf Quality, 10 µL, 960 tips (10 racks × 96 tips)               | 0030 014 383 | 109,00 € |
| Eppendorf Quality, 50 µL, 960 tips (10 racks × 96 tips)               | 0030 014 405 | 109,00 € |
| Eppendorf Quality, 300 µL, 960 tips (10 racks × 96 tips)              | 0030 014 448 | 109,00 € |
| Eppendorf Quality, 1,000 µL, 960 tips (10 racks × 96 tips)            | 0030 014 480 | 109,00 € |
| <b>epT.I.P.S.<sup>®</sup> Motion pipette tips, with filter</b>        |              |          |
| PCR clean and sterile, 10 µL, 960 tips (10 racks × 96 tips)           | 0030 015 193 | 240,00 € |
| PCR clean and sterile, 50 µL, 960 tips (10 racks × 96 tips)           | 0030 015 215 | 240,00 € |
| PCR clean and sterile, 300 µL, 960 tips (10 racks × 96 tips)          | 0030 015 231 | 240,00 € |
| PCR clean and sterile, 1,000 µL, 960 tips (10 racks × 96 tips)        | 0030 015 258 | 240,00 € |
| PCR clean, 10 µL, 960 tips (10 racks × 96 tips)                       | 0030 014 391 | 214,00 € |
| PCR clean, 50 µL, 960 tips (10 racks × 96 tips)                       | 0030 014 413 | 214,00 € |
| PCR clean, 300 µL, 960 tips (10 racks × 96 tips)                      | 0030 014 456 | 214,00 € |
| PCR clean, 1,000 µL, 960 tips (10 racks × 96 tips)                    | 0030 014 499 | 214,00 € |
| <b>epT.I.P.S.<sup>®</sup> Motion as SafeRack tips, without filter</b> |              |          |
| Eppendorf Quality, 50 µL, 960 tips (10 racks × 96 tips)               | 0030 014 600 | 129,00 € |
| Eppendorf Quality, 300 µL, 960 tips (10 racks × 96 tips)              | 0030 014 626 | 129,00 € |
| Eppendorf Quality, 1,000 µL, 960 tips (10 racks × 96 tips)            | 0030 014 642 | 129,00 € |
| <b>epT.I.P.S.<sup>®</sup> Motion as SafeRack tips, with filter</b>    |              |          |
| PCR clean, 50 µL, 960 tips (10 racks × 96 tips)                       | 0030 014 618 | 224,00 € |
| PCR clean, 300 µL, 960 tips (10 racks × 96 tips)                      | 0030 014 634 | 224,00 € |
| PCR clean, 1,000 µL, 960 tips (10 racks × 96 tips)                    | 0030 014 650 | 224,00 € |

## Automated Pipetting

### Device Hardware Upgrades and Software

Eppendorf Product Catalog p. 90

|  | Order no.    | Price       |
|--|--------------|-------------|
| <b>Thermal module</b> , thermal module for heating or cooling (0 – 110 °C) thermo adapters, thermoblock and thermoracks  | 5075 757 001 | 7.011,00 €  |
| <b>Work surface adapter</b> , for raising the epMotion® worktable by 5.5 cm, 4 feet  | 5070 752 001 | 349,00 €    |
| <b>Device Hardware Upgrades and Software</b> , Upgrade set to MultiCon, changes existing epMotion EasyCon to PC controller. Includes PC, keyboard, mouse, epBlue SW, PCR & Prep assistant, cables  | 5075 000 725 | 4.991,00 €  |
| <b>CleanCap upgrade set</b> , UV lamp and HEPA filter system for decontamination and pipetting under clean air conditions, for upgrading an epMotion 5073 (SN>7000) with epBlue version >40.7, optionally available  | 5073 001 333 | 6.392,00 €  |
| <b>epBlue ID barcode and tracking software incl. manual barcode scanner</b> , modular expansion of epBlue for barcode support. Incl. manual barcode reader with stand. Compatible with epBlue version 40.1 and higher  | 5075 002 701 | 5.606,00 €  |
| <b>epBlue GxP software</b> , for use in supporting process environments (according to GLP, GMP, 21 CFR part 11), as configuration option or for pre-installed epMotion® MultiCon versions with epBlue version 40.1 to 40.7, with epBlue GxP software, corresponding firmware, certificates | 5075 002 728 | 15.640,00 € |
| <b>epMotion® Editor 40</b> , software CD ROM with instructions, used to create, edit and simulate application on a PC, compatible with epBlue version >40.x  | 5075 014 220 | 1.910,00 €  |
| <b>epMotion® Editor 40</b> , additional software license   | 5075 014 300 | 422,00 €    |
| <b>Enhanced Feature Set 1</b> , license for epBlue additional features: Calculated no. of samples; Automatic tool selection; Start at command; Email notification, requires service visit and epBlue >40.6   | 5075 000 964 | 833,00 €    |
| <b>Enhanced Feature Set 2</b> , license for epBlue feature: Normalization, requires Service visit and epBlue 40.8 or higher  | 5075 000 981 | 833,00 €    |
| <b>Modification kit TMX for 5075I</b> , to modify a 5075I with Eppendorf ThermoMixer® to a 5075t (for units SN>4000)   | 5075 000 630 | 9.238,00 €  |
| <b>Extension plate 5075I to TMX and VAC</b> , to be ordered with 5075 000.620 and 5075 000.630 (for units SN>4000)   | 5075 000 640 | 1.209,00 €  |
| <b>Expansion board</b> , for supporting the cleanbench work surface  | 5070 751 005 | 517,00 €    |
| <b>Waste bin UV-shield 45 mm</b> , adapter to install at waste position. Allows using UV light by an increased waste capacity. Waste bag holder not included.  | 5075 751 976 | 205,00 €    |
| <b>Waste bin UV-shield 100 mm, incl. adapter</b> , adapter to install at waste position. Allows using UV light by an increased waste capacity for epMotions with raised worktable (5.5 cm, work surface adapters included). Waste bag holder not included.                                 | 5075 751 992 | 758,00 €    |

### epMotion® Dispensing tools

Eppendorf Product Catalog p. 91

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>TS 10 single channel dispensing tool</b> , 1-channel, 0,2 – 10 µL  | 5280 000 100 | 2.940,00 € |
| <b>TS 50 single-channel dispensing tool</b> , 1-channel, 1,0 – 50 µL  | 5280 000 010 | 2.558,00 € |
| <b>TS 300 single-channel dispensing tool</b> , 1-channel, 20 – 300 µL   | 5280 000 037 | 2.558,00 € |
| <b>TS 1000 single-channel dispensing tool</b> , 1-channel, 40 – 1,000 µL  | 5280 000 053 | 2.558,00 € |
| <b>TM 10-8 eight-channel dispensing tool</b> , 8-channel, 0,2 – 10 µL   | 5280 000 304 | 5.486,00 € |
| <b>TM 50-8 eight-channel dispensing tool</b> , 8-channel, 1,0 – 50 µL   | 5280 000 215 | 4.770,00 € |
| <b>TM 300-8 eight-channel dispensing tool</b> , 8-channel, 20 – 300 µL  | 5280 000 231 | 4.770,00 € |
| <b>TM 1000-8 eight-channel dispensing tool</b> , 8-channel, 40 – 1,000 µL   | 5280 000 258 | 4.770,00 € |
| <b>Holder for 6 dispensing tools</b>  | 5075 774 003 | 317,00 €   |
| <b>Holder for gripper</b>   | 5075 759 004 | 517,00 €   |
| <b>Gripper</b> , incl. holder, for transporting plates on the worktable and automatic operation of the vacuum manifold, additionally order Gripper Tower 5075 751 895 for 5073 with SN > 6000 | 5282 000 018 | 8.232,00 € |

### epMotion® 96 Accessories

Eppendorf Product Catalog p. 91

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Loading frame for epMotion® 96</b>   | 5069 072 005 | 334,00 €   |
| <b>Adapter 384-well plate for epMotion® 96</b>  | 5069 073 001 | 126,00 €   |
| <b>Upgrade set 2-position slider for epMotion® 96</b>   | 5069 074 008 | Inquire    |
| <b>O-ring tool set</b> , with separator, blade, hook, mounting cone, mounting cylinder  | 5069 858 505 | 341,00 €   |
| <b>Modification kit for 2-position slider (4.5 mm increments)</b> , including constant-force spring, shielding plates (epMotion 96 and 96xl), screws, guiding aid and holding devices (epMotion 96 as of serial number 2008, all epMotion 96xl) | 5069 878 506 | 2.638,00 € |

## epMotion® Reservoir Racks and Modules

Eppendorf Product Catalog p. 92

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>ReservoirRack</b> , for presenting 10 mL, 30 mL and 100 mL reagent reservoirs. Up to 7 reservoirs or Reservoir Rack modules can be presented.   | 5075 754 002 | 173,00 € |
| <b>Adapter sleeves</b> , to retrofit the 1.5/2.0 mL thermoracks for use with 0.5 mL vessels, 25 pcs.   | 5075 772 000 | 285,00 € |
| <b>Gripper Tower</b> , enables positioning epT.I.P.S.® Motion Reload trays in the parking position of the gripper and holds the gripper when it is not in use; only for use with epMotion® 5073 as of serial number 6000   | 5075 751 895 | 561,00 € |
| <b>Tip Tool</b> , tool for exclusively transferring epT.I.P.S.® Motion from racks and trays to the ReservoirRack modules for tips, incl. wall mount for safe storage of the tool   | 5285 000 000 | 403,00 € |
| <b>Reservoir Rack Module NGS</b> , for use with epMotion® ReservoirRacks, for distributing samples and reagents for NGS-specific applications, space for 20 x ILMN tubes, 4 x Eppendorf Tubes® 5.0 mL, 4 x Eppendorf Tubes 1.5 mL  | 5075 751 917 | 816,00 € |
| <b>Reservoir Rack Module PCR</b> , for use with epMotion® ReservoirRacks, for distributing samples and reagents for PCR, space for 48 x PCR tubes 0.2 mL incl. 2 x strip adapters 0.2 mL, (alternatively, space sufficient for 16 index tubes and up to 32 x PCR tubes 0.2 mL), 12 x Eppendorf Tubes 1.5/2.0 mL, 1 x tube 5.0 mL | 5075 751 933 | 875,00 € |
| <b>Reservoir Rack Module Tips</b> , for use with epMotion® ReservoirRacks, for providing up to 16 epT.I.P.S.® Motion per module and up to 7 modules per ReservoirRack (112 tips in total), set consisting of 7 modules   | 5075 751 950 | 524,00 € |
| <b>TipHolder</b> , for epT.I.P.S.® Motion Reloads  | 5075 751 399 | 457,00 € |
| <b>TipHolder 73</b> , enables positioning epT.I.P.S.® Motion Reload trays above a rack of epT.I.P.S.® Motion and thus doubles the tip capacity in an ANSI/SLAS position, reusable and autoclavable, epBlue 40.7 or higher required   | 5075 751 879 | 514,00 € |
| <b>ReservoirRack Module TC</b> , for use in epMotion® ReservoirRacks, temperable, 8 x PCR tubes 0.2 mL   | 5075 799 049 | 360,00 € |
| <b>ReservoirRack Module TC</b> , for use in epMotion® ReservoirRacks, temperable, 4 x Safe-Lock tubes 0.5/1.5/2.0 mL   | 5075 799 081 | 382,00 € |
| <b>ReservoirRack Module TC</b> , for use in epMotion® ReservoirRacks, temperable, 4 x Eppendorf Tubes® 5.0 mL  | 5075 799 340 | 375,00 € |
| <b>ReservoirRack Module TC</b> , for use in epMotion® ReservoirRacks, temperable, 4 x reaction vessels Ø 12 mm   | 5075 799 103 | 382,00 € |
| <b>ReservoirRack Module TC</b> , for use in epMotion® ReservoirRacks, temperable, 4 x reaction vessels Ø 16 mm   | 5075 799 120 | 382,00 € |
| <b>ReservoirRack Module TC</b> , for use in epMotion® ReservoirRacks, temperable, 4 x reaction vessels Ø 17 mm   | 5075 799 162 | 382,00 € |
| <b>ReservoirRack Module TC</b> , for use in epMotion® ReservoirRacks, temperable, 2 x reaction vessels Ø 29 mm   | 5075 799 189 | 382,00 € |
| <b>ReservoirRack Module TC</b> , for use in epMotion® ReservoirRacks, temperable, 1 x reservoir 10 mL  | 5075 799 421 | 156,00 € |
| <b>ReservoirRack Module TC</b> , for use in epMotion® ReservoirRacks, temperable, 1 x reservoir 30 mL  | 5075 799 146 | 156,00 € |
| <b>ReservoirRack Module TC</b> , for use in epMotion® ReservoirRacks, temperable, 1 x reservoir 100 mL   | 5075 799 260 | 156,00 € |
| <b>epMotion® reservoir 10 mL</b> , small-volume container for reagent presentation on the epMotion®, can only be inserted with the Reservoir Rack; production batch testing, certified, 10 mL, PCR clean, 10 x 5 reservoirs in the bag, polypropylene (PP)   | 0030 126 521 | 181,00 € |
| <b>epMotion® reservoir 30 mL</b> , large-volume container for reagent presentation on the epMotion®, can only be inserted with the Reservoir Rack; production batch testing, certified, 30 mL, PCR clean, 10 x 5 reservoirs in the bag, polypropylene (PP)   | 0030 126 505 | 181,00 € |
| <b>epMotion® reservoir 100 mL</b> , large-volume container for reagent presentation on the epMotion®, can only be inserted with the Reservoir Rack; production batch testing, certified, 100 mL, PCR clean, 10 x 5 reservoirs in the bag, polypropylene (PP)   | 0030 126 513 | 181,00 € |
| <b>epMotion® reservoir 400 mL</b> , large-volume container for reagent presentation on the epMotion®; production batch testing, certified, 400 mL, PCR clean, 400 mL, 10 reservoirs, PCR clean, autoclavable, polypropylene, polypropylene (PP)  | 5075 751 364 | 241,00 € |
| <b>epMotion® reservoir 195 mL</b> , low dead volume reservoir, single well with 96 cavities, PCR clean, 4 x 5 reservoirs, 195 mL, polypropylene (PP)   | 0030 126 556 | 360,00 € |

## Automated Pipetting

### epMotion® Thermoracks and Plate Accessories

Eppendorf Product Catalog p. 93

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Thermorack CB</b> , for up to 384 × 0.1 mL strip tubes in the Qiagen® Rotor-Gene®   | 5075 767 031 | 1.232,00 € |
| <b>Eppendorf Magnum FLX® Magnet Adapter</b> , universal ring magnet plate for rapid magnetic bead separation. Supports most 96-well plates (standard, PCR, deepwell) | 5075 751 836 | 1.904,00 € |
| <b>Thermorack Rotor/Tubes</b> , for temperature controlled pipetting of Qiagen® Rotor-Disc® 72/100 and 20 tubes 1.5/2.0 mL   | 5075 751 526 | 986,00 €   |
| <b>Thermoblock PCR 96 OC</b> , for temperature control and safe piercing of sealed, semi-skirted PCR plates with orientation control                                 | 5075 751 666 | 801,00 €   |
| <b>Thermoadapter for PCR</b> , for temperature control of PCR plates, 96 wells, skirted  | 5075 787 008 | 745,00 €   |
| <b>Thermoadapter for PCR</b> , for temperature control of PCR plates, 384 wells, skirted   | 5075 788 004 | 816,00 €   |
| <b>Thermoadapter Frosty</b> , combination of height adapter and PCR-Cooler, for cooling skirted PCR plates   | 5075 789 000 | 678,00 €   |
| <b>Lid for Thermorack Rotor/Tubes</b> , with adapter for placing the lid on the epMotion® worktable  | 5075 751 640 | 308,00 €   |
| <b>Thermoblock DWP 2000</b> , for 2 mL Eppendorf Deepwell Plates 96/2000   | 5075 751 330 | 574,00 €   |
| <b>Thermoadapter DWP 96</b> , for Deepwell Plates 96/1000  | 5075 751 054 | 1.436,00 € |
| <b>Thermoadapter LC Sample</b> , for MagNA Pure LC sample cartridge  | 5075 751 305 | 737,00 €   |
| <b>Thermorack</b> , for 24 Safe-Lock tubes, temperable, 0.5/1.5/2.0 mL   | 5075 769 000 | 1.006,00 € |
| <b>Thermorack</b> , for 24 Safe-Lock tubes, temperable, 1.5/2.0 mL   | 5075 771 004 | 957,00 €   |
| <b>Thermorack TMX</b> , for 24 Safe-Lock micro test tubes, 0.5 mL/1.5 mL/2.0 mL  | 5075 751 160 | 965,00 €   |
| <b>Thermorack TMX</b> , for 24 Safe-Lock tubes, temperable, 1.5/2.0 mL   | 5075 751 186 | 965,00 €   |
| <b>Thermorack</b> , for 24 Cryo tubes  | 5075 777 055 | 960,00 €   |
| <b>Thermoblock for PCR</b> , for temperature control of 96 vessels 0.2 mL or one PCR plate 96 wells  | 5075 766 000 | 745,00 €   |
| <b>Thermoblock for PCR</b> , for temperature control of PCR plate 384 wells  | 5075 767 007 | 816,00 €   |
| <b>Thermoadapter for Microplate 96/V/U</b> , for holding and temperature control of Eppendorf Microplates with V or U bottom   | 5075 751 577 | 579,00 €   |

### epMotion® Tuberacks, Adapters and Vacuum Modules

Eppendorf Product Catalog p. 94

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>PrepRack for 24 Eppendorf Safe-Lock Tubes 2 mL</b> , for nucleic acid preparation                                  | 5073 751 006 | 1.035,00 € |
| <b>Rack for single test tubes</b> , for presenting Eppendorf reaction vessels, glass or plastic tubes, not temperable |              |            |
| Ø 17 mm × 100 mm max. length  | 5075 761 009 | 317,00 €   |
| Ø 17 mm × 60 mm max. length   | 5075 775 000 | 317,00 €   |
| Ø 16 mm × 100 mm max. length  | 5075 760 002 | 317,00 €   |
| Ø 16 mm × 60 mm max. length   | 5075 776 006 | 317,00 €   |
| Ø 15 mm × 100 mm max. length  | 5075 792 028 | 317,00 €   |
| Ø 15 mm × 60 mm max. length   | 5075 792 044 | 317,00 €   |
| Ø 14 mm × 100 mm max. length  | 5075 792 001 | 317,00 €   |
| Ø 14 mm × 60 mm max. length   | 5075 792 060 | 317,00 €   |
| Ø 13 mm × 100 mm max. length  | 5075 762 005 | 317,00 €   |
| Ø 13 mm × 60 mm max. length   | 5075 792 087 | 317,00 €   |
| Ø 12 mm × 100 mm max. length  | 5075 763 001 | 317,00 €   |
| Ø 12 mm × 60 mm max. length   | 5075 792 109 | 317,00 €   |
| <b>Rack</b> , for 24 HPLC tubes, Ø 12 mm × 40 mm max. length  | 5075 792 125 | 317,00 €   |
| <b>Rack</b> , for 96 screw cap tubes, 1.5/2.0 mL  | 5075 791 005 | 808,00 €   |
| <b>Rack for 24 tubes</b> , for 24 Safe-Lock tubes, no tempering   |              |            |
| 0.5/1.5/2.0 mL  | 5075 751 453 | 570,00 €   |
| 1.5/2.0 mL  | 5075 751 275 | 544,00 €   |



## epMotion® Tuberacks, Adapters and Vacuum Modules

Eppendorf Product Catalog p. 94

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Vac holder</b> , adapter for Vac Frame   | 5075 778 009 | 444,00 €   |
| <b>Vacuum lid</b> , lid for vacuum unit   | 5075 779 005 | 308,00 €   |
| <b>Silicone mat for Vac Lid</b>   | 5075 793 008 | 209,00 €   |
| <b>Vac Frame 1</b> , 16 mm height, recommended for filter plates from 5Prime, Invitex   | 5075 784 009 | 931,00 €   |
| <b>Vac Frame 2</b> , 26 mm height, recommended for filter plates from MACHERY NAGEL®, Promega® or Qiagen®   | 5075 785 005 | 931,00 €   |
| <b>Anticontamination plate adapter</b> , for vacuum processes with foaming solutions in multiwell plates, set of 10   | 5075 794 004 | 49,90 €    |
| <b>Vac Thermo lid</b>   | 5075 796 007 | 1.583,00 € |
| <b>Collection plate adapter MN</b> , type for VAC station   | 5075 785 064 | 106,00 €   |
| <b>Elution plate adapter</b> , for processing micro tubes in racks  | 5075 785 030 | 108,00 €   |
| <b>Rack Smart</b> , for holding one tube rack of SmartCycler® reaction tubes  | 5075 790 009 | 291,00 €   |
| <b>Rack LC 20 µL/100 µL</b> , for up to 96 × 20 µL or 100 µL LightCycler capillaries for use with Centrifuge 5804/5804 R or 5810/5810 R and other devices, set of 2 | 5075 795 000 | 2.329,00 € |
| <b>Adapter for Covaris Rack 12 x 8 Strips</b> , for E-series rack 12 place 8 microTUBEStrip V2 from Covaris PN 500444   | 5075 751 801 | 180,00 €   |
| <b>Rack ILMN tubes</b> , Rack for single tubes, 40 × ø 8.4 mm and 12 × ø 11.2 mm, transportable, autoclavable   | 5075 751 747 | 580,00 €   |

## epMotion® Waste Accessories

Eppendorf Product Catalog p. 95

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Waste bag holder</b> , for epMotion® 5070/5073/5075 waste position  | 5075 753 103 | 128,00 € |
| <b>Waste container with ring</b> , for holding used pipette tips   | 5075 753 006 | 134,00 € |
| <b>Liquid Waste Tub 400 mL</b> , to mount in waste bag holder and waste container of epMotion®, with lid, autoclavable, working volume 400 mL  | 5075 751 720 | 130,00 € |
| <b>Waste bags bio-hazardous waste</b> , for epMotion® 5070/5073/5075 waste position, up to 7 L volume, ideal for disposing of biologically hazardous waste, autoclavable, PP material, thickness 50 µm, transparent, 50 bags | 5075 751 763 | 33,40 €  |
| <b>Waste bags</b> , for epMotion® 5070/5073/5075 waste position, up to 7 L volume, ideal for disposing of non hazardous waste, autoclavable, PP material, thickness 50 µm, transparent, 50 bags                              | 5075 751 780 | 29,70 €  |

## epMotion® Height adapters

Eppendorf Product Catalog p. 95

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>epMotion® height adapter</b> , for leveling labware, enables faster plate processing |              |          |
| 85 mm   | 5075 751 003 | 444,00 € |
| 55 mm   | 5075 752 000 | 444,00 € |
| 40 mm, for pipette tips   | 5075 755 009 | 444,00 € |

## Temperature Control and Mixing

### Eppendorf ThermoMixer® C

Eppendorf Product Catalog p. 100

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Eppendorf ThermoMixer® C</b> , basic device without thermoblock, 220 – 240 V/50 – 60 Hz (EU) | 5382 000 015 | 3.017,00 € |

### Eppendorf ThermoStat C

Eppendorf Product Catalog p. 101

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Eppendorf ThermoStat C</b> , basic device without thermoblock, 220 – 240 V/50 – 60 Hz (EU) | 5383 000 019 | 2.063,00 € |
| <b>ThermoTop®</b> , with <i>condens.protect</i> technology                                    | 5308 000 003 | 528,00 €   |
| <b>SmartExtender</b> , 1.5 mL   | 5322 000 008 | 822,00 €   |

### Eppendorf ThermoMixer® F

Eppendorf Product Catalog p. 102

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Eppendorf ThermoMixer® FP</b> , with thermoblock for Microplates and Deepwell Plates, including lid, 220 – 240 V/50 – 60 Hz (EU) | 5385 000 016 | 3.081,00 € |
| <b>Eppendorf ThermoMixer® F0.5</b> , with thermoblock for 24 reaction vessels 0.5 mL, 200 – 240 V/50 – 60 Hz (EU)                   | 5386 000 010 | 2.958,00 € |
| <b>Eppendorf ThermoMixer® F1.5</b> , with thermoblock for 24 reaction vessels 1.5 mL, 220 – 240 V/50 – 60 Hz (EU)                   | 5384 000 012 | 2.886,00 € |
| <b>Eppendorf ThermoMixer® F2.0</b> , with thermoblock for 24 reaction vessels 2.0 mL, 220 – 240 V/50 – 60 Hz (EU)                   | 5387 000 013 | 2.886,00 € |

### MixMate®

Eppendorf Product Catalog p. 103

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>MixMate®</b> , basic device without tube holders, 230 V/50 – 60 Hz (EU)  | 5353 000 510 | 1.998,00 € |
| <b>Tube Holder 0.5 mL</b> , for 24 × 0.5 mL reaction vessels  | 5353 040 121 | 106,00 €   |
| <b>Tube Holder 1.5/2.0 mL</b> , for 24 × 1.5 mL or 2.0 mL reaction vessels  | 5353 040 130 | 106,00 €   |
| <b>Tube Holder 5/15 mL</b> , for 8 × 5/15 mL conical tubes  | 5353 040 148 | 106,00 €   |
| <b>Tube Holder 25/50 mL</b> , for 4 × 25/50 mL conical tubes  | 5353 040 156 | 106,00 €   |
| <b>Tube Holder PCR 96</b> , for 96 × 0.2 mL PCR tubes, PCR strips or 1 × 96-well PCR plate, semi-skirted or unskirted | 5353 040 113 | 106,00 €   |

### ThermoTop®

Eppendorf Product Catalog p. 104

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>ThermoTop®</b> , with <i>condens.protect</i> technology | 5308 000 003 | 528,00 € |

### SmartExtender

Eppendorf Product Catalog p. 104

|                               | Order no.    | Price    |
|-------------------------------|--------------|----------|
| <b>SmartExtender</b> , 1.5 mL | 5322 000 008 | 822,00 € |

## IsoTherm-System®

| Eppendorf Product Catalog p. 104  | Order no.    | Price    |
|---|--------------|----------|
| <b>IsoTherm-System® starter set</b> , includes IsoSafe, IsoRack, 0 °C IsoPack and -21 °C IsoPack for 1.5/2.0 mL tubes | 3880 001 018 | 173,00 € |
| <b>IsoPack and IsoRack set</b> , includes one IsoPack and one IsoRack for 1.5/2.0 mL vessels, 0 °C                    | 3880 001 166 | 56,90 €  |
| for 1.5/2.0 mL vessels, -21 °C  | 3880 001 174 | 56,90 €  |

## Eppendorf SmartBlock for Eppendorf ThermoMixer® C and Eppendorf ThermoStat C

| Eppendorf Product Catalog p. 105  | Order no.    | Price    |
|---|--------------|----------|
| <b>Eppendorf SmartBlock cryo thaw</b> , thermoblock for 24 cryo tubes, 1.5 mL – 2 mL, max. external diameter 12.7 mm, flat base. Replaces order no. 5367000025. | 5318 000 001 | 626,00 € |
| <b>Eppendorf SmartBlock 0.5 mL</b> , thermoblock for 24 reaction vessels 0.5 mL, incl. Transfer Rack 0.5 mL   | 5361 000 031 | 626,00 € |
| <b>Eppendorf SmartBlock 1.5 mL</b> , thermoblock for 24 reaction vessels 1.5 mL, incl. Transfer Rack 1.5/2.0 mL   | 5360 000 038 | 626,00 € |
| <b>Eppendorf SmartBlock 2.0 mL</b> , thermoblock for 24 reaction vessels 2.0 mL, incl. Transfer Rack 1.5/2.0 mL   | 5362 000 035 | 626,00 € |
| <b>Eppendorf SmartBlock 5 mL</b> , thermoblock for 8 Eppendorf Tubes® 5.0 mL, 8 × 5.0 mL reaction vessels   | 5309 000 007 | 626,00 € |
| <b>Eppendorf SmartBlock 15 mL</b> , thermoblock for 8 conical tubes 15 mL   | 5366 000 021 | 626,00 € |
| <b>Adapter 25 mL for Eppendorf SmartBlock 50 mL</b> , Set consisting of 4 adapters and removal tool   | 5365 000 109 | 164,00 € |
| <b>Eppendorf SmartBlock 50 mL</b> , thermoblock for 4 conical tubes 50 mL   | 5365 000 028 | 626,00 € |
| <b>Eppendorf SmartBlock 12 mm</b> , thermoblock for 24 reaction vessels, diameter up to 12 mm, height 35 mm – 76 mm   | 5364 000 024 | 778,00 € |
| <b>Eppendorf SmartBlock PCR 96</b> , thermoblock for PCR plates 96, incl. lid   | 5306 000 006 | 778,00 € |
| <b>Eppendorf SmartBlock PCR 384</b> , thermoblock for PCR plates 384, incl. lid   | 5307 000 000 | 778,00 € |
| <b>Eppendorf SmartBlock DWP 500</b> , thermoblock for Eppendorf Deepwell Plates 96/500 µL, incl. lid  | 5316 000 004 | 778,00 € |
| <b>Eppendorf SmartBlock DWP 1000</b> , thermoblock for Eppendorf Deepwell Plate 96/1,000 µL, incl. lid  | 5310 000 002 | 778,00 € |
| <b>Eppendorf SmartBlock plates</b> , thermoblock for microplates and Deepwell Plates, incl. lid   | 5363 000 039 | 778,00 € |
| <b>Lid</b> , for Eppendorf ThermoMixer® F0.5/F1.5/F2.0 and FP, for Eppendorf SmartBlock 0.5 – 2.0 mL, plates, PCR 96, PCR 384, and DWP                          | 5363 000 233 | 140,00 € |
| <b>Transfer Rack 0.5 mL</b> , for Eppendorf SmartBlock 0.5 mL   | 3880 000 305 | 49,10 €  |
| <b>Transfer Rack 1.5/2.0 mL</b> , for Eppendorf SmartBlock 1.5 mL and Eppendorf SmartBlock 2.0 mL   | 3880 000 151 | 49,10 €  |

## Tubes and Plates

### Eppendorf Safe-Lock Tubes

| Eppendorf Product Catalog p. 113   | Order no.    | Price   |
|--|--------------|---------|
| <b>Eppendorf Safe-Lock Tubes, 0,5 mL</b>   |              |         |
| Eppendorf Quality, colorless, 500 tubes. Replaces order no. 0030121503.  | 0030 121 023 | 48,00 € |
| Eppendorf Quality, yellow, 500 tubes   | 0030 121 112 | 63,60 € |
| Eppendorf Quality, red, 500 tubes  | 0030 121 120 | 63,60 € |
| Eppendorf Quality, blue, 500 tubes   | 0030 121 139 | 63,60 € |
| Eppendorf Quality, green, 500 tubes  | 0030 121 147 | 63,60 € |
| Eppendorf Quality, assorted colors (per 100), 500 tubes  | 0030 121 708 | 52,00 € |
| PCR clean, amber (light protection), 500 tubes   | 0030 121 155 | 69,30 € |
| PCR clean, colorless, 500 tubes  | 0030 123 301 | 53,10 € |
| Biopur®, colorless, 50 tubes, individually packed  | 0030 121 570 | 32,90 € |
| Forensic DNA Grade, colorless, 500 tubes (10 bags x 50 tubes)  | 0030 123 603 | 57,70 € |
| <b>Eppendorf Safe-Lock Tubes, 1,5 mL</b>   |              |         |
| Eppendorf Quality, colorless, 1.000 tubes. Replaces order no. 0030121872.  | 0030 120 086 | 48,50 € |
| Eppendorf Quality, yellow, 1.000 tubes   | 0030 120 159 | 68,00 € |
| Eppendorf Quality, red, 1.000 tubes  | 0030 120 167 | 68,00 € |
| Eppendorf Quality, blue, 1.000 tubes   | 0030 120 175 | 68,00 € |
| Eppendorf Quality, green, 1.000 tubes  | 0030 120 183 | 68,00 € |
| Eppendorf Quality, assorted colors (per 200), 1.000 tubes  | 0030 121 694 | 52,00 € |
| PCR clean, amber (light protection), 1.000 tubes   | 0030 120 191 | 74,10 € |
| PCR clean, colorless, 1.000 tubes  | 0030 123 328 | 53,10 € |
| Biopur®, colorless, 100 tubes, individually packed   | 0030 121 589 | 33,60 € |
| Forensic DNA Grade, colorless, 500 tubes (10 bags x 50 tubes)  | 0030 123 611 | 30,00 € |
| <b>Eppendorf Safe-Lock Tubes, 2,0 mL</b>   |              |         |
| Eppendorf Quality, colorless, 1.000 tubes. Replaces order no. 0030121880.  | 0030 120 094 | 66,30 € |
| Eppendorf Quality, yellow, 1.000 tubes   | 0030 120 205 | 89,00 € |
| Eppendorf Quality, red, 1.000 tubes  | 0030 120 213 | 89,00 € |
| Eppendorf Quality, blue, 1.000 tubes   | 0030 120 221 | 89,00 € |
| Eppendorf Quality, green, 1.000 tubes  | 0030 120 230 | 89,00 € |
| Eppendorf Quality, assorted colors (per 200), 1.000 tubes  | 0030 121 686 | 71,60 € |
| PCR clean, amber (light protection), 1.000 tubes   | 0030 120 248 | 96,00 € |
| PCR clean, colorless, 1.000 tubes  | 0030 123 344 | 72,30 € |
| Biopur®, colorless, 100 tubes, individually packed   | 0030 121 597 | 44,80 € |
| Forensic DNA Grade, colorless, 500 tubes (10 bags x 50 tubes)  | 0030 123 620 | 41,20 € |
| <b>Storage Box 8 x 8</b> , freezer box, for 64 tubes 1 – 2 mL, 3 pcs., height 53 mm, 2 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code | 0030 140 524 | 30,10 € |
| <b>Eppendorf Tube Rack</b> , 36 positions, for 1.5 mL and 2.0 mL tubes, polypropylene, numbered positions, autoclavable, 2 pcs.  | 0030 119 819 | 48,50 € |

### Eppendorf Tubes® 3810X

| Eppendorf Product Catalog p. 114   | Order no.    | Price   |
|--|--------------|---------|
| <b>Eppendorf Tubes® 3810X, 1,5 mL, g-safe® centrifugation stability</b>  |              |         |
| Eppendorf Quality, colorless, 1.000 pcs.   | 0030 125 150 | 42,60 € |
| Eppendorf Quality, yellow, 1.000 pcs.  | 0030 125 207 | 69,60 € |
| Eppendorf Quality, red, 1.000 pcs.   | 0030 125 193 | 69,60 € |
| Eppendorf Quality, blue, 1.000 pcs.  | 0030 125 177 | 69,60 € |
| Eppendorf Quality, green, 1.000 pcs.   | 0030 125 185 | 69,60 € |
| PCR clean, colorless, 1.000 pcs.   | 0030 125 215 | 45,50 € |
| <b>Storage Box 8 x 8</b> , freezer box, for 64 tubes 1 – 2 mL, 3 pcs., height 53 mm, 2 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code | 0030 140 524 | 30,10 € |
| <b>Eppendorf Tube Rack</b> , 36 positions, for 1.5 mL and 2.0 mL tubes, polypropylene, numbered positions, autoclavable, 2 pcs.  | 0030 119 819 | 48,50 € |

## Eppendorf Tubes® 5.0 mL

| Eppendorf Product Catalog p. 114-115  | Order no.    | Price    |
|---|--------------|----------|
| <b>Eppendorf Tubes® 5.0 mL with snap cap, 5,0 mL</b>  |              |          |
| Eppendorf Quality, colorless, 200 tubes (2 bags × 100 tubes)  | 0030 119 401 | 43,60 €  |
| PCR clean, colorless, 200 tubes (2 bags × 100 tubes)  | 0030 119 460 | 48,80 €  |
| Biopur®, colorless, 50 tubes, individually wrapped  | 0030 119 479 | 53,90 €  |
| Forensic DNA Grade, colorless, 200 tubes (4 bags × 50 tubes)  | 0030 119 606 | 55,40 €  |
| sterile, colorless, 200 tubes (10 bags × 20 tubes)  | 0030 119 487 | 53,90 €  |
| PCR clean, amber (light protection), 200 tubes (2 bags × 100 tubes)   | 0030 119 452 | 55,60 €  |
| <b>DNA LoBind® Tubes</b> , DNA LoBind®, 5,0 mL, with snap cap, PCR clean, colorless, 200 tubes (4 bags × 50 tubes)  | 0030 108 310 | 58,80 €  |
| <b>Protein LoBind® Tubes</b> , Protein LoBind®, 5,0 mL, with snap cap, PCR clean, colorless, 100 tubes (2 bags × 50 tubes)  | 0030 108 302 | 46,40 €  |
| <b>Eppendorf Tubes® 5.0 mL with screw cap, 5,0 mL</b>   |              |          |
| Eppendorf Quality, colorless, 200 tubes (2 bags × 100 tubes). Replaces order no. 0030119517.  | 0030 122 305 | 55,10 €  |
| PCR clean, colorless, 200 tubes (2 bags × 100 tubes)  | 0030 122 313 | 63,60 €  |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 200 tubes (2 bags × 100 tubes)  | 0030 122 321 | 70,80 €  |
| Forensic DNA Grade, colorless, 200 tubes (2 bags × 100 tubes)   | 0030 122 402 | 72,10 €  |
| PCR clean, amber (light protection), 200 tubes (2 bags × 100 tubes)   | 0030 122 330 | 75,50 €  |
| <b>Eppendorf Tubes® 5.0 mL with snap cap starter pack, 5,0 mL, 400 tubes (2 packages with 2 bags of 100 each), 2 racks (with 12 locations each), white, 8 universal adapters for rotors with bore holes for 15 mL conical tubes, PCR clean, colorless</b>   | 0030 119 380 | 155,00 € |
| <b>Eppendorf Tubes® 5.0 mL with screw cap, starter pack, 5,0 mL, 400 tubes (2 packages with 2 bags of 100 each), 2 racks (with 12 locations each), white, 8 universal adapters for rotors with bore holes for 15 mL conical tubes, PCR clean, colorless</b> | 0030 122 364 | 166,00 € |

## Eppendorf Tubes® 5.0 mL – System Accessories

| Eppendorf Product Catalog p. 115  | Order no.    | Price      |
|---|--------------|------------|
| <b>Tube Clip for Eppendorf Tubes® 5.0 mL, 5,0 mL, Eppendorf Quality, 10 pcs.</b>  | 0030 119 509 | 25,70 €    |
| <b>Eppendorf Tube Rack, 12 positions, for 5.0 mL and 15 mL tubes, polypropylene, numbered positions, autoclavable, 2 pcs.</b>   | 0030 119 827 | 48,50 €    |
| <b>Storage Box 5 × 5, freezer box, for 25 tubes 5 mL, 4 pcs., height 64 mm, 2.5 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code</b>         | 0030 140 532 | 42,50 €    |
| <b>Storage Box 5 × 5, freezer box, for 25 tubes 5 mL screw cap, 2 pcs., height 76 mm, 3 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code</b> | 0030 140 613 | 21,50 €    |
| <b>Rotor FA-45-12-17, aerosol-tight, incl. rotor lid, aluminum</b>  | 5409 700 006 | 895,00 €   |
| <b>Rotor FA-45-16-17, aerosol-tight, incl. rotor lid, for 16 × Eppendorf Tubes® 5.0 mL</b>  | 5427 750 002 | 895,00 €   |
| <b>Rotor FA-45-20-17, incl. rotor lid, for 20 Eppendorf Tubes® 5.0 mL</b>   | 5820 765 003 | 1.433,00 € |
| <b>Adapter, for 14 Eppendorf Tubes® 5.0 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.</b>   | 5825 739 000 | 299,00 €   |
| <b>Adapter, for 1 Eppendorf Tubes® 5.0 mL, 17 mm × 60 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs.</b>   | 5702 733 008 | 90,00 €    |
| <b>Adapter, for 1 Eppendorf Tubes® 5.0 mL, for Rotor F-35-6-30, large rotor bore, 2 pcs.</b>  | 5427 747 001 | 96,00 €    |
| <b>Adapter, for 1 Eppendorf Tubes® 5.0 mL, for Rotor F-34-6-38, 2 pcs.</b>  | 5804 777 000 | 83,40 €    |
| <b>Adapter, for 1 Eppendorf Tubes® 5.0 mL, for Rotor F-35-6-30, small rotor bore, 2 pcs.</b>  | 5427 746 005 | 96,00 €    |
| <b>Adapter, for 1 Eppendorf Tubes® 5.0 mL with snap cap, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.</b>   | 5820 730 005 | 83,40 €    |
| <b>Universal adapter, for 1 Eppendorf Tubes® 5.0 mL, for rotors for conical 15 mL tubes, 8 pcs.</b>   | 5820 732 008 | 25,50 €    |
| <b>Adapter, for 8 Eppendorf Tubes® 5.0 mL, for Rotor S-4-72, 2 pcs.</b>   | 5804 793 005 | 194,00 €   |
| <b>Exchangeable thermoblocks, for 8 × Eppendorf Tubes® 5.0 mL</b>   | 5309 000 333 | Inquire    |
| <b>Eppendorf SmartBlock 5 mL, thermoblock for 8 Eppendorf Tubes® 5.0 mL, 8 × 5.0 mL reaction vessels</b>  | 5309 000 007 | 626,00 €   |
| <b>ReservoirRack Module TC, for use in epMotion® ReservoirRacks, temperable, 4 × Eppendorf Tubes® 5.0 mL</b>  | 5075 799 340 | 375,00 €   |

## Tubes and Plates

### Eppendorf Conical Tubes

Eppendorf Product Catalog p. 116-117

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Eppendorf Conical Tubes 25 mL with SnapTec® cap, 25 mL</b>  |              |          |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, amber (light protection), 150 tubes (6 bags × 25 tubes)                                   | 0030 118 430 | 55,30 €  |
| Eppendorf Quality, colorless, 200 tubes (5 bags × 40 tubes)  | 0030 118 405 | 62,70 €  |
| PCR clean, colorless, 200 tubes (5 bags × 40 tubes)  | 0030 118 413 | 67,40 €  |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 150 tubes (6 bags × 25 tubes)  | 0030 118 421 | 55,80 €  |
| <b>Eppendorf Conical Tubes 25 mL with SnapTec® cap, 25 mL, Starter Pack, 4 Single Tube Stands, 6 adapters for rotors with bore for 50 mL conical tubes</b> |              |          |
| PCR clean, colorless, 200 tubes (5 bags × 40 tubes)  | 0030 118 502 | 111,00 € |
| <b>Eppendorf Conical Tubes 25 mL with screw cap, 25 mL</b>   |              |          |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, amber (light protection), 200 tubes (8 bags × 25 tubes)                                   | 0030 122 445 | 75,90 €  |
| Eppendorf Quality, colorless, 200 tubes (4 bags × 50 tubes)  | 0030 122 410 | 67,40 €  |
| PCR clean, colorless, 200 tubes (4 bags × 50 tubes)  | 0030 122 429 | 72,00 €  |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 200 tubes (8 bags × 25 tubes)  | 0030 122 437 | 76,60 €  |
| <b>Eppendorf Conical Tubes 25 mL with screw cap, 25 mL, Starter Pack, 4 Single Tube Stands, 6 adapters for rotors with bore for 50 mL conical tubes</b>    |              |          |
| PCR clean, colorless, 200 tubes (4 bags × 50 tubes)  | 0030 122 500 | 115,00 € |

### Eppendorf Conical Tubes

Eppendorf Product Catalog p. 118-119

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Eppendorf Conical Tubes, 15 mL</b>  |              |          |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 500 tubes (10 bags × 50 tubes)   | 0030 122 151 | 171,00 € |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, amber (light protection), 200 tubes (4 bags × 50 tubes)   | 0030 122 194 | 68,50 €  |
| Forensic DNA Grade, colorless, 100 tubes, individually wrapped   | 0030 122 259 | 84,30 €  |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 500 tubes (20 racks × 25 tubes)  | 0030 122 160 | 187,00 € |
| <b>Eppendorf Conical Tubes, 50 mL</b>  |              |          |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 500 tubes (20 bags × 25 tubes)   | 0030 122 178 | 205,00 € |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, amber (light protection), 200 tubes (8 bags × 25 tubes)   | 0030 122 224 | 82,90 €  |
| Forensic DNA Grade, colorless, 48 tubes, individually wrapped  | 0030 122 267 | 48,10 €  |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 300 tubes (12 racks × 25 tubes)  | 0030 122 186 | 161,00 € |
| <b>Eppendorf Conical Tubes SnapTec® 50, 45 mL</b>  |              |          |
| Eppendorf Quality, colorless, 240 tubes (8 bags × 30 tubes)  | 0030 118 650 | 96,00 €  |
| PCR clean, colorless, 240 tubes (8 bags × 30 tubes)  | 0030 118 669 | 99,40 €  |
| sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 200 tubes (10 bags × 20 tubes)   | 0030 118 677 | 84,80 €  |
| <b>Eppendorf Tube Rack, 12 positions, for 5.0 mL and 15 mL tubes, polypropylene, numbered positions, autoclavable, 2 pcs.</b>  | 0030 119 827 | 48,50 €  |
| <b>Eppendorf Tube Rack, 12 positions, 6 for 5.0 mL and 15 mL tubes and 6 for 25 mL and 50 mL tubes, polypropylene, numbered positions, autoclavable, 2 pcs.</b>                              | 0030 119 835 | 48,50 €  |
| <b>Storage Box 5 × 5, freezer box, for 25 tubes 15 mL, 2 pcs., height 127 mm, 5 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code</b>                  | 0030 140 583 | 22,10 €  |
| <b>Storage Box 3 × 3, freezer box, for 9 tubes 50 mL and 4 tubes 15 mL, 2 pcs., height 127 mm, 5 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code</b> | 0030 140 591 | 22,10 €  |

## Eppendorf Tubes® BioBased

New

| Eppendorf Product Catalog p. 121  | Order no.    | Price    |
|---|--------------|----------|
| <b>Eppendorf Tubes® BioBased</b>  |              |          |
| 15 mL, screw cap, sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 500 tubes (10 bags × 50 tubes)  | 0030 122 526 | 185,15 € |
| 25 mL, screw cap, sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 200 tubes (8 bags × 25 tubes)   | 0030 122 534 | 82,92 €  |
| 50 mL, screw cap, sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 500 tubes (20 bags × 25 tubes)  | 0030 122 542 | 221,95 € |
| 5.0 mL, screw cap, sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 200 tubes (2 bags × 100 tubes) | 0030 122 518 | 76,59 €  |

## Protein LoBind® Tubes

| Eppendorf Product Catalog p. 122  | Order no.    | Price    |
|---|--------------|----------|
| <b>Protein LoBind® Tubes</b> , Protein LoBind®, 25 mL, with SnapTec® cap, PCR clean, colorless, 200 tubes (5 bags × 40 tubes) | 0030 108 531 | Inquire  |
| <b>Protein LoBind® Tubes</b> , Protein LoBind®, 25 mL, with screw cap, PCR clean, colorless, 200 tubes (4 bags × 50 tubes)    | 0030 122 283 | 103,00 € |
| <b>Protein LoBind® Tubes</b> , Protein LoBind®, 0,5 mL, PCR clean, colorless, 100 tubes (2 bags × 50 tubes)                   | 0030 108 094 | 15,30 €  |
| <b>Protein LoBind® Tubes</b> , Protein LoBind®, 1,5 mL, PCR clean, colorless, 100 tubes (2 bags × 50 tubes)                   | 0030 108 116 | 11,50 €  |
| <b>Protein LoBind® Tubes</b> , Protein LoBind®, 2,0 mL, PCR clean, colorless, 100 tubes (2 bags × 50 tubes)                   | 0030 108 132 | 18,30 €  |
| <b>Protein LoBind® Tubes</b> , Protein LoBind®, 5,0 mL, with snap cap, PCR clean, colorless, 100 tubes (2 bags × 50 tubes)    | 0030 108 302 | 46,40 €  |
| <b>Protein LoBind® Tubes</b> , Protein LoBind®, 5,0 mL, with screw cap, PCR clean, colorless, 200 tubes (2 bags × 100 tubes)  | 0030 122 356 | 96,40 €  |
| <b>Protein LoBind® Tubes</b> , Protein LoBind®, 15 mL, conical tube, PCR clean, colorless, 200 tubes (4 bags × 50 tubes)      | 0030 122 216 | 97,70 €  |
| <b>Protein LoBind® Tubes</b> , Protein LoBind®, 50 mL, conical tube, PCR clean, colorless, 200 tubes (8 bags × 25 tubes)      | 0030 122 240 | 114,00 € |

## DNA LoBind® Tubes

| Eppendorf Product Catalog p. 123   | Order no.    | Price   |
|--|--------------|---------|
| <b>DNA LoBind® Tubes</b> , DNA LoBind®, 25 mL, conical tubes, PCR clean, colorless, 200 tubes (5 bags × 40 tubes)    | 0030 108 523 | 75,90 € |
| <b>DNA LoBind® Tubes</b> , DNA LoBind®, 25 mL, with screw cap, PCR clean, colorless, 200 tubes (4 bags × 50 tubes)   | 0030 122 275 | 78,10 € |
| <b>DNA LoBind® Tubes</b> , DNA LoBind®, 0,5 mL, PCR clean, colorless, 1 bag × 500 tubes                              | 0030 108 400 | 45,50 € |
| <b>DNA LoBind® Tubes</b> , DNA LoBind®, 1,5 mL, PCR clean, colorless, 1 bag × 500 tubes                              | 0030 108 418 | 25,00 € |
| <b>DNA LoBind® Tubes</b> , DNA LoBind®, 2,0 mL, PCR clean, colorless, 1 bag × 500 tubes                              | 0030 108 426 | 29,60 € |
| <b>DNA LoBind® Tubes</b> , DNA LoBind®, 0,5 mL, PCR clean, colorless, 250 tubes (5 bags × 50 tubes)                  | 0030 108 035 | 28,70 € |
| <b>DNA LoBind® Tubes</b> , DNA LoBind®, 1,5 mL, PCR clean, colorless, 250 tubes (5 bags × 50 tubes)                  | 0030 108 051 | 15,80 € |
| <b>DNA LoBind® Tubes</b> , DNA LoBind®, 2,0 mL, PCR clean, colorless, 250 tubes (5 bags × 50 tubes)                  | 0030 108 078 | 18,90 € |
| <b>DNA LoBind® Tubes</b> , DNA LoBind®, 5,0 mL, with snap cap, PCR clean, colorless, 200 tubes (4 bags × 50 tubes)   | 0030 108 310 | 58,80 € |
| <b>DNA LoBind® Tubes</b> , DNA LoBind®, 5,0 mL, with screw cap, PCR clean, colorless, 200 tubes (2 bags × 100 tubes) | 0030 122 348 | 70,80 € |
| <b>DNA LoBind® Tubes</b> , DNA LoBind®, 15 mL, conical tubes, PCR clean, colorless, 200 tubes (4 bags × 50 tubes)    | 0030 122 208 | 73,20 € |
| <b>DNA LoBind® Tubes</b> , DNA LoBind®, 50 mL, conical tubes, PCR clean, colorless, 200 tubes (8 bags × 25 tubes)    | 0030 122 232 | 87,80 € |

## Tubes and Plates

### Eppendorf Microplates 96 and 384

| Eppendorf Product Catalog p. 124                               | Order no.    | Price    |
|--|--------------|----------|
| <b>Microplate 96/F</b> , wells clear, RecoverMax® well design  |              |          |
| PCR clean, white, 80 plates (5 bags × 16 plates)               | 0030 601 106 | 235,00 € |
| sterile, white, 80 plates (5 bags × 16 plates)                 | 0030 602 102 | 269,00 € |
| <b>Microplate 96/U</b> , wells clear, RecoverMax® well design  |              |          |
| PCR clean, white, 80 plates (5 bags × 16 plates)               | 0030 601 203 | 235,00 € |
| sterile, white, 80 plates (5 bags × 16 plates)                 | 0030 602 200 | 269,00 € |
| <b>Microplate 96/V</b> , wells clear, RecoverMax® well design  |              |          |
| PCR clean, white, 80 plates (5 bags × 16 plates)               | 0030 601 300 | 235,00 € |
| sterile, white, 80 plates (5 bags × 16 plates)                 | 0030 602 307 | 269,00 € |
| <b>Microplate 384/F</b> , wells clear, RecoverMax® well design |              |          |
| PCR clean, white, 80 plates (5 bags × 16 plates)               | 0030 621 107 | 596,00 € |
| sterile, white, 80 plates (5 bags × 16 plates)                 | 0030 622 103 | 689,00 € |
| <b>Microplate 384/V</b> , wells clear, RecoverMax® well design |              |          |
| PCR clean, white, 80 plates (5 bags × 16 plates)               | 0030 621 301 | 596,00 € |
| sterile, white, 80 plates (5 bags × 16 plates)                 | 0030 622 308 | 689,00 € |

### Assay/Reader Microplates

| Eppendorf Product Catalog p. 125                 | Order no.    | Price    |
|--|--------------|----------|
| <b>Microplate 96/F</b> , wells white             |              |          |
| PCR clean, gray, 80 plates (5 bags × 16 plates)  | 0030 601 475 | 328,00 € |
| <b>Microplate 96/U</b> , wells white             |              |          |
| PCR clean, gray, 80 plates (5 bags × 16 plates)  | 0030 601 572 | 328,00 € |
| <b>Microplate 96/V</b> , wells white             |              |          |
| PCR clean, gray, 80 plates (5 bags × 16 plates)  | 0030 601 670 | 328,00 € |
| <b>Microplate 384/V</b> , wells white            |              |          |
| PCR clean, gray, 80 plates (5 bags × 16 plates)  | 0030 621 670 | 839,00 € |
| <b>Microplate 96/F</b> , wells black             |              |          |
| PCR clean, white, 80 plates (5 bags × 16 plates) | 0030 601 700 | 328,00 € |
| <b>Microplate 96/U</b> , wells black             |              |          |
| PCR clean, white, 80 plates (5 bags × 16 plates) | 0030 601 807 | 328,00 € |
| <b>Microplate 96/V</b> , wells black             |              |          |
| PCR clean, white, 80 plates (5 bags × 16 plates) | 0030 601 904 | 328,00 € |
| <b>Microplate 384/V</b> , wells black            |              |          |
| PCR clean, white, 80 plates (5 bags × 16 plates) | 0030 621 905 | 839,00 € |



## Eppendorf Deepwell Plates 96 and 384

| Eppendorf Product Catalog p. 126                         | Order no.    | Price      |
|--|--------------|------------|
| <b>Deepwell Plate 96/500 µL</b> , wells clear, 500 µL    |              |            |
| PCR clean, white, 40 plates (5 bags × 8 plates)          | 0030 501 101 | 361,00 €   |
| PCR clean, yellow, 40 plates (5 bags × 8 plates)         | 0030 501 110 | 361,00 €   |
| PCR clean, green, 40 plates (5 bags × 8 plates)          | 0030 501 136 | 361,00 €   |
| PCR clean, blue, 40 plates (5 bags × 8 plates)           | 0030 501 144 | 361,00 €   |
| PCR clean, white, 120 plates (10 bags × 12 plates)       | 0030 505 107 | 982,00 €   |
| sterile, white, 40 plates (5 bags × 8 plates)            | 0030 502 108 | 419,00 €   |
| sterile, yellow, 40 plates (5 bags × 8 plates)           | 0030 502 116 | 419,00 €   |
| sterile, green, 40 plates (5 bags × 8 plates)            | 0030 502 132 | 419,00 €   |
| sterile, blue, 40 plates (5 bags × 8 plates)             | 0030 502 140 | 419,00 €   |
| sterile, white, 120 plates (10 bags × 12 plates)         | 0030 506 103 | 1.129,00 € |
| <b>Deepwell Plate 96/1000 µL</b> , wells clear, 1.000 µL |              |            |
| PCR clean, white, 20 plates (5 bags × 4 plates)          | 0030 501 209 | 166,00 €   |
| PCR clean, yellow, 20 plates (5 bags × 4 plates)         | 0030 501 217 | 166,00 €   |
| PCR clean, green, 20 plates (5 bags × 4 plates)          | 0030 501 233 | 166,00 €   |
| PCR clean, blue, 20 plates (5 bags × 4 plates)           | 0030 501 241 | 166,00 €   |
| PCR clean, white, 80 plates (10 bags × 8 plates)         | 0030 505 204 | 596,00 €   |
| sterile, white, 20 plates (5 bags × 4 plates)            | 0030 502 205 | 190,00 €   |
| sterile, yellow, 20 plates (5 bags × 4 plates)           | 0030 502 213 | 190,00 €   |
| sterile, green, 20 plates (5 bags × 4 plates)            | 0030 502 230 | 190,00 €   |
| sterile, blue, 20 plates (5 bags × 4 plates)             | 0030 502 248 | 190,00 €   |
| sterile, white, 80 plates (10 bags × 8 plates)           | 0030 506 200 | 689,00 €   |
| <b>Deepwell Plate 96/2000 µL</b> , wells clear, 2.000 µL |              |            |
| PCR clean, white, 20 plates (5 bags × 4 plates)          | 0030 501 306 | 193,00 €   |
| PCR clean, yellow, 20 plates (5 bags × 4 plates)         | 0030 501 314 | 193,00 €   |
| PCR clean, green, 20 plates (5 bags × 4 plates)          | 0030 501 330 | 193,00 €   |
| PCR clean, blue, 20 plates (5 bags × 4 plates)           | 0030 501 349 | 193,00 €   |
| PCR clean, white, 80 plates (10 bags × 8 plates)         | 0030 505 301 | 707,00 €   |
| sterile, white, 20 plates (5 bags × 4 plates)            | 0030 502 302 | 227,00 €   |
| sterile, yellow, 20 plates (5 bags × 4 plates)           | 0030 502 310 | 227,00 €   |
| sterile, green, 20 plates (5 bags × 4 plates)            | 0030 502 337 | 227,00 €   |
| sterile, blue, 20 plates (5 bags × 4 plates)             | 0030 502 345 | 227,00 €   |
| sterile, white, 80 plates (10 bags × 8 plates)           | 0030 506 308 | 819,00 €   |
| <b>Deepwell Plate 384/200 µL</b> , wells clear, 200 µL   |              |            |
| PCR clean, white, 40 plates (5 bags × 8 plates)          | 0030 521 102 | 409,00 €   |
| sterile, white, 40 plates (5 bags × 8 plates)            | 0030 522 109 | 471,00 €   |
| PCR clean, white, 120 plates (10 bags × 12 plates)       | 0030 525 108 | 1.107,00 € |
| sterile, white, 120 plates (10 bags × 12 plates)         | 0030 526 104 | 1.271,00 € |

## Protein LoBind® Plates

| Eppendorf Product Catalog p. 127  | Order no.    | Price      |
|---|--------------|------------|
| <b>Deepwell Plate 96/500 µL</b> , Protein LoBind®, wells colorless, 500 µL    |              |            |
| PCR clean, white, 40 plates (5 bags × 8 plates)                               | 0030 504 100 | 729,00 €   |
| PCR clean, yellow, 40 plates (5 bags × 8 plates)                              | 0030 504 119 | 729,00 €   |
| PCR clean, white, 120 plates (10 bags × 12 plates)                            | 0030 508 106 | 1.968,00 € |
| <b>Deepwell Plate 96/1000 µL</b> , Protein LoBind®, wells colorless, 1.000 µL |              |            |
| PCR clean, white, 20 plates (5 bags × 4 plates)                               | 0030 504 208 | 335,00 €   |
| PCR clean, yellow, 20 plates (5 bags × 4 plates)                              | 0030 504 216 | 335,00 €   |
| PCR clean, white, 80 plates (10 bags × 8 plates)                              | 0030 508 203 | 1.199,00 € |
| <b>Deepwell Plate 96/2000 µL</b> , Protein LoBind®, wells colorless, 2.000 µL |              |            |
| PCR clean, white, 20 plates (5 bags × 4 plates)                               | 0030 504 305 | 395,00 €   |
| <b>Deepwell Plate 384/200 µL</b> , Protein LoBind®, wells colorless, 200 µL   |              |            |
| PCR clean, white, 40 plates (5 bags × 8 plates)                               | 0030 524 101 | 819,00 €   |
| PCR clean, white, 120 plates (10 bags × 12 plates)                            | 0030 528 107 | 2.209,00 € |
| <b>Microplate 384/V-PP</b> , Protein LoBind®, wells clear, 140 µL             |              |            |
| PCR clean, white, 80 plates (5 bags × 16 plates)                              | 0030 624 300 | 1.079,00 € |
| PCR clean, white, 240 plates (10 bags × 24 plates)                            | 0030 628 306 | 2.428,00 € |

## Tubes and Plates

### DNA LoBind® Plates

| Eppendorf Product Catalog p. 128  | Order no.    | Price      |
|---|--------------|------------|
| <b>Deepwell Plate 96/500µL</b> , DNA LoBind®, wells clear, 500 µL, LoBind®    |              |            |
| PCR clean, white, 40 plates (5 bags × 8 plates)                               | 0030 503 104 | 402,00 €   |
| PCR clean, blue, 40 plates (5 bags × 8 plates)                                | 0030 503 147 | 402,00 €   |
| PCR clean, white, 120 plates (10 bags × 12 plates)                            | 0030 507 100 | 1.085,00 € |
| <b>Deepwell Plate 96/1000µL</b> , DNA LoBind®, wells clear, 1.000 µL, LoBind® |              |            |
| PCR clean, white, 20 plates (5 bags × 4 plates)                               | 0030 503 201 | 183,00 €   |
| PCR clean, blue, 20 plates (5 bags × 4 plates)                                | 0030 503 244 | 183,00 €   |
| PCR clean, white, 80 plates (10 bags × 8 plates)                              | 0030 507 207 | 657,00 €   |
| <b>Deepwell Plate 384/200µL</b> , DNA LoBind®, wells clear, 200 µL, LoBind®   |              |            |
| PCR clean, white, 40 plates (5 bags × 8 plates)                               | 0030 523 105 | 452,00 €   |
| PCR clean, white, 120 plates (10 bags × 12 plates)                            | 0030 527 100 | 1.217,00 € |
| <b>Microplate 96/V-PP</b> , DNA LoBind®, wells clear, 350 µL, LoBind®         |              |            |
| PCR clean, white, 80 plates (5 bags × 16 plates)                              | 0030 603 303 | 252,00 €   |
| <b>Microplate 384/V-PP</b> , DNA LoBind®, wells clear, 140 µL, LoBind®        |              |            |
| PCR clean, white, 80 plates (5 bags × 16 plates)                              | 0030 623 304 | 689,00 €   |
| PCR clean, white, 240 plates (10 bags × 24 plates)                            | 0030 627 300 | 1.550,00 € |

### Sealing options for Sample Preparation

| Eppendorf Product Catalog p. 129   | Order no.    | Price    |
|--|--------------|----------|
| <b>Eppendorf Storage Film</b> , self-adhesive, PCR clean, 100 pcs. (2 bags × 50 pcs.)                | 0030 127 870 | 147,00 € |
| <b>Eppendorf Storage Foil</b> , self-adhesive, PCR clean, 100 pcs.                                   | 0030 127 889 | 132,00 € |
| <b>Eppendorf Sealing Mat</b> , for DWP 96/2000, PCR clean, 50 pcs. (5 bags × 10 pcs.)                | 0030 127 960 | 166,00 € |
| <b>Eppendorf Sealing Mat</b> , for DWP 96/1000 and DWP 96/500, PCR clean, 50 pcs. (5 bags × 10 pcs.) | 0030 127 978 | 166,00 € |
| <b>Eppendorf Plates® Lid</b> , for MTP and DWP, PCR clean, 80 pcs. (5 bags × 16 pcs.)                | 0030 131 517 | 129,00 € |
| <b>Eppendorf Plates® Lid</b> , for MTP and DWP, sterile, 80 pcs. (5 bags × 16 pcs.)                  | 0030 131 525 | 161,00 € |
| <b>Eppendorf Heat Sealing Film</b> , PCR clean, 100 pcs.   | 0030 127 838 | 193,00 € |
| <b>Eppendorf Heat Sealing Foil</b> , PCR clean, 100 pcs.   | 0030 127 854 | 136,00 € |

### Eppendorf Tube and Cuvette Racks

| Eppendorf Product Catalog p. 129   | Order no.    | Price   |
|--|--------------|---------|
| <b>Eppendorf Tube Rack</b> , 48 positions, for 0.5 mL tubes, polypropylene, numbered positions, autoclavable, 2 pcs.   | 0030 119 800 | 48,50 € |
| <b>Eppendorf Tube Rack</b> , 12 positions, for 5.0 mL and 15 mL tubes, polypropylene, numbered positions, autoclavable, 2 pcs.                                   | 0030 119 827 | 48,50 € |
| <b>Eppendorf Tube Rack</b> , 12 positions, 6 for 5.0 mL and 15 mL tubes and 6 for 25 mL and 50 mL tubes, polypropylene, numbered positions, autoclavable, 2 pcs. | 0030 119 835 | 48,50 € |
| <b>Eppendorf Cuvette Rack</b> , 30 positions, for glass and plastic cuvettes, polypropylene, numbered positions, autoclavable, 2 pcs.                            | 0030 119 851 | 48,50 € |
| <b>Eppendorf Tube Rack</b> , 36 positions, for cryogenic tubes, polypropylene, numbered positions, autoclavable, 2 pcs.  | 0030 119 843 | 48,50 € |
| <b>Eppendorf Tube Rack</b> , 36 positions, for 1.5 mL and 2.0 mL tubes, polypropylene, numbered positions, autoclavable, 2 pcs.                                  | 0030 119 819 | 48,50 € |

### Micropestle

| Eppendorf Product Catalog p. 130   | Order no.    | Price   |
|--|--------------|---------|
| <b>Micropestle</b> , for resuspending pellets in 1.5/2.0 mL vessels, 10 pcs. | 0030 120 973 | 46,40 € |

## Eppendorf Storage Boxes

Eppendorf Product Catalog p. 130

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Storage Box 10 × 10</b> , freezer box, for 100 cryogenic tubes w. int. thread, 3 pcs., height 53 mm, 2 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code  | 0030 140 508 | 29,30 € |
| <b>Storage Box 9 × 9</b> , freezer box, for 81 screw cap (cryog.) tubes 1 – 2 mL, 3 pcs., height 53 mm, 2 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code  | 0030 140 516 | 30,10 € |
| <b>Storage Box 9 × 9</b> , freezer box, for 81 screw cap (cryog.) tubes 3 mL, 2 pcs., height 76 mm, 3 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code      | 0030 140 540 | 20,90 € |
| <b>Storage Box 9 × 9</b> , freezer box, for 81 screw cap (cryog.) tubes 4 – 5 mL, 2 pcs., height 102 mm, 4 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code | 0030 140 567 | 21,50 € |
| <b>Storage Box 8 × 8</b> , freezer box, for 64 tubes 1 – 2 mL, 3 pcs., height 53 mm, 2 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code                     | 0030 140 524 | 30,10 € |
| <b>Storage Box 5 × 5</b> , freezer box, for 25 tubes 5 mL, 4 pcs., height 64 mm, 2.5 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code                       | 0030 140 532 | 42,50 € |
| <b>Storage Box 5 × 5</b> , freezer box, for 25 tubes 5 mL screw cap, 2 pcs., height 76 mm, 3 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code               | 0030 140 613 | 21,50 € |
| <b>Storage Box 5 × 5</b> , freezer box, for 25 tubes 15 mL, 2 pcs., height 127 mm, 5 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code                       | 0030 140 583 | 22,10 € |
| <b>Storage Box 3 × 3</b> , freezer box, for 9 tubes 25 mL, 2 pcs., height 89 mm, 3.5 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code                       | 0030 140 729 | 22,40 € |
| <b>Storage Box 3 × 3</b> , freezer box, for 9 tubes 50 mL and 4 tubes 15 mL, 2 pcs., height 127 mm, 5 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code      | 0030 140 591 | 22,10 € |

## Eppendorf µCuvette® G1.0

Eppendorf Product Catalog p. 132

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Eppendorf µCuvette® G1.0</b> , Microvolume measuring cell for Eppendorf BioPhotometer® and BioSpectrometer® | 6138 000 018 | 2.165,00 € |

## UVette®

Eppendorf Product Catalog p. 132

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>UVette® 220 nm – 1,600 nm</b> , original Eppendorf plastic cuvette, individually packed, (including RNase-free, DNase-free) and protein-free, PCR clean, 50 – 2,000 µL, 80 pcs. | 0030 106 300 | 65,90 €  |
| <b>UVette® starter set</b> , 80 UVettes and 1 universal adapter for light path height 15 mm, convertible to 8.5 mm, 50 – 2,000 µL  | 4099 100 007 | 98,90 €  |
| <b>UVette® routine pack 220 nm – 1,600 nm</b> , reclosable box, Eppendorf Quality, 50 – 2,000 µL, 200 pcs.   | 0030 106 318 | 138,00 € |
| <b>Adapter for UVette®</b> , for photometer/spectrophotometer with light path height   |              |          |
| 8.5 mm   | 4099 001 009 | 75,90 €  |
| 10 mm  | 4099 002 005 | 75,90 €  |
| 15 mm  | 4099 003 001 | 75,90 €  |
| 20 mm  | 4099 005 004 | 75,90 €  |
| <b>Eppendorf Cuvette Rack</b> , 30 positions, for glass and plastic cuvettes, polypropylene, numbered positions, autoclavable, 2 pcs.  | 0030 119 851 | 48,50 €  |

## Eppendorf Vis Cuvettes

Eppendorf Product Catalog p. 133

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Eppendorf macro Vis Cuvette</b> , plastic cuvette for measurements in the Vis range, max. filling volume 4,500 µL                  |              |          |
| 1.000 cuvettes (10 boxes × 100 cuvettes)  | 0030 079 345 | 136,00 € |
| <b>Eppendorf semi-micro Vis Cuvette</b> , plastic cuvette for measurements in the Vis range   |              |          |
| max. filling volume 3,000 µL  | 0030 079 353 | 144,00 € |
| <b>Eppendorf Cuvette Rack</b> , 30 positions, for glass and plastic cuvettes, polypropylene, numbered positions, autoclavable, 2 pcs. | 0030 119 851 | 48,50 €  |

## Centrifuges and Rotors

### MiniSpin®/MiniSpin® plus

Eppendorf Product Catalog p. 142

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Centrifuge MiniSpin®</b> , keypad, non-refrigerated, with Rotor F-45-12-11, 230 V/50 – 60 Hz (EU)                                | 5452 000 010 | 1.088,00 € |
| <b>Centrifuge MiniSpin® plus</b> , keypad, non-refrigerated, with Rotor F-45-12-11, 230 V/50 – 60 Hz (EU)                           | 5453 000 015 | 1.414,00 € |
| <b>Rotor F-45-12-11</b> , angle 45°, 12 places, max. tube diameter 11 mm, incl. rotor lid and rotor nut, MiniSpin®/MiniSpin® plus   | 5452 720 008 | 266,00 €   |
| <b>Rotor F-55-16-5 PCR</b> , angle 55°, 16 places, max. tube diameter 5 mm, incl. rotor lid (aluminum), MiniSpin®/MiniSpin® plus    | 5452 727 007 | 298,00 €   |
| <b>Adapter</b> , for 1 PCR tube (0.2 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.  | 5425 715 005 | 15,70 €    |
| <b>Adapter</b> , for 1 sample tube (0.5 mL, max. Ø 8 mm) or 1 Microtainer® (0.6 mL, max. Ø 8 mm), for all 1.5/2.0 mL rotors, 6 pcs. | 5425 716 001 | 15,70 €    |
| <b>Adapter</b> , for 1 micro test tube (0.4 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.                                     | 5425 717 008 | 15,70 €    |

### Centrifuge 5418 R

Eppendorf Product Catalog p. 143

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Centrifuge 5418 R</b> , rotary knobs, refrigerated, with Rotor FA-45-18-11 and rotor lid, 230 V/50 – 60 Hz (EU)                  | 5401 000 010 | 4.955,00 € |
| <b>Rotor FA-45-18-11</b> , aerosol-tight, 45° angle, 18 places, max. tube Ø 11 mm, incl. rotor lid (aluminum)                       | 5418 707 005 | 594,00 €   |
| <b>Adapter</b> , for 1 PCR tube (0.2 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.  | 5425 715 005 | 15,70 €    |
| <b>Adapter</b> , for 1 sample tube (0.5 mL, max. Ø 8 mm) or 1 Microtainer® (0.6 mL, max. Ø 8 mm), for all 1.5/2.0 mL rotors, 6 pcs. | 5425 716 001 | 15,70 €    |
| <b>Adapter</b> , for 1 micro test tube (0.4 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.                                     | 5425 717 008 | 15,70 €    |

### Centrifuge 5420

Eppendorf Product Catalog p. 144

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Centrifuge 5420</b> , keypad, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU)   | 5420 000 318 | 1.764,00 € |
| <b>Centrifuge 5420</b> , keypad, non-refrigerated, with aerosol-tight Rotor FA-24x2, 230 V/50 – 60 Hz (EU)                          | 5420 000 113 | 2.175,00 € |
| <b>Rotor F-32x0.2-PCR</b> , incl. rotor lid, for 32 × 0.2 mL PCR tubes or 4 × 8 PCR tube strips, for 5420                           | 5495 404 007 | 512,00 €   |
| <b>Rotor FA-24x2</b> , aerosol-tight, incl. rotor lid, aerosol-tight, for 24 × 1.5/2.0 mL tubes, for 5420                           | 5495 400 001 | 647,00 €   |
| <b>Adapter</b> , for 1 PCR tube (0.2 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.  | 5425 715 005 | 15,70 €    |
| <b>Adapter</b> , for 1 sample tube (0.5 mL, max. Ø 8 mm) or 1 Microtainer® (0.6 mL, max. Ø 8 mm), for all 1.5/2.0 mL rotors, 6 pcs. | 5425 716 001 | 15,70 €    |
| <b>Adapter</b> , for 1 micro test tube (0.4 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.                                     | 5425 717 008 | 15,70 €    |

### Centrifuge 5425/5425 R

Eppendorf Product Catalog p. 145

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Centrifuge 5425</b> , keypad, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU)  | 5405 000 310 | 2.222,00 € |
| <b>Centrifuge 5425</b> , rotary knobs, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU)                                      | 5405 000 719 | 2.222,00 € |
| <b>Centrifuge 5425</b> , keypad, non-refrigerated, with Rotor FA-24x2, 230 V/50 – 60 Hz (EU)                                       | 5405 000 115 | 2.762,00 € |
| <b>Centrifuge 5425</b> , rotary knobs, non-refrigerated, with Rotor FA-24x2, 230 V/50 – 60 Hz (EU)                                 | 5405 000 514 | 2.762,00 € |
| <b>Centrifuge 5425</b> , rotary knobs, non-refrigerated, Molecular Biology Solution, with Rotor FA-18x2-Kit, 230 V/50 – 60 Hz (EU) | 5405 000 174 | 2.639,00 € |
| <b>Centrifuge 5425 R</b> , keypad, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)  | 5406 000 313 | 5.224,00 € |
| <b>Centrifuge 5425 R</b> , rotary knobs, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)  | 5406 000 712 | 5.224,00 € |
| <b>Centrifuge 5425 R</b> , keypad, refrigerated, with Rotor FA-24x2, 230 V/50 – 60 Hz (EU)   | 5406 000 119 | 5.617,00 € |
| <b>Centrifuge 5425 R</b> , rotary knobs, refrigerated, with Rotor FA-24x2, 230 V/50 – 60 Hz (EU)                                   | 5406 000 518 | 5.617,00 € |

## Rotors for Centrifuge 5425/5425 R

Eppendorf Product Catalog p. 146

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Rotor FA-24x2</b> , aerosol-tight, incl. rotor lid, for 24 × 1.5/2.0 mL tubes  | 5495 500 006 | 702,00 € |
| <b>Rotor FA-24x2-PTFE</b> , aerosol-tight, incl. rotor lid, PTFE coated, for 24 × 1.5/2.0 mL tubes                                  | 5495 503 005 | 894,00 € |
| <b>Rotor FA-10x5</b> , aerosol-tight, incl. rotor lid, for 10 × 5.0 mL tubes (snap and screw cap)                                   | 5495 505 008 | 894,00 € |
| <b>Rotor FA-18x2-Kit</b> , aerosol-tight, incl. rotor lid, for 18 × 1.5/2.0 mL tubes or spin columns                                | 5495 508 007 | 557,00 € |
| <b>Rotor F-32x0.2-PCR</b> , incl. rotor lid, for 32 × 0.2 mL PCR tubes or 4 × 8 PCR strips  | 5495 510 001 | 511,00 € |
| <b>Rotor S-96x0.2-PCR</b> , for 96 × 0.2 mL PCR tubes or 12 × 8 PCR strips, incl. buckets   | 5495 512 004 | 702,00 € |
| <b>Adapter</b> , for 1 PCR tube (0.2 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.  | 5425 715 005 | 15,70 €  |
| <b>Adapter</b> , for 1 sample tube (0.5 mL, max. Ø 8 mm) or 1 Microtainer® (0.6 mL, max. Ø 8 mm), for all 1.5/2.0 mL rotors, 6 pcs. | 5425 716 001 | 15,70 €  |
| <b>Adapter</b> , for 1 vessel 1.5 – 2.0 mL, for all 5.0 mL rotors, 4 pcs.   | 5820 768 002 | 37,00 €  |
| <b>Adapter</b> , for 1 cryo tube, for all 5.0 mL rotors, 4 pcs.   | 5820 769 009 | 37,00 €  |
| <b>Adapter</b> , for 1 HPLC tube, for all 5.0 mL rotors, 4 pcs.   | 5820 770 007 | 37,00 €  |
| <b>Adapter</b> , for 1 micro test tube (0.4 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.                                     | 5425 717 008 | 15,70 €  |

## Centrifuge 5427 R

New

Eppendorf Product Catalog p. 145

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Centrifuge 5427 R</b> , rotary knobs, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)          | 5429 000 010 | Inquire |
| <b>Centrifuge 5427 R</b> , rotary knobs, refrigerated, with Rotor FA-45-30-11, 230 V/50 – 60 Hz (EU) | 5429 000 310 | Inquire |
| <b>Centrifuge 5427 R</b> , rotary knobs, refrigerated, with Rotor FA-45-48-11, 230 V/50 – 60 Hz (EU) | 5429 000 210 | Inquire |

## Rotors for Centrifuge 5427 R

Eppendorf Product Catalog p. 148

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor FA-45-48-11</b> , aerosol-tight, incl. rotor lid, for 5427R  | 5409 710 001 | 780,00 €   |
| <b>Rotor F-45-48-11</b> , incl. rotor lid, aluminum   | 5409 712 004 | 565,00 €   |
| <b>Rotor FA-45-30-11</b> , aerosol-tight, incl. rotor lid, aluminum   | 5409 706 004 | 668,00 €   |
| <b>Rotor F-45-30-11</b> , incl. rotor lid, aluminum   | 5409 708 007 | 483,00 €   |
| <b>Rotor FA-45-24-11</b> , aerosol-tight, incl. rotor lid, aluminum   | 5409 702 009 | 724,00 €   |
| <b>Rotor FA-45-24-11-Kit</b> , aerosol-tight, incl. rotor lid, aluminum   | 5409 704 001 | 858,00 €   |
| <b>Rotor S-24-11-AT</b> , aerosol-tight, stainless steel  | 5409 715 003 | 1.151,00 € |
| <b>Rotor F-45-48-5-PCR</b> , aluminum   | 5409 714 007 | 585,00 €   |
| <b>Rotor FA-45-12-17</b> , aerosol-tight, incl. rotor lid, aluminum   | 5409 700 006 | 895,00 €   |
| <b>Adapter</b> , for 1 PCR tube (0.2 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.  | 5425 715 005 | 15,70 €    |
| <b>Adapter</b> , for 1 sample tube (0.5 mL, max. Ø 8 mm) or 1 Microtainer® (0.6 mL, max. Ø 8 mm), for all 1.5/2.0 mL rotors, 6 pcs. | 5425 716 001 | 15,70 €    |
| <b>Adapter</b> , for 1 vessel 1.5 – 2.0 mL, for all 5.0 mL rotors, 4 pcs.   | 5820 768 002 | 37,00 €    |
| <b>Adapter</b> , for 1 cryo tube, for all 5.0 mL rotors, 4 pcs.   | 5820 769 009 | 37,00 €    |
| <b>Adapter</b> , for 1 HPLC tube, for all 5.0 mL rotors, 4 pcs.   | 5820 770 007 | 37,00 €    |
| <b>Adapter</b> , for 1 micro test tube (0.4 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.                                     | 5425 717 008 | 15,70 €    |

## Centrifuges and Rotors

### Centrifuge 5430 and Centrifuge 5430 R

Eppendorf Product Catalog p. 149

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Centrifuge 5430</b> , keypad, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU)   | 5427 000 015 | 3.508,00 € |
| <b>Centrifuge 5430</b> , rotary knobs, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU)                                     | 5427 000 610 | 3.508,00 € |
| <b>Centrifuge 5430</b> , keypad, non-refrigerated, with Rotor FA-45-30-11 incl. rotor lid, 230 V/50 – 60 Hz (EU)                  | 5427 000 210 | 3.956,00 € |
| <b>Centrifuge 5430</b> , rotary knobs, non-refrigerated, with Rotor FA-45-30-11 incl. rotor lid, 230 V/50 – 60 Hz (EU)            | 5427 000 415 | 3.956,00 € |
| <b>Centrifuge 5430</b> , rotary knobs, non-refrigerated, Low-Throughput Plate Solution, with Rotor A-2-MTP, 230 V/50 – 60 Hz (EU) | 5427 000 178 | 4.610,00 € |
| <b>Centrifuge 5430 R</b> , keypad, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)   | 5428 000 205 | 7.134,00 € |
| <b>Centrifuge 5430 R</b> , rotary knobs, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)                                       | 5428 000 610 | 7.134,00 € |
| <b>Centrifuge 5430 R</b> , keypad, refrigerated, with Rotor FA-45-30-11 incl. rotor lid, 230 V/50 – 60 Hz (EU)                    | 5428 000 010 | 7.402,00 € |
| <b>Centrifuge 5430 R</b> , rotary knobs, refrigerated, with Rotor FA-45-30-11 incl. rotor lid, 230 V/50 – 60 Hz (EU)              | 5428 000 410 | 7.402,00 € |

### Rotors for Centrifuge 5430/5430 R

Eppendorf Product Catalog p. 152-153

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor FA-45-16-17</b> , aerosol-tight, incl. rotor lid, for 16 x Eppendorf Tubes® 5.0 mL   | 5427 750 002 | 895,00 €   |
| <b>Adapter</b> , for 1 PCR tube (0.2 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.  | 5425 715 005 | 15,70 €    |
| <b>Adapter</b> , for 1 sample tube (0.5 mL, max. Ø 8 mm) or 1 Microtainer® (0.6 mL, max. Ø 8 mm), for all 1.5/2.0 mL rotors, 6 pcs.                 | 5425 716 001 | 15,70 €    |
| <b>Rotor FA-45-48-11</b> , aerosol-tight, incl. rotor lid, for 5430/5430R   | 5427 754 008 | 780,00 €   |
| <b>Rotor F-45-48-11</b> , incl. rotor lid   | 5427 755 004 | 565,00 €   |
| <b>Rotor S-24-11-AT</b> , aerosol-tight, incl. rotor lid, stainless steel   | 5427 757 007 | 1.151,00 € |
| <b>Rotor FA-45-30-11</b> , aerosol-tight, incl. rotor lid, for 5430/5430R   | 5427 753 001 | 668,00 €   |
| <b>Rotor FA-45-24-11-Kit</b> , aerosol-tight, incl. rotor lid   | 5427 752 005 | 858,00 €   |
| <b>Rotor F-45-30-11</b> , incl. rotor lid, PTFE coated  | 5427 712 003 | 496,00 €   |
| <b>Rotor FA-45-24-11-HS</b> , aerosol-tight, incl. rotor lid, PTFE coated   | 5427 710 000 | 693,00 €   |
| <b>Rotor F-35-6-30</b> , incl. rotor lid, with 6 adapters each for 15 mL/50 mL conical tubes  | 5427 716 009 | 1.346,00 € |
| <b>Rotor A-2-MTP</b> , incl. wind shield and rotor lid, with 2 buckets  | 5427 700 005 | 1.346,00 € |
| <b>Rotor F-45-64-5-PCR</b> , incl. rotor lid and adapters   | 5427 714 006 | 819,00 €   |
| <b>Rotor F-45-18-17-Cryo</b> , incl. rotor lid and cryovial adapters  | 5427 705 007 | 489,00 €   |
| <b>Adapter</b> , for 1 cryo tubes (max. Ø 13 mm) or sealable centrifugation tubes (Ø 12.2 mm), max. length 50 mm, for Rotor F-45-18-17-Cryo, 6 pcs. | 5702 752 002 | 72,60 €    |
| <b>Adapter</b> , for 1 vessel 1.5 – 2.0 mL, for all 5.0 mL rotors, 4 pcs.   | 5820 768 002 | 37,00 €    |
| <b>Adapter</b> , for 1 cryo tube, for all 5.0 mL rotors, 4 pcs.   | 5820 769 009 | 37,00 €    |
| <b>Adapter</b> , for 1 HPLC tube, for all 5.0 mL rotors, 4 pcs.   | 5820 770 007 | 37,00 €    |
| <b>Adapter</b> , for 1 micro test tube (0.4 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.   | 5425 717 008 | 15,70 €    |
| <b>Adapter</b> , for 1 HPLC vessel 1.5 mL, for Rotor F-45-18-17-Cryo, 18 pcs.   | 5427 708 006 | 136,00 €   |

### Adapter options 6 x 15/50 mL conical tube rotor

Eppendorf Product Catalog p. 152-153

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Adapter</b> , for 1 tube 65 – 89 mm, for Rotor F-35-6-30, small rotor bore, 2 pcs.                                       | 5427 740 007 | 96,00 €  |
| <b>Adapter</b> , for 1 tube 65 – 89 mm, for Rotor F-35-6-30, large rotor bore, 2 pcs.                                       | 5427 742 000 | 96,00 €  |
| <b>Adapter</b> , for 1 tube 90 – 110 mm, for Rotor F-35-6-30, small rotor bore, 2 pcs.                                      | 5427 741 003 | 96,00 €  |
| <b>Adapter</b> , for 1 tube 90 – 110 mm, for Rotor F-35-6-30, large rotor bore, 2 pcs.                                      | 5427 743 006 | 96,00 €  |
| <b>Adapter</b> , for 1 round-bottom tube and blood collection tube 7 – 15 mL, for Rotor F-35-6-30, small rotor bore, 2 pcs. | 5427 732 004 | 97,10 €  |
| <b>Adapter</b> , for 1 round-bottom tube and blood collection tube 7 – 15 mL, for Rotor F-35-6-30, large rotor bore, 2 pcs. | 5427 734 007 | 97,10 €  |
| <b>Adapter</b> , for 1 round-bottom tube and blood collection tube 9 – 15 mL, for Rotor F-35-6-30, small rotor bore, 2 pcs. | 5427 735 003 | 97,10 €  |
| <b>Adapter</b> , for 1 round-bottom tube and blood collection tube 9 – 15 mL, for Rotor F-35-6-30, large rotor bore, 2 pcs. | 5427 738 002 | 97,10 €  |
| <b>Adapter</b> , for 1 Eppendorf Tubes® 5.0 mL, for Rotor F-35-6-30, small rotor bore, 2 pcs.                               | 5427 746 005 | 96,00 €  |
| <b>Adapter</b> , for 1 Eppendorf Tubes® 5.0 mL, for Rotor F-35-6-30, large rotor bore, 2 pcs.                               | 5427 747 001 | 96,00 €  |
| <b>Adapter</b> , for 1 conical tube 15 mL, for Rotor F-35-6-30, small rotor bore, 2 pcs.                                    | 5427 726 004 | 88,70 €  |
| <b>Adapter</b> , for 1 round-bottom tube 30 mL, for Rotor F-35-6-30, large rotor bore, 2 pcs.                               | 5427 736 000 | 78,70 €  |
| <b>Adapter</b> , for 1 round-bottom tube 50 mL, for Rotor F-35-6-30, large rotor bore, 2 pcs.                               | 5427 737 006 | 78,70 €  |
| <b>Adapter</b> , for 1 conical tube 50 mL, for Rotor F-35-6-30, large rotor bore, 2 pcs.                                    | 5427 727 000 | 52,60 €  |
| <b>Adapter</b> , for 1 Centriprep® 100 – 125 mm, for Rotor F-35-6-30, large rotor bore, 6 pcs.                              | 5427 723 005 | 144,00 € |

## Centrifuge 5702/5702 R/5702 RH

Eppendorf Product Catalog p. 152-153

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Centrifuge 5702</b> , rotary knobs, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU)  | 5702 000 010 | 1.940,00 € |
| <b>Centrifuge 5702</b> , rotary knobs, non-refrigerated, with Rotor A-4-38 incl. adapters for 15/50 mL conical tubes, 2 sets of 2 adapters, 230 V/50 – 60 Hz (EU)  | 5702 000 320 | 3.189,00 € |
| <b>Centrifuge 5702 R</b> , rotary knobs, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)  | 5703 000 010 | 4.601,00 € |
| <b>Centrifuge 5702 R</b> , rotary knobs, refrigerated, with Rotor A-4-38 incl. adapters for 13/16 mm blood collection tubes, 2 sets of 2 adapters, 230 V/50 – 60 Hz (EU)   | 5703 000 320 | 5.725,00 € |
| <b>Centrifuge 5702 R</b> , rotary knobs, refrigerated, with Rotor A-4-38 incl. adapters for 15/50 mL conical tubes, 2 sets of 2 adapters, 230 V/50 – 60 Hz (EU)  | 5703 000 425 | 5.584,00 € |
| <b>Centrifuge 5702 RH</b> , rotary knobs, refrigerated/heated, without rotor, 230 V/50 – 60 Hz (EU)  | 5704 000 010 | 5.199,00 € |
| <b>Centrifuge 5702</b> , rotary knobs, non-refrigerated, with Rotor A-4-38 incl. adapters for round-bottom tubes 4 – 10 mL and adapters for round-bottom tubes 2.6 – 7 mL, 2 sets of 2 adapters, 230 V/50 – 60 Hz (EU) | 5702 000 515 | 3.354,00 € |

## Rotor A-4-38 with round buckets

Eppendorf Product Catalog p. 152-153

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor A-4-38 with round buckets</b> , 8 positions, max. Ø 38 mm, incl. 4 round buckets 100 mL            | 5702 720 003 | 1.068,00 € |
| <b>16 x 125 mm bucket</b> , for 16 x 125 mm round bottom tubes, non aerosol-tight, for rotor A-4-38, 4 pcs. | 5702 764 000 | 452,00 €   |
| <b>Rotor A-4-38</b> , 8 positions, max. Ø 38 mm, without buckets  | 5702 740 004 | 583,00 €   |
| <b>Round bucket 100 mL</b> , for Rotor A-4-38, 4 pcs.   | 5702 722 006 | 603,00 €   |
| <b>Caps for 100 mL round buckets</b> , aerosol-tight, 2 pcs.  | 5702 721 000 | 64,30 €    |

## Adapter for Rotor A-4-38 with round buckets

Eppendorf Product Catalog p. 152-153

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Adapter</b> , for 1 round-bottom tubes 85/100 mL, 38 mm x 106 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs.               | 5702 718 009 | 100,00 € |
| <b>Adapter</b> , for 1 conical tubes 15 mL, 17.2 mm x 121 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs.                      | 5702 732 001 | 100,00 € |
| <b>Adapter</b> , for 1 conical tube 25 mL, 25 mm x 100 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs.                         | 5702 717 002 | 100,00 € |
| <b>Adapter</b> , for 1 conical tube 50 mL, 30 mm x 115 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs.                         | 5702 734 004 | 100,00 € |
| <b>Adapter</b> , for 2 conical tubes 15 mL, 17.2 mm x 121 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs.                      | 5702 723 002 | 180,00 € |
| <b>Adapter</b> , for 4 tubes 1.5 – 2.0 mL, 11 mm x 43 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs.                          | 5702 745 006 | 186,00 € |
| <b>Adapter</b> , for 4 round-bottom tubes 4 – 10 mL, 16 mm x 100 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs.               | 5702 735 000 | 180,00 € |
| <b>Adapter</b> , for 4 round-bottom tubes 9 – 15 mL, 17.5 mm x 100 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs.             | 5702 724 009 | 180,00 € |
| <b>Adapter</b> , for 4 round-bottom tubes 2.6 – 7 mL, 13.5 mm x 100 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs.            | 5702 719 005 | 180,00 € |
| <b>Adapter</b> , for 1 Eppendorf Tubes® 5.0 mL, 17 mm x 60 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs.                     | 5702 733 008 | 90,00 €  |
| <b>Adapter</b> , for 5 round-bottom tubes 1 – 1.4 mL, 8.5 mm x 100 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs.             | 5702 736 007 | 180,00 € |
| <b>Adapter</b> , for 5 round-bottom tubes 2 – 7 mL, 12.5 mm x 100 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs.              | 5702 737 003 | 180,00 € |
| <b>Adapter</b> , for 1 conical tube 25 mL screw cap <sup>1)</sup> , 30 mm x 83 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs. | 5702 734 004 | 100,00 € |
| <b>Adapter</b> , for 1 conical tube 25 mL snap cap <sup>1)</sup> , 30 mm x 83 mm, for 100 mL round bucket in Rotor A-4-38, 2 pcs.  | 5702 734 004 | 100,00 € |

<sup>1)</sup> Please order adapter for 25 mL conical tubes snap cap (5820 733 004) or adapter for 25 mL conical tubes screw cap (5820 734 000).



## Centrifuges and Rotors

### Rotor A-4-38 with rectangular buckets

Eppendorf Product Catalog p. 154-155

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Rotor A-4-38</b> , 8 positions, max. Ø 38 mm, without buckets | 5702 740 004 | 583,00 € |
| <b>Rectangular bucket 90 mL</b> , for Rotor A-4-38, 4 pcs.       | 5702 709 000 | 363,00 € |

### Adapter for Rotor A-4-38 with rectangular buckets

Eppendorf Product Catalog p. 154-155

|   | Order no.    | Price   |
|---|--------------|---------|
| <b>Adapter</b> , for 10 round-bottom tubes 5 – 7 mL, for 90 mL rectangular bucket in Rotor A-4-38, 13 mm × 100 mm, 2 pcs. | 5702 710 008 | 91,50 € |
| <b>Adapter</b> , for 2 round-bottom tubes 25 mL, for 90 mL rectangular bucket in Rotor A-4-38, 25 mm × 100 mm, 2 pcs.     | 5702 716 006 | 91,50 € |
| <b>Adapter</b> , for 4 round-bottom tubes 20 mL, for 90 mL rectangular bucket in Rotor A-4-38, 22 mm × 100 mm, 2 pcs.     | 5702 713 007 | 91,50 € |
| <b>Adapter</b> , for 6 round-bottom tubes 15 mL, for 90 mL rectangular bucket in Rotor A-4-38, 17.5 mm × 100 mm, 2 pcs.   | 5702 712 000 | 64,30 € |
| <b>Adapter</b> , for 8 round-bottom tubes 9 mL, for 90 mL rectangular bucket in Rotor A-4-38, 14.5 mm × 100 mm, 2 pcs.    | 5702 711 004 | 91,50 € |
| <b>Adapter</b> , for 9 round-bottom tubes 7.5 mm, for 90 mL rectangular bucket in Rotor A-4-38, Ø12 mm, 2 pcs.            | 5702 763 004 | 90,00 € |

### Rotors for Centrifuge 5702/5702 R/5702 RH

Eppendorf Product Catalog p. 154-155

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Rotor A-8-17</b> , 8 × 15 mL, without adapter, without rubber mats   | 5702 700 002 | 288,00 € |
| <b>Rotor F-35-30-17</b> , incl. 30 steel sleeves for 15 mL vessels, with 20 adapters for conical tubes and 30 rubber mats                           | 5702 704 008 | 706,00 € |
| <b>Rotor F-35-30-17</b> , incl. 10 steel sleeves for 15 mL vessels, with 10 adapters for conical tubes and 10 rubber mats                           | 5702 705 004 | 470,00 € |
| <b>Rotor F-45-24-11</b> , 45° angle, 24 places, max. Ø 11 mm, without lid   | 5702 746 002 | 428,00 € |
| <b>Rotor F-45-18-17-Cryo</b> , without rotor lid, without adapters  | 5702 747 009 | 492,00 € |
| <b>Adapter</b> , for 1 PCR tube (0.2 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.  | 5425 715 005 | 15,70 €  |
| <b>Adapter</b> , for 1 sample tube (0.5 mL, max. Ø 8 mm) or 1 Microtainer® (0.6 mL, max. Ø 8 mm), for all 1.5/2.0 mL rotors, 6 pcs.                 | 5425 716 001 | 15,70 €  |
| <b>Adapter</b> , for 1 cryo tubes (max. Ø 13 mm) or sealable centrifugation tubes (Ø 12.2 mm), max. length 50 mm, for Rotor F-45-18-17-Cryo, 6 pcs. | 5702 752 002 | 72,60 €  |
| <b>Adapter</b> , for 1 micro test tube (0.4 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.   | 5425 717 008 | 15,70 €  |
| <b>Adapter</b> , for 1 HPLC vessel 1.5 mL, for Rotor F-45-18-17-Cryo, 18 pcs.   | 5427 708 006 | 136,00 € |

### Centrifuge 5804/5804 R and 5810/5810 R

Eppendorf Product Catalog p. 156-157

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Centrifuge 5804</b> , keypad, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU) | 5804 000 010 | 4.089,00 € |
| <b>Centrifuge 5804 R</b> , keypad, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)   | 5805 000 010 | 7.266,00 € |
| <b>Centrifuge 5810</b> , keypad, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU) | 5810 000 010 | 5.369,00 € |
| <b>Centrifuge 5810 R</b> , keypad, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)   | 5811 000 015 | 8.999,00 € |

### Centrifuge Packages

Eppendorf Product Catalog p. 156-157

|   | Order no.    | Price       |
|---|--------------|-------------|
| <b>Centrifuge 5804</b> , keypad, non-refrigerated, with Rotor A-4-44 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU) | 5804 000 320 | 5.466,00 €  |
| <b>Centrifuge 5804 R</b> , keypad, refrigerated, with Rotor A-4-44 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU)   | 5805 000 320 | 8.327,00 €  |
| <b>Centrifuge 5804</b> , keypad, non-refrigerated, with Rotor S-4-72 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU) | 5804 000 520 | 6.569,00 €  |
| <b>Centrifuge 5804 R</b> , keypad, refrigerated, with Rotor S-4-72 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU)   | 5805 000 620 | 9.151,00 €  |
| <b>Centrifuge 5810</b> , keypad, non-refrigerated, with Rotor A-4-62 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU) | 5810 000 320 | 7.161,00 €  |
| <b>Centrifuge 5810 R</b> , keypad, refrigerated, with Rotor A-4-62 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU)   | 5811 000 325 | 10.430,00 € |
| <b>Centrifuge 5810</b> , keypad, non-refrigerated, with Rotor A-4-81 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU) | 5810 000 420 | 8.051,00 €  |



## Centrifuge Packages

Eppendorf Product Catalog p. 156-157

|   | Order no.    | Price       |
|---|--------------|-------------|
| <b>Centrifuge 5810 R</b> , keypad, refrigerated, with Rotor A-4-81 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU)   | 5811 000 420 | 11.320,00 € |
| <b>Centrifuge 5810</b> , keypad, non-refrigerated, Classic Cell Culture Solution, with Rotor S-4-104 including 4 adapters for 5/15 mL tubes and 4 adapters for 50 mL conical tubes, 230 V/50 – 60 Hz (EU) | 5810 000 525 | 8.622,00 €  |
| <b>Centrifuge 5810 R</b> , keypad, refrigerated, with Rotor S-4-104 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU)  | 5811 000 620 | 11.890,00 € |
| <b>Centrifuge 5810</b> , keypad, non-refrigerated, with Rotor S-4-104 incl. adapters for 13/16 mm round-bottom tubes, 230 V/50 – 60 Hz (EU)   | 5810 000 620 | 8.622,00 €  |
| <b>Centrifuge 5810 R</b> , keypad, refrigerated, with Rotor S-4-104 incl. adapters for 13/16 mm round-bottom tubes, 230 V/50 – 60 Hz (EU)   | 5811 000 820 | 11.890,00 € |

## Rotors for Centrifuge 5804/5804 R/5810/5810 R

Eppendorf Product Catalog p. 160-161

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor FA-45-20-17</b> , incl. rotor lid, for 20 Eppendorf Tubes® 5.0 mL                      | 5820 765 003 | 1.433,00 € |
| <b>Adapter</b> , for 1 vessel 1.5 – 2.0 mL, for all 5.0 mL rotors, 4 pcs.                       | 5820 768 002 | 37,00 €    |
| <b>Adapter</b> , for 1 cryo tube, for all 5.0 mL rotors, 4 pcs.                                 | 5820 769 009 | 37,00 €    |
| <b>Adapter</b> , for 1 HPLC tube, for all 5.0 mL rotors, 4 pcs.                                 | 5820 770 007 | 37,00 €    |
| <b>Rotor FA-45-48-11</b> , aerosol-tight, incl. rotor lid, for 48 × 1.5/2.0 mL tubes            | 5820 760 001 | 1.433,00 € |
| <b>Rotor F-45-30-11</b> , incl. rotor lid, for 30 × 1.5/2.0 mL reaction vessels                 | 5804 715 004 | 696,00 €   |
| <b>Rotor FA-45-30-11</b> , aerosol-tight, incl. rotor lid, for 30 × 1.5/2.0 mL reaction vessels | 5804 726 006 | 863,00 €   |

## Rotors for Centrifuge 5804/5804 R/5810/5810 R

Eppendorf Product Catalog p. 160-161

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor F-35-48-17</b> , for 24 × 15 mL conical tubes, incl. 24 steel sleeves and adapters   | 5820 771 003 | 706,00 €   |
| <b>Rotor F-35-48-17</b> , for 40 × 15 mL conical tubes, incl. 48 steel sleeves and adapters   | 5820 772 000 | 988,00 €   |
| <b>Rotor T-60-11</b> , incl. rotor lid, for 1.5/2.0 mL tubes, without adapter                 | 5804 730 003 | 1.582,00 € |
| <b>Adapter</b> , for 10 reaction vessels 1.5/2.0 mL, for Rotor T-60-11, 6 pcs.                | 5804 731 000 | 220,00 €   |
| <b>Adapter</b> , for 20 reaction vessels 0.4 mL, for Rotor T-60-11, 6 pcs.                    | 5804 732 006 | 220,00 €   |
| <b>Rotor F-45-48-PCR</b> , for 6 × 8, 6 × 5 PCR tube strips, 48 × 0.2 mL individual PCR tubes | 5804 735 005 | 651,00 €   |

## Rotors for Centrifuge 5804/5804 R/5810/5810 R

Eppendorf Product Catalog p. 160-161

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Adapter</b> , for 1 cell culture flask, for Rotors A-4-81-MTP/Flex   | 5825 719 000 | 190,00 € |
| <b>Adapter</b> , for 1 PCR plate (96 wells), for 96-well PCR plates, used in A-4-81-MTP/Flex, A-4-62-MTP, A-2-DWP-AT and A-2-DWP, 2 pcs.                          | 5825 711 009 | 159,00 € |
| <b>Adapter</b> , for 1 PCR plate (384 wells), for 384-well PCR plates, used in A-4-81-MTP/Flex, A-4-62-MTP, A-2-DWP-AT and A-2-DWP, 2 pcs.                        | 5825 713 001 | 126,00 € |
| <b>IsoRack adapter</b> , for 24 vessels in the IsoRack 0.5 mL, for Rotor A-4-81 (MTP), for 24 × 0.5 mL tubes in the IsoRack, 2 pcs.                               | 5825 708 008 | 154,00 € |
| <b>IsoRack adapter</b> , for 24 vessels in the IsoRack 1.5/2.0 mL, for Rotor A-4-81 (MTP), for Rotor A-2-DWP-AT, for 24 × 1.5/2.0 mL tubes in the IsoRack, 2 pcs. | 5825 709 004 | 154,00 € |
| <b>CombiSlide Adapter</b> , for 12 slides Combi Slide, for Rotors A-2-MTP and A-4-81 (MTP), 2 pcs.  | 5825 706 005 | 314,00 € |

## Rotor A-2-DWP

Eppendorf Product Catalog p. 161

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor A-2-DWP</b> , for Deepwell plates, incl. 2 buckets | 5804 740 009 | 1.782,00 € |

## Rotor A-2-DWP-AT

Eppendorf Product Catalog p. 161

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Rotor A-2-DWP-AT</b> , incl. 2 buckets, 2 aerosol-tight caps and 2 plate carriers | 5820 710 004 | 2.433,00 € |

## Centrifuges and Rotors

### Rotor S-4-104

Eppendorf Product Catalog p. 162-163

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor S-4-104</b>  |              |            |
| incl. 4 x 750 mL round buckets  | 5820 740 000 | 3.074,00 € |
| aerosol-tight capable, incl. 4 plate buckets  | 5820 754 001 | 3.074,00 € |
| without buckets   | 5820 755 008 | 1.435,00 € |
| <b>Inlay for microfluidic card</b> , for Rotor S-4-104, 4 pcs.                                | 5820 751 002 | 76,50 €    |
| <b>Plate carrier</b> , Rotoren S-4-104, S-4x750, 2 pcs.                                       | 5820 756 004 | 303,00 €   |
| <b>Plate bucket</b> , for Rotor S-4-104, aerosol-tight capable, incl. plate carrier, 2 pcs.   | 5820 744 006 | 791,00 €   |
| <b>Plate bucket</b> , for Rotor S-4-104, aerosol-tight capable, incl. plate carrier, 4 pcs.   | 5820 743 000 | 1.460,00 € |
| <b>Plate bucket (open)</b> , for Rotor S-4-104, incl. plate carrier, 4 pcs.                   | 5820 757 000 | 1.694,00 € |
| <b>Round bucket 750 mL</b> , for Rotor S-4-104, 2 pcs.  | 5820 742 003 | 887,00 €   |
| <b>Round bucket 750 mL</b> , for Rotor S-4-104, aerosol-tight capable, 4 pcs.                 | 5820 741 007 | 1.694,00 € |
| <b>Caps</b> , rotors S-4-104, S-4x750, plate bucket, aerosol-tight, 2 pcs.                    | 5820 748 001 | 213,00 €   |
| <b>Caps for 750 mL/1000 mL</b> , for rotors S-4-104, S-4x750, S-4x1000, aerosol-tight, 2 pcs. | 5820 747 005 | 195,00 €   |
| <b>Wide-neck bottle 750 mL</b> , for Rotors S-4-104, S-4x750, S-4x1000, 2 pcs.                | 5820 708 000 | 198,00 €   |

### Adapter for Rotor S-4-104

Eppendorf Product Catalog p. 162-163

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Adapter</b> , for 7 conical tubes 25 mL screw cap <sup>1)</sup> , for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.     | 5825 733 002 | 299,00 € |
| <b>Adapter</b> , for 7 conical tubes 25 mL snap cap <sup>1)</sup> , for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.      | 5825 733 002 | 299,00 € |
| <b>Adapter</b> , for 7 conical tubes Conical SnapTec® 50, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                | 5825 733 002 | 299,00 € |
| <b>Adapter</b> , for 14 tubes 14 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                                     | 5825 748 000 | 299,00 € |
| <b>Adapter</b> , for tubes 30 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.  | 5825 755 006 | 295,00 € |
| <b>Adapter</b> , for 27 round-bottom tubes 5 mL FACS, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                    | 5825 747 003 | 299,00 € |
| <b>Adapter</b> , for bottle 500 mL Corning®, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                             | 5825 745 000 | 299,00 € |
| <b>Adapter</b> , for 7 conical tubes 50 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                              | 5825 733 002 | 299,00 € |
| <b>Adapter</b> , for 14 Eppendorf Tubes® 5.0 mL and 15 mL conical tubes, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs. | 5825 734 009 | 299,00 € |
| <b>Adapter</b> , for 20 round-bottom tubes 5.5 – 12 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                  | 5825 736 001 | 299,00 € |
| <b>Adapter</b> , for 20 round-bottom tubes 9 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                         | 5825 743 008 | 299,00 € |
| <b>Adapter</b> , for 23 round-bottom tubes 2.6 – 8 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                   | 5825 738 004 | 299,00 € |
| <b>Adapter</b> , for 1 bottle 250 mL flat, 175 – 225/ 250 mL conical, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.    | 5825 741 005 | 299,00 € |
| <b>Adapter</b> , for 50 reaction vessels 1.5/2.0 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                     | 5825 740 009 | 299,00 € |
| <b>Adapter</b> , for 5 tubes 50 mL conical, skirted, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                     | 5825 732 006 | 299,00 € |
| <b>Adapter</b> , for 1 Eppendorf bottle 750 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                          | 5825 744 004 | 159,00 € |

<sup>1)</sup> Please order adapter for 25 mL conical tubes snap cap (5820 733 004) or adapter for 25 mL conical tubes screw cap (5820 734 000).

## Rotor A-4-81

Eppendorf Product Catalog p. 164-165

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor A-4-81</b> , for 500 mL rectangular buckets or MTP/Flex buckets, incl. 4 x 500 mL rectangular buckets                            | 5810 718 007 | 2.841,00 € |
| <b>Rotor A-4-81</b> , without buckets   | 5810 743 001 | 1.336,00 € |
| <b>Rotor A-4-81-MTP/Flex</b> , Swing-bucket rotor, incl. 4 MTP/Flex buckets   | 5810 725 003 | 2.841,00 € |
| <b>Bucket for 7 x 50 mL conical tubes</b> , for Rotor A-4-81, 4 pcs.  | 5825 730 003 | 1.721,00 € |
| <b>Rectangular bucket 500 mL</b> , for Rotor A-4-81, 4 pcs.   | 5810 730 007 | 1.840,00 € |
| <b>Cap</b> , aerosol-tight, for 500 mL rectangular bucket, 2 pcs.   | 5810 724 007 | 180,00 €   |
| <b>Wide-neck bottle 500 mL</b> , for Rotors A-4-81, S-4x500, rectangular, 2 pcs.  | 5820 707 003 | 200,00 €   |
| <b>Wide-neck bottle 400 mL</b> , for Rotors A-4-81, S-4x500, 400 mL, lid blue, 2 pcs.   | 5810 729 009 | 200,00 €   |
| <b>MTP/Flex buckets</b> , for Rotors A-4-81, S-4x500, for use with IsoRack and cell culture flask adapters as well as MTP and DWP, 2 pcs. | 5810 742 005 | 791,00 €   |
| <b>MTP/Flex buckets</b> , for Rotors A-4-81, S-4x500, for use with IsoRack and cell culture flask adapters as well as MTP and DWP, 4 pcs. | 5810 741 009 | 1.460,00 € |

## Adapter for Rotor A-4-81

Eppendorf Product Catalog p. 164-165

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Bucket for 7 x 50 mL conical tubes</b> , for Rotor A-4-81, 4 pcs.   | 5825 730 003 | 1.721,00 € |
| <b>MTP/Flex buckets</b> , for Rotors A-4-81, S-4x500, for use with IsoRack and cell culture flask adapters as well as MTP and DWP, 2 pcs.  | 5810 742 005 | 791,00 €   |
| <b>Adapter</b> , for 5 conical tubes 25 mL screw cap, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.  | 5810 723 000 | 196,00 €   |
| <b>Adapter</b> , for 5 conical tubes 25 mL snap cap, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.   | 5810 723 000 | 196,00 €   |
| <b>Adapter</b> , for 5 conical tubes 50 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.  | 5810 723 000 | 196,00 €   |
| <b>Bucket for 7 x 25 mL conical tubes, screw cap<sup>1)</sup></b> , for Rotor A-4-81, 4 pcs.   | 5825 730 003 | 1.721,00 € |
| <b>Bucket for 7 x 25 mL snap cap conical tubes<sup>1)</sup></b> , for Rotor A-4-81, 4 pcs.   | 5825 730 003 | 1.721,00 € |
| <b>Bucket for 7 x 25 mL snap cap conical tubes<sup>1)</sup></b> , for Rotor A-4-81, 4 pcs.   | 5825 730 003 | 1.721,00 € |
| <b>Adapter</b> , for 16 round-bottom tubes 3 – 15 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.  | 5810 748 003 | 196,00 €   |
| <b>Adapter</b> , for 16 round-bottom tubes 7 – 17 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.  | 5810 721 008 | 196,00 €   |
| <b>Adapter</b> , for 18 round-bottom tubes 2.6 – 7 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.   | 5810 747 007 | 196,00 €   |
| <b>Adapter</b> , for 1 bottle 180 – 250 mL flat/ 175 – 225 mL conical <sup>2)</sup> , for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, Please note that these conical bottles need to be used together with a conical form insert provided by the manufacturer, 2 pcs. | 5825 722 000 | 197,00 €   |
| <b>Adapter</b> , for 1 bottle 400 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.  | 5810 728 002 | 164,00 €   |
| <b>Adapter</b> , for 20 reaction vessels 1.5/2.0 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.   | 5810 745 004 | 196,00 €   |
| <b>Adapter</b> , for 20 round-bottom tubes 1.2 – 5 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.   | 5810 746 000 | 196,00 €   |
| <b>Adapter</b> , for 25 round-bottom tubes 2.6 – 5 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.   | 5810 720 001 | 196,00 €   |
| <b>Adapter</b> , for 12 Eppendorf Tubes® 5.0 mL and 15 mL conical tubes, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.   | 5810 722 004 | 196,00 €   |
| <b>Adapter</b> , for 5 conical tubes 50 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.  | 5810 723 000 | 196,00 €   |

<sup>1)</sup> Please order adapter for 25 mL conical tubes snap cap (5820 733 004) or adapter for 25 mL conical tubes screw cap (5820 734 000). <sup>2)</sup> Please note that these conical bottles need to be used together with a conical form insert provided by the manufacturer.

## Centrifuges and Rotors

### Rotor A-4-62

Eppendorf Product Catalog p. 166-167

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Rotor A-4-62</b> , incl. 4 x 250 mL rectangular buckets         | 5810 709 008 | 1.935,00 € |
| <b>Rotor A-4-62-MTP</b> , with 4 MTP buckets                       | 5810 711 002 | 1.853,00 € |
| <b>MTP bucket for Rotor A-4-62</b> , for 4 MTP or 1 DWP, 4 pcs.    | 5810 702 003 | 1.156,00 € |
| <b>Caps</b> , for rectangular bucket 250 mL, aerosol-tight, 2 pcs. | 5810 710 006 | 135,00 €   |
| <b>Rectangular bucket 250 mL</b> , for Rotor A-4-62, 4 pieces      | 5810 716 004 | 1.011,00 € |

### Adapter for Rotor A-4-62

Eppendorf Product Catalog p. 166-167

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Adapter</b> , for 4 conical tubes 50 mL, for 250 mL rectangular bucket in Rotor A-4-62, flat bottom shape, please order conical inlay 5804728009 when using conical tubes   | 5810 763 002 | 162,00 €   |
| <b>Adapter</b> , for 3 conical tubes 25 mL conical, for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.  | 5810 758 009 | 184,00 €   |
| <b>Adapter</b> , for 3 conical tubes 25 mL conical, for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.  | 5810 758 009 | 184,00 €   |
| <b>Adapter</b> , for 3 conical tubes 50 mL, for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.  | 5810 758 009 | 184,00 €   |
| <b>Adapter</b> , for 12 round-bottom tubes 3 – 15 mL, for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.  | 5810 753 007 | 171,00 €   |
| <b>Adapter</b> , for 12 round-bottom tubes 7 – 17 mL, for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.  | 5810 754 003 | 171,00 €   |
| <b>Adapter</b> , for 15 round-bottom tubes 2.6 – 7 mL, for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.   | 5810 752 000 | 171,00 €   |
| <b>Adapter</b> , for 16 reaction vessels 1.5 – 2.0 mL, for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.   | 5810 751 004 | 162,00 €   |
| <b>Adapter</b> , for 1 vessel 80 – 120 mL, for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.   | 5810 761 000 | 162,00 €   |
| <b>Adapter</b> , for 1 bottle 200 – 250 mL flat/ 175 – 225 mL conical, for 250 mL rectangular bucket in Rotor A-4-62, Please note that these conical bottles need to be used together with a conical form insert provided by the manufacturer., 2 pcs. | 5810 770 009 | 189,00 €   |
| <b>Adapter</b> , for 25 round-bottom tubes 1.2 – 5 mL, for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.   | 5810 750 008 | 162,00 €   |
| <b>Adapter</b> , for 2 round-bottom tubes 50 – 75 mL, for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.  | 5810 760 003 | 162,00 €   |
| <b>Adapter</b> , for 3 conical tubes 50 mL, for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.  | 5810 758 009 | 184,00 €   |
| <b>Adapter</b> , for 4 round-bottom tubes 18 – 30 mL, for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.  | 5810 757 002 | 162,00 €   |
| <b>Adapter</b> , for 4 tubes 30 – 50 mL round-bottom/ 50 mL conical, skirted <sup>1)</sup> , for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.   | 5810 759 005 | 162,00 €   |
| <b>Adapter</b> , for 4 conical tubes 50 mL, for 250 mL rectangular bucket in Rotor A-4-62, flat bottom shape, please order conical inlay 5804728009 when using conical tubes   | 5810 763 002 | 162,00 €   |
| <b>Adapter</b> , for 8 round-bottom tubes 7 – 18 mL, for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.   | 5810 756 006 | 162,00 €   |
| <b>Adapter</b> , for 9 Eppendorf Tubes® 5.0 mL and 15 mL conical tubes, for 250 mL rectangular bucket in Rotor A-4-62, 2 pcs.  | 5810 755 000 | 184,00 €   |
| <b>MTP bucket for Rotor A-4-62</b> , for 4 MTP or 1 DWP, 4 pcs.  | 5810 702 003 | 1.156,00 € |

<sup>1)</sup> Please order adapter for 50 mL conical tubes skirted 5804 737 008.

### Rotor S-4-72

Eppendorf Product Catalog p. 168-169

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Rotor S-4-72</b> , incl. 4 x 250 mL round buckets | 5804 746 007 | 2.167,00 € |

### Adapter for Rotor S-4-72

Eppendorf Product Catalog p. 168-169

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Adapter</b> , for 4 conical tubes 25 mL screw cap <sup>1)</sup> , for Rotor S-4-72, 2 pcs.      | 5804 784 006 | 194,00 € |
| <b>Adapter</b> , for 4 conical tubes 25 mL snap cap <sup>1)</sup> , for Rotor S-4-72, 2 pcs.       | 5804 784 006 | 194,00 € |
| <b>Adapter</b> , for 4 conical tubes 50 mL, for Rotor S-4-72, 2 pcs.                               | 5804 784 006 | 194,00 € |
| <b>Adapter</b> , for 13 round-bottom tubes 5.5 – 12 mL, for Rotor S-4-72, 2 pcs.                   | 5804 791 002 | 194,00 € |
| <b>Adapter</b> , for 14 round-bottom tubes 2.6 – 8 mL, for Rotor S-4-72, 2 pcs.                    | 5804 789 008 | 194,00 € |
| <b>Adapter</b> , for 26 tubes 1.5/2.0 mL, for Rotor S-4-72, 2 pcs.                                 | 5804 794 001 | 194,00 € |
| <b>Adapter</b> , for 2 conical tubes, skirted 50 mL, for Rotor S-4-72, 2 pcs.                      | 5804 785 002 | 194,00 € |
| <b>Adapter</b> , for 4 conical tubes 50 mL, for Rotor S-4-72, 2 pcs.                               | 5804 784 006 | 194,00 € |
| <b>Adapter</b> , for 8 conical tubes 15 mL, for Rotor S-4-72, 2 pcs.                               | 5804 783 000 | 194,00 € |
| <b>Adapter</b> , for 12 round-bottom tubes 9 mL, for Rotor S-4-72, 2 pcs.                          | 5804 792 009 | 194,00 € |
| <b>Adapter</b> , for 8 Eppendorf Tubes® 5.0 mL, for Rotor S-4-72, 2 pcs.                           | 5804 793 005 | 194,00 € |
| <b>Adapter</b> , for 1 bottle 250 mL flat/ 175 mL Falcon® <sup>2)</sup> , for Rotor S-4-72, 2 pcs. | 5804 787 005 | 194,00 € |

<sup>1)</sup> Please order adapter for 25 mL conical tubes snap cap (5820 733 004) or adapter for 25 mL conical tubes screw cap (5820 734 000). <sup>2)</sup> Please note that these conical bottles need to be used together with a conical form insert provided by the manufacturer.

## Rotor A-4-44

Eppendorf Product Catalog p. 170-171

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Rotor A-4-44</b> , incl. 4 × 100 mL rectangular buckets                         | 5804 709 004 | 1.401,00 € |
| <b>Bucket for 2 × 50 mL conical tubes</b> , for Rotor A-4-44, 4 pcs.               | 5804 706 005 | 606,00 €   |
| <b>Rotor A-4-44</b> , incl. 4 buckets for 2 × 50 mL conical tubes                  | 5804 719 000 | 1.690,00 € |
| <b>Rectangular bucket 100 mL</b> , for Rotor A-4-44, 4 pcs.                        | 5804 741 005 | 864,00 €   |
| <b>Caps</b> , for 100 mL rectangular bucket, aerosol-tight, 2 pcs.                 | 5804 712 005 | 117,00 €   |
| <b>Adapter Insert</b> , insert for 50 mL conical tubes with skirted bottom, 8 pcs. | 5804 737 008 | 102,00 €   |

## Adapter for Rotor A-4-44

Eppendorf Product Catalog p. 170-171

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Bucket for 2 × 50 mL conical tubes</b> , for Rotor A-4-44, 4 pcs.   | 5804 706 005 | 606,00 € |
| <b>Bucket for 2 × 50 mL conical tubes</b> , for Rotor A-4-44, 4 pcs.   | 5804 706 005 | 606,00 € |
| <b>Adapter</b> , for 1 conical tubes 25 mL screw cap <sup>1)</sup> , for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.         | 5804 718 003 | 155,00 € |
| <b>Adapter</b> , for 1 conical tube 50 mL, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                                   | 5804 718 003 | 155,00 € |
| <b>Adapter</b> , for 1 conical tubes 25 mL snap cap <sup>1)</sup> , for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.          | 5804 718 003 | 155,00 € |
| <b>Adapter</b> , for 1 conical tube 25 mL screw cap, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                         | 5804 758 005 | 155,00 € |
| <b>Adapter</b> , for 1 conical tube 50 mL, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                                   | 5804 758 005 | 155,00 € |
| <b>Adapter</b> , for 1 conical tube 25 mL snap cap, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                          | 5804 758 005 | 155,00 € |
| <b>Adapter</b> , for 12 reaction vessels 1.5 – 2.0 mL, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                       | 5804 751 000 | 144,00 € |
| <b>Adapter</b> , for 14 round-bottom tubes 1.2 – 5 mL, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                       | 5804 750 004 | 144,00 € |
| <b>Adapter</b> , for 1 round-bottom tube 80 – 100 mL, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                        | 5804 761 006 | 144,00 € |
| <b>Adapter</b> , for 1 tube 30 – 50 mL round-bottom/ 50 mL conical, skirted, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs. | 5804 759 001 | 144,00 € |
| <b>Adapter</b> , for 1 conical tube 50 mL, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                                   | 5804 758 005 | 155,00 € |
| <b>Adapter</b> , for 1 round-bottom tube 50 – 75 mL, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                         | 5804 760 000 | 144,00 € |
| <b>Adapter</b> , for 2 round-bottom tubes 18 – 30 mL, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                        | 5804 757 009 | 144,00 € |
| <b>Adapter</b> , for 4 Eppendorf Tubes® 5.0 mL and 15 mL conical tubes, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.      | 5804 755 006 | 155,00 € |
| <b>Adapter</b> , for 4 round-bottom tubes 7 – 18 mL, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                         | 5804 756 002 | 144,00 € |
| <b>Adapter</b> , for 6 round-bottom tubes 7 – 17 mL, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                         | 5804 754 000 | 150,00 € |
| <b>Adapter</b> , for 7 round-bottom tubes 3 – 15 mL, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                         | 5804 753 003 | 150,00 € |
| <b>Adapter</b> , for 9 round-bottom tubes 2.6 – 7 mL, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                        | 5804 752 007 | 150,00 € |
| <b>Adapter</b> , for 1 conical tubes 50 mL, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                                  | 5804 718 003 | 155,00 € |
| <b>Adapter</b> , for 2 conical tubes 15 mL conical, for 100 mL rectangular bucket in Rotor A-4-44, 2 pcs.                          | 5804 717 007 | 155,00 € |

<sup>1)</sup> Please order adapter for 25 mL conical tubes snap cap (5820 733 004) or adapter for 25 mL conical tubes screw cap (5820 734 000).

## Centrifuges and Rotors

### Rotor F-34-6-38

Eppendorf Product Catalog p. 172-173

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor F-34-6-38</b> , for 6 × 85 mL tubes, incl. rotor lid | 5804 727 002 | 1.966,00 € |
| <b>Rotor lid</b> , for F-34-6-38, for Centrifuge 58xx         | 5804 727 509 | 185,00 €   |

### Adapter for Rotor F-34-6-38

Eppendorf Product Catalog p. 172-173

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Adapter</b> , for 1 conical tube 50 mL, for Rotor F-34-6-38, 2 pcs.                      | 5804 775 007 | 89,20 €  |
| <b>Adapter</b> , for 1 conical tube 25 mL conical, for Rotor F-34-6-38, 2 pcs.              | 5804 775 007 | 89,20 €  |
| <b>Adapter</b> , for 1 conical tube 15 mL conical, for Rotor F-34-6-38, 2 pcs.              | 5804 776 003 | 89,20 €  |
| <b>Adapter</b> , for 1 round-bottom tube 15 – 18 mL, for Rotor F-34-6-38, 2 pcs.            | 5804 772 008 | 104,00 € |
| <b>Adapter</b> , for 1 round-bottom tube 30 mL, for Rotor F-34-6-38, 2 pcs.                 | 5804 773 004 | 84,10 €  |
| <b>Adapter</b> , for 1 conical tube 50 mL conical, for Rotor F-34-6-38, 2 pcs.              | 5804 775 007 | 89,20 €  |
| <b>Adapter</b> , for 1 round-bottom tube 50 mL, for Rotor F-34-6-38, 2 pcs.                 | 5804 774 000 | 84,10 €  |
| <b>Adapter</b> , for 1 Eppendorf Tubes® 5.0 mL, for Rotor F-34-6-38, 2 pcs.                 | 5804 777 000 | 83,40 €  |
| <b>Adapter</b> , for 2 round-bottom tubes 7 – 15 mL, for Rotor F-34-6-38, 2 pcs.            | 5804 771 001 | 104,00 € |
| <b>Adapter</b> , for 3 round-bottom tubes 4 – 7 mL (13x100 mm), for Rotor F-34-6-38, 2 pcs. | 5804 739 000 | 104,00 € |
| <b>Adapter</b> , for 3 round-bottom tubes 2 – 5 mL (13x75 mm), for Rotor F-34-6-38, 2 pcs.  | 5804 738 004 | 104,00 € |
| <b>Adapter</b> , for 4 reaction vessels 1.5 – 2.0 mL, for Rotor F-34-6-38, 2 pcs.           | 5804 770 005 | 138,00 € |

### Rotor FA-45-6-30

Eppendorf Product Catalog p. 174-175

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Rotor FA-45-6-30</b> , aerosol-tight, incl. rotor lid, for 6 × 15/50 mL conical tubes | 5820 715 006 | 2.382,00 € |

### Adapter for Rotor FA-45-6-30

Eppendorf Product Catalog p. 174-175

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Adapter</b> , for 1 conical tube 15 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.  | 5820 717 009 | 84,70 € |
| <b>Adapter</b> , for 1 round-bottom tube 16 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                                       | 5820 720 000 | 84,70 € |
| <b>Adapter</b> , for 1 round-bottom tube 4 – 8 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                                    | 5820 725 001 | 83,40 € |
| <b>Adapter</b> , for 1 round-bottom tube 2.6 – 5 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                                  | 5820 726 008 | 83,40 € |
| <b>Adapter</b> , for 1 round-bottom tube 7.5 – 12 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                                 | 5820 727 004 | 83,40 € |
| <b>Adapter</b> , for 1 round-bottom tube 5.5 – 10 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                                 | 5820 728 000 | 83,40 € |
| <b>Adapter</b> , for 1 round-bottom tube 7 – 17 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                                   | 5820 729 007 | 83,40 € |
| <b>Adapter</b> , for 1 round-bottom tube 30 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                                       | 5820 721 006 | 84,70 € |
| <b>Adapter</b> , for 1 Eppendorf Conical Tube 25 mL with screw cap 25 mL screw cap, for rotors for conical 50 mL tubes, 6 pcs. | 5820 734 000 | 33,00 € |
| <b>Adapter</b> , for 1 Eppendorf Conical Tube 25 mL with snap cap 25 mL snap cap, for rotors for conical 50 mL tubes, 6 pcs.   | 5820 733 004 | 32,40 € |
| <b>Adapter</b> , for 1 Eppendorf Tubes® 5.0 mL with snap cap, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                         | 5820 730 005 | 83,40 € |

## Centrifuge 5910 Ri

Eppendorf Product Catalog p. 176

|   | Order no.    | Price       |
|---|--------------|-------------|
| <b>Centrifuge 5910 Ri</b> , touch interface, refrigerated, without rotor, 230 V/50 – 60 Hz (EU) | 5943 000 011 | 10.610,00 € |

## Centrifuge Packages 5910 Ri

Eppendorf Product Catalog p. 176-177

|   | Order no.    | Price       |
|---|--------------|-------------|
| <b>Centrifuge 5910 Ri</b> , touch interface, refrigerated, High Speed Solution, with Rotor FA-6x50, 230 V/50 – 60 Hz (EU)   | 5943 000 149 | 12.300,00 € |
| <b>Centrifuge 5910 Ri</b> , touch interface, refrigerated, with Rotor S-4x750, incl. round buckets and adapters for 5 mL/15 mL/50 mL conical tubes, 230 V/50 – 60 Hz (EU)   | 5943 000 516 | 13.390,00 € |
| <b>Centrifuge 5910 Ri</b> , touch interface, refrigerated, with Rotor S-4xUniversal, incl. universal buckets and adapters for 5 mL/15 mL/50 mL conical tubes and plates, 230 V/50 – 60 Hz (EU)  | 5943 000 316 | 14.020,00 € |
| <b>Centrifuge 5910 Ri</b> , touch interface, refrigerated, Big Volume Harvesting Solution, with Rotor S-4xUniversal, incl. universal buckets and 1 L adapters (4x), Eppendorf 1 L bottles (4x), 230 V/50 – 60 Hz (EU)   | 5943 000 173 | 14.140,00 € |
| <b>Centrifuge 5910 Ri</b> , touch interface, refrigerated, Versatility Solution, with Rotor S-4xUniversal, incl. universal buckets and adapters for 5 mL/15 mL/50 mL conical tubes, 175 – 250 mL bottles and plates, 230 V/50 – 60 Hz (EU)                                | 5943 000 216 | 14.130,00 € |
| <b>Centrifuge 5910 Ri</b> , touch interface, refrigerated, with Rotor S-4xUniversal, incl. universal buckets, adapters for 5 mL/15 mL/50 mL conical tubes, 175 – 250 mL bottles and plates, adapters for 1000 mL bottles incl. 4 x 1000 mL bottles, 230 V/50 – 60 Hz (EU) | 5943 000 904 | 14.930,00 € |

## Centrifuge 5920 R

Eppendorf Product Catalog p. 176

|   | Order no.    | Price       |
|---|--------------|-------------|
| <b>Centrifuge 5920 R</b> , keypad, refrigerated, without rotor, 230 V/50 – 60 Hz (EU) | 5948 000 010 | 11.500,00 € |

## Centrifuge Packages 5920 R

Eppendorf Product Catalog p. 176-177

|  | Order no.    | Price       |
|--|--------------|-------------|
| <b>Centrifuge 5920 R</b> , keypad, refrigerated, with Rotor S-4xUniversal-Large, incl. universal buckets and adapters for 5 mL/15 mL/25 mL/50 mL conical tubes and plates, 230 V/50 – 60 Hz (EU) | 5948 000 911 | 15.200,00 € |
| <b>Centrifuge 5920 R</b> , keypad, refrigerated, with Rotor S-4xUniversal-Large, incl. universal buckets and adapters for round-bottom tubes with diameter 13/16 mm, 230 V/50 – 60 Hz (EU)       | 5948 000 810 | 14.830,00 € |
| <b>Centrifuge 5920 R</b> , keypad, refrigerated, with Rotor S-4x1000, incl. high-capacity buckets and adapters for 15 mL/50 mL conical tubes, 230 V/50 – 60 Hz (EU)                              | 5948 000 210 | 14.490,00 € |
| <b>Centrifuge 5920 R</b> , keypad, refrigerated, with Rotor S-4x1000, incl. plate/tube buckets and adapters for 15 mL/50 mL conical tubes, 230 V/50 – 60 Hz (EU)                                 | 5948 000 410 | 14.000,00 € |
| <b>Centrifuge 5920 R</b> , keypad, refrigerated, High-Throughput Screening Solution, with Rotor S-4x1000, incl. plate/tube bucket and plate carrier, 230 V/50 – 60 Hz (EU)                       | 5948 000 174 | 14.260,00 € |
| <b>Centrifuge 5920 R</b> , keypad, refrigerated, with Rotor S-4x1000, incl. round buckets and adapters for 15 mL/50 mL conical tubes, 230 V/50 – 60 Hz (EU)                                      | 5948 000 315 | 14.000,00 € |

## Centrifuges and Rotors

### Various accessories for centrifugation

Eppendorf Product Catalog p. 178-179

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Wide-neck bottle 400 mL</b> , for Rotors A-4-81, S-4x500, 400 mL, lid blue, 2 pcs. | 5810 729 009 | 200,00 € |
| <b>Wide-neck bottle 500 mL</b> , for Rotors A-4-81, S-4x500, rectangular, 2 pcs.      | 5820 707 003 | 200,00 € |
| <b>Wide-neck bottle 750 mL</b> , for Rotors S-4-104, S-4x750, S-4x1000, 2 pcs.        | 5820 708 000 | 198,00 € |
| <b>Wide-neck bottle 1,000 mL</b> , for Rotor S-4xUniversal, 2 pcs.                    | 5910 770 006 | 212,00 € |
| <b>Wide-neck bottle 400 mL</b> , for Rotor S-4x400, 2 pcs.                            | 5910 707 002 | 212,00 € |

### Rotor S-4xUniversal

Eppendorf Product Catalog p. 180-181

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Rotor S-4xUniversal</b> , incl. universal buckets                           | 5895 200 001 | 3.643,00 € |
| <b>Universal buckets</b> , for Rotor S-4xUniversal, 2 pcs.                     | 5895 202 004 | 1.035,00 € |
| <b>Universal buckets</b> , for Rotor S-4xUniversal, 4 pcs.                     | 5895 203 000 | 2.057,00 € |
| <b>Caps</b> , for Rotor S-4xUniversal, aerosol-tight, 2 pcs.                   | 5910 750 005 | 213,00 €   |
| <b>Wide-neck bottle 1,000 mL</b> , for Rotor S-4xUniversal, 2 pcs.             | 5910 770 006 | 212,00 €   |
| <b>Wide-neck bottle 750 mL</b> , for Rotors S-4-104, S-4x750, S-4x1000, 2 pcs. | 5820 708 000 | 198,00 €   |

### Adapter for Rotor S-4xUniversal

Eppendorf Product Catalog p. 180-181

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Adapter</b> , for 7 conical tubes 50 mL or 1 bottle 250 mL or 1 microplate 250 mL, for Rotor S-4xUniversal, 2 pcs.                           | 5910 751 001 | 420,00 € |
| <b>Adapter</b> , for 7 conical tubes 50 mL or 1 bottle 250 mL or 1 microplate 25 mL screw cap, for Rotor S-4xUniversal, 2 pcs.                  | 5910 751 001 | 420,00 € |
| <b>Adapter</b> , for 7 conical tubes 50 mL or 1 bottle 250 mL or 1 microplate Microplate 96/384 wells, for Rotor S-4xUniversal, 2 pcs.          | 5910 751 001 | 420,00 € |
| <b>Adapter</b> , for 7 conical tubes 50 mL or 1 bottle 250 mL or 1 microplate 50 mL, for Rotor S-4xUniversal, 2 pcs.                            | 5910 751 001 | 420,00 € |
| <b>Adapter</b> , for 7 conical tubes 50 mL or 1 bottle 250 mL or 1 microplate Microplate 96/384 wells, for Rotor S-4xUniversal, 2 pcs.          | 5910 751 001 | 420,00 € |
| <b>Adapter</b> , for 17 conical tubes 15 mL or 16 Eppendorf Tubes 5 mL or 1 microplate MTP, DWP, for Rotor S-4xUniversal, 2 pcs.                | 5910 752 008 | 299,00 € |
| <b>Adapter</b> , for 17 conical tubes 15 mL or 16 Eppendorf Tubes 5 mL or 1 microplate 15 mL conical, for Rotor S-4xUniversal, 2 pcs.           | 5910 752 008 | 299,00 € |
| <b>Adapter</b> , for 9 conical tubes Microplate 96/384 wells, for Rotor S-4xUniversal, 2 pcs.   | 5910 769 008 | 359,00 € |
| <b>Adapter</b> , for 9 conical tubes 25 mL snap cap, for Rotor S-4xUniversal, 2 pcs.  | 5910 769 008 | 359,00 € |
| <b>Adapter</b> , for 9 conical tubes 25 mL screw cap, for Rotor S-4xUniversal, 2 pcs.   | 5910 769 008 | 359,00 € |
| <b>Adapter</b> , for 5 MTP, DWP, PCR plates, for Rotor S-4xUniversal, 2 pcs.  | 5910 753 004 | 299,00 € |
| <b>Adapter</b> , for 1 1,000 mL bottle 1,000 mL, for Centrifuge 5910 R with Rotor S-4xUniversal, 2 pcs.   | 5910 756 003 | 299,00 € |
| <b>Adapter</b> , for 46 tubes 12 x 75 mm, 5 mL FACS, for Rotor S-4xUniversal, 2 pcs.  | 5910 758 006 | 299,00 € |
| <b>Adapter</b> , for 30 round-bottom tubes 4 mL – 8 mL, for Rotor S-4xUniversal, 2 pcs.   | 5910 755 007 | 299,00 € |
| <b>Adapter</b> , for 26 round-bottom tubes 7.5 mL – 12 mL, for Rotor S-4xUniversal, 2 pcs.  | 5910 754 000 | 299,00 € |
| <b>Adapter</b> , for 21 tubes 9 mL, for Rotor S-4xUniversal, 2 pcs.   | 5910 762 003 | 299,00 € |
| <b>Adapter for 500 mL bottle</b> , for 1 bottle 500 mL Corning®, for Rotor S-4xUniversal, 2 pcs.  | 5910 760 000 | 299,00 € |
| <b>Adapter</b> , for 7 conical tubes 50 mL or 1 bottle 250 mL or 1 microplate 50 mL conical, for Rotor S-4xUniversal, 2 pcs.                    | 5910 751 001 | 420,00 € |
| <b>Adapter</b> , for 9 conical tubes 50 mL conical, for Rotor S-4xUniversal, 2 pcs.   | 5910 769 008 | 359,00 € |
| <b>Adapter</b> , for 7 tubes 50 mL conical, skirted, for Rotor S-4xUniversal, 2 pcs.  | 5910 764 006 | 299,00 € |
| <b>Adapter</b> , for 17 conical tubes 15 mL or 16 Eppendorf Tubes 5 mL or 1 microplate Eppendorf Tubes® 5.0 mL, for Rotor S-4xUniversal, 2 pcs. | 5910 752 008 | 299,00 € |
| <b>Adapter</b> , for 1 bottle 750 mL, for Rotor S-4xUniversal, 2 pcs.   | 5910 757 000 | 299,00 € |



## Rotor S-4x500

Eppendorf Product Catalog p. 182-183

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor S-4x500</b> , for 500 mL rectangular buckets or MTP/Flex buckets   |              |            |
| incl. 4 x 500 mL rectangular buckets  | 5895 170 005 | 2.645,00 € |
| without buckets   | 5895 171 001 | 1.311,00 € |
| <b>Rectangular bucket 500 mL</b> , for Rotor A-4-81, 4 pcs.   | 5810 730 007 | 1.840,00 € |
| <b>Cap</b> , aerosol-tight, for 500 mL rectangular bucket, 2 pcs.   | 5810 724 007 | 180,00 €   |
| <b>Wide-neck bottle 500 mL</b> , for Rotors A-4-81, S-4x500, rectangular, 2 pcs.  | 5820 707 003 | 200,00 €   |
| <b>Wide-neck bottle 400 mL</b> , for Rotors A-4-81, S-4x500, 400 mL, lid blue, 2 pcs.   | 5810 729 009 | 200,00 €   |
| <b>MTP/Flex buckets</b> , for Rotors A-4-81, S-4x500, for use with IsoRack and cell culture flask adapters as well as MTP and DWP, 2 pcs. | 5810 742 005 | 791,00 €   |
| <b>MTP/Flex buckets</b> , for Rotors A-4-81, S-4x500, for use with IsoRack and cell culture flask adapters as well as MTP and DWP, 4 pcs. | 5810 741 009 | 1.460,00 € |

## Adapter for Rotor S-4x500

Eppendorf Product Catalog p. 182-183

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>MTP/Flex buckets</b> , for Rotors A-4-81, S-4x500, for use with IsoRack and cell culture flask adapters as well as MTP and DWP, 2 pcs.  | 5810 742 005 | 791,00 € |
| <b>Adapter</b> , for 5 conical tubes 25 mL screw cap, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.  | 5810 723 000 | 196,00 € |
| <b>Adapter</b> , for 5 conical tubes 25 mL snap cap, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.   | 5810 723 000 | 196,00 € |
| <b>Adapter</b> , for 5 conical tubes 50 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.  | 5810 723 000 | 196,00 € |
| <b>Adapter</b> , for 16 round-bottom tubes 3 – 15 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.  | 5810 748 003 | 196,00 € |
| <b>Adapter</b> , for 16 round-bottom tubes 7 – 17 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.  | 5810 721 008 | 196,00 € |
| <b>Adapter</b> , for 18 round-bottom tubes 2.6 – 7 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.   | 5810 747 007 | 196,00 € |
| <b>Adapter</b> , for 1 bottle 180 – 250 mL flat/ 175 – 225 mL conical <sup>1)</sup> , for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, Please note that these conical bottles need to be used together with a conical form insert provided by the manufacturer, 2 pcs. | 5825 722 000 | 197,00 € |
| <b>Adapter</b> , for 1 bottle 400 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.  | 5810 728 002 | 164,00 € |
| <b>Adapter</b> , for 20 reaction vessels 1.5/2.0 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.   | 5810 745 004 | 196,00 € |
| <b>Adapter</b> , for 20 round-bottom tubes 1.2 – 5 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.   | 5810 746 000 | 196,00 € |
| <b>Adapter</b> , for 25 round-bottom tubes 2.6 – 5 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.   | 5810 720 001 | 196,00 € |
| <b>Adapter</b> , for 12 Eppendorf Tubes® 5.0 mL and 15 mL conical tubes, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.   | 5810 722 004 | 196,00 € |
| <b>Adapter</b> , for 5 conical tubes 50 mL, for 500 mL rectangular bucket in Rotor A-4-81 and Rotor S-4x500, 2 pcs.  | 5810 723 000 | 196,00 € |

<sup>1)</sup> Please note that these conical bottles need to be used together with a conical form insert provided by the manufacturer.

## Centrifuges and Rotors

### Rotor S-4x400

| Eppendorf Product Catalog p. 184-185                       | Order no.    | Price      |
|--|--------------|------------|
| <b>Rotor S-4x400</b> , incl. 400 mL round bucket           | 5895 180 000 | 2.704,00 € |
| <b>Caps</b> , for Rotor S-4x400, aerosol-tight, 2 pcs.     | 5910 700 008 | 195,00 €   |
| <b>Wide-neck bottle 400 mL</b> , for Rotor S-4x400, 2 pcs. | 5910 707 002 | 212,00 €   |

### Adapter for Rotor S-4x400

| Eppendorf Product Catalog p. 184-185  | Order no.    | Price    |
|---|--------------|----------|
| <b>Adapter</b> , for 4 tubes 25 mL screw cap <sup>1)</sup> , for Rotor S-4x400, 2 pcs.                | 5910 701 004 | 251,00 € |
| <b>Adapter</b> , for 4 tubes 25 mL snap cap <sup>1)</sup> , for Rotor S-4x400, 2 pcs.                 | 5910 701 004 | 251,00 € |
| <b>Adapter</b> , for 4 tubes 50 mL, for Rotor S-4x400, 2 pcs.   | 5910 701 004 | 251,00 € |
| <b>Adapter</b> , for 26 tubes 1.5 mL/2.0 mL, for Rotor S-4x400, 2 pcs.                                | 5910 708 009 | 251,00 € |
| <b>Adapter</b> , for 17 tubes 12 mm x 75 mm, for Rotor S-4x400, 2 pcs.                                | 5910 711 000 | 251,00 € |
| <b>Adapter</b> , for 7 Eppendorf Tubes® 5.0 mL and 15 mL conical, for Rotor S-4x400, 2 pcs.           | 5910 702 000 | 251,00 € |
| <b>Adapter</b> , for 1 bottle 175 – 225 mL conical, for Rotor S-4x400, 2 pcs.                         | 5910 714 009 | 251,00 € |
| <b>Adapter</b> , for 1 bottle 175 mL – 250 mL, for Rotor S-4x400, 2 pcs.                              | 5910 705 000 | 251,00 € |
| <b>Adapter</b> , for 1 bottle 400 mL, for Centrifuge 5910 R with Rotor S-4x400, 2 pcs.                | 5910 706 006 | 159,00 € |
| <b>Adapter</b> , for 15 tubes 4 mL – 8 mL, for Rotor S-4x400, 2 pcs.                                  | 5910 703 007 | 251,00 € |
| <b>Adapter</b> , for 3 tubes 50 mL skirted, Centriprep and Oak Ridge 30 mL, for Rotor S-4x400, 2 pcs. | 5910 712 006 | 251,00 € |
| <b>Adapter</b> , for 4 tubes 50 mL conical, for Rotor S-4x400, 2 pcs.                                 | 5910 701 004 | 251,00 € |
| <b>Adapter</b> , for 11 tubes 7.5 mL – 12 mL, for Rotor S-4x400, 2 pcs.                               | 5910 704 003 | 251,00 € |
| <b>Adapter</b> , for 8 tubes 9 mL, for Rotor S-4x400, 2 pcs.  | 5910 709 005 | 251,00 € |

<sup>1)</sup> Please order adapter for 25 mL conical tubes snap cap (5820 733 004) or adapter for 25 mL conical tubes screw cap (5820 734 000).

### Rotor S-4x750

| Eppendorf Product Catalog p. 186-187  | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor S-4x750</b> , aerosol-tight, without buckets   | 5895 121 004 | 1.406,00 € |
| <b>Rotor S-4x750</b> , aerosol-tight, incl. round buckets                                     | 5895 120 008 | 3.068,00 € |
| <b>Plate bucket</b> , for Rotor S-4x750, aerosol-tight capable, incl. plate carrier, 2 pcs.   | 5895 125 000 | 862,00 €   |
| <b>Round bucket</b> , for Rotor S-4x750, 2 pcs.   | 5895 123 007 | 862,00 €   |
| <b>Round bucket</b> , for Rotor S-4x750, 4 pcs.   | 5895 122 000 | 1.667,00 € |
| <b>Plate bucket</b> , for Rotor S-4x750, aerosol-tight capable, incl. plate carrier, 4 pcs.   | 5895 124 003 | 1.667,00 € |
| <b>Caps</b> , rotors S-4-104, S-4x750, plate bucket, aerosol-tight, 2 pcs.                    | 5820 748 001 | 213,00 €   |
| <b>Caps for 750 mL/1000 mL</b> , for rotors S-4-104, S-4x750, S-4x1000, aerosol-tight, 2 pcs. | 5820 747 005 | 195,00 €   |
| <b>Wide-neck bottle 750 mL</b> , for Rotors S-4-104, S-4x750, S-4x1000, 2 pcs.                | 5820 708 000 | 198,00 €   |

### Adapter for Rotor S-4x750

| Eppendorf Product Catalog p. 186-187   | Order no.    | Price    |
|--|--------------|----------|
| <b>Plate bucket</b> , for Rotor S-4x750, aerosol-tight capable, incl. plate carrier, 2 pcs.  | 5895 125 000 | 862,00 € |
| <b>Adapter</b> , for 7 conical tubes 25 mL screw cap <sup>1)</sup> , for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.     | 5825 733 002 | 299,00 € |
| <b>Adapter</b> , for 7 conical tubes 25 mL snap cap <sup>1)</sup> , for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.      | 5825 733 002 | 299,00 € |
| <b>Adapter</b> , for 7 conical tubes Conical SnapTec® 50, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                | 5825 733 002 | 299,00 € |
| <b>Adapter</b> , for 14 tubes 14 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                                     | 5825 748 000 | 299,00 € |
| <b>Adapter</b> , for tubes 30 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.  | 5825 755 006 | 295,00 € |
| <b>Adapter</b> , for 27 round-bottom tubes 5 mL FACS, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                    | 5825 747 003 | 299,00 € |
| <b>Adapter</b> , for bottle 500 mL Corning®, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                             | 5825 745 000 | 299,00 € |
| <b>Adapter</b> , for 7 conical tubes 50 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                              | 5825 733 002 | 299,00 € |
| <b>Adapter</b> , for 14 Eppendorf Tubes® 5.0 mL and 15 mL conical tubes, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs. | 5825 734 009 | 299,00 € |
| <b>Adapter</b> , for 20 round-bottom tubes 5.5 – 12 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                  | 5825 736 001 | 299,00 € |

<sup>1)</sup> Please order adapter for 25 mL conical tubes snap cap (5820 733 004) or adapter for 25 mL conical tubes screw cap (5820 734 000).

## Adapter for Rotor S-4x750

| Eppendorf Product Catalog p. 186-187  | Order no.    | Price    |
|---|--------------|----------|
| <b>Adapter</b> , for 20 round-bottom tubes 9 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                      | 5825 743 008 | 299,00 € |
| <b>Adapter</b> , for 23 round-bottom tubes 2.6 – 8 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                | 5825 738 004 | 299,00 € |
| <b>Adapter</b> , for 1 bottle 250 mL flat, 175 – 225/ 250 mL conical, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs. | 5825 741 005 | 299,00 € |
| <b>Adapter</b> , for 50 reaction vessels 1.5/2.0 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                  | 5825 740 009 | 299,00 € |
| <b>Adapter</b> , for 5 tubes 50 mL conical, skirted, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                  | 5825 732 006 | 299,00 € |
| <b>Adapter</b> , for 1 Eppendorf bottle 750 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                       | 5825 744 004 | 159,00 € |

<sup>1)</sup> Please order adapter for 25 mL conical tubes snap cap (5820 733 004) or adapter for 25 mL conical tubes screw cap (5820 734 000).

## Rotor S-4xUniversal-Large

| Eppendorf Product Catalog p. 188-191  | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor S-4xUniversal-Large</b> , incl. universal buckets                        | 5895 190 006 | 3.494,00 € |
| <b>Buckets</b> , for Rotor S-4xUniversal-Large, 4 pcs.                            | 5895 192 009 | 2.252,00 € |
| <b>Caps</b> , Rotor S-4xUniversal-Large, universal buckets, aerosol-tight, 2 pcs. | 5920 752 006 | 213,00 €   |

## Adapter for Rotor S-4xUniversal-Large

| Eppendorf Product Catalog p. 188-191   | Order no.    | Price    |
|--|--------------|----------|
| <b>Plate carrier</b> , for Rotor S-4xUniversal-Large, universal buckets, 2 pcs.  | 5920 737 007 | 299,00 € |
| <b>Adapters for Rotors</b> , for conical tubes 25/50 mL and plates, for Rotor S-4xUniversal-Large, 2 pcs.                          | 5920 756 001 | 353,00 € |
| <b>Adapters for Rotors</b> , for conical tubes 25/50 mL and plates, for Rotor S-4xUniversal-Large, 2 pcs.                          | 5920 756 001 | 353,00 € |
| <b>Adapters for Rotors</b> , for conical tubes 25/50 mL and plates, for Rotor S-4xUniversal-Large, 2 pcs.                          | 5920 756 001 | 353,00 € |
| <b>Adapters for Rotors</b> , for conical tubes 50 mL, 175 – 250 mL bottles and plates, for Rotor S-4xUniversal-Large, 2 pcs.       | 5920 755 005 | 411,00 € |
| <b>Adapters for Rotors</b> , for conical tubes 25 mL, 175 – 250 mL bottles and plates, for Rotor S-4xUniversal-Large, 2 pcs.       | 5920 755 005 | 411,00 € |
| <b>Adapters for Rotors</b> , for conical tubes 50 mL, 175 – 250 mL bottles and plates, for Rotor S-4xUniversal-Large, 2 pcs.       | 5920 755 005 | 411,00 € |
| <b>Adapters for Rotors</b> , for conical tubes 50 mL, 175 – 250 mL bottles and plates, for Rotor S-4xUniversal-Large, 2 pcs.       | 5920 755 005 | 411,00 € |
| <b>Adapters for Rotors</b> , for conical tubes 50 mL, 175 – 250 mL bottles and plates, for Rotor S-4xUniversal-Large, 2 pcs.       | 5920 755 005 | 411,00 € |
| <b>Adapters for Rotors</b> , for conical tubes 25/50 mL and plates, for Rotor S-4xUniversal-Large, 2 pcs.                          | 5920 756 001 | 353,00 € |
| <b>Adapters for Rotors</b> , for conical tubes 5/15 mL and plates, for Rotor S-4xUniversal-Large, 2 pcs.                           | 5920 757 008 | 294,00 € |
| <b>Adapters for Rotors</b> , for conical tubes 5/15 mL and plates, for Rotor S-4xUniversal-Large, 2 pcs.                           | 5920 757 008 | 294,00 € |
| <b>Adapters for Rotors</b> , for 92 tubes 1.5/2.0 mL, for Rotor S-4xUniversal-Large, 2 pcs.  | 5920 747 002 | 299,00 € |
| <b>Adapters for Rotors</b> , for 24 tubes 14 mL, for Rotor S-4xUniversal-Large, 2 pcs.   | 5920 751 000 | 299,00 € |
| <b>Adapters for Rotors</b> , for 28 tubes 9 mL, for Rotor S-4xUniversal-Large, 2 pcs.  | 5920 746 006 | 299,00 € |
| <b>Adapters for Rotors</b> , for 1 bottle 250 mL conical (Corning®, Nunc®), for Rotor S-4xUniversal-Large, 2 pcs.                  | 5920 750 003 | 299,00 € |
| <b>Adapters for Rotors</b> , for 2 bottles 250 mL flat, 175 – 225 mL conical <sup>1)</sup> , for Rotor S-4xUniversal-Large, 2 pcs. | 5920 740 008 | 299,00 € |
| <b>Adapters for Rotors</b> , for 12 tubes 30 mL, for Rotor S-4xUniversal-Large, 2 pcs.   | 5920 743 007 | 299,00 € |
| <b>Adapters for Rotors</b> , for 1 bottle 500 mL conical Corning®, for Rotor S-4xUniversal-Large, 2 pcs.                           | 5920 744 003 | 299,00 € |
| <b>Adapters for Rotors</b> , for 1 bottle 500 mL flat Nalgene®, for Rotor S-4xUniversal-Large, 2 pcs.                              | 5920 745 000 | 299,00 € |
| <b>Adapters for Rotors</b> , for 12 tubes 50 mL Nalgene® Oak Ridge, for Rotor S-4xUniversal-Large, 2 pcs.                          | 5920 753 002 | 299,00 € |
| <b>Adapters for Rotors</b> , for 7 tubes 50 mL conical, skirted, for Rotor S-4xUniversal-Large, 2 pcs.                             | 5920 748 009 | 299,00 € |
| <b>Adapters for Rotors</b> , for 51 tubes 5 mL FACS, for Rotor S-4xUniversal-Large, 2 pcs.   | 5920 742 000 | 299,00 € |
| <b>Adapters for Rotors</b> , for 1 bottle 750 mL, for Rotor S-4xUniversal-Large, 2 pcs.  | 5920 741 004 | 299,00 € |
| <b>Adapters for Rotors</b> , for 49 tubes 2.6 – 8 mL, for Rotor S-4xUniversal-Large, 2 pcs.  | 5920 739 000 | 299,00 € |
| <b>Adapters for Rotors</b> , for 37 tubes 5.5 – 12 mL, for Rotor S-4xUniversal-Large, 2 pcs.                                       | 5920 738 003 | 299,00 € |
| <b>Adapters for Rotors</b> , for 3 Micro Fluidic Cards, for Rotor S-4xUniversal-Large, 2 pcs.                                      | 5920 749 005 | 299,00 € |

<sup>1)</sup> Please note that these conical bottles need to be used together with a conical form insert provided by the manufacturer.

## Centrifuges and Rotors

### Rotor S-4x1000 high-capacity buckets

Eppendorf Product Catalog p. 192-193

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Rotor S-4x1000</b> , incl. high-capacity buckets      | 5895 118 003 | 3.563,00 € |
| <b>High-Capacity bucket</b> , Rotor S-4x1000, 2 pcs.     | 5895 107 001 | 1.340,00 € |
| <b>High-Capacity bucket</b> , for Rotor S-4x1000, 4 pcs. | 5895 106 005 | 2.682,00 € |

### Adapter for Rotor S-4x1000 high-capacity buckets

Eppendorf Product Catalog p. 192-193

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Adapters for Rotors</b> , for 13 conical tubes 25 mL, Rotor S-4x1000, for high-capacity buckets, 2 pcs.  | 5920 715 003 | 299,00 € |
| <b>Adapters for Rotors</b> , for 13 conical tubes 25 mL, Rotor S-4x1000, for high-capacity buckets, 2 pcs.  | 5920 715 003 | 299,00 € |
| <b>Adapters for Rotors</b> , for 49 tubes 2.6 – 8 mL, Rotor S-4x1000, for high-capacity buckets, 2 pcs.   | 5920 718 002 | 299,00 € |
| <b>Adapters for Rotors</b> , for Eppendorf Tubes® 5.0 mL and 15 mL conical tubes, Rotor S-4x1000, for high-capacity buckets, 2 pcs.                                     | 5920 716 000 | 299,00 € |
| <b>Adapters for Rotors</b> , for 36 tubes 5.5 – 12 mL, Rotor S-4x1000, for high-capacity buckets, 2 pcs.  | 5920 720 007 | 299,00 € |
| <b>Adapters for Rotors</b> , for 29 tubes 9 mL, Rotor S-4x1000, for high-capacity buckets, 2 pcs.   | 5920 722 000 | 299,00 € |
| <b>Adapters for Rotors</b> , for 2 bottles 250 mL flat, 175 – 225 mL conical <sup>1)</sup> and Nunc® 250 mL conical , Rotor S-4x1000, for high-capacity buckets, 2 pcs. | 5920 717 006 | 299,00 € |
| <b>Adapters for Rotors</b> , for 13 conical tubes 50 mL, Rotor S-4x1000, for high-capacity buckets, 2 pcs.  | 5920 715 003 | 299,00 € |
| <b>Adapters for Rotors</b> , for 52 tubes 5 mL FACS, Rotor S-4x1000, for high-capacity buckets, 2 pcs.  | 5920 724 002 | 303,00 € |
| <b>Plate carrier</b> , Rotor S-4x1000, for high-capacity buckets, 2 pcs.  | 5920 729 004 | 299,00 € |

<sup>1)</sup> Please note that these conical bottles need to be used together with a conical form insert provided by the manufacturer.

### Rotor S-4x1000 plate/tube buckets

Eppendorf Product Catalog p. 194-195

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor S-4x1000</b> , aerosol-tight, incl. plate/tube buckets                                     | 5895 117 007 | 2.961,00 € |
| <b>Caps</b> , Rotor S-4x1000: plate/tube bucket; Rotor S-4x750: plate bucket, aerosol-tight, 2 pcs. | 5895 111 009 | 213,00 €   |
| <b>Plate carrier</b> , for Rotor S-4x1000 with plate/tube buckets, 2 pcs.                           | 5920 705 008 | 299,00 €   |
| <b>Plate/tube bucket</b> , for Rotor S-4x1000, 4 pcs.   | 5895 104 002 | 1.781,00 € |
| <b>Plate/tube bucket</b> , for Rotor S-4x1000, 2 pcs.   | 5895 105 009 | 921,00 €   |

### Adapter for Rotor S-4x1000 plate/tube buckets

Eppendorf Product Catalog p. 194-195

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Plate carrier</b> , for Rotor S-4x1000 with plate/tube buckets, 2 pcs.   | 5920 705 008 | 299,00 € |
| <b>Adapters for Rotors</b> , for 10 conical tubes 25 mL, for Rotor S-4x1000 with plate/tube buckets, 2 pcs.                             | 5920 709 003 | 299,00 € |
| <b>Adapters for Rotors</b> , for 10 conical tubes 25 mL, for Rotor S-4x1000 with plate/tube buckets, 2 pcs.                             | 5920 709 003 | 299,00 € |
| <b>Adapters for Rotors</b> , for 10 conical tubes 50 mL, for Rotor S-4x1000 with plate/tube buckets, 2 pcs.                             | 5920 709 003 | 299,00 € |
| <b>Adapters for Rotors</b> , for 6 conical tubes 25 mL, for Rotor S-4x1000 with plate/tube buckets, 2 pcs.                              | 5920 711 008 | 299,00 € |
| <b>Adapters for Rotors</b> , for 6 conical tubes 25 mL, for Rotor S-4x1000 with plate/tube buckets, 2 pcs.                              | 5920 711 008 | 299,00 € |
| <b>Adapters for Rotors</b> , for 6 conical tubes 50 mL, for Rotor S-4x1000 with plate/tube buckets, 2 pcs.                              | 5920 711 008 | 299,00 € |
| <b>Adapters for Rotors</b> , for 35 tubes 2.6 – 8 mL , for Rotor S-4x1000 with plate/tube buckets, 2 pcs.                               | 5920 706 004 | 299,00 € |
| <b>Adapters for Rotors</b> , for 22 Eppendorf Tubes® 5.0 mL and 15 mL conical tubes, for Rotor S-4x1000 with plate/tube buckets, 2 pcs. | 5920 710 001 | 299,00 € |
| <b>Adapters for Rotors</b> , for 16 Eppendorf Tubes® 5.0 mL and 15 mL conical tubes, for Rotor S-4x1000 with plate/tube buckets, 2 pcs. | 5920 712 004 | 299,00 € |
| <b>Adapters for Rotors</b> , for 33 tubes 5.5 – 12 mL, for Rotor S-4x1000 with plate/tube buckets, 2 pcs.                               | 5920 707 000 | 299,00 € |
| <b>Adapters for Rotors</b> , for 28 tubes 9 mL, for Rotor S-4x1000 with plate/tube buckets, 2 pcs.                                      | 5920 708 007 | 299,00 € |
| <b>Adapters for Rotors</b> , for 10 conical tubes 50 mL, for Rotor S-4x1000 with plate/tube buckets, 2 pcs.                             | 5920 709 003 | 299,00 € |
| <b>Adapters for Rotors</b> , for 6 conical tubes 50 mL, for Rotor S-4x1000 with plate/tube buckets, 2 pcs.                              | 5920 711 008 | 299,00 € |

## Rotor S-4x1000 round buckets

Eppendorf Product Catalog p. 196-197

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Round bucket</b> , for Rotor S-4x1000, 2 pcs.  | 5895 103 006 | 921,00 €   |
| <b>Round bucket</b> , for Rotor S-4x1000, 4 pcs.  | 5895 102 000 | 1.781,00 € |
| <b>Rotor S-4x1000 round buckets</b> , aerosol-tight, incl. 4 × 1,000 mL round buckets         | 5895 100 007 | 2.961,00 € |
| <b>Caps for 750 mL/1000 mL</b> , for rotors S-4-104, S-4x750, S-4x1000, aerosol-tight, 2 pcs. | 5820 747 005 | 195,00 €   |
| <b>Wide-neck bottle 750 mL</b> , for Rotors S-4-104, S-4x750, S-4x1000, 2 pcs.                | 5820 708 000 | 198,00 €   |

## Adapter for Rotor S-4x1000 round buckets

Eppendorf Product Catalog p. 196-197

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Plate/tube bucket</b> , for Rotor S-4x1000, 2 pcs.  | 5895 105 009 | 921,00 € |
| <b>Adapter</b> , for 7 conical tubes 25 mL screw cap <sup>1)</sup> , for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.     | 5825 733 002 | 299,00 € |
| <b>Adapter</b> , for 6 conical tubes 25 mL snap cap <sup>1)</sup> , for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.      | 5825 733 002 | 299,00 € |
| <b>Adapter</b> , for 7 conical tubes Conical SnapTec® 50, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                | 5825 733 002 | 299,00 € |
| <b>Adapter</b> , for 14 tubes 14 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                                     | 5825 748 000 | 299,00 € |
| <b>Adapter</b> , for 1 bottle 1 Liter Nalgene®, for Rotor S-4x1000 round buckets, 2 pcs.   | 5920 700 006 | 299,00 € |
| <b>Adapter</b> , for tubes 30 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.  | 5825 755 006 | 295,00 € |
| <b>Adapter</b> , for 27 round-bottom tubes 5 mL FACS, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                    | 5825 747 003 | 299,00 € |
| <b>Adapter</b> , for bottle 500 mL Corning®, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                             | 5825 745 000 | 299,00 € |
| <b>Adapter</b> , for 1 bottle 500 mL Nalgene®, for Rotor S-4x1000 round buckets, 2 pcs.  | 5920 703 005 | 295,00 € |
| <b>Adapter</b> , for 1 bottle 600 mL TPP, for Rotor S-4x1000 round buckets, 2 pcs.   | 5920 701 002 | 295,00 € |
| <b>Adapter</b> , for 7 conical tubes 50 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                              | 5825 733 002 | 299,00 € |
| <b>Adapter</b> , for 14 Eppendorf Tubes® 5.0 mL and 15 mL conical tubes, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs. | 5825 734 009 | 299,00 € |
| <b>Adapter</b> , for 20 round-bottom tubes 5.5 – 12 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                  | 5825 736 001 | 299,00 € |
| <b>Adapter</b> , for 20 round-bottom tubes 9 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                         | 5825 743 008 | 299,00 € |
| <b>Adapter</b> , for 23 round-bottom tubes 2.6 – 8 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                   | 5825 738 004 | 299,00 € |
| <b>Adapter</b> , for 1 bottle 250 mL flat, 175 – 225/ 250 mL conical, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.    | 5825 741 005 | 299,00 € |
| <b>Adapter</b> , for 50 reaction vessels 1.5/2.0 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                     | 5825 740 009 | 299,00 € |
| <b>Adapter</b> , for 5 tubes 50 mL conical, skirted, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                     | 5825 732 006 | 299,00 € |
| <b>Adapter</b> , for 1 Eppendorf bottle 750 mL, for Rotor S-4-104, Rotor S-4x1000 round buckets and Rotor S-4x750, 2 pcs.                          | 5825 744 004 | 159,00 € |

<sup>1)</sup> Please order adapter for 25 mL conical tubes snap cap (5820 733 004) or adapter for 25 mL conical tubes screw cap (5820 734 000).

## Centrifuges and Rotors

### Rotor FA-6x250

Eppendorf Product Catalog p. 200-201

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor FA-6x250</b> , aerosol-tight, for 6 x 250 mL bottles, incl. Eppendorf QuickLock® rotor lid | 5895 175 007 | 4.455,00 € |

### Adapter for Rotor FA-6x250

Eppendorf Product Catalog p. 200-201

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Adapter</b> , for 1 conical tube 25 mL, for Rotor FA-6x250, 2 pcs.                               | 5920 760 009 | 299,00 € |
| <b>Adapter</b> , for 1 conical tube 50 mL, for Rotor FA-6x250, 2 pcs.                               | 5920 760 009 | 299,00 € |
| <b>Adapters for Rotors</b> , for 50 mL conical skirted tubes 50 mL, for Rotor FA-6x250, 2 pcs.      | 5920 766 007 | 299,00 € |
| <b>Adapters for Rotors</b> , for 9 FACS round tubes 5 mL, for Rotor FA-6x250, 2 pcs.                | 5920 765 000 | 299,00 € |
| <b>Adapters for Rotors</b> , for 8 blood collection tubes 2.6 – 8.0 mL, for Rotor FA-6x250, 2 pcs.  | 5920 763 008 | 299,00 € |
| <b>Adapters for Rotors</b> , for 7 blood collection tubes 5.5 – 12.0 mL, for Rotor FA-6x250, 2 pcs. | 5920 762 001 | 299,00 € |
| <b>Adapters for Rotors</b> , for 7 blood collection tubes 9 mL, for Rotor FA-6x250, 2 pcs.          | 5920 764 004 | 299,00 € |
| <b>Adapters for Rotors</b> , for 7 Oak Ridge tubes 10 mL, for Rotor FA-6x250, 2 pcs.                | 5920 769 006 | 299,00 € |
| <b>Adapters for Rotors</b> , for 1 Oak Ridge tubes 85 mL, for Rotor FA-6x250, 2 pcs.                | 5920 768 000 | 299,00 € |
| <b>Adapters for Rotors</b> , for 7 Oak Ridge tubes 16 mL, for Rotor FA-6x250, 2 pcs.                | 5920 770 004 | 299,00 € |
| <b>Adapters for Rotors</b> , for 3 Oak Ridge tubes 30 mL, for Rotor FA-6x250, 2 pcs.                | 5920 767 003 | 299,00 € |
| <b>Adapters for Rotors</b> , for 1 Oak Ridge tubes 50 mL, for Rotor FA-6x250, 2 pcs.                | 5920 771 000 | 299,00 € |
| <b>Adapter</b> , for 1 conical tube 50 mL, for Rotor FA-6x250, 2 pcs.                               | 5920 760 009 | 299,00 € |

### Rotor FA-6x50

Eppendorf Product Catalog p. 202-203

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Rotor FA-6x50</b> , aerosol-tight, incl. rotor lid, 6 x 50 mL conical tubes | 5895 150 004 | 2.336,00 € |

### Adapter for Rotor FA-6x50

Eppendorf Product Catalog p. 202-203

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Adapter</b> , for 1 conical tube 15 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.  | 5820 717 009 | 84,70 € |
| <b>Adapter</b> , for 1 round-bottom tube 16 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                                       | 5820 720 000 | 84,70 € |
| <b>Adapter</b> , for 1 round-bottom tube 4 – 8 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                                    | 5820 725 001 | 83,40 € |
| <b>Adapter</b> , for 1 round-bottom tube 2.6 – 5 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                                  | 5820 726 008 | 83,40 € |
| <b>Adapter</b> , for 1 round-bottom tube 7.5 – 12 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                                 | 5820 727 004 | 83,40 € |
| <b>Adapter</b> , for 1 round-bottom tube 5.5 – 10 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                                 | 5820 728 000 | 83,40 € |
| <b>Adapter</b> , for 1 round-bottom tube 7 – 17 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                                   | 5820 729 007 | 83,40 € |
| <b>Adapter</b> , for 1 round-bottom tube 30 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                                       | 5820 721 006 | 84,70 € |
| <b>Adapter</b> , for 1 tube 35 mL, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.  | 5820 722 002 | 84,70 € |
| <b>Adapter</b> , for 1 Eppendorf Conical Tube 25 mL with screw cap 25 mL screw cap, for rotors for conical 50 mL tubes, 6 pcs. | 5820 734 000 | 33,00 € |
| <b>Adapter</b> , for 1 Eppendorf Conical Tube 25 mL with snap cap 25 mL snap cap, for rotors for conical 50 mL tubes, 6 pcs.   | 5820 733 004 | 32,40 € |
| <b>Adapter</b> , for 1 Eppendorf Tubes® 5.0 mL with snap cap, for Rotor FA-45-6-30 and FA-6x50, 2 pcs.                         | 5820 730 005 | 83,40 € |

## Rotor FA-20x5

Eppendorf Product Catalog p. 204

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Rotor FA-20x5</b> , aerosol-tight, for 20 Eppendorf Tubes® 5.0 mL, incl. Eppendorf QuickLock® rotor lid | 5895 130 003 | 1.433,00 € |
| <b>Adapter</b> , for 1 vessel 1.5 – 2.0 mL, for all 5.0 mL rotors, 4 pcs.                                  | 5820 768 002 | 37,00 €    |
| <b>Adapter</b> , for 1 cryo tube, for all 5.0 mL rotors, 4 pcs.  | 5820 769 009 | 37,00 €    |
| <b>Adapter</b> , for 1 HPLC tube, for all 5.0 mL rotors, 4 pcs.  | 5820 770 007 | 37,00 €    |

## Rotor FA-48x2

Eppendorf Product Catalog p. 204

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Adapter</b> , for 1 PCR tube (0.2 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.  | 5425 715 005 | 15,70 €    |
| <b>Adapter</b> , for 1 sample tube (0.5 mL, max. Ø 8 mm) or 1 Microtainer® (0.6 mL, max. Ø 8 mm), for all 1.5/2.0 mL rotors, 6 pcs. | 5425 716 001 | 15,70 €    |
| <b>Rotor FA-48x2</b> , 48 × 1.5/2 mL tubes, incl. aerosol-tight Eppendorf QuickLock® rotor lid                                      | 5895 135 005 | 1.433,00 € |
| <b>Adapter</b> , for 1 micro test tube (0.4 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.                                     | 5425 717 008 | 15,70 €    |

## Rotor F-48x15

Eppendorf Product Catalog p. 204

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Rotor F-48x15</b> , for 40 × 15 mL conical tubes, incl. 48 steel sleeves and adapters | 5895 160 000 | 1.019,00 € |

## Rotor FA-30x2

Eppendorf Product Catalog p. 205

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor FA-30x2</b> , incl. aerosol-tight Eppendorf QuickLock® rotor lid | 5895 155 006 | 1.019,00 € |

## Mobile Tables

Eppendorf Product Catalog p. 205

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Mobile table for under bench storage (low design)</b> , fits under standard 90 cm lab bench with centrifuge lid closed for Centrifuge 5804/5804 R, 5810/5810 R, 5910 R, 5910 Ri and 5920 R | 5921 000 010 | 1.613,00 € |
| <b>Mobile table for ergonomic loading and unloading (high design)</b> , fits in any empty floor space for Centrifuge 5804/5804 R, 5810/5810 R, 5910 R, 5910 Ri and 5920 R                     | 5921 000 012 | 2.016,00 € |

## Centrifuges and Rotors

### Concentrator plus and accessories

Eppendorf Product Catalog p. 206-207

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Concentrator plus basic device</b> , w/o pump, with Rotor F-45-48-11<br>230 V/50 – 60 Hz (EU)   | 5305 000 100 | 4.652,00 € |
| <b>Concentrator plus complete system</b> , with integrated diaphragm vacuum pump<br>with Rotor F-45-48-11, 230 V/50 – 60 Hz (EU)                                 | 5305 000 304 | 8.290,00 € |
| without rotor, 230 V/50 – 60 Hz (EU)   | 5305 000 509 | 8.079,00 € |
| <b>Concentrator plus complete system</b> , with integrated diaphragm vacuum pump, with connection, e.g., for a gel dryer<br>without rotor, 230 V/50 – 60 Hz (EU) | 5305 000 703 | 9.619,00 € |
| <b>Adapter</b> , for 1 PCR tube (0.2 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.   | 5425 715 005 | 15,70 €    |
| <b>Adapter</b> , for 1 sample tube (0.5 mL, max. Ø 8 mm) or 1 Microtainer® (0.6 mL, max. Ø 8 mm), for all 1.5/2.0 mL rotors, 6 pcs.                              | 5425 716 001 | 15,70 €    |
| <b>CombiSlide Adapter</b> , for 12 slides Combi Slide, for Rotors A-2-MTP and A-4-81 (MTP), 2 pcs.   | 5825 706 005 | 314,00 €   |
| <b>Adapter</b> , for 1 PCR tube 0.2 mL, for Rotor F-45-72-8, 6 pcs.  | 5425 723 008 | 15,70 €    |
| <b>Rotor spacer</b> , required for stacking rotors F-45-72-8 and F-45-48-11  | 5301 316 005 | 29,30 €    |
| <b>Work tray</b> , for semi-/unskirted PCR plates, PCR strips and 0.2 mL PCR tubes, for Rotor A-2-VC, 10 pcs.  | 0030 124 235 | Inquire    |
| <b>Frame for work tray</b> , for Rotor A-2-VC, 5 pcs.  | 0030 124 243 | 54,30 €    |
| <b>Adapter</b> , for 1 micro test tube (0.4 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.  | 5425 717 008 | 15,70 €    |

### Rotors for Concentrator plus

Eppendorf Product Catalog p. 206-207

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Rotor F-45-72-8</b> , 72 places for 0.5 mL reaction vessels  | 5490 034 007 | 415,00 €   |
| <b>Rotor F-45-70-11</b> , 70 places for 1.5/2.0 mL reaction vessels   | 5490 032 004 | 415,00 €   |
| <b>Rotor F-45-48-11</b> , 48 places for 1.5/2.0 mL reaction vessels   | 5490 030 001 | 415,00 €   |
| <b>Rotor F-45-22-17</b> , 22 places for Eppendorf Tubes® 5.0 mL (17 x 59.4 – 67.4 mm)                                   | 5490 048 008 | 525,00 €   |
| <b>Rotor F-40-36-12</b> , 36 places for 1.5 mL flat-bottom tubes (12 x 32 mm)   | 5490 040 007 | 525,00 €   |
| <b>Rotor F-45-36-15</b> , 36 places for 3.0/5.0 mL flat-bottom tubes (15 x 48 mm)                                       | 5490 035 003 | 415,00 €   |
| <b>Rotor F-45-16-20</b> , 16 places for 6.5/10.0 mL flat-bottom tubes (20 x 42 – 55 mm)                                 | 5490 043 006 | 525,00 €   |
| <b>Rotor F-40-18-19</b> , 18 places for 10.0 mL flat-bottom tubes (19 x 66 mm)  | 5490 037 006 | 525,00 €   |
| <b>Rotor F-45-12-31</b> , 12 places for 20.0 mL flat-bottom tubes (31 x 55 mm)  | 5490 044 002 | Inquire    |
| <b>Rotor F-35-8-24</b> , 8 places for 25.0 mL flat-bottom tubes (24 x 86-90 mm)   | 5490 039 009 | 525,00 €   |
| <b>Rotor F-45-24-12</b> , 24 places for 6.0/8.0 mL round-bottom tubes (12 x 67-100 mm)                                  | 5490 036 000 | 415,00 €   |
| <b>Rotor F-50-8-16</b> , 8 places for 8 mL – 12 mL round-bottom tubes (16 x 105 – 120 mm)                               | 5490 041 003 | 525,00 €   |
| <b>Rotor F-50-8-18</b> , 8 places for 10 mL – 18 mL round-bottom tubes (18 x 105-128 mm)                                | 5490 042 000 | 525,00 €   |
| <b>Rotor F-45-8-17</b> , 8 places for 15 mL conical tubes (17 x 118 – 123 mm)   | 5490 038 002 | 525,00 €   |
| <b>Rotor F-35-6-30</b> , 6 places for 15 mL conical tubes (17 x 116-123 mm) and 50 mL conical tubes (29.5 x 116-123 mm) | 5490 047 001 | 525,00 €   |
| <b>Rotor A-2-VC</b> , incl. 2 buckets   | 5490 045 009 | 1.116,00 € |



## High-Speed Centrifuge CR22N

Eppendorf Product Catalog p. 208

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>High-Speed Centrifuge CR22N</b> , refrigerated, 230 V/50/60 Hz (no plug), 30 A  | 5721 220 012 | Inquire |
| <b>CR22N Harvesting Bundle (4 x 1.5 L)</b> , refrigerated, High-Speed Centrifuge CR22N incl. rotor R9A2 and 4 x 1.5 L triangular bottles, 230 V/50/60 Hz (no plug), 30 A | 5721 261 412 | Inquire |

## Rotors for CR22N

Eppendorf Product Catalog p. 209

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Rotor R9A2</b> , fixed-angle rotor, incl. 4 x 1.5 L triangular bottles, for Centrifuge CR22N  | 5721 221 014 | Inquire |
| <b>Rotor R9A2</b> , fixed-angle rotor, incl. 4 x 1.5 L triangular bottles, for Centrifuge CR30NX | 5721 302 115 | Inquire |

## Rotors for CR22/21N

**New**

Eppendorf Product Catalog p. 210-211

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Rotor R10A3</b> , fixed-angle rotor, Harvesting Kit 6 x 0.5 L, for Centrifuge CR22N incl. 6 x 500 mL PA bottles and caps  | 5721 221 013 | Inquire |
| <b>Rotor R10A5</b> , fixed-angle rotor, suitable for 6 x 0.5 L Nalgene bottles, for Centrifuge CR22N   | 5721 221 012 | Inquire |
| <b>Rotor R10H</b> , horizontal rotor, for Centrifuge CR22N and for Centrifuge CR21N  | 5721 227 000 | Inquire |
| <b>Rotor R11D2</b> , soil dehydration rotor, for Centrifuge CR22N and for Centrifuge CR21N   | 5721 228 000 | Inquire |
| <b>Rotor R12A3</b> , fixed-angle rotor, incl. 6 x 250 mL PP bottles, for Centrifuge CR22N  | 5721 221 011 | Inquire |
| <b>Rotor R12A6</b> , fixed-angle rotor, 4 x 500 mL incl. 4 x 500PA bottles, for Centrifuge CR22N   | 5721 221 010 | Inquire |
| <b>Rotor R13A</b> , fixed-angle rotor, incl. 6 x 250 mL PP bottles, for Centrifuge CR22N   | 5721 221 009 | Inquire |
| <b>Rotor R14A</b> , fixed-angle rotor, incl. 4 x 250 mL bottles, for Centrifuge CR22N  | 5721 221 008 | Inquire |
| <b>Rotor R15A</b> , fixed-angle rotor, High-Speed Pelleting Kit 10 x 15/50 mL for Centrifuge CR22N incl. 100 x 15/50 high-speed conical tubes                                      | 5721 221 007 | Inquire |
| <b>Rotor R18A</b> , fixed-angle rotor, incl. 100 x 50 mL high-speed conical tubes, for Centrifuge CR22N  | 5721 221 006 | Inquire |
| <b>Rotor R20A2</b> , fixed-angle rotor, High-Speed Pelleting Kit 8 x 50 mL for Centrifuge CR22N incl. 8 x 50 mL PA tube and M-PP lid   | 5721 302 106 | Inquire |
| <b>Rotor R21A2</b> , fixed-angle rotor, incl. 6 x 30 mL PA bottles, for Centrifuge CR22N   | 5721 221 004 | Inquire |
| <b>Rotor R22A4</b> , fixed-angle rotor, for 30 x 2 mL mico tubes, for Centrifuge CR22N   | 5721 221 003 | Inquire |
| <b>Rotor R22A6</b> , fixed-angle rotor, incl. 12 x 10 mL PA bottles, for Centrifuge CR22N  | 5721 221 002 | Inquire |
| <b>Rotor R27A</b> , fixed-angle rotor, incl. 8 x 50 mL PP bottles, for Centrifuge CR22N  | 5721 202 101 | Inquire |
| <b>Rotor R3S</b> , swing-bucket rotor, incl. 4 buckets, for 500 mL bottle, for Centrifuge CR22N  | 5721 224 003 | Inquire |
| <b>Rotor R4SS</b> , swing-bucket rotor, Low-speed Pelleting Kit 15/50 mL incl. adapter for 15/50 mL Eppendorf Conical Tubes (4 pcs. each), no tubes included, for Centrifuge CR22N | 5721 264 002 | Inquire |
| <b>Rotor R5S4</b> , swing-bucket rotor, for up to 12 Microplates, for Centrifuge CR22N   | 5721 224 001 | Inquire |
| <b>Rotor R8S</b> , swing-bucket rotor, incl. 4 buckets, for 15/50 mL Eppendorf Conical Tubes, for Centrifuge CR22N   | 5721 202 400 | Inquire |
| <b>Rotor R9A</b> , fixed-angle rotor, incl. 4 x 1 L PP bottles, for Centrifuge CR22N   | 5721 221 015 | Inquire |
| <b>Rotor R9A2</b> , fixed-angle rotor, incl. 4 x 1.5 L triangular bottles, for Centrifuge CR22N  | 5721 221 014 | Inquire |

## Centrifuges and Rotors

### High-Speed Centrifuge CR30NX

Eppendorf Product Catalog p. 212

|   | Order no.    | Price   |
|---|--------------|---------|
| <b>High-Speed Centrifuge CR30NX</b> , refrigerated, 230 V/50/60 Hz (no plug), 30 A, 1-phase   | 5721 300 012 | Inquire |
| <b>CR30NX Harvesting Bundle (4 x 1.5 L)</b> , refrigerated, High Speed Centrifuge CR30NX incl. rotor R9A2 and 4 x 1.5 L triangular bottles, 230 V/50/60 Hz (no plug), 30 A, 1-phase | 5721 351 512 | Inquire |

### Rotors for CR30NX

Eppendorf Product Catalog p. 213

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Rotor R30AT</b> , fixed-angle rotor, Super-Speed Pelleting Kit 8 x 50 mL for Centrifuge CR30NX, incl. 8 x 50PP bottle | 5721 302 101 | Inquire |
| <b>Rotor R25ST</b> , swing-bucket rotor, Super-Speed Pelleting Kit 6 x 40 mL for Centrifuge CR30NX, incl. 50 x 40PA tube | 5721 302 400 | Inquire |

### Rotors for CR30NX

**New**

Eppendorf Product Catalog p. 214-215

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Rotor R10A3</b> , fixed-angle rotor, Harvesting Kit 6 x 0.5 L for Centrifuge CR30NX, incl. 6 x 500 mL PA bottles and caps   | 5721 302 112 | Inquire |
| <b>Rotor R10A5</b> , fixed-angle rotor, suitable for 6 x 0.5 L Nalgene bottles, for Centrifuge CR30NX  | 5721 302 113 | Inquire |
| <b>Rotor R10C2</b> , continuous flow rotor, Continuous Flow Kit for Centrifuge CR30NX, rotor nominal capacity: 3,200 mL  | 5721 352 601 | Inquire |
| <b>Rotor R10H</b> , horizontal rotor, for plates, for Centrifuge CR30NX  | 5721 302 700 | Inquire |
| <b>Rotor R11D2</b> , soil dehydration rotor, 4 x 100 mL, for Centrifuge CR30NX   | 5721 302 800 | Inquire |
| <b>Rotor R16A2</b> , fixed-angle rotor, incl. 100 x 15/50 mL high-speed conical tubes, for Centrifuge CR30NX   | 5721 302 110 | Inquire |
| <b>Rotor R16A3</b> , fixed-angle rotor, for 6 x 250 mL bottles, including 6 x 250 mL PP bottles, for Centrifuge CR30NX   | 5721 302 111 | Inquire |
| <b>Rotor R17A</b> , fixed-angle rotor, incl. 14 x 50 mL high-speed conical tubes, for Centrifuge CR30NX  | 5721 302 109 | Inquire |
| <b>Rotor R18C2</b> , continuous flow rotor, batch unit, for Centrifuge CR30NX  | 5721 302 600 | Inquire |
| <b>Rotor R19A</b> , fixed-angle rotor, incl. 6 x 80 mL PA bottles, for Centrifuge CR30NX   | 5721 302 108 | Inquire |
| <b>Rotor R19A2</b> , fixed-angle rotor, High-Speed Pelleting Kit 8 x 50 mL conical tubes for Centrifuge CR30NX, incl. 100 x 50 mL high-speed conical tubes plus caps | 5721 302 107 | Inquire |
| <b>Rotor R20A2</b> , fixed-angle rotor, High-Speed Pelleting Kit 8 x 50 mL, conical tubes for Centrifuge CR30NX incl. 8 x 50 mL PA tube and M-PP lid                 | 5721 221 005 | Inquire |
| <b>Rotor R21A2</b> , fixed-angle rotor, incl. 6 x 30 mL PA bottles, for Centrifuge CR30NX  | 5721 302 105 | Inquire |
| <b>Rotor R22A4</b> , fixed-angle rotor, for 30 x 2 mL micro tubes, for Centrifuge CR30NX   | 5721 302 104 | Inquire |
| <b>Rotor R22A6</b> , fixed-angle rotor, incl. 12 x 10 mL PA bottles, for Centrifuge CR30NX   | 5721 302 103 | Inquire |
| <b>Rotor R25ST</b> , swing-bucket rotor, Super-Speed Pelleting Kit 6 x 40 mL for Centrifuge CR30NX, incl. 50 x 40PA tube   | 5721 302 400 | Inquire |
| <b>Rotor R27A</b> , fixed-angle rotor, incl. 8 x 50 mL PP bottles, for Centrifuge CR30NX   | 5721 302 102 | Inquire |
| <b>Rotor R30AT</b> , fixed-angle rotor, Super-Speed Pelleting Kit 8 x 50 mL for Centrifuge CR30NX, incl. 8 x 50PP bottle   | 5721 302 101 | Inquire |
| <b>Rotor R4SS</b> , swing-bucket rotor, incl. 4 buckets, for 15/50 mL Eppendorf Conical Tubes, for Centrifuge CR30NX   | 5721 352 402 | Inquire |
| <b>Rotor R5S4</b> , swing-bucket rotor, for up to 12 microplates, for Centrifuge CR30NX  | 5721 302 401 | Inquire |
| <b>Rotor R7A</b> , fixed-angle rotor, incl. 6 x 1 L PP bottles, for Centrifuge CR30NX  | 5721 352 100 | Inquire |
| <b>Rotor R9A</b> , fixed-angle rotor, incl. 4 x 1 L bottles, for Centrifuge CR30NX   | 5721 302 114 | Inquire |
| <b>Rotor R9A2</b> , fixed-angle rotor, incl. 4 x 1.5 L triangular bottles, for Centrifuge CR30NX   | 5721 302 115 | Inquire |

## Ultracentrifuge CP-NX Series

Eppendorf Product Catalog p. 216

|   | Order no.    | Price   |
|---|--------------|---------|
| <b>Centrifuge CP100NX</b> , refrigerated, touch interface, 230 V/50/60 Hz (no plug), 16 A | 5720 110 012 | Inquire |
| <b>Centrifuge CP90NX</b> , refrigerated, touch interface, 230 V/50/60 Hz (no plug), 16 A  | 5720 109 012 | Inquire |
| <b>Centrifuge CP80NX</b> , refrigerated, touch interface, 230 V/50/60 Hz (no plug), 16 A  | 5720 108 012 | Inquire |

## Rotors for CP-NX

**New**

Eppendorf Product Catalog p. 216-217

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Rotor P21A2</b> , fixed-angle rotor, incl. 6 x 230 mL PA bottle, for CP-NX Series         | 5720 211 205 | Inquire |
| <b>Rotor P28S</b> , swing-bucket rotor, incl. 100 x 40 mL PA tubes, for CP-NX Series         | 5720 214 005 | Inquire |
| <b>Rotor P32ST</b> , swing-bucket rotor, for CP-NX Series                                    | 5720 214 003 | Inquire |
| <b>Rotor P40ST</b> , swing-bucket rotor, incl. 100 x 13 mL PA tubes, for CP-NX Series        | 5720 214 002 | Inquire |
| <b>Rotor P45AT</b> , fixed-angle rotor, for 6 x 94 mL tubes, for CP-NX Series                | 5720 211 109 | Inquire |
| <b>Rotor P50A3</b> , fixed-angle rotor, incl. 300 x 1.5 mL micro tubes (C), for CP-NX Series | 5720 211 108 | Inquire |
| <b>Rotor P50AT2</b> , fixed-angle rotor, for 12 x 40 mL tubes, for CP-NX Series              | 5720 211 107 | Inquire |
| <b>Rotor P55ST2</b> , swing-bucket rotor, incl. 100 x 5 mL PA tubes, for CP-NX Series        | 5720 214 001 | Inquire |
| <b>Rotor P70AT</b> , fixed-angle rotor, for 8 x 40 mL tubes, for CP-NX Series                | 5720 211 105 | Inquire |
| <b>Rotor P70AT2</b> , fixed-angle rotor, for 12 x 12 mL tubes, for CP-NX Series              | 5720 211 104 | Inquire |
| <b>Rotor P80AT</b> , fixed-angle rotor, for 8 x 12 mL tubes                                  | 5720 211 103 | Inquire |
| <b>Rotor P90AT</b> , fixed-angle rotor, for 8 x 12 mL tubes                                  | 5720 211 102 | Inquire |
| <b>Rotor P100AT2</b> , fixed-angle rotor, for 8 x 6.5 mL tubes, for Centrifuge CP-NX         | 5720 211 101 | Inquire |
| <b>Rotor P100VT</b> , vertical rotor, for 8 x 5 mL tubes, for CP-NX Series                   | 5720 213 001 | Inquire |

## Micro-Ultracentrifuges - CS150FNX and CS150NX

Eppendorf Product Catalog p. 217

|   | Order no.    | Price   |
|---|--------------|---------|
| <b>Centrifuge CS150FNX</b> , refrigerated, floorstanding, touch interface, 230 V/50/60 Hz (EU), 8 A | 5720 121 511 | Inquire |
| <b>Centrifuge CS150NX</b> , refrigerated, benchtop, touch interface, 230 V/50/60 Hz (EU), 8 A       | 5720 131 511 | Inquire |

## Rotors for CS-FNX/NX

**New**

Eppendorf Product Catalog p. 218-219

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Rotor S50A</b> , fixed-angle rotor, incl. 20 x 25 mL thick-walled tubes, for CS-FNX/NX                          | 5720 221 012 | Inquire |
| <b>Rotor S50ST</b> , swing-bucket rotor, incl. 50 x 7 mL PA tubes, for CS-FNX/NX                                   | 5720 224 002 | Inquire |
| <b>Rotor S52ST</b> , swing-bucket rotor, incl. 100 x 5 mL PA tube, for CS-FNX/NX                                   | 5720 224 001 | Inquire |
| <b>Rotor S55A2</b> , fixed-angle rotor, incl. 300 x 1.5 mL micro tubes, for CS-FNX/NX                              | 5720 221 011 | Inquire |
| <b>Rotor S55S</b> , swing-bucket rotor, incl. 100 x 2.2 mL PA tubes, for CS-FNX/NX                                 | 5720 224 000 | Inquire |
| <b>Rotor S58A</b> , fixed-angle rotor, incl. 20 x 10 mL PC thick-walled tubes, for CS-FNX/NX                       | 5720 221 010 | Inquire |
| <b>Rotor S70AT</b> , fixed-angle rotor, incl. 100 x 0.5 mL PC tubes, for CS-FNX/NX                                 | 5720 221 009 | Inquire |
| <b>Rotor S80AT2</b> , fixed-angle rotor, incl. 100 x 0.5 mL PC tubes, for CS-FNX/NX                                | 5720 221 007 | Inquire |
| <b>Rotor S80AT3</b> , fixed-angle rotor, incl. 100 x 6 mL thick-walled tubes, for CS-FNX/NX Series, for Japan only | 5720 221 008 | Inquire |
| <b>Rotor S100AT3</b> , fixed-angle rotor, incl. 100 x 0.23 mL PC tubes, for CS-FNX/NX                              | 5720 221 005 | Inquire |
| <b>Rotor S100AT4</b> , fixed-angle rotor, incl. 100 x 3 mL PC tubes, for CS-FNX/NX                                 | 5720 221 006 | Inquire |
| <b>Rotor S100NT</b> , neo-angle rotor, incl. 100 x 4 mL PA seal tube, for CS-FNX/NX                                | 5720 222 002 | Inquire |
| <b>Rotor S110AT</b> , fixed-angle rotor, incl. 100 x 4 mL PC tubes, for CS-FNX/NX                                  | 5720 221 004 | Inquire |
| <b>Rotor S120AT2</b> , fixed-angle rotor, incl. 100 x 1 mL PC tubes, for CS-FNX/NX                                 | 5720 221 002 | Inquire |
| <b>Rotor S120AT3</b> , fixed-angle rotor, incl. 100 x 0.5 mL PC tubes, for CS-FNX/NX                               | 5720 221 003 | Inquire |
| <b>Rotor S120NT</b> , neo-angle rotor, incl. 100 x 2 mL PA seal tubes, for CS-FNX/NX                               | 5720 222 001 | Inquire |
| <b>Rotor S120VT</b> , vertical rotor, incl. 100 x 2 mL PA seal tubes, for CS-FNX/NX                                | 5720 223 000 | Inquire |
| <b>Rotor S140AT</b> , fixed-angle rotor, incl. 100 x 1 mL PC tubes, for CS-FNX/NX                                  | 5720 221 001 | Inquire |
| <b>Rotor S140NT</b> , neo-angle rotor, incl. 100 x 2 mL PA Seal tubes, for CS-FNX/NX                               | 5720 222 000 | Inquire |
| <b>Rotor S150AT</b> , fixed-angle rotor, incl. 100 x 1 mL PC tubes, for CS-FNX/NX                                  | 5720 221 000 | Inquire |

## PCR

## Mastercycler® X50

Eppendorf Product Catalog p. 226-227

|  | Order no.    | Price       |
|--|--------------|-------------|
| <b>Mastercycler® X50s</b> , 110 – 230 V/50 – 60 Hz (EU/JP/TW/US/South America), silver block, 96-well or 0.1/0.2 mL tubes, with touchscreen interface    | 6311 000 010 | 10.530,00 € |
| <b>Mastercycler® X50a</b> , 110 – 230 V/50 – 60 Hz (EU/JP/TW/US/South America), aluminium block, 96-well or 0.1/0.2 mL tubes, with touchscreen interface | 6313 000 018 | 9.603,00 €  |
| <b>Mastercycler® X50h</b> , 110 – 230 V/50 – 60 Hz (EU/JP/TW/US/South America), aluminium block, 384-well, high pressure lid, with touchscreen interface | 6316 000 019 | 10.530,00 € |
| <b>Mastercycler® X50i</b> , 110 – 230 V/50 – 60 Hz (EU/JP/TW/US/South America), silver block, 96-well or 0.1/0.2 mL tubes, eco unit                      | 6301 000 012 | 8.655,00 €  |
| <b>Mastercycler® X50l</b> , 110 – 230 V/50 – 60 Hz (EU/JP/TW/US/South America), aluminium block, 96-well or 0.1/0.2 mL tubes, eco unit                   | 6303 000 010 | 8.071,00 €  |
| <b>Mastercycler® X50t</b> , 110 – 230 V/50 – 60 Hz (EU/JP/TW/US/South America), aluminium block, 384-well, high pressure lid, eco unit                   | 6306 000 010 | 9.241,00 €  |

## CycleManager X50 Software

Eppendorf Product Catalog p. 227

|                                  | Order no.    | Price    |
|----------------------------------|--------------|----------|
| <b>CycleManager X50 Software</b> | 6349 000 014 | 268,00 € |

## Mastercycler® nexus

Eppendorf Product Catalog p. 228

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Mastercycler® nexus gradient</b> , 230 V/50 – 60 Hz (EU)            | 6331 000 017 | 9.104,00 € |
| <b>Mastercycler® nexus</b> , 230 V/50 – 60 Hz (EU)                     | 6333 000 014 | 7.606,00 € |
| <b>Mastercycler® nexus gradient eco</b> , 230 V/50 – 60 Hz (EU)        | 6334 000 018 | 7.406,00 € |
| <b>Mastercycler® nexus eco</b> , 230 V/50 – 60 Hz (EU)                 | 6332 000 010 | 5.872,00 € |
| <b>CAN bus connecting cable</b> , 150 cm                               | 5341 611 000 | 75,20 €    |
| <b>CAN bus cable</b> , connecting cable between 2 thermomodules, 50 cm | 5341 612 006 | 58,50 €    |

## Mastercycler® nexus X2

Eppendorf Product Catalog p. 229

|   | Order no.    | Price       |
|---|--------------|-------------|
| <b>Mastercycler® nexus GX2</b> , 230 V/50 – 60 Hz (EU)  | 6336 000 015 | 10.910,00 € |
| <b>Mastercycler® nexus GX2e</b> , 230 V/50 – 60 Hz (EU) | 6338 000 012 | 8.273,00 €  |
| <b>Mastercycler® nexus X2</b> , 230 V/50 – 60 Hz (EU)   | 6337 000 019 | 9.112,00 €  |
| <b>Mastercycler® nexus X2e</b> , 230 V/50 – 60 Hz (EU)  | 6339 000 016 | 6.954,00 €  |

## Eppendorf twin.tec® PCR Plates

Eppendorf Product Catalog p. 230

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Eppendorf twin.tec® Trace PCR Plate 96</b> , skirted, 150 µL |              |          |
| <b>New</b> PCR clean, colorless, 25 plates                      | 0030 129 768 | 142,00 € |
| <b>New</b> PCR clean, crystal blue, 25 plates                   | 0030 129 776 | 142,00 € |
| <b>New</b> PCR clean, fuchsia, 25 plates                        | 0030 129 784 | 142,00 € |
| <b>Eppendorf twin.tec® PCR Plate 96</b> , skirted, 150 µL       |              |          |
| PCR clean, colorless, 25 plates                                 | 0030 128 648 | 142,00 € |
| PCR clean, yellow, 25 plates                                    | 0030 128 656 | 142,00 € |
| PCR clean, green, 25 plates                                     | 0030 128 664 | 142,00 € |
| PCR clean, blue, 25 plates                                      | 0030 128 672 | 142,00 € |
| PCR clean, red, 25 plates                                       | 0030 128 680 | 142,00 € |
| Forensic DNA Grade, colorless, 10 plates, individually wrapped  | 0030 129 601 | 70,70 €  |

## Eppendorf twin.tec® PCR Plates

Eppendorf Product Catalog p. 230

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Eppendorf twin.tec® Trace PCR Plate 96</b> , skirted, 150 µL                   |              |          |
| <b>Eppendorf twin.tec® PCR Plate 96</b> , semi-skirted, 250 µL                    |              |          |
| PCR clean, colorless, 25 plates   | 0030 128 575 | 142,00 € |
| PCR clean, yellow, 25 plates  | 0030 128 583 | 142,00 € |
| PCR clean, green, 25 plates   | 0030 128 591 | 142,00 € |
| PCR clean, blue, 25 plates  | 0030 128 605 | 142,00 € |
| PCR clean, red, 25 plates   | 0030 128 613 | 142,00 € |
| Forensic DNA Grade, colorless, 10 plates, individually wrapped                    | 0030 129 610 | 70,70 €  |
| <b>Eppendorf twin.tec® PCR Plate 96, low profile</b> , unskirted, 150 µL          |              |          |
| PCR clean, colorless, 20 plates   | 0030 133 307 | 90,10 €  |
| PCR clean, blue, 20 plates  | 0030 133 331 | 90,10 €  |
| <b>Eppendorf twin.tec® PCR Plate 96</b> , unskirted, 250 µL                       |              |          |
| PCR clean, colorless, 20 plates   | 0030 133 366 | 90,10 €  |
| PCR clean, blue, 20 plates  | 0030 133 390 | 90,10 €  |
| <b>Eppendorf twin.tec® PCR Plate 96, divisible</b> , unskirted, divisible, 150 µL |              |          |
| PCR clean, colorless, 20 plates   | 0030 133 358 | 163,00 € |
| PCR clean, blue, 20 plates  | 0030 133 382 | 163,00 € |
| <b>Eppendorf twin.tec® PCR Plate 96, divisible</b> , unskirted, divisible, 250 µL |              |          |
| PCR clean, colorless, 20 plates   | 0030 133 374 | 163,00 € |
| PCR clean, blue, 20 plates  | 0030 133 404 | 163,00 € |
| <b>Eppendorf twin.tec® microbiology PCR Plate 96</b> , skirted, 150 µL            |              |          |
| PCR clean, colorless, 10 plates   | 0030 129 300 | 85,40 €  |
| PCR clean, blue, 10 plates  | 0030 129 318 | 85,40 €  |
| <b>Eppendorf twin.tec® microbiology PCR Plate 96</b> , semi-skirted, 250 µL       |              |          |
| PCR clean, colorless, 10 plates   | 0030 129 326 | 85,40 €  |
| PCR clean, blue, 10 plates  | 0030 129 334 | 85,40 €  |
| <b>Eppendorf twin.tec® microbiology PCR Plate 384</b> , skirted, 45 µL            |              |          |
| PCR clean, colorless, 10 plates   | 0030 129 342 | 127,00 € |
| PCR clean, blue, 10 plates  | 0030 129 350 | 127,00 € |
| <b>Eppendorf twin.tec® PCR Plate 384</b> , skirted, 45 µL                         |              |          |
| PCR clean, colorless, 25 plates   | 0030 128 508 | 219,00 € |
| PCR clean, yellow, 25 plates  | 0030 128 516 | 219,00 € |
| PCR clean, green, 25 plates   | 0030 128 524 | 219,00 € |
| PCR clean, blue, 25 plates  | 0030 128 532 | 219,00 € |
| PCR clean, red, 25 plates   | 0030 128 540 | 219,00 € |
| Forensic DNA Grade, colorless, 10 plates, individually wrapped                    | 0030 129 628 | 70,70 €  |

## Eppendorf twin.tec® PCR Plates

Eppendorf Product Catalog p. 231

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Eppendorf twin.tec® PCR Plate 96</b> , skirted, 150 µL      |              |            |
| PCR clean, colorless, 300 plates                               | 0030 128 770 | 1.514,00 € |
| PCR clean, green, 300 plates                                   | 0030 128 796 | 1.514,00 € |
| PCR clean, yellow, 300 plates                                  | 0030 128 788 | 1.514,00 € |
| PCR clean, blue, 300 plates                                    | 0030 128 842 | 1.514,00 € |
| PCR clean, red, 300 plates                                     | 0030 128 850 | 1.514,00 € |
| <b>Eppendorf twin.tec® PCR Plate 96</b> , semi-skirted, 250 µL |              |            |
| PCR clean, colorless, 300 plates                               | 0030 128 869 | 1.514,00 € |
| <b>Eppendorf twin.tec® PCR Plate 384</b> , 40 µL               |              |            |
| PCR clean, colorless, 300 plates                               | 0030 128 931 | 2.245,00 € |
| PCR clean, yellow, 300 plates                                  | 0030 128 940 | 2.245,00 € |
| PCR clean, green, 300 plates                                   | 0030 128 958 | 2.245,00 € |
| PCR clean, blue, 300 plates                                    | 0030 128 966 | 2.245,00 € |
| PCR clean, red, 300 plates                                     | 0030 128 974 | 2.245,00 € |

## PCR

## Eppendorf twin.tec® PCR Plates LoBind®

Eppendorf Product Catalog p. 231

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Eppendorf twin.tec® PCR Plate 96 LoBind®, semi-skirted, 250 µL</b> |              |            |
| PCR clean, colorless, 25 plates                                       | 0030 129 504 | 172,00 €   |
| <b>Eppendorf twin.tec® PCR Plate 96 LoBind®, skirted, 150 µL</b>      |              |            |
| PCR clean, colorless, 25 plates                                       | 0030 129 512 | 172,00 €   |
| PCR clean, red, 25 plates   | 0030 129 598 | 187,00 €   |
| PCR clean, blue, 25 plates  | 0030 129 580 | 187,00 €   |
| PCR clean, yellow, 25 plates  | 0030 129 679 | 187,00 €   |
| PCR clean, green, 25 plates   | 0030 129 660 | 187,00 €   |
| <b>Eppendorf twin.tec® PCR Plate 96 LoBind®, skirted</b>              |              |            |
| PCR clean, yellow, 300 plates   | 0030 129 563 | 1.749,00 € |
| PCR clean, orange, 300 plates   | 0030 129 571 | 1.749,00 € |
| PCR clean, green, 300 plates  | 0030 129 555 | 1.749,00 € |
| <b>Eppendorf twin.tec® PCR Plate 384 LoBind®, skirted, 45 µL</b>      |              |            |
| PCR clean, colorless, 25 plates                                       | 0030 129 547 | 264,00 €   |

Eppendorf twin.tec® *real-time* PCR plates

Eppendorf Product Catalog p. 232

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Eppendorf twin.tec® 96 <i>real-time</i> PCR Plate, skirted, 150 µL</b>                |              |          |
| PCR clean, blue, wells white, 25 plates  | 0030 132 505 | 181,00 € |
| PCR clean, white, wells white, 25 plates   | 0030 132 513 | 181,00 € |
| <b>Eppendorf twin.tec® 96 <i>real-time</i> PCR Plate, semi-skirted, 250 µL</b>           |              |          |
| PCR clean, blue, wells white, 25 plates  | 0030 132 530 | 181,00 € |
| PCR clean, white, wells white, 25 plates   | 0030 132 548 | 181,00 € |
| <b>Eppendorf twin.tec® 96 <i>real-time</i> PCR Plate, unskirted, low profile, 150 µL</b> |              |          |
| PCR clean, white, wells white, 20 plates   | 0030 132 700 | 115,00 € |
| PCR clean, blue, wells white, 20 plates  | 0030 132 718 | 115,00 € |
| <b>Eppendorf twin.tec® 96 <i>real-time</i> PCR Plate, skirted, 150 µL</b>                |              |          |
| Forensic DNA Grade, white, wells white, 10 plates, individually wrapped                  | 0030 129 636 | 70,70 €  |
| <b>Eppendorf twin.tec® 96 <i>real-time</i> PCR Plate, semi-skirted, 250 µL</b>           |              |          |
| Forensic DNA Grade, white, wells white, 10 plates, individually wrapped                  | 0030 129 644 | 70,70 €  |

## Eppendorf PCR Tubes

Eppendorf Product Catalog p. 233

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Eppendorf PCR Tube Strips, 0,1 mL</b>  |              |          |
| PCR clean, without lids (10 × 12 strips)  | 0030 124 804 | 142,00 € |
| PCR clean, with Cap Strips, domed (10 × 12 strips)                              | 0030 124 812 | 192,00 € |
| PCR clean, with Cap Strips, flat (10 × 12 strips)                               | 0030 124 820 | 209,00 € |
| <b>Eppendorf PCR Tubes, 0,2 mL</b>  |              |          |
| PCR clean, colorless, 1.000 pcs.  | 0030 124 332 | 153,00 € |
| Forensic DNA Grade, colorless, 500 pcs. (5 bags × 100 pcs.)                     | 0030 124 707 | 86,60 €  |
| <b>Eppendorf PCR Tube Strips, 0,2 mL</b>  |              |          |
| PCR clean, 8-tube strips, 120 pcs. (960 tubes)                                  | 0030 124 359 | 186,00 € |
| <b>Eppendorf PCR Tubes, 0,5 mL, PCR clean, with hinged lid, clear, 500 pcs.</b> | 0030 124 537 | 52,40 €  |

## Eppendorf Fast PCR Tube Strips

Eppendorf Product Catalog p. 233

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Eppendorf Fast PCR Tube Strips</b>                              |              |          |
| 0,1 mL, without lids, PCR clean, 120 pcs. (10 × 12 pcs.)           | 0030 124 901 | 170,00 € |
| 0,1 mL, with cap strips, flat, PCR clean, 120 pcs. (10 × 12 pcs.)  | 0030 124 910 | 246,00 € |
| 0,1 mL, with cap strips, domed, PCR clean, 120 pcs. (10 × 12 pcs.) | 0030 124 928 | 229,00 € |

## PCR cap strips

Eppendorf Product Catalog p. 233

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Cap Strips</b> , 8-strips, domed, 120 pcs. (10 bags × 12 pcs.)                      | 0030 124 839 | 62,20 €  |
| <b>Cap Strips</b> , 8-strips, flat, 120 pcs. (10 × 12 pcs.)                            | 0030 124 847 | 78,60 €  |
| <i>real-time</i> <b>PCR Tube Strips</b> , without caps, 120 pcs. (10 × 12 pcs.)        | 0030 132 882 | 173,00 € |
| <b>Masterclear® Cap Strips and real-time PCR Tube Strips</b> , 120 pcs. (10 × 12 pcs.) | 0030 132 890 | 274,00 € |
| <b>Masterclear® Cap Strip</b> , 120 pcs. (10 × 12 pcs.)                                | 0030 132 874 | 105,00 € |

## PCR-Cooler

Eppendorf Product Catalog p. 234

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>PCR-Cooler 0.2 mL Starter Set</b> , 1 pink, 1 blue | 3881 000 015 | 114,00 € |
| <b>PCR-Cooler 0.2 mL</b> , Pink                       | 3881 000 023 | 62,30 €  |
| <b>PCR-Cooler 0.2 mL</b> , Blue                       | 3881 000 031 | 62,30 €  |

## PCR Rack

Eppendorf Product Catalog p. 234

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>PCR-Rack</b> , assorted colors, 10 pcs. | 0030 124 545 | 158,00 € |

## PCR films and foils

Eppendorf Product Catalog p. 234

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Eppendorf Heat Sealing Film</b> , PCR clean, 100 pcs.                    | 0030 127 838 | 193,00 € |
| <b>Eppendorf Heat Sealing Foil</b> , PCR clean, 100 pcs.                    | 0030 127 854 | 136,00 € |
| <b>Eppendorf PCR Film</b> , self-adhesive, PCR clean, 100 pcs.              | 0030 127 781 | 159,00 € |
| <b>Eppendorf PCR Foil</b> , self-adhesive, PCR clean, 100 pcs.              | 0030 127 790 | 159,00 € |
| <b>Masterclear® real-time PCR Film</b> , self-adhesive, PCR clean, 100 pcs. | 0030 132 947 | 290,00 € |

## HeatSealer

Eppendorf Product Catalog p. 235

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>HeatSealer S100</b> , incl. PCR Plate Adapter, 100 – 230 V/50 – 60 Hz (EU) | 5391 000 001 | 2.847,00 € |
| <b>HeatSealer S200</b> , incl. PCR Plate Adapter, 100 – 230 V/50 – 60 Hz (EU) | 5392 000 005 | 4.087,00 € |
| <b>Plate Adapter low profile</b> , for HeatSealer S100                        | 5391 070 018 | 346,00 €   |
| <b>Plate Adapter high profile</b> , for HeatSealer S100                       | 5391 070 034 | 346,00 €   |
| <b>96-well-PCR Plate Adapter</b> , for HeatSealer S100                        | 5391 070 026 | 346,00 €   |
| <b>Plate Adapter high profile</b> , for HeatSealer S200                       | 5392 070 011 | 346,00 €   |
| <b>Plate Adapter low profile</b> , for HeatSealer S200                        | 5392 070 020 | 346,00 €   |
| <b>96-well Plate Adapter</b> , for HeatSealer S200                            | 5392 070 038 | 346,00 €   |
| <b>Eppendorf Heat Sealing Film</b> , PCR clean, 100 pcs.                      | 0030 127 838 | 193,00 €   |
| <b>Eppendorf Heat Sealing Foil</b> , PCR clean, 100 pcs.                      | 0030 127 854 | 136,00 €   |

## Freezer

## Eppendorf ULT Upright Freezers

Eppendorf Product Catalog p. 246

|  | Order no.  | Price       |
|--|------------|-------------|
| <b>CryoCube® F740h</b> , 740 L, classic interface, VIP, green cooling liquids, air cooling, handle left side, 3 shelves, 230 V/50 Hz (EU)  | F740340011 | 14.900,00 € |
| <b>CryoCube® F740h</b> , 740 L, classic interface, VIP, green cooling liquids, air cooling, handle right side, 3 shelves, 230 V/50 Hz (EU)   | F740340021 | 15.650,00 € |
| <b>CryoCube® F740h</b> , 740 L, classic interface, VIP, green cooling liquids, air cooling, handle left side, 5 shelves, 230 V/50 Hz (EU)  | F740340031 | 14.900,00 € |
| <b>CryoCube® F740h</b> , 740 L, classic interface, VIP, green cooling liquids, air cooling, handle right side, 5 shelves, 230 V/50 Hz (EU)   | F740340041 | 15.650,00 € |
| <b>CryoCube® F740hi</b> , 740 L, touchscreen interface, VIP, green cooling liquids, air cooling, handle left side, 3 shelves, 230 V/50 Hz (EU)   | F740320011 | 16.490,00 € |
| <b>CryoCube® F740hi</b> , 740 L, touchscreen interface, VIP, green cooling liquids, air cooling, handle right side, 3 shelves, 230 V/50 Hz (EU)  | F740320021 | 17.270,00 € |
| <b>CryoCube® F740hi</b> , 740 L, touchscreen interface, VIP, green cooling liquids, air cooling, handle left side, 5 shelves, 230 V/50 Hz (EU)   | F740320031 | 16.490,00 € |
| <b>CryoCube® F740hi</b> , 740 L, touchscreen interface, VIP, green cooling liquids, air cooling, handle right side, 5 shelves, 230 V/50 Hz (EU)  | F740320041 | 17.270,00 € |
| <b>CryoCube® F740hiw</b> , 740 L, touchscreen interface, VIP, green cooling liquids, water cooling, handle left side, 3 shelves, 230 V/50 Hz (EU)  | F740320111 | 16.960,00 € |
| <b>CryoCube® F740hiw</b> , 740 L, touchscreen interface, VIP, green cooling liquids, water cooling, handle left side, 5 shelves, 230 V/50 Hz (EU)  | F740320131 | 16.960,00 € |
| <b>CryoCube® F570n</b> , 570 L, classic interface, green cooling liquids, air cooling, handle left side, 5 shelves, 230 V/50 Hz (EU/GB)  | F571300031 | 12.040,00 € |
| <b>CryoCube® F570h</b> , 570 L, classic interface, VIP, green cooling liquids, air cooling, handle left side, 5 shelves, 230 V/50 Hz (EU/GB)   | F571340031 | 13.060,00 € |
| <b>CryoCube® F570hw</b> , 570 L, classic interface, VIP, green cooling liquids, water cooling, handle left side, 5 shelves, 230 V/50 Hz (EU/GB)  | F571340131 | 14.130,00 € |
| <b>Innova® U535</b> , 535 L, classic interface, VIP, classic cooling liquids, air cooling, handle left side, 3 shelves, 230 V/50 Hz (EU/GB)  | U9430-0001 | 12.860,00 € |
| <b>CryoCube® F440h</b> , 440 L, classic interface, VIP, green cooling liquids, air cooling, handle left side, 5 shelves, 230 V/50 Hz (EU/GB)   | F440340031 | 12.090,00 € |
| <b>CryoCube® F440n</b> , 440 L, classic interface, green cooling liquids, air cooling, handle left side, 5 shelves, 230 V/50 Hz (EU/GB)  | F440300031 | 11.090,00 € |
| <b>CryoCube® F101h</b> , 101 L, classic interface, VIP, green cooling liquids, air cooling, handle right side, 2 shelves, 230 V/50 Hz (EU), chart recorder, CO <sub>2</sub> backup system, LN <sub>2</sub> backup system | F101340061 | 8.490,00 €  |

## Eppendorf ULT Chest Freezers

Eppendorf Product Catalog p. 246

|  | Order no.  | Price       |
|--|------------|-------------|
| <b>CryoCube® FC660h</b> , 660 L, classic interface, VIP, green cooling liquids, air cooling, 230 V/50 Hz (EU/GB) | F660320001 | 14.350,00 € |

CO<sub>2</sub> and LN<sub>2</sub> Back-up System for Innova® U101, U360, U535, U725/ -G, C585, C760, CryoCube® F101h and F740 series

Eppendorf Product Catalog p. 247

|  | Order no.  | Price      |
|--|------------|------------|
| <b>CO<sub>2</sub> back-up system</b> , 100 V/50 – 60 Hz (JP)   | F652999005 | Inquire    |
| <b>CO<sub>2</sub> back-up system</b> , 120 – 220 V/60 Hz, Innova® and Green 'G', ('narrow')                    | U9043-0002 | Inquire    |
| <b>LN<sub>2</sub> back-up system</b> , 100 V/50 – 60 Hz (JP), Innova® Green 'G' / F740, LN <sub>2</sub> Backup | F652999006 | Inquire    |
| <b>LN<sub>2</sub> back-up system</b> , 115 – 220 V/60 Hz, Innova® and Green 'G', ('narrow')                    | U9044-0002 | 1.445,00 € |
| <b>CO<sub>2</sub> back-up system</b> , 230 V/50 Hz, Innova®, Green 'G', and CryoCube® 740 series, ('narrow')   | U9043-0004 | 1.712,00 € |
| <b>LN<sub>2</sub> back-up system</b> , 230 V/50 Hz, Innova®, Green 'G', and CryoCube® 740 series, ('narrow')   | U9044-0004 | 1.712,00 € |



## CO<sub>2</sub> and LN<sub>2</sub> Back-up Systems for CryoCube® F440/ F570/ FC660 series, Premium U410, HEF® U410

| Eppendorf Product Catalog p. 247  | Order no.  | Price      |
|---|------------|------------|
| <b>CO<sub>2</sub> back-up system</b> , 120 – 220 V/60 Hz, CryoCube® F440/ F570/ FC660 series, Premium, and HEF® | U9043-0006 | Inquire    |
| <b>LN<sub>2</sub> back-up system</b> , 120 – 220 V/60 Hz, CryoCube® F440/F570/FC660 series, Premium, and HEF®   | U9044-0006 | Inquire    |
| <b>CO<sub>2</sub> back-up system</b> , 230 V/50 Hz, CryoCube® F440/F570/FC660 series, Premium, and HEF®         | U9043-0008 | 1.712,00 € |
| <b>LN<sub>2</sub> back-up system</b> , 230 V/50 Hz, CryoCube® F440/ F570/FC660 series, Premium, and HEF®        | U9044-0008 | 1.712,00 € |

## Chart Recorder

| Eppendorf Product Catalog p. 247                               | Order no.  | Price      |
|--|------------|------------|
| <b>Chart Recorder (Type I)</b>                                 | P0625-2100 | 1.074,00 € |
| <b>Chart Recorder pens (Type I)</b> , 3 pcs.                   | K0660-0051 | 60,70 €    |
| <b>Chart Recorder paper (Type I)</b> , 0 to -50 °C, 60 pcs.    | P0625-2111 | 85,90 €    |
| <b>Chart Recorder paper (Type I)</b> , -50 to -100 °C, 60 pcs. | P0625-2110 | 85,90 €    |

## Freezer Racks

| Eppendorf Product Catalog p. 248-249   | Order no.    | Price    |
|--|--------------|----------|
| <b>Freezer Rack: Innova® U101, CryoCube® F101h</b>                                       |              |          |
| 2 in/53 mm, drawer, stainless steel  | 6001 062 210 | 170,00 € |
| 2.5 in/64 mm, drawer, stainless steel  | 6001 062 910 | 170,00 € |
| 3 in/76 mm, drawer, stainless steel  | 6001 062 310 | 160,00 € |
| 4 in/102 mm, drawer, stainless steel   | 6001 062 410 | 144,00 € |
| 5 in/127 mm, drawer, stainless steel   | 6001 062 510 | 144,00 € |
| 2 in/53 mm, side access, stainless steel   | 6001 061 210 | 127,00 € |
| 2.5 in/64 mm, side access, stainless steel   | 6001 061 910 | 127,00 € |
| 3 in/76 mm, side access, stainless steel   | 6001 061 310 | 121,00 € |
| 4 in/102 mm, side access, stainless steel  | 6001 061 410 | 106,00 € |
| 5 in/127 mm, side access, stainless steel  | 6001 061 510 | 121,00 € |
| <b>Freezer Rack: CryoCube® F740 series (3-compartment, Standard), Innova® U360, U535</b> |              |          |
| 2 in/53 mm, drawer, stainless steel  | 6001 012 210 | 202,00 € |
| 2.5 in/64 mm, drawer, stainless steel  | 6001 012 910 | 187,00 € |
| 3 in/76 mm, drawer, stainless steel  | 6001 012 310 | 182,00 € |
| 4 in/102 mm, drawer, stainless steel   | 6001 012 410 | 170,00 € |
| 5 in/127 mm, drawer, stainless steel   | 6001 012 510 | 170,00 € |
| 2 in/53 mm, side access, stainless steel   | 6001 011 210 | 178,00 € |
| 2.5 in/64 mm, side access, stainless steel   | 6001 011 910 | 179,00 € |
| 3 in/76 mm, side access, stainless steel   | 6001 011 310 | 176,00 € |
| 4 in/102 mm, side access, stainless steel  | 6001 011 410 | 136,00 € |
| 5 in/127 mm, side access, stainless steel  | 6001 011 510 | 136,00 € |
| DWP (53 mm), side access, stainless steel  | 6001 011 110 | 209,00 € |

## Freezer

## Freezer Racks

Eppendorf Product Catalog p. 248-249

Order no.

Price

**Freezer Rack: Innova® U101, CryoCube® F101h****Freezer Rack: CryoCube® F440/ F570 series, CryoCube® F740 series (5-compartment, compartment 1-4), Premium U410, HEF® U410**

|  |              |          |
|--|--------------|----------|
| 2 in/53 mm, drawer, stainless steel  | 6001 022 210 | 172,00 € |
| 2.5 in/64 mm, drawer, stainless steel  | 6001 022 910 | 170,00 € |
| 3 in/76 mm, drawer, stainless steel  | 6001 022 310 | 167,00 € |
| 4 in/102 mm, drawer, stainless steel   | 6001 022 410 | 167,00 € |
| 2 in/53 mm, side access, stainless steel   | 6001 021 210 | 138,00 € |
| 2.5 in/64 mm, side access, stainless steel   | 6001 021 910 | 136,00 € |
| 3 in/76 mm, side access, stainless steel   | 6001 021 310 | 127,00 € |
| 4 in/102 mm, side access, stainless steel  | 6001 021 410 | 128,00 € |
| DWP (53 mm), side access, stainless steel  | 6001 021 110 | 184,00 € |
| <b>Freezer Rack: CryoCube® F740 series (5-compartment, MAX for 5<sup>th</sup> comp.)</b> |              |          |
| 2 in/53 mm, side access, stainless steel   | 6001 081 210 | 193,00 € |
| 2.5 in/64 mm, side access, stainless steel   | 6001 081 910 | 181,00 € |
| 3 in/76 mm, side access, stainless steel   | 6001 081 310 | 170,00 € |
| 4 in/102 mm, side access, stainless steel  | 6001 081 410 | 159,00 € |
| 5 in/127 mm, side access, stainless steel  | 6001 081 510 | 144,00 € |
| DWP (53 mm), side access, stainless steel  | 6001 081 110 | 220,00 € |
| 2 in/53 mm, drawer, stainless steel  | 6001 082 210 | 234,00 € |
| 2.5 in/64 mm, drawer, stainless steel  | 6001 082 910 | 195,00 € |
| 3 in/76 mm, drawer, stainless steel  | 6001 082 310 | 200,00 € |
| 4 in/102 mm, drawer, stainless steel   | 6001 082 410 | 182,00 € |
| 5 in/127 mm, drawer, stainless steel   | 6001 082 510 | 177,00 € |

## Freezer Rack: CryoCube® F740 series (3-compartment, MAX)

Eppendorf Product Catalog p. 250

Order no.

Price

**Freezer Rack: CryoCube® F740 series (3-compartment, MAX)**

|  |              |          |
|--|--------------|----------|
| 2 in/53 mm, side access, stainless steel   | 6001 071 210 | 206,00 € |
| 2.5 in/64 mm, side access, stainless steel | 6001 071 910 | 179,00 € |
| 3 in/76 mm, side access, stainless steel   | 6001 071 310 | 176,00 € |
| 4 in/102 mm, side access, stainless steel  | 6001 071 410 | 181,00 € |
| 5 in/127 mm, side access, stainless steel  | 6001 071 510 | 136,00 € |
| DWP (53 mm), side access, stainless steel  | 6001 071 110 | 243,00 € |
| 2 in/53 mm, drawer, stainless steel        | 6001 072 210 | 234,00 € |
| 2.5 in/64 mm, drawer, stainless steel      | 6001 072 910 | 187,00 € |
| 3 in/76 mm, drawer, stainless steel        | 6001 072 310 | 182,00 € |
| 4 in/102 mm, drawer, stainless steel       | 6001 072 410 | 210,00 € |
| 5 in/127 mm, drawer, stainless steel       | 6001 072 510 | 170,00 € |

## Freezer Rack: Innova® C585, C760

Eppendorf Product Catalog p. 250

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Freezer Rack: Innova® C585, C760</b>                         |              |          |
| 2 in/53 mm, side access, stainless steel                        | 6001 040 210 | 114,00 € |
| <b>Freezer Rack: Innova® C585, C760, CryoCube® FC660 series</b> |              |          |
| 2.5 in/64 mm, side access, stainless steel                      | 6001 000 910 | 96,80 €  |
| 3 in/76 mm, side access, stainless steel                        | 6001 000 310 | 92,00 €  |
| 4 in/102 mm, side access, stainless steel                       | 6001 000 410 | 104,00 € |
| 5 in/127 mm, side access, stainless steel                       | 6001 000 510 | 114,00 € |
| <b>Freezer Rack: Innova® C585, C760</b>                         |              |          |
| DWP (53 mm), side access, stainless steel                       | 6001 040 110 | 125,00 € |
| MTP, side access, stainless steel                               | 6001 040 010 | 194,00 € |
| locking rod, 2 in/53 mm, side access, stainless steel           | 6001 040 211 | 125,00 € |
| <b>Freezer Rack: Innova® C585, C760, CryoCube® FC660 series</b> |              |          |
| locking rod, 2.5 in/64 mm, side access, stainless steel         | 6001 000 911 | 109,00 € |
| locking rod, 3 in/76 mm, side access, stainless steel           | 6001 000 311 | 111,00 € |
| locking rod, 4 in/102 mm, side access, stainless steel          | 6001 000 411 | 114,00 € |
| <b>Freezer Rack: Innova® C585, C760</b>                         |              |          |
| locking rod, DWP (53 mm), side access, stainless steel          | 6001 040 111 | 137,00 € |
| locking rod, MTP, side access, stainless steel                  | 6001 040 011 | 209,00 € |

## Freezer Rack: Innova® and CryoCube® FC660 series

Eppendorf Product Catalog p. 251

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Freezer Rack: CryoCube® FC660 series</b>                     |              |          |
| 2 in/53 mm, side access, stainless steel                        | 6001 050 210 | 101,00 € |
| <b>Freezer Rack: Innova® C585, C760, CryoCube® FC660 series</b> |              |          |
| 2.5 in/64 mm, side access, stainless steel                      | 6001 000 910 | 96,80 €  |
| 3 in/76 mm, side access, stainless steel                        | 6001 000 310 | 92,00 €  |
| 4 in/102 mm, side access, stainless steel                       | 6001 000 410 | 104,00 € |
| 5 in/127 mm, side access, stainless steel                       | 6001 000 510 | 114,00 € |
| <b>Freezer Rack: CryoCube® FC660 series</b>                     |              |          |
| DWP (53 mm), side access, stainless steel                       | 6001 050 110 | 113,00 € |
| MTP, side access, stainless steel                               | 6001 050 010 | 190,00 € |
| locking rod, 2 in/53 mm, side access, stainless steel           | 6001 050 211 | 114,00 € |
| <b>Freezer Rack: Innova® C585, C760, CryoCube® FC660 series</b> |              |          |
| locking rod, 2.5 in/64 mm, side access, stainless steel         | 6001 000 911 | 109,00 € |
| locking rod, 3 in/76 mm, side access, stainless steel           | 6001 000 311 | 111,00 € |
| locking rod, 4 in/102 mm, side access, stainless steel          | 6001 000 411 | 114,00 € |
| <b>Freezer Rack: CryoCube® FC660 series</b>                     |              |          |
| locking rod, DWP (53 mm), side access, stainless steel          | 6001 050 111 | 125,00 € |
| locking rod, MTP, side access, stainless steel                  | 6001 050 011 | 202,00 € |

## Aluminum racks

Eppendorf Product Catalog p. 251

|  | Order no.  | Price    |
|--|------------|----------|
| <b>Upright freezer rack: Innova® U101, CryoCube® F101h</b>   |            |          |
| 10 boxes per rack, 2 in/53 mm, drawer, aluminum  | K0641-3003 | 213,00 € |
| 6 boxes per rack, 3 in/76 mm, drawer, aluminum   | K0641-3004 | 165,00 € |
| 4 boxes per rack, 4 in/102 mm, drawer, aluminum  | K0641-3005 | 117,00 € |
| <b>Chest freezer rack: Innova® C585, C760, CryoCube® FC660 series</b>  |            |          |
| 13 boxes per rack, 2 in/53 mm, side access, aluminum   | K0641-1690 | 154,00 € |
| 9 boxes per rack, 3 in/76 mm, side access, aluminum  | K0641-1700 | 145,00 € |
| 6 boxes per rack, 4 in/102 mm, side access, aluminum   | K0641-1750 | 140,00 € |
| <b>Upright freezer rack: CryoCube® F740 series (3-compartment), Innova® U360, U535</b>                       |            |          |
| 28 boxes per rack, 2 in/53 mm, drawer, aluminum  | K0641-3000 | 314,00 € |
| 20 boxes per rack, 3 in/76 mm, drawer, aluminum  | K0641-3001 | 245,00 € |
| 12 boxes per rack, 4 in/102 mm, drawer, aluminum   | K0641-3002 | 174,00 € |
| <b>Upright freezer rack: U410 series, CryoCube® F440/ F570 series, CryoCube® F740 series (5-compartment)</b> |            |          |
| 16 boxes per rack, 2 in/53 mm, drawer, aluminum  | K0641-1900 | 186,00 € |
| 12 boxes per rack, 3 in/76 mm, drawer, aluminum  | K0641-1890 | 171,00 € |
| 8 boxes per rack, 4 in/102 mm, drawer, aluminum  | K0641-1880 | 154,00 € |

## Freezer

## Eppendorf Storage Boxes

Eppendorf Product Catalog p. 252

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Storage Box 10 × 10</b> , freezer box, for 100 cryogenic tubes w. int. thread, 3 pcs., height 53 mm, 2 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code  | 0030 140 508 | 29,30 € |
| <b>Storage Box 9 × 9</b> , freezer box, for 81 screw cap (cryog.) tubes 1 – 2 mL, 3 pcs., height 53 mm, 2 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code  | 0030 140 516 | 30,10 € |
| <b>Storage Box 9 × 9</b> , freezer box, for 81 screw cap (cryog.) tubes 3 mL, 2 pcs., height 76 mm, 3 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code      | 0030 140 540 | 20,90 € |
| <b>Storage Box 9 × 9</b> , freezer box, for 81 screw cap (cryog.) tubes 4 – 5 mL, 2 pcs., height 102 mm, 4 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code | 0030 140 567 | 21,50 € |
| <b>Storage Box 8 × 8</b> , freezer box, for 64 tubes 1 – 2 mL, 3 pcs., height 53 mm, 2 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code                     | 0030 140 524 | 30,10 € |
| <b>Storage Box 5 × 5</b> , freezer box, for 25 tubes 5 mL, 4 pcs., height 64 mm, 2.5 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code                       | 0030 140 532 | 42,50 € |
| <b>Storage Box 5 × 5</b> , freezer box, for 25 tubes 5 mL screw cap, 2 pcs., height 76 mm, 3 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code               | 0030 140 613 | 21,50 € |
| <b>Storage Box 5 × 5</b> , freezer box, for 25 tubes 15 mL, 2 pcs., height 127 mm, 5 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code                       | 0030 140 583 | 22,10 € |
| <b>Storage Box 3 × 3</b> , freezer box, for 9 tubes 25 mL, 2 pcs., height 89 mm, 3.5 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code                       | 0030 140 729 | 22,40 € |
| <b>Storage Box 3 × 3</b> , freezer box, for 9 tubes 50 mL and 4 tubes 15 mL, 2 pcs., height 127 mm, 5 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code      | 0030 140 591 | 22,10 € |

## Freezer Cardboard Storage Boxes and dividers

Eppendorf Product Catalog p. 253

|  | Order no. | Price  |
|--|-----------|--------|
| <b>Freezer Cardboard Storage Boxes, white</b>        |           |        |
| 53 mm tall   | B50-SQ    | 4,40 € |
| 76 mm tall   | B75-SQ    | 5,20 € |
| 102 mm tall  | B95-SQ    | 5,20 € |
| <b>Freezer Boxes Dividers</b>                        |           |        |
| 7 × 7, 17.4 mm max. tube diameter, holds 49 tubes    | D49       | 4,40 € |
| 8 × 8, 15 mm max. tube diameter, holds 64 tubes      | D64       | 4,40 € |
| 9 × 9, 13 mm max. tube diameter, holds 81 tubes      | D81       | 4,40 € |
| 10 × 10, 11.8 mm max. tube diameter, holds 100 tubes | D100      | 4,40 € |

## CryoStorage Vials

Eppendorf Product Catalog p. 254-255

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Eppendorf CryoStorage Vial</b> , 0,5 mL, free of DNase, RNase, human DNA, endotoxin; external thread, pre-capped gray, with 2D SafeCode (DataMatrix/ECC200), linear and human readable code, 960 vials (10 racks × 96 vials), individually packed racks | 0030 079 400 | 1.327,00 € |
| <b>Eppendorf CryoStorage Vial</b> , 1,0 mL, free of DNase, RNase, human DNA, endotoxin; external thread, pre-capped gray, with 2D SafeCode (DataMatrix/ECC200), linear and human readable code, 960 vials (10 racks × 96 vials), individually packed racks | 0030 079 418 | 1.371,00 € |
| <b>Eppendorf CryoStorage Vial</b> , 1,5 mL, free of DNase, RNase, human DNA, endotoxin; external thread, pre-capped gray, with 2D SafeCode (DataMatrix/ECC200), linear and human readable code, 480 vials (10 racks × 48 vials), individually packed racks | 0030 079 426 | 708,00 €   |

## CryoStorage Vials

Eppendorf Product Catalog p. 254-255

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Eppendorf CryoStorage Vial</b> , 2,0 mL, free of DNase, RNase, human DNA, endotoxin; external thread, pre-capped gray, with 2D SafeCode (DataMatrix/ECC200), linear and human readable code, 480 vials (10 racks × 48 vials), individually packed racks          | 0030 079 434 | 731,00 €   |
| <b>Eppendorf CryoStorage Vial</b> , 4,0 mL, free of DNase, RNase, human DNA, endotoxin; external thread, pre-capped gray, with 2D SafeCode (DataMatrix/ECC200), linear and human readable code, 480 vials (10 racks × 48 vials), individually packed racks          | 0030 079 442 | 797,00 €   |
| <b>Eppendorf CryoStorage Vial</b> , sterile, 0,5 mL, free of DNase, RNase, human DNA, endotoxin; external thread, pre-capped gray, with 2D SafeCode (DataMatrix/ECC200), linear and human readable code, 960 vials (10 racks × 96 vials), individually packed racks | 0030 079 450 | 1.371,00 € |
| <b>Eppendorf CryoStorage Vial</b> , sterile, 1,0 mL, free of DNase, RNase, human DNA, endotoxin; external thread, pre-capped gray, with 2D SafeCode (DataMatrix/ECC200), linear and human readable code, 960 vials (10 racks × 96 vials), individually packed racks | 0030 079 469 | 1.416,00 € |
| <b>Eppendorf CryoStorage Vial</b> , sterile, 1,5 mL, free of DNase, RNase, human DNA, endotoxin; external thread, pre-capped gray, with 2D SafeCode (DataMatrix/ECC200), linear and human readable code, 480 vials (10 racks × 48 vials), individually packed racks | 0030 079 477 | 731,00 €   |
| <b>Eppendorf CryoStorage Vial</b> , sterile, 2,0 mL, free of DNase, RNase, human DNA, endotoxin; external thread, pre-capped gray, with 2D SafeCode (DataMatrix/ECC200), linear and human readable code, 480 vials (10 racks × 48 vials), individually packed racks | 0030 079 485 | 753,00 €   |
| <b>Eppendorf CryoStorage Vial</b> , sterile, 4,0 mL, free of DNase, RNase, human DNA, endotoxin; external thread, pre-capped gray, with 2D SafeCode (DataMatrix/ECC200), linear and human readable code, 480 vials (10 racks × 48 vials), individually packed racks | 0030 079 493 | 819,00 €   |

## SafeCode Tubes

Eppendorf Product Catalog p. 256

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Eppendorf Conical Tubes</b> , 15 mL, 2D SafeCode, sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 500 tubes (10 bags × 50 tubes)    | 0030 113 489 | 296,00 € |
| <b>Eppendorf Conical Tubes</b> , 50 mL, 2D SafeCode, sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 500 tubes (20 bags × 25 tubes)    | 0030 113 497 | 328,00 € |
| <b>Eppendorf Tubes® 5.0 mL with screw cap</b> , 5 mL, 2D SafeCode, sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, 200 tubes (2 bags × 100 tubes) | 0030 113 500 | 120,60 € |

## Eppendorf RackScan

Eppendorf Product Catalog p. 257

|   | Order no.    | Price       |
|---|--------------|-------------|
| <b>Eppendorf RackScan s</b> , code scanner for side reading of SBS racks and cryostorage boxes, accessory for the RackScan product family, /50 – 60 Hz                            | 6002 000 014 | 1.487,00 €  |
| <b>Eppendorf RackScan b</b> , scanner (bottom reading) for 2D-coded tubes in SBS racks and cryostorage boxes, incl. software RackScan Analyzer, 100 – 240 V ±10 %/50 – 60 Hz (EU) | 6002 000 006 | 10.150,00 € |

## Micromanipulation

### CellTram® 4r Air/Oil

Eppendorf Product Catalog p. 260-261

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>CellTram® 4r Air</b> , manual pneumatic microinjector, with gears 1:1 and 1:10, for holding and injection                     | 5196 000 013 | 2.066,00 € |
| <b>CellTram® 4r Oil</b> , manual hydraulic microinjector, with gears 1:1 and 1:10, for holding and injection                     | 5196 000 030 | 3.042,00 € |
| <b>Filling and Cleaning set</b> , for CellTram® 4 Oil and CellTram® Oil/vario, incl. filling tube, Luer lock adapter, 2 syringes | 5196 088 000 | 78,00 €    |
| <b>Injection tube Air</b> , white ring mark, I.D. 0.5 mm, length 1.3 m   | 5196 061 004 | 71,90 €    |
| <b>Injection tube Oil</b> , blue ring mark, I.D. 1.0 mm, length 1.3 m  | 5196 089 006 | 71,90 €    |
| <b>Capillary holder 4</b> , incl. 1x Grip head set 4 size 0 (for microcapillaries with outer diameter 1.0 – 1.1 mm)              | 5196 081 005 | 233,00 €   |
| <b>Capillary holder 4 slim shape</b> , for flat angle injections, for microcapillaries with outer diameter 1.0 mm                | 5196 062 000 | 238,00 €   |
| <b>O-ring set 4</b> , incl. 10 o-rings large, 10 o-rings small, 2 distance sleeves, o-ring removal tool for grip head set 4      | 5196 086 007 | 47,50 €    |
| <b>Grip head set 4 size 0</b> , 1 pc, for microcapillaries with outer diameter 1.0 to 1.1 mm                                     | 5196 082 001 | 47,50 €    |
| <b>Grip head set 4 size 1</b> , 1 pc, for microcapillaries with outer diameter 1.2 to 1.3 mm                                     | 5196 083 008 | 47,50 €    |
| <b>Grip head set 4 size 2</b> , 1 pc, for microcapillaries with outer diameter 1.4 to 1.5 mm                                     | 5196 084 004 | 47,50 €    |
| <b>Grip head set 4 size 3</b> , 1 pc, for microcapillaries with outer diameter 0.7 to 0.9 mm                                     | 5196 085 000 | 47,50 €    |

### TransferMan® 4r

Eppendorf Product Catalog p. 262

|  | Order no.    | Price       |
|--|--------------|-------------|
| <b>TransferMan® 4r</b> , micromanipulator with DualSpeed™ joystick for direct and dynamic movement control, 100 – 240 V/50 Hz (EU), For research use only. | 5193 000 012 | 17.470,00 € |

### InjectMan® 4

Eppendorf Product Catalog p. 263

|   | Order no.    | Price       |
|---|--------------|-------------|
| <b>InjectMan® 4</b> , micromanipulator with dynamic movement control, 100 – 240 V/50 – 60 Hz (EU) | 5192 000 019 | 16.940,00 € |

### InjectMan® 4 and TransferMan® 4r Accessories

Eppendorf Product Catalog p. 263

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Positioning aid</b> , for universal capillary holder, capillary holder 4, 2 pcs.   | 5192 072 001 | 111,00 € |
| <b>Spare parts kit</b> , 1 slider (complete), 2 machine screws M2.5x6 (swivel joint), 2 set screws (angle head), 2 compression springs (knurled screw angle head) | 5192 071 005 | 114,00 € |
| <b>Headstage holder</b> , for preamplifier  | 5192 073 008 | 121,00 € |
| <b>Foot control</b> , for FemtoJet® 4i/4x   | 5252 070 020 | 189,00 € |
| <b>Connection cable</b> , for connecting TransferMan® 4r/InjectMan® 4 to FemtoJet® 4i/4x  | 5192 082 007 | 82,10 €  |
| <b>Y-cable FJ4</b> , Y-cable for connecting FemtoJet® 4i/4x with a PC and TransferMan® 4r or InjectMan® 4   | 5192 080 004 | 137,00 € |
| <b>Y-cable PX</b> , Y-cable for connecting Eppendorf PiezoXpert® or a PC with TransferMan® 4r or InjectMan® 4 and FemtoJet® 4i/4x                                 | 5192 081 000 | 137,00 € |

## Adapters for the InjectMan® 4, TransferMan® 4r

Eppendorf Product Catalog p. 264

|  | Order no.    | Price    |
|--|--------------|----------|
| <b>Adapter Leica® 1</b> , for Eppendorf micromanipulation systems, for Leica® DMi8, DMi3000 B, 3000 M, 4000 B, 5000 B, 5000 M, 6000 B, DM IRB E, HC and DM IRE 2 microscopes | 5192 301 000 | 755,00 € |
| <b>Adapter Leica® 2</b> , for Eppendorf micromanipulation systems, for Leica® DM IL LED and HC microscopes   | 5192 302 007 | 755,00 € |
| <b>Adapter Nikon® 1</b> , for Eppendorf micromanipulation systems, for Nikon® Eclipse® Diaphot 200, 300 and Eclipse® Ti-E, Ti-U, Ti-S, TE200, TE300, TE2000 microscopes      | 5192 316 008 | 883,00 € |
| <b>Adapter Nikon® 2</b> , for Eppendorf micromanipulation systems, for Nikon® Eclipse® Ts2R microscope   | 5192 317 004 | 866,00 € |
| <b>Adapter Nikon® 3</b> , for Eppendorf micromanipulation systems, for Nikon® Eclipse Ti2-E/-A/-U microscopes  | 5192 318 000 | 755,00 € |
| <b>Adapter Olympus® 1</b> , for Eppendorf micromanipulation systems, for Olympus® IX50, IX51, IX70, IX71, IX80 and IX81 microscopes  | 5192 306 002 | 755,00 € |
| <b>Adapter Olympus® 2</b> , for Eppendorf micromanipulation systems, for microscopes of Olympus® IX3 series (IX53, IX73, IX83) with standard illumination column             | 5192 307 009 | 680,00 € |
| <b>Adapter Olympus® 3</b> , for Eppendorf micromanipulation systems, for microscopes of Olympus® IX3 series (IX53, IX73, IX83) with special illumination column IX2-ILL30    | 5192 308 005 | 680,00 € |
| <b>Adapter Zeiss® 1</b> , for Eppendorf micromanipulation systems, for Zeiss® Axiovert® 200, Axio Observer A1, D1, Z1 and Axio Observer 3, 5 and 7 microscopes               | 5192 311 006 | 755,00 € |
| <b>Adapter Zeiss® 2</b> , for Eppendorf micromanipulation systems, for Zeiss® Axio Vert.A1 microscope  | 5192 312 002 | 883,00 € |
| <b>Universal stand</b> , for mounting Eppendorf micromanipulators on upright microscopes and stereo microscopes independent of the microscope tripod used                    | 5192 325 007 | 923,00 € |
| <b>Adapter bridge</b> , for mounting the TransferMan® 4 and InjectMan® 4 on the microscope adapters of the predecessors TransferMan NK2, InjectMan NI2 and PatchMan NP2.     | 5192 321 001 | 513,00 € |

## Antivibration pad

Eppendorf Product Catalog p. 265

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Antivibration Pad XS</b> , weight range 4.5 – 6.0 kg   | 5181 301 009 | 551,00 € |
| <b>Antivibration Pad S</b> , weight range 6.0 – 8.0 kg    | 5181 303 001 | 551,00 € |
| <b>Antivibration Pad M</b> , weight range 8.0 – 10.0 kg   | 5181 305 004 | 551,00 € |
| <b>Antivibration Pad L</b> , weight range 10.0 – 12.5 kg  | 5181 307 007 | 551,00 € |
| <b>Antivibration Pad XL</b> , weight range 12.5 – 16.5 kg | 5181 309 000 | 551,00 € |

## PiezoXpert®

Eppendorf Product Catalog p. 266

|  | Order no.    | Price       |
|--|--------------|-------------|
| <b>PiezoXpert®</b> , for piezo-assisted micromanipulation, incl. actuator 2, foot control and spacer plate and grip head 4 size 0, 100 – 240 V/50 – 60 Hz (EU)     | 5194 000 016 | 13.470,00 € |
| <b>Piezo Drill Tip ES</b> , for piezo-assisted mouse ES cell transfer (for research use only), 25 ° tip angle, 15 µm inner diameter, 6 mm flange, sterile, 25 pcs. | 5195 000 095 | 489,00 €    |
| <b>Piezo Drill Tip ICSI</b> , for piezo-assisted mouse ICSI (for research use only), 25 ° tip angle, 6 µm inner diameter, 6 mm flange, sterile, 25 pcs.            | 5195 000 087 | 489,00 €    |
| <b>Actuator 2</b> , for piezo-assisted micromanipulation, as accessory for the PiezoXpert®   | 5194 075 156 | 1.627,00 €  |
| <b>Grip head set 4 size 0</b> , 1 pc, for microcapillaries with outer diameter 1.0 to 1.1 mm   | 5196 082 001 | 47,50 €     |
| <b>Grip head set 4 size 1</b> , 1 pc, for microcapillaries with outer diameter 1.2 to 1.3 mm   | 5196 083 008 | 47,50 €     |
| <b>Grip head set 4 size 2</b> , 1 pc, for microcapillaries with outer diameter 1.4 to 1.5 mm   | 5196 084 004 | 47,50 €     |
| <b>Grip head set 4 size 3</b> , 1 pc, for microcapillaries with outer diameter 0.7 to 0.9 mm   | 5196 085 000 | 47,50 €     |
| <b>Tube adapter</b> , pack of 2, for connecting tubes with an outer diameter of 2 mm or 3 mm with the capillary holder and the capillary holder 4                  | 5194 075 407 | 64,10 €     |

## Micromanipulation

### FemtoJet® 4i and FemtoJet® 4x

Eppendorf Product Catalog p. 267

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>FemtoJet® 4i</b> , programmable microinjector with integrated pressure supply, 100 – 240 V ±10 %/50 – 60 Hz (EU), including foot control                | 5252 000 013 | 7.866,00 € |
| <b>FemtoJet® 4x</b> , programmable microinjector with external pressure supply (not included), 100 – 240 V ±10 %/50 – 60 Hz (EU), including foot control   | 5253 000 017 | 6.409,00 € |
| <b>Capillary holder 4</b> , incl. 1x Grip head set 4 size 0 (for microcapillaries with outer diameter 1.0 – 1.1 mm)  | 5196 081 005 | 233,00 €   |
| <b>Capillary holder 4 slim shape</b> , for flat angle injections, for microcapillaries with outer diameter 1.0 mm  | 5196 062 000 | 238,00 €   |
| <b>Injection tube</b> , for universal capillary holder and capillary holder 4, 2 m   | 5252 070 054 | 63,60 €    |
| <b>O-ring set 4</b> , incl. 10 o-rings large, 10 o-rings small, 2 distance sleeves, o-ring removal tool for grip head set 4                                | 5196 086 007 | 47,50 €    |
| <b>Hand control</b> , for triggering inject function   | 5252 070 011 | 63,00 €    |
| <b>Grip head set 4 size 0</b> , 1 pc, for microcapillaries with outer diameter 1.0 to 1.1 mm   | 5196 082 001 | 47,50 €    |
| <b>Grip head set 4 size 1</b> , 1 pc, for microcapillaries with outer diameter 1.2 to 1.3 mm   | 5196 083 008 | 47,50 €    |
| <b>Grip head set 4 size 2</b> , 1 pc, for microcapillaries with outer diameter 1.4 to 1.5 mm   | 5196 084 004 | 47,50 €    |
| <b>Grip head set 4 size 3</b> , 1 pc, for microcapillaries with outer diameter 0.7 to 0.9 mm   | 5196 085 000 | 47,50 €    |
| <b>Grip head set 4 (slim shape)</b> , for Capillary holder 4 slim shape incl. 6 O-rings and 2 distance sleeves, for capillaries with outer diameter 1.0 mm | 5196 063 007 | 47,50 €    |
| <b>Y-cable FJ4</b> , Y-cable for connecting FemtoJet® 4i/4x with a PC and TransferMan® 4r or InjectMan® 4  | 5192 080 004 | 137,00 €   |
| <b>Y-cable PX</b> , Y-cable for connecting Eppendorf PiezoXpert® or a PC with TransferMan® 4r or InjectMan® 4 and FemtoJet® 4i/4x                          | 5192 081 000 | 137,00 €   |

### Eppendorf Microcapillaries

Eppendorf Product Catalog p. 268-271

|   | Order no.    | Price    |
|---|--------------|----------|
| <b>Biopsy Tip I</b> , capillary for laser-assisted biopsy of cells and organelles (for research use only), 35 ° tip angle, 19 µm inner diameter, 1,9 mm flange, sterile, 25 pcs., sterile, 25 pcs.  | 5195 000 052 | 489,00 € |
| <b>Biopsy Tip II</b> , capillary for laser-assisted biopsy of cells and organelles (for research use only), 15 ° tip angle, 20 µm inner diameter, 0,5 mm flange, sterile, 25 pcs., sterile, 25 pcs. | 5195 000 060 | 489,00 € |
| <b>Femtotips®</b> , injection capillary (for research use only), sterile, 20 pcs., sterile, 20 pcs.   | 5242 952 008 | 137,00 € |
| <b>Femtotip II</b> , injection capillary (for research use only), sterile, 20 pcs., sterile, 20 pcs.  | 5242 957 000 | 137,00 € |
| <b>VacuTip I</b> , holding capillary (for research use only), 35 ° tip angle, 15 µm inner diameter, 100 µm outer diameter, 1 mm flange, sterile, 25 pcs., sterile, 25 pcs.                          | 5195 000 036 | 477,00 € |
| <b>VacuTip II</b> , holding capillary (for research use only), 25 ° tip angle, 60 µm inner diameter, 110 µm outer diameter, 0,5 mm flange, sterile, 25 pcs., sterile, 25 pcs.                       | 5195 000 044 | 489,00 € |
| <b>TransferTip® RP (ICSI)</b> , for sperm injection using the ICSI technique (for research use only), 35 ° tip angle, 4 µm inner diameter, 0,5 mm flange, sterile, 25 pcs., sterile, 25 pcs.        | 5195 000 010 | 626,00 € |
| <b>TransferTip® R (ICSI)</b> , for sperm injection using the ICSI technique (for research use only), 35 ° tip angle, 4 µm inner diameter, 1 mm flange, sterile, 25 pcs., sterile, 25 pcs.           | 5195 000 028 | 626,00 € |
| <b>TransferTip® F (ICSI)</b> , for sperm injection using the ICSI technique (for research use only), 35 ° tip angle, 4 µm inner diameter, 0,4 mm flange, sterile, 25 pcs., sterile, 25 pcs.         | 5195 000 001 | 627,00 € |
| <b>TransferTip® (ES)</b> , for ES cell transfer (for research use only), 20 ° tip angle, 15 µm inner diameter, 1 mm flange, sterile, 25 pcs., sterile, 25 pcs.                                      | 5195 000 079 | 627,00 € |
| <b>Piezo Drill Tip ES</b> , for piezo-assisted mouse ES cell transfer (for research use only), 25 ° tip angle, 15 µm inner diameter, 6 mm flange, sterile, 25 pcs.                                  | 5195 000 095 | 489,00 € |
| <b>Piezo Drill Tip ICSI</b> , for piezo-assisted mouse ICSI (for research use only), 25 ° tip angle, 6 µm inner diameter, 6 mm flange, sterile, 25 pcs.   | 5195 000 087 | 489,00 € |
| <b>Microloader</b> , tip for filling Femtotips® and other glass microcapillaries (for research use only), 0,5 – 20 µL, 100 mm, light gray, 192 pcs. (2 racks × 96 pcs.)                             | 5242 956 003 | 102,00 € |



## CellXpert® C170i

Eppendorf Product Catalog p. 276-278

|  | Order no.    | Price       |
|--|--------------|-------------|
| <b>CellXpert® C170i</b>  |              |             |
| non-segmented inner door, handle right side, 220 – 240 V/50 – 60 Hz (EU)   | 6731 000 011 | 10.340,00 € |
| non-segmented inner door, handle right side, copper package, 220 – 240 V/50 – 60 Hz (EU)   | 6731 000 511 | 14.570,00 € |
| non-segmented inner door, handle right side, humidity monitoring, 220 – 240 V/50 – 60 Hz (EU)  | 6731 000 111 | 11.730,00 € |
| non-segmented inner door, handle right side, water-level monitoring, 220 – 240 V/50 – 60 Hz (EU)   | 6731 000 211 | 11.040,00 € |
| non-segmented inner door, handle right side, water and humidity monitoring package, 220 – 240 V/50 – 60 Hz (EU)  | 6731 000 311 | 12.420,00 € |
| non-segmented inner door, handle right side, 1 – 20 % oxygen control, 220 – 240 V/50 – 60 Hz (EU)  | 6731 001 011 | 12.580,00 € |
| non-segmented inner door, handle right side, 1 – 20 % oxygen control, water and humidity monitoring package, copper package, 220 – 240 V/50 – 60 Hz (EU)                           | 6731 001 811 | 18.920,00 € |
| non-segmented inner door, handle right side, 1 – 20 % oxygen control, copper package, 220 – 240 V/50 – 60 Hz (EU)  | 6731 001 511 | 16.530,00 € |
| non-segmented inner door, handle left side, 220 – 240 V/50 – 60 Hz (EU)  | 6731 000 021 | 10.340,00 € |
| non-segmented inner door, handle left side, 1 – 20 % oxygen control, copper package, 220 – 240 V/50 – 60 Hz (EU)   | 6731 001 521 | 16.530,00 € |
| non-segmented inner door, handle left side, copper package, 220 – 240 V/50 – 60 Hz (EU)  | 6731 000 521 | 14.570,00 € |
| non-segmented inner door, handle left side, humidity monitoring, 220 – 240 V/50 – 60 Hz (EU)   | 6731 000 121 | 11.730,00 € |
| non-segmented inner door, handle left side, water-level monitoring, 220 – 240 V/50 – 60 Hz (EU)  | 6731 000 221 | 11.040,00 € |
| non-segmented inner door, handle left side, water and humidity monitoring package, 220 – 240 V/50 – 60 Hz (EU)   | 6731 000 321 | 12.420,00 € |
| non-segmented inner door, handle left side, 1 – 20 % oxygen control, 220 – 240 V/50 – 60 Hz (EU)   | 6731 001 021 | 12.580,00 € |
| non-segmented inner door, handle left side, 1 – 20 % oxygen control, water and humidity monitoring package, copper package, 220 – 240 V/50 – 60 Hz (EU)                            | 6731 001 821 | 18.920,00 € |
| inner door with 4 door segments, handle right side, 220 – 240 V/50 – 60 Hz (EU)  | 6731 000 041 | 12.220,00 € |
| <b>mAb-Discovery Bundle: Cell Culture Screening Solution</b>   |              |             |
| inner door with 4 door segments, handle right side, water and humidity monitoring package, 220 – 240 V/50 – 60 Hz (EU), For screening assays. Bundle consists of: CellXpert® C170i | 6731 000 341 | 14.300,00 € |
| <b>CellXpert® C170i</b>  |              |             |
| inner door with 4 door segments, handle right side, water and humidity monitoring package, copper package, 220 – 240 V/50 – 60 Hz (EU)   | 6731 000 841 | 18.560,00 € |
| inner door with 4 door segments, handle right side, 1 – 20 % oxygen control, 220 – 240 V/50 – 60 Hz (EU)   | 6731 001 041 | 14.460,00 € |
| inner door with 4 door segments, handle right side, 1 – 20 % oxygen control, water and humidity monitoring package, 220 – 240 V/50 – 60 Hz (EU)                                    | 6731 001 341 | 16.530,00 € |
| inner door with 4 door segments, handle right side, 0.1 – 20 % oxygen control, 220 – 240 V/50 – 60 Hz (EU)   | 6731 002 041 | 15.480,00 € |
| inner door with 4 door segments, handle left side, 220 – 240 V/50 – 60 Hz (EU)   | 6731 000 051 | 12.220,00 € |
| inner door with 4 door segments, handle left side, water and humidity monitoring package, 220 – 240 V/50 – 60 Hz (EU)  | 6731 000 351 | 14.300,00 € |
| inner door with 4 door segments, handle left side, water and humidity monitoring package, copper package, 220 – 240 V/50 – 60 Hz (EU)  | 6731 000 851 | 18.560,00 € |
| inner door with 4 door segments, handle left side, 1 – 20 % oxygen control, water and humidity monitoring package, 220 – 240 V/50 – 60 Hz (EU)                                     | 6731 001 351 | 16.530,00 € |
| inner door with 4 door segments, handle left side, 1 – 20 % oxygen control, 220 – 240 V/50 – 60 Hz (EU)  | 6731 001 051 | 14.460,00 € |
| inner door with 8 door segments, handle right side, 1 – 20 % oxygen control, 220 – 240 V/50 – 60 Hz (EU)   | 6731 001 081 | 14.920,00 € |
| inner door with 8 door segments, handle right side, 0.1 – 20 % oxygen control, 220 – 240 V/50 – 60 Hz (EU)   | 6731 002 081 | 15.960,00 € |
| inner door with 8 door segments, handle left side, 1 – 20 % oxygen control, 220 – 240 V/50 – 60 Hz (EU)  | 6731 001 091 | 14.920,00 € |

## CO<sub>2</sub> Incubators

### CellXpert® C170/C170i shelves, water trays and stacking stand

Eppendorf Product Catalog p. 278

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Stacking stand, lower frame</b> , with castors  | 6731 070 093 | 1.058,00 € |
| <b>Stacking stand, upper frame</b>   | 6731 070 085 | 353,00 €   |
| <b>Universal stacking stand, upper frame</b> , CellXpert with other 170 L CO <sub>2</sub> incubators   | 6731 070 158 | 338,00 €   |
| <b>Set of copper water tray and 4 copper shelves for 170 L incubators</b>                              | 6731 080 013 | 2.798,00 € |
| <b>Copper water tray</b> , for 170 L incubators  | P0628-6260   | 802,00 €   |
| <b>Copper shelves</b> , for 170 L incubators, 2 pcs.   | 6710 859 106 | 1.567,00 € |
| <b>Shelves</b> , for 170 L incubators, 2 pcs.  | 6710 859 009 | 398,00 €   |
| <b>Shelf</b> , for 170 L CO <sub>2</sub> incubators, stainless steel, 1 piece                          | 6731 070 123 | 214,00 €   |
| <b>Shelf rack for 8 shelves</b> , for 170 L incubators   | P0628-6390   | 279,00 €   |
| <b>Water tray</b> , for 170 L incubators   | P0628-6140   | 136,00 €   |
| <b>Door lock</b> , for CellXpert® C170i CO <sub>2</sub> incubators, retrofittable, no service required | 6731 070 166 | 210,00 €   |

### CellXpert® C170

Eppendorf Product Catalog p. 279

|  | Order no.    | Price       |
|--|--------------|-------------|
| <b>CellXpert® C170</b> , non-segmented inner door, handle right side, 220 – 240 V/50 – 60 Hz (EU)        | 6734 000 011 | 8.858,00 €  |
| <b>CellXpert® C170</b> , inner door with 4 door segments, handle right side, 220 – 240 V/50 – 60 Hz (EU) | 6734 000 041 | 10.760,00 € |

### New Brunswick S41i

Eppendorf Product Catalog p. 280

|   | Order no.  | Price       |
|---|------------|-------------|
| <b>New Brunswick S41i</b> , 170 L, CO <sub>2</sub> incubator shaker with inner shelf and touch screen control, 1 shelf, stackable up to 2-fold, 1 (2. optional) shelves, 230 V/50/60 Hz (EU), orbit 2.5 cm (1 in)   | S411230011 | 16.970,00 € |
| <b>Universal platform for New Brunswick S41i</b> , 36 × 61 cm (14 × 24 in), compatible with sticky pads: 5 necessary to cover platform (sold separately, order no: M1250-9700, cutting required), anodized aluminum | M1334-9920 | 736,00 €    |
| <b>Interchangeable platform for New Brunswick S41i</b> , 36 × 61 cm (14 × 24 in)  |            |             |
| 125 mL Erlenmeyer flask dedicated platform, comes complete with 32 clamps installed   | M1334-9921 | 845,00 €    |
| 250 mL Erlenmeyer flask dedicated platform, comes complete with 24 clamps installed   | M1334-9922 | 792,00 €    |
| 500 mL Erlenmeyer flask dedicated platform, comes complete with 15 clamps installed   | M1334-9923 | 779,00 €    |
| 1 L Erlenmeyer flask dedicated platform, comes complete with 11 clamps installed  | M1334-9924 | 769,00 €    |
| 2 L Erlenmeyer flask dedicated platform, comes complete with 6 clamps installed   | M1334-9925 | 769,00 €    |
| 2.8 L Erlenmeyer flask dedicated platform, comes complete with 3 clamps installed   | M1334-9926 | 769,00 €    |
| 4 L Erlenmeyer flask dedicated platform, comes complete with 3 clamps installed   | M1334-9927 | 769,00 €    |
| <b>Stacking kit</b> , allows up to 2 each New Brunswick S41i shakers to be stacked one on top of the other  | P0628-6502 | 1.350,00 €  |
| <b>Additional shelf for New Brunswick S41i</b> , perforated, can easily be added for static cell culture of adherent cells  | P0628-6181 | 278,00 €    |
| <b>Sticky pad</b> , 20 × 20 cm (8 × 8 in), Compatible platform or adapter kit necessary   | M1250-9700 | 85,50 €     |

### CO<sub>2</sub> Incubators Accessories

Eppendorf Product Catalog p. 281

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>HEPA filter 25 mm</b>                                     | P0628-5060   | 71,70 €    |
| <b>In-line gas filter</b> , 2 pcs.                           | 6710 070 251 | Inquire    |
| <b>Automatic gas cylinder CO<sub>2</sub> changeover unit</b> | P0628-5000   | 1.398,00 € |
| <b>Temperature probe</b> , 100 mm tip for gas analyzer       | P0628-7880   | 247,00 €   |
| <b>Temperature probe</b> , 5 mm tip for gas analyzer         | P0628-7881   | 104,00 €   |

## Innova® S44i

Eppendorf Product Catalog p. 286-287

|   | Order no.  | Price       |
|---|------------|-------------|
| <b>Innova® S44i</b> , non-refrigerated, stackable up to 3-fold, 230 V/50/60 Hz (EU), orbit 2.5 cm (1 in)  | S44I300001 | 19.900,00 € |
| <b>Innova® S44i</b> , refrigerated, stackable up to 3-fold, 230 V/50/60 Hz (EU), orbit 2.5 cm (1 in)  | S44I310001 | 22.680,00 € |
| <b>Innova® S44i</b> , refrigerated, stackable up to 3-fold, 230 V/50/60 Hz (EU), orbit 2.5 cm (1 in), photosynthetic LED light bank   | S44I311001 | 25.860,00 € |
| <b>Innova® S44i</b> , non-refrigerated, stackable up to 3-fold, 230 V/50/60 Hz (EU), orbit 5.1 cm (2 in)  | S44I320001 | 19.900,00 € |
| <b>Innova® S44i</b> , refrigerated, stackable up to 3-fold, 230 V/50/60 Hz (EU), orbit 5.1 cm (2 in)  | S44I330001 | 22.680,00 € |
| <b>Innova® S44i</b> , refrigerated, stackable up to 3-fold, 230 V/50/60 Hz (EU), orbit 5.1 cm (2 in), photosynthetic LED light bank   | S44I331001 | 25.860,00 € |
| <b>Positioning base</b> , 10.5 cm (4.1 in), Innova® S44i  | S44I041001 | 731,00 €    |
| <b>Positioning base</b> , 30.7 cm (12.1 in), Innova® S44i   | S44I041002 | 1.339,00 €  |
| <b>Stacking kit</b> , for double stack of Innova® S44i  | S44I041005 | 401,00 €    |
| <b>Stacking kit</b> , for Innova® S44i to Innova® 44/44R  | S44I041006 | 401,00 €    |
| <b>Stacking kit</b> , for triple stack Innova® S44i   | S44I041003 | 1.201,00 €  |
| <b>Static incubation shelf</b> , Innova® S44i   | S44I042001 | 386,00 €    |
| <b>Darkening window shade</b> , removable, Innova® S44i   | S44I042002 | 293,00 €    |
| <b>Universal platform for Innova® S44i</b> , 76 × 62 cm (30 × 24 in), compatible with sticky pads: 9 necessary to cover platform (sold separately, order no: M1250-9700, cutting required), anodized aluminum | S44I040001 | 1.215,00 €  |
| <b>125 mL Erlenmeyer flask dedicated platform</b> , 86 clamps   | S44I040002 | 1.587,00 €  |
| <b>250 mL Erlenmeyer flask dedicated platform</b> , 53 clamps   | S44I040003 | 1.467,00 €  |
| <b>500 mL Erlenmeyer flask dedicated platform</b> , 39 clamps   | S44I040004 | 1.376,00 €  |
| <b>2 L Lauber Thomson flask dedicated platform</b> , 15 clamps, aluminum, for 2 L narrow neck Erlenmeyer, Thomson, or Corning flasks (DIN-ISO 1773) bottom diameter ~166 mm                                   | S44I040010 | 1.512,00 €  |
| <b>1 L Erlenmeyer flask dedicated platform</b> , 23 clamps  | S44I040005 | 1.303,00 €  |
| <b>2 L Erlenmeyer flask dedicated platform</b> , 15 clamps, aluminum, for 2 L wide neck Erlenmeyer flasks (DIN-ISO 24450), bottom diameter ~153 mm  | S44I040006 | 1.230,00 €  |

## New Brunswick Innova® 44/44R

Eppendorf Product Catalog p. 288

|   | Order no.  | Price       |
|---|------------|-------------|
| <b>New Brunswick Innova® 44</b>   |            |             |
| stackable up to 3-fold, 230 V/50/60 Hz (EU/CN), orbit 2.5 cm (1 in)   | M1282-0002 | 15.800,00 € |
| stackable up to 3-fold, 230 V/50/60 Hz (EU/CN), orbit 5.1 cm (2 in)   | M1282-0012 | 15.800,00 € |
| <b>New Brunswick Innova® 44R</b>  |            |             |
| refrigerated, stackable up to 3-fold, 230 V/50/60 Hz (EU/CN), orbit 2.5 cm (1 in)   | M1282-0006 | 19.430,00 € |
| refrigerated, stackable up to 3-fold, 230 V/50/60 Hz (EU/CN), orbit 5.1 cm (2 in)   | M1282-0016 | 19.430,00 € |
| <b>Stacking kit</b> , for stacking one additional 44 or 44R (3 × Innova® 44/44R can be placed in a single stack)  | M1282-0500 | 392,00 €    |
| <b>Remote alarm option</b> , for Innova® 44/44, Innova® 43/43, Innova® 42/42 or Innova® 40/40   | M1320-8029 | 420,00 €    |
| <b>Universal platform for Innova® 44/44R</b> , 76 × 46 cm (30 × 18 in), not interchangeable with other 76 × 46 cm platforms   | M1282-9904 | 1.245,00 €  |
| <b>Aluminum platform for Innova® 44/44R</b> , 76 × 46 cm (30 × 18 in), not interchangeable with other 76 × 46 cm platforms  |            |             |
| 125 mL Erlenmeyer flask dedicated platform  | M1282-9905 | 1.626,00 €  |
| 250 mL Erlenmeyer flask dedicated platform  | M1282-9906 | 1.503,00 €  |
| 500 mL Erlenmeyer flask dedicated platform  | M1282-9907 | 1.410,00 €  |
| 1 L Erlenmeyer flask dedicated platform   | M1282-9908 | 1.334,00 €  |
| <b>2 L Erlenmeyer flask dedicated platform for Innova® 44</b> , 76 × 46 cm (30 × 18 in), 12 clamps, aluminum, for 2 L wide neck Erlenmeyer flasks (DIN-ISO 24450), bottom diameter ~153 mm. For 2 L narrow neck Erlenmeyer, Thomson, or Corning flasks (DIN-ISO 1773) bottom diameter ~166 mm, please refer to universal platform and clamps M1190-9004 | M1282-9909 | 1.260,00 €  |
| <b>Sticky pad platform for Innova® 44(R)</b> , 76 × 46 cm (30 × 18 in), 7 sticky pads necessary to cover platform (sold separately, order no: M1250-9700, cutting required), anodized aluminum  | M1282-9913 | 881,00 €    |

## Shaker

## New Brunswick Innova® 44/44R

Eppendorf Product Catalog p. 288

Order no.

Price

**New Brunswick Innova® 44****Heavy duty 0.95 cm-thick aluminum plate platform for Innova® 44/44R only**, 76 × 46 cm (30 × 18 in), not interchangeable with other 76 × 46 cm platforms.

Heavy-duty platforms are recommended for shaking full loads of flasks, 1 L and larger, containing more than 20 % volume or with baffles.

|  | Order no.  | Price      |
|--|------------|------------|
| Heavy duty universal platform  | M1282-9915 | 1.495,00 € |
| 1 L Erlenmeyer flask dedicated platform, heavy duty  | M1282-9925 | 1.604,00 € |
| 2 L Erlenmeyer flask dedicated platform, heavy duty, for 2 L wide neck Erlenmeyer flasks (DIN-ISO 24450)       | M1282-9926 | 1.519,00 € |
| <b>Positioning bases</b> , increases the height of the shaker for easier access, mandatory for stacked shakers |            |            |
| for Innova® 44/44R, 10.2 cm, short   | M1282-0600 | 713,00 €   |
| for Innova® 44/44R, 30.5 cm, medium  | M1282-0800 | 1.307,00 € |
| for Innova® 44/44R, 40.6 cm, tall  | M1282-0700 | 1.417,00 € |
| <b>12-port gassing manifold</b> , for Innova® 44R - contact Eppendorf for installation options                 | M1282-0291 | 798,00 €   |
| <b>Humidity monitor</b> , for Innova® 44R - contact Eppendorf for installation options                         | M1282-5011 | 1.080,00 € |
| <b>UV germicidal light</b> , for Innova® 44R - contact Eppendorf for installation options                      | M1282-9921 | 2.156,00 € |
| <b>Photosynthetic light bank</b> , for Innova® 44R - contact Eppendorf for installation options                | M1282-9922 | 2.156,00 € |

## New Brunswick Innova® 43/43R

Eppendorf Product Catalog p. 289

Order no.

Price

**New Brunswick Innova® 43**

|  |            |             |
|--|------------|-------------|
| non-stackable, 230 V/50/60 Hz (CN/EU), orbit 2.5 cm (1 in) | M1320-0002 | 16.190,00 € |
| non-stackable, 230 V/50/60 Hz (CN/EU), orbit 5.1 cm (2 in) | M1320-0012 | 16.190,00 € |

**New Brunswick Innova® 43R**

|  |            |             |
|--|------------|-------------|
| refrigerated, non-stackable, 230 V/50/60 Hz (CN/EU), orbit 2.5 cm (1 in)   | M1320-0006 | 18.370,00 € |
| refrigerated, non-stackable, 230 V/50/60 Hz (CN/EU), orbit 5.1 cm (2 in)   | M1320-0016 | 18.370,00 € |
| <b>Remote alarm option</b> , for Innova® 44/44, Innova® 43/43, Innova® 42/42 or Innova® 40/40  | M1320-8029 | 420,00 €    |
| <b>Platforms for Innova® 2300/43/43R, Excella® E10/E25/E25R</b> , 76 × 46 cm (30 × 18 in), aluminum, Universal platform  | M1250-9920 | 1.223,00 €  |
| <b>Platform for Innova® 2300/43/43R, Excella® E10/E25/E25R</b> , 76 × 46 cm (30 × 18 in), aluminum   |            |             |
| 125 mL Erlenmeyer flask dedicated platform   | M1191-9909 | 1.584,00 €  |
| 250 mL Erlenmeyer flask dedicated platform   | M1191-9910 | 1.467,00 €  |
| 500 mL Erlenmeyer flask dedicated platform   | M1191-9911 | 1.374,00 €  |
| 1 L Erlenmeyer flask dedicated platform  | AG-1       | 1.364,00 €  |
| 2 L Erlenmeyer flask dedicated platform, for 2 L wide neck Erlenmeyer flasks (DIN-ISO 24450), bottom diameter ~153 mm. For 2 L narrow neck Erlenmeyer, Thomson, or Corning flasks (DIN-ISO 1773) bottom diameter ~166 mm please refer to universal platform and clamps M1190-9004. | AG-2       | 1.343,00 €  |
| <b>Sticky pad platform for Innova® 43(R)/2300</b> , 76 × 46 cm (30 × 18 in), 7 sticky pads necessary to cover platform (sold separately, order no: M1250-9700, cutting required), anodized aluminum  | M1250-9904 | 839,00 €    |
| <b>Utility tray with non-skid rubber surface for Innova® 2300</b> , usable as a stand-alone – no other platform required   | AG-00      | 986,00 €    |
| <b>12-port gassing manifold</b> , permits direct gassing into flasks   | M1320-0500 | 798,00 €    |
| <b>Humidity monitor</b> , for use with spill pan as water reservoir, for Innova® 43/43R  | M1320-0600 | 1.080,00 €  |
| <b>UV germicidal light</b> , in the airflow path outside the chamber reduces the bio-burden of circulated air, Innova® 43R only  | M1320-0400 | 2.156,00 €  |
| <b>Photosynthetic light bank</b> , for Innova® 43R only  | M1320-0300 | 2.156,00 €  |

## New Brunswick Innova® 42/42R

Eppendorf Product Catalog p. 290-291

|   | Order no.  | Price       |
|---|------------|-------------|
| <b>New Brunswick Innova® 42</b>   |            |             |
| stackable up to 2-fold, 230 V/50/60 Hz (CN/EU), orbit 1.9 cm (0.75 in)  | M1335-0002 | 11.230,00 € |
| stackable up to 2-fold, 230 V/50/60 Hz (CN/EU), orbit 2.5 cm (1 in)   | M1335-0012 | 11.230,00 € |
| <b>New Brunswick Innova® 42R</b>  |            |             |
| refrigerated, stackable up to 2-fold, 230 V/50/60 Hz (CN/EU), orbit 1.9 cm (0.75 in)  | M1335-0006 | 12.750,00 € |
| refrigerated, stackable up to 2-fold, 230 V/50/60 Hz (CN/EU), orbit 2.5 cm (1 in)   | M1335-0016 | 12.750,00 € |
| <b>Stacking kit</b> , for stacking Innova® 42/42R units   | M1335-0800 | 391,00 €    |
| <b>Stacking kit</b> , for stacking Innova® 4200/4230 units on top of Innova® 42/42R shakers   | M1335-0810 | 560,00 €    |
| <b>Remote alarm option</b> , for Innova® 44/44, Innova® 43/43, Innova® 42/42 or Innova® 40/40   | M1320-8029 | 420,00 €    |
| <b>Additional perforated shelf and brackets</b> , for Innova® 42/42R  | M1335-0080 | 260,00 €    |
| <b>Platforms for Innova® 2100/40/40R/42/42R, Excella® E24/E24R, 46 × 46 cm (18 × 18 in), phenolic resin (unless otherwise stated)</b>   |            |             |
| Universal platform, aluminum  | M1250-9902 | 888,00 €    |
| 125 mL Erlenmeyer flask dedicated platform  | M1194-9904 | 947,00 €    |
| 250 mL Erlenmeyer flask dedicated platform  | M1194-9905 | 981,00 €    |
| 500 mL Erlenmeyer flask dedicated platform  | M1194-9906 | 947,00 €    |
| 1 L Erlenmeyer flask dedicated platform   | M1194-9907 | 907,00 €    |
| 2 L Erlenmeyer flask dedicated platform, for 2 L wide neck Erlenmeyer flasks (DIN-ISO 24450), bottom diameter ~153 mm. For 2 L narrow neck Erlenmeyer, Thomson, or Corning flasks (DIN-ISO 1773) bottom diameter ~166 mm refer to universal platform and clamps M1190-9004. | M1194-9908 | 907,00 €    |
| <b>Sticky pad platform for Innova® 40(R)/42(R)/2100, Excella® E24(R)</b> , 46 × 46 cm (18 × 18 in), 4 sticky pads necessary to cover platform (sold separately, order no: M1250-9700), anodized aluminum  | M1250-9903 | 637,00 €    |
| <b>Quick-change platform kit</b> , enables platform attachment without tools, for Innova® 2100/4000/4080/4230   | M1192-9901 | 511,00 €    |
| <b>Culture drawer</b> , (height adjustable) prevents plates from drying   | M1335-0501 | 909,00 €    |
| <b>Utility carrier</b> , holds racks and vessels between x2 rubber-cushioned cross bars, for Innova® 2100/Innova® 40 / Innova® 42/Excella® 24   | M1194-9909 | 799,00 €    |
| <b>Rod kit</b> , for utility carrier (one additional rod), of Innova 40/42/2100 and Excella 24 lab shakers  | M1194-9923 | 249,00 €    |
| <b>Utility tray with non-skid rubber surface for Innova® 2100, 40, 42</b> , usable as a stand-alone – no other platform required  | M1194-9910 | 727,00 €    |
| <b>Drip pan</b> , to catch accidental spills, for universal platform M1250-9902   | M1250-9906 | 229,00 €    |
| <b>Photosynthetic light bank</b> , programmable, Innova® 42R only - contact Eppendorf for installation options  | M1335-0300 | 2.156,00 €  |
| <b>Photosynthetic light bank</b> , programmable, Innova® 42R only - contact Eppendorf for installation options  | M1335-0300 | 2.156,00 €  |
| <b>UV germicidal light</b> , programmable, Innova® 42R only - contact Eppendorf for installation options  | M1335-0400 | 2.156,00 €  |
| <b>UV germicidal light</b> , programmable, Innova® 42R only - contact Eppendorf for installation options  | M1335-0400 | 2.156,00 €  |
| <b>Humidity monitor</b> , for Innova® 42/42R - contact Eppendorf for installation options   | M1335-0500 | 1.080,00 €  |
| <b>Humidity monitor</b> , for Innova® 42/42R - contact Eppendorf for installation options   | M1335-0500 | 1.080,00 €  |
| <b>12-port gassing manifold</b> , for Innova® 42/42R - contact Eppendorf for installation options   | M1335-0600 | 798,00 €    |
| <b>12-port gassing manifold</b> , for Innova® 42/42R - contact Eppendorf for installation options   | M1335-0600 | 798,00 €    |

## Shaker

## New Brunswick Innova® 40/40R

Eppendorf Product Catalog p. 292-293

|   | Order no.  | Price       |
|---|------------|-------------|
| <b>New Brunswick Innova® 40</b>   |            |             |
| non-stackable, 230 V/50/60 Hz (CN/EU), orbit 1.9 cm (0.75 in)   | M1299-0082 | 9.667,00 €  |
| non-stackable, 230 V/50/60 Hz (CN/EU), orbit 2.5 cm (1 in)  | M1299-0092 | 9.667,00 €  |
| <b>New Brunswick Innova® 40R</b>  |            |             |
| refrigerated, non-stackable, 230 V/50/60 Hz (CN/EU), orbit 1.9 cm (0.75 in)   | M1299-0086 | 11.340,00 € |
| refrigerated, non-stackable, 230 V/50/60 Hz (CN/EU), orbit 2.5 cm (1 in)  | M1299-0096 | 11.340,00 € |
| <b>Remote alarm option</b> , for Innova® 44/44, Innova® 43/43, Innova® 42/42 or Innova® 40/40   | M1320-8029 | 420,00 €    |
| <b>Platforms for Innova® 2100/40/40R/42/42R, Excella® E24/E24R, 46 × 46 cm (18 × 18 in), phenolic resin (unless otherwise stated)</b>   |            |             |
| Universal platform, aluminum  | M1250-9902 | 888,00 €    |
| 125 mL Erlenmeyer flask dedicated platform  | M1194-9904 | 947,00 €    |
| 250 mL Erlenmeyer flask dedicated platform  | M1194-9905 | 981,00 €    |
| 500 mL Erlenmeyer flask dedicated platform  | M1194-9906 | 947,00 €    |
| 1 L Erlenmeyer flask dedicated platform   | M1194-9907 | 907,00 €    |
| 2 L Erlenmeyer flask dedicated platform, for 2 L wide neck Erlenmeyer flasks (DIN-ISO 24450), bottom diameter ~153 mm. For 2 L narrow neck Erlenmeyer, Thomson, or Corning flasks (DIN-ISO 1773) bottom diameter ~166 mm refer to universal platform and clamps M1190-9004. | M1194-9908 | 907,00 €    |
| <b>Sticky pad platform for Innova® 40(R)/42(R)/2100, Excella® E24(R)</b> , 46 × 46 cm (18 × 18 in), 4 sticky pads necessary to cover platform (sold separately, order no: M1250-9700), anodized aluminum  | M1250-9903 | 637,00 €    |
| <b>Quick-change platform kit</b> , enables platform attachment without tools, for Innova® 2100/4000/4080/4230   | M1192-9901 | 511,00 €    |
| <b>Utility carrier</b> , holds racks and vessels between x2 rubber-cushioned cross bars, for Innova® 2100/Innova® 40 / Innova® 42/Excella® 24   | M1194-9909 | 799,00 €    |
| <b>Rod kit</b> , for utility carrier (one additional rod), of Innova 40/42/2100 and Excella 24 lab shakers  | M1194-9923 | 249,00 €    |
| <b>Utility tray with non-skid rubber surface for Innova® 2100, 40, 42</b> , usable as a stand-alone – no other platform required  | M1194-9910 | 727,00 €    |
| <b>Drip pan</b> , to catch accidental spills, for universal platform M1250-9902   | M1250-9906 | 229,00 €    |

## Clamps and Accessories for Universal Platforms

Eppendorf Product Catalog p. 294

|   | Order no.  | Price    |
|---|------------|----------|
| <b>Test tube rack, small, 30 × 8 × 12.9 cm</b>  |            |          |
| diameter 8 – 11 mm, 48 tubes  | M1289-0001 | 345,00 € |
| diameter 12 – 15 mm, 34 tubes   | M1289-0002 | 345,00 € |
| diameter 15 – 18 mm, 24 tubes   | M1289-0003 | 345,00 € |
| diameter 18 – 21 mm, 18 tubes   | M1289-0004 | 345,00 € |
| diameter 22 – 26 mm, 13 tubes   | M1289-0005 | 345,00 € |
| diameter 26 – 30 mm, 12 tubes   | M1289-0006 | 345,00 € |
| <b>Test tube rack, medium, 36 × 8 × 12.9 cm</b>   |            |          |
| diameter 8 – 11 mm, 60 tubes  | M1289-0010 | 345,00 € |
| diameter 12 – 15 mm, 44 tubes   | M1289-0020 | 345,00 € |
| diameter 15 – 18 mm, 31 tubes   | M1289-0030 | 345,00 € |
| diameter 18 – 21 mm, 23 tubes   | M1289-0040 | 345,00 € |
| diameter 22 – 26 mm, 16 tubes   | M1289-0050 | 345,00 € |
| diameter 26 – 30 mm, 16 tubes   | M1289-0060 | 345,00 € |
| <b>Test tube rack, large, 45 × 8 × 12.9 cm</b>  |            |          |
| diameter 8 – 11 mm, 80 tubes  | M1289-0100 | 345,00 € |
| diameter 12 – 15 mm, 60 tubes   | M1289-0200 | 345,00 € |
| diameter 15 – 18 mm, 42 tubes   | M1289-0300 | 345,00 € |
| diameter 18 – 21 mm, 30 tubes   | M1289-0400 | 345,00 € |
| diameter 22 – 26 mm, 22 tubes   | M1289-0500 | 345,00 € |
| diameter 26 – 30 mm, 20 tubes   | M1289-0600 | 345,00 € |
| <b>Clamp, for 1 L media bottles</b>   | ACSB-1000S | 110,00 € |
| <b>Clamp, for 500 mL media bottles and TPP Bioreactors 450 mL or 600 mL</b>             | ACSB-500S  | 110,00 € |
| <b>Angled test tube rack holder, for user-supplied test tube racks and microplates</b>  | TTR-210    | 360,00 € |
| <b>Angled Test Tube Rack Spacer Bar, for use with TTR-210</b>                           |            |          |
| to accommodate test tubes racks that are less than 13 mm (5 in) wide                    | TTR-215    | 121,00 € |
| <b>Microplate rack, single layer, for single layer of five plates</b>                   | TTR-221    | 360,00 € |
| <b>Microplate rack, stack, stack three deepwell or nine standard plates</b>             | M1289-0700 | 345,00 € |
| <b>Sticky tape, 23 m × 3.8 cm (75 feet × 1.5 in) roll</b>                               | M1250-9600 | 300,00 € |
| <b>Sticky pad, 20 × 20 cm (8 × 8 in), Compatible platform or adapter kit necessary</b>  | M1250-9700 | 85,50 €  |
| <b>Drip pan, to catch accidental spills, for use with universal platform M1250-9920</b> | M1250-9921 | 249,00 € |
| <b>Drip pan, to catch accidental spills, for universal platform M1250-9902</b>          | M1250-9906 | 229,00 € |

## Clamps for universal platforms

Eppendorf Product Catalog p. 296-297

|   | Order no.  | Price    |
|---|------------|----------|
| <b>Clamp</b>  |            |          |
| for 10 mL Erlenmeyer flasks   | ACE-10S    | 21,10 €  |
| for 25 mL Erlenmeyer flasks   | M1190-9004 | 21,20 €  |
| for 50 mL Erlenmeyer flasks   | M1190-9000 | 22,30 €  |
| for 125 mL Erlenmeyer flasks  | M1190-9001 | 23,50 €  |
| for 250 mL Erlenmeyer flasks  | M1190-9002 | 33,50 €  |
| for 500 mL Erlenmeyer flasks  | M1190-9003 | 39,60 €  |
| for 1 L Erlenmeyer flasks   | ACE-1000S  | 52,10 €  |
| for 2 L wide neck Erlenmeyer flasks (DIN-ISO 24450)                         | ACE-2000S  | 81,80 €  |
| for 3 L Erlenmeyer flasks   | ACE-3000S  | 95,00 €  |
| for 4 L Erlenmeyer flasks   | ACE-4000S  | 93,40 €  |
| for 5 L Erlenmeyer flasks   | ACE-5000S  | 101,00 € |
| for 6 L Erlenmeyer flasks   | ACE-6000S  | 112,00 € |
| for 2 L narrow neck Erlenmeyer (DIN-ISO 1773), Thomson, or Corning flasks   | M1190-9005 | 119,00 € |
| for 2.8 L Fernbach flasks   | ACFE-2800S | 104,00 € |
| <b>Clamp, for 1 L media bottles</b>   | ACSB-1000S | 110,00 € |
| <b>Clamp, for 500 mL media bottles and TPP Bioreactors 450 mL or 600 mL</b> | ACSB-500S  | 110,00 € |

## Centrifuges and Rotors

### MiniSpin®/MiniSpin® plus

Eppendorf Product Catalog p. 336

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Centrifuge MiniSpin® (EU-IVD)</b> , keypad, non-refrigerated, with Rotor F-45-12-11, 230 V/50 – 60 Hz (EU)                       | 5452 000 018 | 1.142,00 € |
| <b>Centrifuge MiniSpin® plus (EU-IVD)</b> , keypad, non-refrigerated, with Rotor F-45-12-11, 230 V/50 – 60 Hz (EU)                  | 5453 000 011 | 1.484,00 € |
| <b>Rotor F-45-12-11</b> , angle 45°, 12 places, max. tube diameter 11 mm, incl. rotor lid and rotor nut, MiniSpin®/MiniSpin® plus   | 5452 720 008 | 266,00 €   |
| <b>Rotor F-55-16-5 PCR</b> , angle 55°, 16 places, max. tube diameter 5 mm, incl. rotor lid (aluminum), MiniSpin®/MiniSpin® plus    | 5452 727 007 | 298,00 €   |
| <b>Adapter</b> , for 1 PCR tube (0.2 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.  | 5425 715 005 | 15,70 €    |
| <b>Adapter</b> , for 1 sample tube (0.5 mL, max. Ø 8 mm) or 1 Microtainer® (0.6 mL, max. Ø 8 mm), for all 1.5/2.0 mL rotors, 6 pcs. | 5425 716 001 | 15,70 €    |
| <b>Adapter</b> , for 1 micro test tube (0.4 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.                                     | 5425 717 008 | 15,70 €    |

### Centrifuge 5418 and Centrifuge 5418 R

Eppendorf Product Catalog p. 337

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Centrifuge 5418 R (EU-IVD)</b> , rotary knobs, refrigerated, with Rotor FA-45-18-11 and rotor lid, 230 V/50 – 60 Hz (EU)         | 5401 000 013 | 5.204,00 € |
| <b>Adapter</b> , for 1 PCR tube (0.2 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.  | 5425 715 005 | 15,70 €    |
| <b>Adapter</b> , for 1 sample tube (0.5 mL, max. Ø 8 mm) or 1 Microtainer® (0.6 mL, max. Ø 8 mm), for all 1.5/2.0 mL rotors, 6 pcs. | 5425 716 001 | 15,70 €    |
| <b>Adapter</b> , for 1 micro test tube (0.4 mL, max. Ø 6 mm), for all 1.5/2.0 mL rotors, 6 pcs.                                     | 5425 717 008 | 15,70 €    |

### Centrifuge 5420

Eppendorf Product Catalog p. 337

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Centrifuge 5420 (EU-IVD)</b> , keypad, non-refrigerated, without rotor, 230 V/50/60 Hz (EU)                      | 5420 000 210 | 1.881,00 € |
| <b>Centrifuge 5420 (EU-IVD)</b> , keypad, non-refrigerated, with aerosol-tight Rotor FA-24x2, 230 V/50 – 60 Hz (EU) | 5420 000 016 | 2.292,00 € |

### Centrifuge 5425/5425 R

Eppendorf Product Catalog p. 338

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Centrifuge 5425 (EU-IVD)</b> , keypad, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU)  | 5405 000 212 | 2.340,00 € |
| <b>Centrifuge 5425 (EU-IVD)</b> , rotary knobs, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU)                                      | 5405 000 611 | 2.340,00 € |
| <b>Centrifuge 5425 (EU-IVD)</b> , keypad, non-refrigerated, with Rotor FA-24x2, 230 V/50 – 60 Hz (EU)                                       | 5405 000 018 | 2.902,00 € |
| <b>Centrifuge 5425 (EU-IVD)</b> , rotary knobs, non-refrigerated, with Rotor FA-24x2, 230 V/50 – 60 Hz (EU)                                 | 5405 000 417 | 2.902,00 € |
| <b>Centrifuge 5425 (EU-IVD)</b> , rotary knobs, non-refrigerated, Molecular Biology Solution, with Rotor FA-18x2-Kit, 230 V/50 – 60 Hz (EU) | 5405 000 158 | 2.639,00 € |
| <b>Centrifuge 5425 R (EU-IVD)</b> , keypad, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)  | 5406 000 216 | 5.485,00 € |
| <b>Centrifuge 5425 R (EU-IVD)</b> , rotary knobs, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)  | 5406 000 615 | 5.485,00 € |
| <b>Centrifuge 5425 R (EU-IVD)</b> , keypad, refrigerated, with Rotor FA-24x2, 230 V/50 – 60 Hz (EU)   | 5406 000 011 | 5.895,00 € |
| <b>Centrifuge 5425 R (EU-IVD)</b> , rotary knobs, refrigerated, with Rotor FA-24x2, 230 V/50 – 60 Hz (EU)                                   | 5406 000 410 | 5.895,00 € |

### Centrifuge 5427 R

**New**

Eppendorf Product Catalog p. 339

|  | Order no.    | Price   |
|--|--------------|---------|
| <b>Centrifuge 5427 R</b> , rotary knobs, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)          | 5429 000 015 | Inquire |
| <b>Centrifuge 5427 R</b> , rotary knobs, refrigerated, with Rotor FA-45-30-11, 230 V/50 – 60 Hz (EU) | 5429 000 315 | Inquire |
| <b>Centrifuge 5427 R</b> , rotary knobs, refrigerated, with Rotor FA-45-48-11, 230 V/50 – 60 Hz (EU) | 5429 000 215 | Inquire |



## Centrifuge 5430 and Centrifuge 5430 R

Eppendorf Product Catalog p. 340

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Centrifuge 5430 (EU-IVD)</b> , keypad, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU)                                | 5427 000 011 | 3.682,00 € |
| <b>Centrifuge 5430 (EU-IVD)</b> , rotary knobs, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU)                          | 5427 000 615 | 3.682,00 € |
| <b>Centrifuge 5430 (EU-IVD)</b> , keypad, non-refrigerated, with Rotor FA-45-30-11 incl. rotor lid, 230 V/50 – 60 Hz (EU)       | 5427 000 216 | 4.154,00 € |
| <b>Centrifuge 5430 (EU-IVD)</b> , rotary knobs, non-refrigerated, with Rotor FA-45-30-11 incl. rotor lid, 230 V/50 – 60 Hz (EU) | 5427 000 410 | 4.154,00 € |
| <b>Centrifuge 5430 R (EU-IVD)</b> , keypad, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)                                  | 5428 000 210 | 7.491,00 € |
| <b>Centrifuge 5430 R (EU-IVD)</b> , rotary knobs, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)                            | 5428 000 619 | 7.491,00 € |
| <b>Centrifuge 5430 R (EU-IVD)</b> , keypad, refrigerated, with Rotor FA-45-30-11 incl. rotor lid, 230 V/50 – 60 Hz (EU)         | 5428 000 015 | 7.773,00 € |
| <b>Centrifuge 5430 R (EU-IVD)</b> , rotary knobs, refrigerated, with Rotor FA-45-30-11 incl. rotor lid, 230 V/50 – 60 Hz (EU)   | 5428 000 414 | 7.773,00 € |

## Centrifuge 5702/5702 R/5702 RH

Eppendorf Product Catalog p. 341

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>Centrifuge 5702 (EU-IVD)</b> , rotary knobs, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU)  | 5702 000 019 | 2.037,00 € |
| <b>Centrifuge 5702 (EU-IVD)</b> , rotary knobs, non-refrigerated, with Rotor A-4-38 incl. adapters for 15/50 mL conical tubes, 2 sets of 2 adapters, 230 V/50 – 60 Hz (EU)  | 5702 000 329 | 3.350,00 € |
| <b>Centrifuge 5702 R (EU-IVD)</b> , rotary knobs, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)  | 5703 000 012 | 4.831,00 € |
| <b>Centrifuge 5702 R (EU-IVD)</b> , rotary knobs, refrigerated, with Rotor A-4-38 incl. adapters for 13/16 mm blood collection tubes, 2 sets of 2 adapters, 230 V/50 – 60 Hz (EU)   | 5703 000 322 | 6.012,00 € |
| <b>Centrifuge 5702 R (EU-IVD)</b> , rotary knobs, refrigerated, with Rotor A-4-38 incl. adapters for 15/50 mL conical tubes, 2 sets of 2 adapters, 230 V/50 – 60 Hz (EU)  | 5703 000 420 | 5.865,00 € |
| <b>Centrifuge 5702 RH (EU-IVD)</b> , rotary knobs, refrigerated/heated, without rotor, 230 V/50 – 60 Hz (EU)  | 5704 000 016 | 5.461,00 € |
| <b>Centrifuge 5702 (EU-IVD)</b> , rotary knobs, non-refrigerated, with Rotor A-4-38 incl. adapters for round-bottom tubes 4 – 10 mL and adapters for round-bottom tubes 2.6 – 7 mL, 2 sets of 2 adapters, 230 V/50 – 60 Hz (EU) | 5702 000 510 | 3.440,00 € |

## Centrifuge 5804 and Centrifuge 5804 R

Eppendorf Product Catalog p. 342

|  | Order no.    | Price      |
|--|--------------|------------|
| <b>Centrifuge 5804 (EU-IVD)</b> , keypad, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU)   | 5804 000 013 | 4.294,00 € |
| <b>Centrifuge 5804 (EU-IVD)</b> , keypad, non-refrigerated, with Rotor A-4-44 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU) | 5804 000 323 | 5.740,00 € |
| <b>Centrifuge 5804 (EU-IVD)</b> , keypad, non-refrigerated, with Rotor S-4-72 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU) | 5804 000 528 | 6.896,00 € |
| <b>Centrifuge 5804 R (EU-IVD)</b> , keypad, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)   | 5805 000 017 | 7.629,00 € |
| <b>Centrifuge 5804 R (EU-IVD)</b> , keypad, refrigerated, with Rotor A-4-44 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU)   | 5805 000 327 | 8.743,00 € |
| <b>Centrifuge 5804 R (EU-IVD)</b> , keypad, refrigerated, with Rotor S-4-72 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU)   | 5805 000 629 | 9.608,00 € |

## Centrifuges and Rotors

### Centrifuge 5810 and Centrifuge 5810 R

Eppendorf Product Catalog p. 342

|  | Order no.    | Price       |
|--|--------------|-------------|
| <b>Centrifuge 5810 (EU-IVD)</b> , keypad, non-refrigerated, without rotor, 230 V/50 – 60 Hz (EU)   | 5810 000 017 | 5.638,00 €  |
| <b>Centrifuge 5810 (EU-IVD)</b> , keypad, non-refrigerated, with Rotor A-4-62 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU)       | 5810 000 327 | 7.520,00 €  |
| <b>Centrifuge 5810 (EU-IVD)</b> , keypad, non-refrigerated, with Rotor A-4-81 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU)       | 5810 000 424 | 8.454,00 €  |
| <b>Centrifuge 5810 (EU-IVD)</b> , keypad, non-refrigerated, with Rotor S-4-104 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU)      | 5810 000 521 | 9.052,00 €  |
| <b>Centrifuge 5810 (EU-IVD)</b> , keypad, non-refrigerated, with Rotor S-4-104 incl. adapters for 13/16 mm round-bottom tubes, 230 V/50 – 60 Hz (EU) | 5810 000 629 | 9.052,00 €  |
| <b>Centrifuge 5810 R (EU-IVD)</b> , keypad, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)   | 5811 000 010 | 9.449,00 €  |
| <b>Centrifuge 5810 R (EU-IVD)</b> , keypad, refrigerated, with Rotor A-4-62 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU)         | 5811 000 320 | 10.950,00 € |
| <b>Centrifuge 5810 R (EU-IVD)</b> , keypad, refrigerated, with Rotor A-4-81 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU)         | 5811 000 428 | 11.890,00 € |
| <b>Centrifuge 5810 R (EU-IVD)</b> , keypad, refrigerated, with Rotor S-4-104 incl. adapters for 15/50 mL conical tubes, 230 V/50 – 60 Hz (EU)        | 5811 000 622 | 12.480,00 € |
| <b>Centrifuge 5810 R (EU-IVD)</b> , keypad, refrigerated, with Rotor S-4-104 incl. adapters for 13/16 mm round-bottom tubes, 230 V/50 – 60 Hz (EU)   | 5811 000 827 | 12.480,00 € |

### Centrifuge 5910 Ri

Eppendorf Product Catalog p. 343

|   | Order no.    | Price       |
|---|--------------|-------------|
| <b>Centrifuge 5910 Ri (EU-IVD)</b> , touch interface, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)  | 5943 000 017 | 11.150,00 € |
| <b>Centrifuge 5910 Ri (EU-IVD)</b> , touch interface, refrigerated, with Rotor S-4x750, incl. round buckets and adapters for 5 mL/15 mL/50 mL conical tubes, 230 V/50 – 60 Hz (EU)  | 5943 000 511 | 13.870,00 € |
| <b>Centrifuge 5910 Ri (EU-IVD)</b> , touch interface, refrigerated, with Rotor S-4xUniversal, incl. universal buckets and adapters for 5 mL/15 mL/50 mL conical tubes and plates, 230 V/50 – 60 Hz (EU)                       | 5943 000 319 | 14.490,00 € |
| <b>Centrifuge 5910 Ri (EU-IVD)</b> , touch interface, refrigerated, with Rotor S-4xUniversal, incl. universal buckets and adapters for 5 mL/15 mL/50 mL conical tubes, 175 – 250 mL bottles and plates, 230 V/50 – 60 Hz (EU) | 5943 000 211 | 14.600,00 € |

### Centrifuge 5920 R

Eppendorf Product Catalog p. 344

|   | Order no.    | Price       |
|---|--------------|-------------|
| <b>Centrifuge 5920 R (EU-IVD)</b> , keypad, refrigerated, without rotor, 230 V/50 – 60 Hz (EU)  | 5948 000 018 | 12.070,00 € |
| <b>Centrifuge 5920 R (EU-IVD)</b> , keypad, refrigerated, with Rotor S-4x1000, incl. high-capacity buckets and adapters for 15 mL/50 mL conical tubes, 230 V/50 – 60 Hz (EU)                              | 5948 000 212 | 15.210,00 € |
| <b>Centrifuge 5920 R (EU-IVD)</b> , keypad, refrigerated, with Rotor S-4x1000, incl. plate/tube buckets and adapters for 15 mL/50 mL conical tubes, 230 V/50 – 60 Hz (EU)                                 | 5948 000 417 | 14.710,00 € |
| <b>Centrifuge 5920 R (EU-IVD)</b> , keypad, refrigerated, with Rotor S-4x1000, incl. round buckets and adapters for 15 mL/50 mL conical tubes, 230 V/50 – 60 Hz (EU)                                      | 5948 000 310 | 14.710,00 € |
| <b>Centrifuge 5920 R (EU-IVD)</b> , keypad, refrigerated, with Rotor S-4xUniversal-Large, incl. universal buckets and adapters for round-bottom tubes with diameter 13/16 mm, 230 V/50 – 60 Hz (EU)       | 5948 000 816 | 15.560,00 € |
| <b>Centrifuge 5920 R (EU-IVD)</b> , keypad, refrigerated, with Rotor S-4xUniversal-Large, incl. universal buckets and adapters for 5 mL/15 mL/25 mL/50 mL conical tubes and plates, 230 V/50 – 60 Hz (EU) | 5948 000 914 | 15.900,00 € |

## TransferMan® 4m

Eppendorf Product Catalog p. 345-346

|  | Order no.    | Price       |
|--|--------------|-------------|
| <b>TransferMan® 4m</b> , micromanipulator with DualSpeed™ joystick for direct and dynamic movement control, 100 – 240 V/50 – 60 Hz (EU), The TransferMan® 4m has been designed and manufactured to assist in the positioning of micro tools (e.g., microcapillaries or micro glass pipettes) under the microscope, during the performance of assisted reproductive technologies (ART) such as intracytoplasmic sperm injection (ICSI), or the collection of genetic material for the purpose of genetic preimplantation testing (PGT). | 5191 000 015 | 18.270,00 € |

## CellTram® 4m Air/Oil

Eppendorf Product Catalog p. 347-348

|   | Order no.    | Price      |
|---|--------------|------------|
| <b>CellTram® 4m Air</b> , manual pneumatic microinjector, with gears 1:1 and 1:10, for holding and injection, The CellTram® 4m Air has been designed and manufactured to support the injection and aspiration of cells or cell components during the performance of assisted reproductive technologies (ART) such as intracytoplasmic sperm injection (ICSI), or the collection of genetic material for the purpose of genetic preimplantation testing (PGT), using micro tools (e.g., microcapillaries or micro glass pipettes). | 5196 000 021 | 2.275,00 € |
| <b>CellTram® 4m Oil</b> , manual hydraulic microinjector, with gears 1:1 and 1:10, for holding and injection, The CellTram® 4m Oil has been designed and manufactured to support the injection and aspiration of cells or cell components during the performance of assisted reproductive technologies (ART) such as intracytoplasmic sperm injection (ICSI), or the collection of genetic material for the purpose of genetic preimplantation testing (PGT), using micro tools (e.g., microcapillaries or micro glass pipettes). | 5196 000 048 | 3.345,00 € |
| <b>Filling and Cleaning set</b> , for CellTram® 4 Oil and CellTram® Oil/vario, incl. filling tube, Luer lock adapter, 2 syringes  | 5196 088 000 | 78,00 €    |
| <b>Injection tube Air</b> , white ring mark, I.D. 0.5 mm, length 1.3 m  | 5196 061 004 | 71,90 €    |
| <b>Injection tube Oil</b> , blue ring mark, I.D. 1.0 mm, length 1.3 m   | 5196 089 006 | 71,90 €    |
| <b>Capillary holder 4</b> , incl. 1x Grip head set 4 size 0 (for microcapillaries with outer diameter 1.0 – 1.1 mm)   | 5196 081 005 | 233,00 €   |
| <b>O-ring set 4</b> , incl. 10 o-rings large, 10 o-rings small, 2 distance sleeves, o-ring removal tool for grip head set 4   | 5196 086 007 | 47,50 €    |
| <b>Grip head set 4 size 0</b> , 1 pc, for microcapillaries with outer diameter 1.0 to 1.1 mm  | 5196 082 001 | 47,50 €    |
| <b>Grip head set 4 size 1</b> , 1 pc, for microcapillaries with outer diameter 1.2 to 1.3 mm  | 5196 083 008 | 47,50 €    |
| <b>Grip head set 4 size 2</b> , 1 pc, for microcapillaries with outer diameter 1.4 to 1.5 mm  | 5196 084 004 | 47,50 €    |
| <b>Grip head set 4 size 3</b> , 1 pc, for microcapillaries with outer diameter 0.7 to 0.9 mm  | 5196 085 000 | 47,50 €    |

## Application and Technical support



Your local distributor:  
[www.eppendorf.com/contact](http://www.eppendorf.com/contact)



Eppendorf Online Shop  
[www.eppendorf.com/eshop](http://www.eppendorf.com/eshop)

Eppendorf Bioprocessing-Catalog. More info at [www.eppendorf.com/catalog](http://www.eppendorf.com/catalog)

