

# Position of the Barcode Label

For all plate types, the label is located in the middle of the side(s) defined. It is not possible to position a label on the top of the plate. The height of the label is stated in the table below:

Plate	Total height of plate	Position of label *
twin.tec PCR Plates 96 skirted	6328006500	2.5 mm (top)
twin.tec PCR Plates 96 semi-skirted	6328007000	1 mm (top)
twin.tec PCR Plates 384	6328007507	1 mm (top)
Microplates	6328008007	7 mm (base)
Deepwell Plates 96/500 µL	6328008503	8.5 mm (base)
Deepwell Plates 96/1000 µL	6328009003	8.5 mm (base)
Deepwell Plates 96/2000 µL	6328009500	8.5 mm (base)
Deepwell Plates 384/200 µL	0056000003	8.5 mm (base)

\*Distance from the top of the plate to the top of the barcode label (top) from the base of the plate to the base of the barcode label (base)

Technical drawings of the plates with all relevant dimensions can be downloaded as a PDF file from the relevant product page on [www.eppendorf.com](http://www.eppendorf.com).

**Your local distributor:** [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)

Eppendorf® and the Eppendorf Brand Design are registered trademarks of Eppendorf SE, Germany. All rights reserved, including graphics and images. Copyright © 2022 by Eppendorf SE.