

Technical Data Sheet

Eppendorf_Microplate UV-VIS, 96/F
English (EN)




Plate dimensions [mm]

Length		Width		Height		
l	l1	w	w1	h	h1	fh
127.8	123.3	85.5	81.2	14.4	0,4	6.5

Well dimensions [mm]

Diameter		Thickness	Depth	Positions of well center [mm]			
d	d1	bt	h2	P1	P2	P3	P4
6.5	6.0	0.025	14.0	9.0	11.2	14.4	9.0

Stack height [mm]

without lid (sh)

28.4

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · 22331 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com