

Technical Data Sheet

Eppendorf Deepwell Plate 384/200 μ L

English (EN)




Plate dimensions [mm]

Length			Width			Height		Stack height
l	l1	l2	w	w1	w2	h	fh	sh (without lid)
127.8	124.5	123.2	85.5	82.2	81.0	25.1	6.1	48.8

Well dimensions [mm]

Depth				Well width	Bottom [°]	Positions of well center [mm]				Bottom space [mm]	
h1	h2	d	d1	a	Angle	P1	P2	P3	P4	b	b1
21.3	0,8	3.8	3.4	120	120	12.1	4.5	8.99	4.5	1.8	1.4

Your local distributor: www.eppendorf.com/contact

Eppendorf AG · 22331 Hamburg · Germany

eppendorf@eppendorf.com · www.eppendorf.com