

epServices  
for premium performance



# epMotion®

## epMotion Performance Plans

epMotion and its dispensing tools are precision instruments. As such, they should be regularly maintained and calibrated, which will help ensure continuing high accuracy, precision and reproducibility.

With our Performance Plans options we help you maintain premium performance of your epMotion over an extended period of time. Our qualified service technicians will take care of your instrument.

Our Certification Services include all function tests, calibration and adjustment services needed for Installation and Operational Qualification (IQ/OQ) of your instrument.

## Dispensing Tools Calibration

Calibration of the dispensing tools are based on ISO 8655, supporting your GLP documentation.

## Value of service

- > Confidence in high quality results and instrument reliability
- > Simple service process
- > Professional maintenance

**Keep your process in motion!**

eppendorf

Performance tested on

Model: \_\_\_\_\_

Serial no.: \_\_\_\_\_

Serviced by: \_\_\_\_\_

Service no.: \_\_\_\_\_

Next service: \_\_\_\_\_

Date: \_\_\_\_\_

epServices  
for premium performance

# epMotion® Performance Plans

Service Operation	ESSENTIAL CHECK	ADVANCED MAINTENANCE	PREMIUM SERVICE	Installation Qualification (IQ)	Operational Qualification (OQ)
<b>Order Number for epMotion® 5070</b>	0082 030 013	0082 030 014	0082 030 015	0082 030 017	0082 030 018
<b>Order Number for epMotion® 5073</b>	0082 030 033	0082 030 034	0082 030 035	0082 030 037	0082 030 038
<b>Order Number for epMotion® 5075</b>	0082 030 023	0082 030 024	0082 030 025	0082 030 027	0082 030 028
<b>Cleaning and Maintenance</b>					
Cleaning of housing/accessories	■	■	■		
Checking epMotion® for external visible damages	■	■	■	■	■
Exchange of CAN bus cables/O-rings of dispensing tools		■	■		
Software update*	■	■	■	■	■
<b>Functional Check</b>					
Instrument set-up, initial operation, instruction of user				■	
Tip ejector, level sensor and carrier	■	■	■		■
Barcode Reader with epBlue™ ID*	■	■	■		
Thermal modules* and vacuum chamber*	■	■	■	■	■
TMX module*/TMX module w. Magnetic Finger module*	■	■	■	■	■
UV light/CleanCap (HEPA Filter)*	■	■	■	■	■
<b>Memory</b>					
Evaluation of stored data		■	■		
<b>Verification</b>					
Calibration/adjustment of level sensor and carrier		■	■	■	■
Calibration/adjustment of vacuum chamber*		■	■		
Verification of operating parameters acc. to specifications					■***
<b>Documentation</b>					
Check List provided	■	■	■		
Dated service sticker to confirm Eppendorf service	■	■	■	■	■
IQ Report and Certificate				■	
OQ Report and Certificate					■
<b>Supporting Information</b>					
Contract period	one year	one year	one year	n/a	n/a
Number of preventive services included	one	one	one	n/a	n/a
Cost of repairs/parts replacement outside scope of preventive maintenance visit (where Eppendorf product warranty has expired)	not included	not included	discount on parts, labor, travel time**	n/a	n/a

\* Where feature is present.

\*\* Please contact your local Eppendorf Service representative for details of discounts and/or any potential additional charges applying. Additionally, if you require a Customized Service, exchanges of spare parts, we will be very pleased to discuss your requirements.

\*\*\* The detailed certification service contents are listed in the OQ Report.

Performance Plans are available in selected countries only and service offer may differ.

## Order number for combined Installation and Operational Qualification (IQ/OQ) services:

- > For epMotion® 5070: 0082 030 019
- > For epMotion® 5073: 0082 030 039
- > For epMotion® 5075: 0082 030 029

Please visit [www.eppendorf.com/epServices](http://www.eppendorf.com/epServices) for information on services for epMotion GxP and epMotion 96.

## Order numbers for additional services:

- > 1-chan. dispensing tool calibration (< epBlue™ 40.0): 5070 005 036
- > 8-chan. dispensing tool calibration (< epBlue™ 40.0): 5070 005 044
- > 1-chan. dispensing tool calibration (≥ epBlue™ 40.0): 0050 600 255
- > 8-chan. dispensing tool calibration (≥ epBlue™ 40.0): 0050 600 259
- > Thermal module calibration: 5075 005 018
- > TMX Module calibration: 5075 005 174
- > HEPA filter exchange: 5073 920 061
- > UV lamp exchange: 5073 920 070

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)

Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com/epServices](http://www.eppendorf.com/epServices)