

Hamburg, July 1st, 2021

Statement on the Relaunch of epT.I.P.S.[®] Packaging Option “Box 2.0”

Product Change Notification Letter (GLP & EU-IVD) – Applies to Europe

Dear valued customer,

This letter is to inform you about the design update of the epT.I.P.S.[®] Box by Eppendorf. The newly introduced Box is named “epT.I.P.S.[®] Box 2.0”.

The design update affects all three Box sizes available for the epT.I.P.S. pipette tips products. New item numbers will be introduced for the epT.I.P.S. Box 2.0 replacing all former box item numbers (Boxes and Sets). A list of predecessor and new item numbers can be found in **Table 1** at the end of this notification letter.

The following design and functionality features have been optimized:

New design features:

- > Fully enclosed containment on the back of boxes (size 1 and 2) and a narrow slot on the back of box 3 only to prevent hidden contamination
- > Lid with recesses for non-slip stackability of boxes during transport
- > New color of box bottom part: winter white
- > New wave design of box button

Additional functionalities:

- > Light touch button for ergonomic handling
- > The new epT.I.P.S. Box 2.0 ensures compatibility with all existing epT.I.P.S. Reloads and future product line extensions

Due to the design update the outer dimensions of the boxes will vary in comparison to the predecessor boxes:

Box sizes	epT.I.P.S. [®] Box, Predecessor (LxWxH)	epT.I.P.S. [®] Box 2.0, New (LxWxH)
Box size 1	138.8 × 106.3 × 65.5 mm	139.9 × 111.3 × 73 mm
Box size 2	140.2 × 109.8 × 85.3 mm	141.8 × 113.9 × 91 mm
Box size 3	144.1 × 119.4 × 140.1 mm	147.1 × 118.9 × 142 mm

There are no changes to the base material of boxes (Polycarbonate of identical type).

However, the design and functionality of the included pipette tips remain unchanged as the tips do not differ in handling nor in any other specification or property.

The products containing the new Box 2.0 will be available starting Q3, 2021 with no changes to the existing pricing. During transitioning period towards the new Box 2.0 design we will work closely with you to make the process as smooth as possible. We greatly appreciate your support and loyalty.

The statement applies to the following products in Europe:

	Volumes of included Tips in μ L	Current Item No. GLP predecessor	New Item No. GLP	Current Item No. IVD predecessor	New Item No. IVD
epT.I.P.S.[®] Box > epT.I.P.S.[®] Box 2.0	0.1-20 / 384	N/A	N/A	0030 076.087	0030 076 370
	0.5-100 / 384	N/A	N/A	0030 076.095	0030 076 389
	0.1-10	0030 073.010	0030 076 125	0030 073.002	0030 076 397
	0.1-20	0030 073.037	0030 076 133	0030 073.029	0030 076 400
	0.5-20	0030 073.053	0030 076 141	0030 073.045	0030 076 419
	200	0030 073.070	0030 076 150	0030 073.061	0030 076 427
	300	0030 073.096	0030 076 168	0030 073.088	0030 076 435
	1,000	0030 073.118	0030 076 176	0030 073.100	0030 076 443
	1,250	0030 073.134	0030 076 184	0030 073.126	0030 076 451
	1,250L	0030 073.657	0030 076 192	0030 073.622	0030 076 460
	2.5	0030 073.150	0030 076 206	0030 073.142	0030 076 478
	5	0030 073.177	0030 076 214	0030 073.169	0030 076 486
epT.I.P.S.[®] Set	0.1-20 / 384	0030 076.028	0030 076 273	N/A	N/A
	0.5-100 / 384	0030 076.036	0030 076 281	N/A	N/A
	0.1-10	0030 073.215	0030 076 290	0030 073.207	0030 076 494
	0.1-20	0030 073.231	0030 076 303	0030 073.223	0030 076 508
	0.5-20	0030 073.258	0030 076 311	0030 073.240	0030 076 516
	200	0030 073.274	0030 076 320	0030 073.266	0030 076 524
	300	0030 073.290	0030 076 338	0030 073.282	0030 076 532
	1,000	0030 073.312	0030 076 346	0030 073.304	0030 076 540
	1,250	0030 073.339	0030 076 354	0030 073.320	0030 076 559
	1,250L	N/A	N/A	N/A	N/A
	2,5	0030 073.355	0030 076 362	0030 073.347	0030 076 567

Table 1: List of predecessor and new item numbers.

Your local distributor: www.eppendorf.com/contact
 Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com