

Technical Data Sheet

Eppendorf Deepwell Plate 96/1000 μ L

English (EN)




Plate dimensions [mm]

Length			Width			Height		Stack height
l	l1	l2	w	w1	w2	h	fh	sh (without lid)
127.8	124.5	123.2	85.5	82.2	81.0	44.1	6.1	86.8

Well dimensions [mm]

Well dimensions [mm]				Bottom [°]	Positions of well center [mm]				Bottom space [mm]	
Depth		Diameter		Angle						
h1	h2	d	d1	a	P1	P2	P3	P4	b	b1
38.7	2.4	7.0	6.0	90	14.4	9.0	11.2	9.0	3.2	1.4

Your local distributor: www.eppendorf.com/contact

Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
 eppendorf@eppendorf.com · www.eppendorf.com