


**Technical Data Sheet**  
Eppendorf Microplate 384/F  
English (EN)



**Plate dimensions [mm]**

Length			Width			Height		Stack height
$l$	$l_1$	$l_2$	$w$	$w_1$	$w_2$	$h$	$fh$	$sh$ (without lid)
127.8	124.8	123.4	85.5	83.1	81.5	14.35	6.1	27.2

**Well dimensions [mm]**

Well dimensions [mm]				Positions of well center [mm]				Bottom space [mm]	
Depth		Well width							
$h_1$	$h_3$	$d$	$d_1$	$P_1$	$P_2$	$P_3$	$P_4$	$b$	$b_1$
11.55	2.8	3.8	3.6	12.1	4.5	9.0	4.5	2.2	1.4

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
Eppendorf AG · 22331 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)