


faster
safer
ultra

Eppendorf Ultracentrifuges



Ultra Your Workflow

Discover the Eppendorf way of ultracentrifugation.

Ultracentrifuges are indispensable tools for modern laboratories, offering powerful separation and purification capabilities. They play a vital role in various applications and workflows, enabling efficient isolation of very small particles. Did you know that Eppendorf Ultracentrifuges are specifically designed to enhance your daily lab efficiency?

We believe, as indispensable tools for many laboratories, ultracentrifuges must adhere to the highest standards of both convenience and safety.

Learn how Eppendorf ultracentrifuges accelerate your workflow, maximize user comfort and increase lab safety.

Get in touch!



Defining the Ultra Standard

Smart features for your lab's advantage



Ultra durable –

Extend rotor longevity with Rotor Lifetime Management

Experience the advanced capabilities of Rotor Lifetime Management (RLM), designed to optimize the lifespan of your ultracentrifuge rotors. RLM detects rotor ID and meticulously records actual run time and speed, ensuring your equipment performs at its best for longer.

Ultra convenient –

Accelerate your workflow with Visual Balancing

Safe time with Visual Balancing - a feature designed to simplify your workflow and enhance safety. Say goodbye to the tedious process of sample balancing with a lab scale.



Discover our ultracentrifuge portfolio
www.eppendorf.link/ultracentrifuge-portfolio



Your local distributor: www.eppendorf.com/contact

Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.link/ultra

Eppendorf and the Eppendorf Brand Design are registered trademarks of Eppendorf SE, Germany. All rights reserved, including graphics and images. AA0X XXX XX/EN1/PDF/1025/ME. Copyright 2025 by Eppendorf SE.