

Eppendorf Certificate

Certificate of Quality

EPA/HEPA* Standard of ep Dualfilter T.I.P.S.®

For the filters in the following pipette tips we hereby confirm a filter efficiency of class ISO 25 E in accordance with DIN EN ISO 29463-5.

Volume	Order no.		
	old (international)	old (North America)	new
ep Dualfilter T.I.P.S.® Racks PCR clean and sterile			
0.1 - 10 µL S, 34 mm	0030 077.504	022491202	0030 078 500
0.1 - 10 µL M, 40 mm	0030 077.512	022491211	0030 078 519
0.5 - 20 µL L, 46 mm	0030 077.520	022491229	0030 078 527
2 - 20 µL, 53 mm	0030 077.539	022491270	0030 078 535
2 - 100 µL, 53 mm	0030 077.547	022491237	0030 078 543
2 - 200 µL, 55 mm	0030 077.555	022491296	0030 078 551
20 - 300 µL, 55 mm	0030 077.563	022491245	0030 078 560
50 - 1000 µL, 76 mm	0030 077.571	022491253	0030 078 578
50 - 1250 µL L, 103 mm	0030 077.750	022494002	0030 078 594
0.25 - 2500 µL, 115 mm	-	022491881	0030 078 586
0.2 - 5000 µL L, 175 mm	0030 077.725	0030077725	0030 078 624
100 - 5000 µL, 120 mm	0030 077.580	022491261	0030 078 616
0.5 - 10 mL L, 243 mm	0030 077.598	022491288	-
ep Dualfilter T.I.P.S.® Racks Forensic Grade			
0.1 - 10 µL M, 40 mm	0030 077.768	0030077768	0030 078 810
2 - 20 µL, 53 mm	0030 077.776	0030077776	0030 078 829
2 - 200 µL, 55 mm	0030 077.784	0030077784	0030 078 837
50 - 1000 µL, 76 mm	0030 077.792	0030077792	0030 078 845
ep Dualfilter T.I.P.S.® LoRetention® Racks PCR clean and sterile			
0.1 - 10 µL S, 34 mm	0030 077.610	022493000	0030 078 632
0.5 - 20 µL L, 46 mm	0030 077.628	022493002	0030 078 640
2 - 100 µL, 53 mm	0030 077.644	022493006	0030 078 659
2 - 200 µL, 55 mm	0030 077.741	0030077741	0030 078 667
20 - 300 µL, 55 mm	0030 077.636	022493004	0030 078 675
50 - 1000 µL, 76 mm	0030 077.652	022493008	0030 078 683
ep Dualfilter T.I.P.S.® SealMax Racks PCR clean and sterile			
0.1 - 10 µL S, 34 mm	0030 077.806	0030078691	0030 078 691
0.5 - 20 µL L, 46 mm	0030 077.814	0030078705	0030 078 705
2 - 100 µL, 53 mm	0030 077.822	0030078713	0030 078 713
2 - 200 µL, 55 mm	0030 077.830	0030078721	0030 078 721
20 - 300 µL, 55 mm	0030 077.849	0030078730	0030 078 730
50 - 1000 µL, 76 mm	0030 077.857	0030078748	0030 078 748
ep Dualfilter T.I.P.S.® 384 PCR clean and sterile			
0.1 - 20 µL, 42 mm	0030 078.403	0030078403	0030 078 853
0.5 - 100 µL, 52.8 mm	0030 078.411	0030078411	0030 078 861

Eppendorf Certificate

There is no standardized test method available yet to determine the efficiency of filters in pipette tips. Hence, the method established is defined following ISO 29463-5 as below:

Aerosol permeability is measured using a defined polydisperse NaCl aerosol. For distinguished particle diameters, penetrating aerosol particles are counted using a Scanning Mobility Particle Sizer where the ratio of this value to the reference blank value is the permeability. EPA/HEPA* classification according to DIN EN ISO 29463-5 rates the filter efficiency at the particle size of highest penetration.

All tests were performed by an independent accredited laboratory.

* Efficient particulate air filter/High efficiency particulate air filter

Hamburg, September 2020

Page 2 of 2

0015 033.516-02



Joana Tziolis
Expert Liquid Handling
LH-Consumables



Monika Schneider
Head of Global Quality Management &
Regulatory Affairs

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified

Eppendorf®, the Eppendorf Brand Design, ep Dualfilter T.I.P.S.® and LoRetention® are registered trademarks of Eppendorf AG, Germany. All rights reserved, incl. graphics and images. Copyright © 2020 by Eppendorf AG.