



Spacer plate for X-head (TransferMan® NK 2, PatchMan NP 2)





Installing the actuator on the TransferMan NK 2 or PatchMan NP 2

Mount the spacer plate under the pressure plate of the X-head




Mounting example:

Actuator on TransferMan NK 2, on the right side of the microscope




Distanzplatte für X-Kopf (TransferMan® NK 2, PatchMan NP 2)





Aktuator am TransferMan NK 2 oder PatchMan NP 2 installieren

Distanzplatte unter die Andruckplatte des X-Kopfs montieren



Montagebeispiel:

Aktuator an TransferMan NK 2, rechtsseitig am Mikroskop



eppendorf
In touch with life

eppendorf® is a registered trademark
5194.200.058-00/0511

Your local distributor: www.eppendorf.com/worldwide

Eppendorf AG · 22331 Hamburg · Germany · Tel: +49 40 53801-0 · Fax: +49 40 538 01-556 · E-mail: eppendorf@eppendorf.com
Eppendorf North America, Inc. · 102 Motor Parkway · Hauppauge, N.Y. 11788-5178 · USA
Tel: +1 516 334 7500 · Toll free phone: +1 800-645-3050 · Fax: +1 516 334 7506 · E-mail: info@eppendorf.com

Application Support Europe: Tel: +49 1803 666 789 (Preis je nach Tarif im Ausland; 9 ct/min aus dem dt. Festnetz; Mobilfunkhöchstpreis 42 ct/min)

support@eppendorf.com

North America: Tel: +1 800 645 3050 · E-mail: techserv@eppendorf.com

Asia Pacific: Tel: +60 3 8023 6869 · E-mail: support_asiapacific@eppendorf.com