


Easy visual
volume check!



Form Follows Function

Safe and contactless pipetting with epT.I.P.S.[®] 5 mL L tips

Another long distance tip completes the epT.I.P.S. elongated pipette tips. The pipette tip epT.I.P.S. 5 mL L with a length of 175 mm easily reaches the bottom of e.g. 15 mL centrifugation tubes or graduated flasks, avoiding sample loss or cross contamination of pipette and vessel.

The long and slim design of this elongated 5 mL pipette tip gives you free access to your sample with

reduced risk of touching the sides of tubes or the neck of your flask. Even minimal volumes can be recovered safely.

Graduations allow easy visual recognition and control of liquid level when pipetting, dispensing or diluting with this 5 mL L epT.I.P.S. pipette tip.



Product Benefits

- > Extended reach minimizes sample loss
- > No cross-contamination of pipette and sample
- > Easy visual volume check
- > Low ergonomic stress
- > Available in Eppendorf Quality, Eppendorf Biopur® and as ep Dualfilter T.I.P.S.® in PCR clean/Sterile

Contamination free pipetting with epT.I.P.S.®
5 mL L tips is possible with following vessels:

Vessel	Volume
Erlenmeyer flasks	100–500 mL
Schott® glass bottles	25–500 mL
Conical tubes	15–50 mL
Cell culture flasks	25–75 cm ³
Reaction tubes	0.1–5.0 mL
Cryogenic vials	1.0–4.5 mL
Oak ridge flasks	30–85 mL
EDTA tubes	5.0–10.0 mL
Round bottom flasks	5.0–15.0 mL



Description

Color code

Purity

ep Dualfilter T.I.P.S.®, Racks, 5 × 24 = 120 tips

0.2–5 mL L, 175 mm

■ purple

PCR clean / Sterile

epT.I.P.S.®, Bulk, 3 × 100 = 300 tips

0.2–5 mL L, 175 mm

■ purple

Eppendorf Quality

epT.I.P.S.®, Racks, 5 × 24 = 120 tips

0.2–5 mL L, 175 mm

■ purple

Biopur®

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com/consumables