



## Simply Better Pipetting

### With the Eppendorf Xplorer® electronic pipette family

People who give 100 % every day deserve the best laboratory equipment: Eppendorf Xplorer electronic pipettes will introduce you to new levels of simplicity, precision and reproducibility while eliminating many of the potential error risks associated with manual pipetting.

Covering a volume range of 0.5 µL to 10 mL, Xplorer (plus)

is available as single-channel pipette, with 8, 12, 16 or 24 channels, or as one of several Move It® variants with adjustable cone spacing. Look no further if you want to experience precisely adjustable parameters, ergonomic pipetting and accurate results with full control over the motor-operated piston at all times.

### Why upgrade from mechanical pipettes?



#### > Better reproducibility

With digital input of settings and due to their smooth, motor-operated piston movement, Xplorer electronic pipettes guarantee consistent, user-independent results.

#### > Less risk of Repetitive Strain Injuries (RSI)

Settings and the pipetting process are controlled via light-touch buttons and require virtually no operating force. The result: Improved ergonomics and overall better comfort for your thumb and wrist.

#### > Increased efficiency

Featuring various operating modes (e.g. dispensing to fill plates quickly), Xplorer electronic pipettes are versatile precision tools making liquid handling quick and easy – especially in high-throughput applications.

Why Eppendorf Xplorer®? Experience the difference!

- > **Intuitive handling – as easy as with a mechanical pipette:** Innovative multifunction rocker with »up is up« (aspiration) and »down is down« (dispensing) logic
- > **No more getting lost in submenus:** Dial for quick mode selection, well-arranged color display with high viewing angle and guidance in 9 different menu languages
- > **Renowned ergonomics:** Based on the standards of the Eppendorf PhysioCare Concept®, Xplorer pipettes are lightweight, well-balanced and feature spring-loaded tip cones\* to reduce tip attachment forces  
\*in all Xplorer pipettes up to 1 mL (except Move It® 300 µL)
- > **On-target performance – fast or slow:** Repeatable pipetting results are safeguarded via motor-operated piston. 8 speed levels support better results in a number of applications with sensitive samples
- > **Easy adjustment without tools:** Pre-programmed adjustments for altitude, liquid types or elongated tips and individual adjustment in 1, 2 or 3 points additionally increase pipetting accuracy
- > **Options for advanced users and regulated labs:** Up to 10 pipetting modes, programming of workflows, password protection, service interval reminder and more...

Eppendorf Xplorer



- > 1-, 8- & 12-channel pipettes
- > ACT®-labeled 1-channel models help you make eco-conscious purchase decisions
- > Volume range: 0.5 µL to 10 mL (multi-channel up to 1,200 µL)
- > 5 operating modes – ideal for stepping into electronic pipetting
- > Up to 8 hours operating time; work can continue while charging
- > Autoclavable lower parts for easy decontamination

Learn more

Eppendorf Xplorer plus



- > 1-, 8-, 12-, 16- & 24-channel pipettes
- > ACT®-labeled 1-channel models help you make eco-conscious purchase decisions
- > Same feature set as Xplorer with additional modes & benefits:
- > 10 operating modes for advanced users
- > Programmable workflows and fixed volumes
- > Password protection
- > Service interval warning

Learn more

Eppendorf Xplorer plus  
Move It



- > 4-, 6-, 8- & 12-channel pipettes with adjustable cone spacing
- > Same feature set as Xplorer plus (even autoclavable lower parts)
- > Ideal for quick & easy transfer of multiple samples between different vessel formats
- > Unique design without tubing connections increases precision and durability as well as user and sample safety

Learn more

epServices  
for premium performance

Treat your pipettes well!

Regular maintenance, calibration, and adjustment services by Eppendorf assure your pipettes continue to generate reproducible results.

[www.eppendorf.com/epServices](http://www.eppendorf.com/epServices)

Discover your  
service options



Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com/pipettes](http://www.eppendorf.com/pipettes)