


epMotion® 96 and 96xl: iPod® Compatibility List

epMotion 96 and 96xl should be operated via iPod touch®. Many model generations of iPod touch® exist, but only a few can be charged in the epMotion 96 and 96xl docking station. This document shows the connection requirements and iPod touch models. Operation via iPhone® or iPad® is possible without charging the device in the docking station.




Apple® iPod touch needs a lightning connection



Lightning connector and cable

Apple® iPod touch generation and type numbers



Generation 5

- > A1509
- > A1421

Generation 6

- > A1574

Generation 7

- > A2178

iPod® touch type number can be found on back of the device