

Rotor Overview for Centrifuge 59xx and Centrifuge 58xx

Swing-Bucket Rotors	Centrifuge 5920 R	Centrifuge 5910 R	Centrifuge 5810 /R	Centrifuge 5804 /R	AT caps available (Y/N)	Max. speed (rcf)	Max. capacity 50 mL conicals	Max. capacity 5/15 mL conicals	Max. capacity 13 mm blood collection tubes	Max. capacity 16 mm blood collection tubes	Max. capacity 1.5/2.0 mL tubes	Max. capacity MTP	Max. capacity DWP	Bottle usage possible (max size)*	Order no.
S-4x Universal-Large	■	-	-	-	■	4,402	40	96	196	148	368	20	8	(4 x 750 mL)	5895 190 006
S-4x Universal	-	■	-	-	■	4,415	36	68	120	104	-	20	4	(4 x 750 mL)	5895 200 001
S-4x1000 HC	■	-	-	-	-	3,153	52	108	196	144	-	24	8	-	5895 118 003
S-4x1000 RB	■	-	-	-	■	3,413	28	56	92	80	200	28	8	(4 x 1000 mL)	5895 100 007
S-4x1000 PTB	■	-	-	-	■	3,076	40	88	140	132	-	28	8	-	5895 117 007
S-4x750	■	■	-	-	■	4,800	28	56	92	80	200	16	4	(4 x 750 mL)	5895 120 008
S-4-104	-	-	■	-	■	3,214	28	56	92	80	200	16	4	(4 x 750 mL)	5820 740 000
S-4x500	-	■	-	-	■	3,220	20	48	100	64	96	16	4	(4 x 500 mL)	5895 170 005
A-4-81	-	-	■	-	■	3,220	20	48	100	64	96	16	4	(4 x 500 mL)	5810 718 007
S-4x400	-	■	-	-	■	5,212	16	28	60	44	104	-	-	(4 x 400 mL)	5895 180 000
A-4-62	-	-	■	-	■	3,220	16	36	60	48	64	16	4	(4 x 250 mL)	5810 709 008
S-4-72	-	-	■	■	-	3,234	16	32	56	52	104	-	-	(4 x 250 mL)	5804 746 007
A-4-44	-	-	■	■	■	4,500	4/8	16	36	28	48	-	-	-	5804 709 004
A-2-DWP	-	-	■	■	-	2,250	-	-	-	-	-	10	2	-	5804 740 009
A-2-DWP-AT	-	-	■	-	■	3,486	-	-	-	-	-	12	2	-	5820 710 004
Fixed-Angle Rotors															
FA-6x250	■	■	-	-	-	15,054	6	-	-	-	n.a.	-	-	(6 x 250 mL)	5895 175 007
F-34-6-38	-	-	■	■	-	18,500	6	6	18	12	24	-	-	(6 x 85 mL)	5804 727 002
FA-6x50	■	■	-	-	-	20,130	6	6	6	6	6	-	-	-	5895 150 004
FA-45-6-30	-	-	■	■	-	20,130	6	6	6	6	6	-	-	-	5820 715 006
F-48x15	-	■	-	-	-	5,005	-	48	-	-	n.a.	-	-	-	5895 160 000
F-35-48-17	-	-	■	■	-	5,005	-	48	-	-	n.a.	-	-	-	5820 772 000
FA-20x5	■	■	-	-	-	20,913	-	20 (5 mL only)	-	-	6	-	-	-	5895 130 003
FA-45-20-17	-	-	■	■	-	20,913	-	20 (5 mL only)	-	-	6	-	-	-	5820 765 003
FA-48x2	■	■	-	-	-	21,194	-	-	-	-	48	-	-	-	5895 135 005
FA-45-48-11	-	-	■	■	-	19,083	-	-	-	-	48	-	-	-	5820 760 001
FA-30x2	-	■	-	-	-	20,984	-	-	-	-	30	-	-	-	5895 155 006
FA-45-30-11/ F-45-30-11	-	-	■	■	-	20,800	-	-	-	-	30	-	-	-	5804 726 006/ 5804 715 004
T-60-11	-	-	■	■	-	16,400	-	-	-	-	60	-	-	-	5804 730 003
F-45-48-PCR	-	-	■	■	-	15,350	-	-	-	-	n.a.	-	-	-	5804 735 005

■ Yes * Please refer to bottle compatibility chart for detailed information

www.eppendorf.com/centrifugation

Your local distributor: www.eppendorf.com/contact
 Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
 E-mail: eppendorf@eppendorf.com · www.eppendorf.com