

# Combination of epT.I.P.S.<sup>®</sup> Pipette Tips with Eppendorf Research<sup>®</sup> plus Pipettes

Eppendorf Research <sup>®</sup> plus	0.1 – 10 µL dark gray	0.1 – 10 µL M medium gray	0.1 – 20 µL medium gray	0.5 – 20 µL L light gray	0.1 – 20 µL pearl white	5 – 100 µL light yellow	2 – 20 µL yellow	2 – 100 µL yellow	2 – 200 µL yellow	20 – 300 µL orange	50 – 1,000 µL blue	50 – 1,250 µL green	50 – 1,250 µL L dark green	0.25 – 2.5 mL red	0.1 – 5 mL violet	0.2 – 5 mL L violet	0.5 – 10 mL turquoise	0.5 – 10 mL L turquoise	GEloader™/ Microloader™ 0.5-20 µL
epT.I.P.S. <sup>®</sup>	■	■	■	■	□	■	■	■	■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S. <sup>®</sup> LoRetention	■			■					■										
ep Dualfilter T.I.P.S. <sup>®</sup>	■	■		■	■	■	■	■	■	■	■				■	■	■	■	
ep Dualfilter T.I.P.S. <sup>®</sup> LoRetention	■			■						■									
ep Dualfilter T.I.P.S. <sup>®</sup> SealMax	■			■				■	■	■	■								
Fixed volume, single channel pipettes																			
10 µL	■	■	■	■															■
20 µL	■	■	■	■															■
10 µL	■						■	■	■										
20 µL	■						■	■	■	■									
25 µL/50 µL/100 µL	■							■	■	■									
200 µL	■							■	■	■									
200/250/500/1,000 µL	■										■	■	■						
Variable volume, single-, 8- and 12-channel pipettes																			
0.1 – 2.5 µL	■	■	■	■															
0.5 – 10 µL	■	■	■	■															■
2 – 20 µL	■	◆	◆	■															■
2 – 20 µL	■						■	■	■	◆									
10 – 100 µL	■							◆	■	■									
20 – 200 µL	■							◆	◆	■									
30 – 300 µL	■							◆	◆	◆									
100 – 1,000 µL	■										■	■	■						
120 – 1,200 µL	■											■	■						
0.25 – 2.5 mL	■													■					
0.5 – 5 mL	■														■	■			
1 – 10 mL	■																■	■	
Variable volume, 16- and 24-channel pipettes																			
1 – 20 µL	□				■														
5 – 100 µL	■					■													

■ normal volume ◆ limited volume



For more information go to [www.eppendorf.com/eptips](http://www.eppendorf.com/eptips)

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)

Eppendorf®, the Eppendorf Brand Design, epT.I.P.S.<sup>®</sup>, ep Dualfilter T.I.P.S.<sup>®</sup>, Eppendorf Research<sup>®</sup> and Move It<sup>®</sup> are registered trademarks of Eppendorf SE, Germany. GEloader™ and Microloader™ are trademarks of Eppendorf SE, Germany. All rights reserved, including graphics and images. Copyright © 2023 by Eppendorf SE.

# Combination of epT.I.P.S.<sup>®</sup> Pipette Tips with Eppendorf Reference<sup>®</sup> 2 Pipettes

Eppendorf Reference <sup>®</sup> 2	0.1 – 10 µL S dark gray	0.1 – 10 µL M medium gray	0.1 – 20 µL medium gray	0.5 – 20 µL L light gray	2 – 20 µL yellow	2 – 100 µL yellow	2 – 200 µL yellow	20 – 300 µL orange	50 – 1,000 µL blue	50 – 1,250 µL green	50 – 1,250 µL L dark green	0.25 – 2.5 mL red	0.1 – 5 mL violet	0.2 – 5 mL L violet	0.5 – 10 mL turquoise	0.5 – 10 mL L turquoise	GEloader™/ Microloader™ 0.5-20 µL
epT.I.P.S. <sup>®</sup>	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S. <sup>®</sup> LoRetention	■			■			■		■								
ep Dualfilter T.I.P.S. <sup>®</sup>	■	■		■	■	■	■	■	■		■		■	■	■		
ep Dualfilter T.I.P.S. <sup>®</sup> LoRetention	■			■				■	■								
ep Dualfilter T.I.P.S. <sup>®</sup> SealMax	■			■			■	■	■								
Fixed volume																	
1 µL/2 µL	■	■	■														■
5 µL/10 µL	■	■	■	■													■
20 µL	■		■	■													
10 µL	■				■	■	■										
20 µL	■				■	■	■	■									
25 µL/50 µL/100 µL	■					■	■	■									
200 µL	■						■	■									
200/250/500/1,000 µL	■								■	■	■						
2 mL/2.5 mL	■											■					
Adjustable volume																	
0.1 – 2.5 µL	■	■	■	■													
0.5 – 10 µL	■	■	■	■													■
2 – 20 µL	■	◆	◆	■													■
2 – 20 µL	■				■	■	■	◆									
10 – 100 µL	■				◆	■	■	■									
20 – 200 µL	■				◆	◆	■	■									
30 – 300 µL	■				◆	◆	◆	■									
100 – 1,000 µL	■								■	■	■						
0.25 – 2.5 mL	■											■					
0.5 – 5 mL	■												■	■			
1 – 10 mL	■														■	■	

■ normal volume ◆ limited volume



For more information go to [www.eppendorf.com/eptips](http://www.eppendorf.com/eptips)

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany  
 eppendorf@eppendorf.com · www.eppendorf.com

[www.eppendorf.com](http://www.eppendorf.com)

Eppendorf®, the Eppendorf Brand Design, epT.I.P.S.<sup>®</sup>, ep Dualfilter T.I.P.S.<sup>®</sup> and Eppendorf Reference<sup>®</sup> are registered trademarks of Eppendorf SE, Germany. GEloader™ and Microloader™ are trademarks of Eppendorf SE, Germany. All rights reserved, including graphics and images. Copyright © 2023 by Eppendorf SE.

# Combination of epT.I.P.S.<sup>®</sup> Pipette Tips with Eppendorf Xplorer<sup>®</sup>/Xplorer<sup>®</sup> plus Pipettes

Eppendorf Xplorer <sup>®</sup> / Eppendorf Xplorer <sup>®</sup> plus	0.1 – 10 µL dark gray	0.1 – 10 µL M medium gray	0.1 – 20 µL medium gray	0.5 – 20 µL L light gray	0.1 – 20 µL pearl white	5 – 100 µL light yellow	2 – 20 µL yellow	2 – 100 µL yellow	2 – 200 µL yellow	20 – 300 µL orange	50 – 1,000 µL blue	50 – 1,250 µL green	50 – 1,250 µL L dark green	0.25 – 2.5 mL red	0.1 – 5 mL violet	0.2 – 5 mL L violet	0.5 – 10 mL turquoise	0.5 – 10 mL L turquoise	GEloader™/ Microloader™ 0.5 – 20 µL	
epT.I.P.S. <sup>®</sup>	■	■	■	■	□	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S. <sup>®</sup> LoRetention	■			■																
ep Dualfilter T.I.P.S. <sup>®</sup>	■	■		■	■	■	■	■	■	■	■		■		■	■	■	■		
ep Dualfilter T.I.P.S. <sup>®</sup> LoRetention	■			■																
ep Dualfilter T.I.P.S. <sup>®</sup> SealMax	■			■				■	■	■	■									
Variable volume, single-, 8- and 12-channel pipettes																				
0.5 – 10 µL	■	■	■	■																■
1 – 20 µL	■	◆	◆	■	■															
5 – 100 µL	■							◆	■	■										
10 – 200 µL	■							◆	◆	■										
15 – 300 µL	■							◆	◆	◆										
50 – 1,000 µL	■										■		■							
50 – 1,200 µL	■											■	■							
0.125 – 2.5 mL	■													■						
0.2 – 5 mL	■														■	■				
0.5 – 10 mL	■																■	■		
Variable volume, 16- and 24-channel pipettes																				
1 – 20 µL	□				■															
5 – 100 µL	■					■														

■ normal volume ◆ limited volume



For more information go to  
[www.eppendorf.com/eptips](http://www.eppendorf.com/eptips)

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)