

Combination of epT.I.P.S.[®] Pipette Tips with Eppendorf Research[®] plus Pipettes

Eppendorf Research [®] plus Eppendorf Research [®] plus Move It [®]	0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	0.1 – 20 µL, light pink	5 – 100 µL, light yellow	2 – 20 µL, yellow	2 – 100 µL, yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.1 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise
epT.I.P.S. [®] BioBased Standard	■		■	■					■	■	■	■	■	■	■	■		
epT.I.P.S. [®] Standard	■		■	■	■				■	■	■	■	■	■	■	■	■	■
epT.I.P.S. [®] BioBased Reloads			■	■					■	■	■	■	■					
ep Dualfilter T.I.P.S. [®]	■	■		■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
ep Dualfilter T.I.P.S. [®] BioBased Reloads		■		■					■	■	■	■	■					
ep Dualfilter T.I.P.S. [®] SealMax [®]	■			■					■	■	■	■	■					
ep Dualfilter T.I.P.S. [®] SealMax [®] BioBased Reloads				■					■	■	■	■	■					
Fixed volume, single-channel pipettes																		
10 µL	●	■	■	●														
20 µL			●	■														
10 µL	■						■	■	■	●								
20 µL	■						■	■	■	●								
25 µL/50 µL/100 µL	■							■	■	●								
200 µL	■								■	●								
200 µL/250 µL/500 µL/1,000 µL	■										■	●	●					
Variable volume, 1-, 8- and 12-channel pipettes																		
0.1 – 2.5 µL	■	●	●															
0.5 – 10 µL	■	●	■	●														
2 – 20 µL	■	●	●	■														
2 – 20 µL	■						■	■	■	●								
10 – 100 µL	■						●	■	■	●								
20 – 200 µL	■						●	●	■	●								
30 – 300 µL	■						●	●	●	■								
100 – 1,000 µL	■										■	●	●					
120 – 1,200 µL	■											●	■					
0.25 – 2.5 mL	■													■				
0.5 – 5 mL	■														■	●		
1 – 10 mL	■																■	●
Variable volume, 16- and 24-channel pipettes																		
1 – 20 µL	■				■													
5 – 100 µL	■					■												
Move It[®], variable volume, 4-, 6-, 8- and 12-channel pipettes																		
1 – 20 µL	■				■													
5 – 100 µL	■					■												
30 – 300 µL	■						●	●	●	■								
120 – 1,200 µL	■											●	■					

■ Test tip (other tips may require adjustment) ● Full volume range ◆ Limited volume range

Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf[®], the Eppendorf Brand Design, epT.I.P.S.[®], ep Dualfilter T.I.P.S.[®], SealMax[®], Eppendorf Research[®], and Move It[®] are registered trademarks of Eppendorf SE, Germany. All rights reserved, including graphics and images. Copyright © 2025 by Eppendorf SE. Order No: AA01 050 120 – EN1/0425/BE.



For more information go to
www.eppendorf.link/eptips

Combination of epT.I.P.S.[®] Pipette Tips with Eppendorf Reference[®] 2 Pipettes

Eppendorf Reference [®] 2	0.1 – 10 µL dark gray	0.1 – 10 µL M medium gray	0.1 – 20 µL medium gray	0.5 – 20 µL L light gray	2 – 20 µL yellow	2 – 100 µL yellow	2 – 200 µL yellow	20 – 300 µL orange	50 – 1,000 µL blue	50 – 1,250 µL green	50 – 1,250 µL L dark green	0.1 – 2.5 mL red	0.1 – 5 mL violet	0.2 – 5 mL L violet	0.5 – 10 mL turquoise	0.5 – 10 mL L turquoise
epT.I.P.S. [®] BioBased Standard	■		■	■			■	■	■	■	■	■	■	■	■	■
epT.I.P.S. [®]	■		■	■			■	■	■	■	■	■	■	■	■	■
epT.I.P.S. [®] BioBased Reloads			■	■	◆		■	■	■	■	■	■	■	■	■	■
ep Dualfilter T.I.P.S. [®]	■	■		■	■		■	■	■	■	■	■	■	■	■	■
ep Dualfilter T.I.P.S. [®] BioBased Reloads		■		■	■		■	■	■	■	■					
ep Dualfilter T.I.P.S. [®] SealMax [®]	■			■	■		■	■	■	■						
ep Dualfilter T.I.P.S. [®] SealMax [®] BioBased Reloads				■	■		■	■	■	■						
Fixed volume, single-channel pipettes																
1 µL/2 µL	■	●	●													
5 µL/10 µL	■	●	■	●												
20 µL	■		●	■												
10 µL	■				◆	■	■	●								
20 µL	■				◆	■	■	●								
25 µL/50 µL/100 µL	■					■	■	●								
200 µL	■						■	●								
200 µL/250 µL/500 µL/1,000 µL	■								■	●	●					
2 mL/2.5 mL	■											■				
Variable volume, 1-, 8- and 12-channel pipettes																
0.1 – 2.5 µL	■	●	●													
0.5 – 10 µL	■	●	■	●												
2 – 20 µL	■	◆	◆	●	■											
2 – 20 µL	■				◆	◆	■	●								
10 – 100 µL	■				◆	■	■	●								
20 – 200 µL	■				◆	■	■	●								
30 – 300 µL	■				◆	◆	◆	■								
100 – 1,000 µL	■								■	■	■					
0.25 – 2.5 mL	■											■				
0.5 – 5 mL	■												■	■		
1 – 10 mL	■														■	■

■ Test tip (other tips may require adjustment) ● Full volume range ◆ Limited volume range

Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf[®], the Eppendorf Brand Design, epT.I.P.S.[®], ep Dualfilter T.I.P.S.[®], SealMax[®] and Eppendorf Reference[®] are registered trademarks of Eppendorf SE, Germany. All rights reserved, including graphics and images. Copyright © 2025 by Eppendorf SE. Order No: AA01 050 120 – EN1/0425/BE



For more information go to www.eppendorf.link/eptips

Combination of epT.I.P.S.[®] Pipette Tips with Eppendorf Xplorer[®]/Xplorer[®] plus Pipettes

Eppendorf Xplorer [®] Eppendorf Xplorer [®] plus Eppendorf Xplorer [®] plus Move It [®]	0.1 – 10 µL dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL medium gray	0.5 – 20 µL L light gray	0.1 – 20 µL light pink	0.5 – 100 µL light yellow	2 – 20 µL yellow	2 – 100 µL yellow	2 – 200 µL yellow	20 – 300 µL orange	50 – 1,000 µL blue	50 – 1,250 µL green	50 – 1,250 µL L dark green	0.1 – 2.5 mL red	0.1 – 5 mL violet	0.2 – 5 mL L violet	0.5 – 10 mL turquoise	0.5 – 10 mL L turquoise
epT.I.P.S. [®] BioBased Standard	■		■	■					■	■	■	■	■	■	■	■	■	■
epT.I.P.S. [®]	■		■	■	■				■	■	■	■	■	■	■	■	■	■
epT.I.P.S. [®] BioBased Reloads			■						■	■	■	■	■					
ep Dualfilter T.I.P.S. [®]	■	■		■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
ep Dualfilter T.I.P.S. [®] BioBased Reloads		■		■			■	■	■	■	■	■	■					
ep Dualfilter T.I.P.S. [®] SealMax [®]	■			■				■	■	■	■	■	■					
ep Dualfilter T.I.P.S. [®] SealMax [®] BioBased Reloads				■				■	■	■	■	■	■					
Variable volume, 1-, 8- and 12-channel pipettes																		
0.5 – 10 µL	■	●	■	■	●													
1 – 20 µL	■	◆	◆	●	■													
5 – 100 µL	■						◆	■	■	●								
10 – 200 µL	■						◆	◆	■	●								
15 – 300 µL	■						◆	◆	◆	■								
50 – 1,000 µL	■										■	●	●					
50 – 1,200 µL	■										◆	■	●					
0.1 – 2.5 mL	■													■				
0.2 – 5 mL	■														■	●		
0.5 – 10 mL	■																■	●
Variable volume, 16- and 24-channel pipettes																		
1 – 20 µL	■				■													
5 – 100 µL	■					■												
Move It[®], variable volume, 4-, 6-, 8-, and 12-channel pipettes																		
1 – 20 µL	■				■													
5 – 100 µL	■					■												
15 – 300 µL	■						◆	◆	◆	■								
50 – 1,200 µL	■											■	●					

■ Test tip (other tips may require adjustment) ● Full volume range ◆ Limited volume range



Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf[®], the Eppendorf Brand Design, epT.I.P.S.[®], ep Dualfilter T.I.P.S.[®], Eppendorf Xplorer[®] and Move It[®] are registered trademarks of Eppendorf SE, Germany. All rights reserved, including graphics and images. Copyright © 2025 by Eppendorf SE. Order No: AA01 050 120 – EN1/0425

For more information go to
www.eppendorf.link/eptips