

# Centrifuge 5702 / 5702 R / 5702 RH



## General Benefits

- > **Small footprint:** Medium capacity (40 x round bottom) at a very compact footprint
- > **High versatility:** Spin a large variety of different round bottom tubes
- > **Optimal conditions:** Optimum pellet formation due to swing-out rotor for contamination-free removal of the supernatant
- > **Strong sample safety:** Timer function and program buttons (Centrifuge 5702 R) for reproducible results

Centrifuge 5702/R/RH		Rotor A-4-38	Round bucket					Rectangular bucket			16 x 125 mm bucket		
			5702 737 003	5702 719 005	5702 735 000	5702 724 009	5702 723 002	5702 763 004	5702 710 008	5702 711 004	5702 712 000	5702 764 000	
Tube type	BD Vacutainer® 	Microtainer®	16*	16*	–	–	–	36*	40*	–	–	–	
		1.8 – 4.8 mL (13 x 75 mm)	16	16	–	–	–	36	–	–	–	–	
		4.0 – 7.0 mL (13 x 100 mm)	16	16	–	–	–	–	40	32	–	–	
		8.0 – 11.0 mL (16 x 100 mm)	–	–	16	16	–	–	–	–	–	24	–
		CPT™ 4.0 mL (13 x 100 mm)	16	16	–	–	–	–	40	–	–	–	–
		CPT™ 8.0 mL (16 x 125 mm)	–	–	–	–	–	–	–	–	–	–	4
	Greiner® Vacuette 	MiniCollect®	16	16	–	–	–	36*	40*	–	–	–	–
		0.9 – 4.0 mL (13 x 75 mm)	16	16	–	–	–	36	40	–	–	–	–
		3.5 – 6.0 mL (13 x 100 mm)	16	16	–	–	–	–	40	–	–	–	–
		6.0 – 10.0 mL (16 x 100 mm)	–	–	16	16	–	–	–	–	–	24	–
	Sarstedt Monovette® 	1.6 – 4.0 mL (ø 11/13 mm)	16	–	–	–	8	36	–	–	–	–	–
		4.0 – 5.0 mL (ø 11 mm)	16	–	–	–	–	36	–	–	–	–	–
		4.0 – 10.0 mL (ø 13 mm)	–	16	–	–	–	–	20	32	–	–	–
		4.0 – 10.0 mL (ø 15 mm)	–	–	16	16	–	–	–	–	–	24	–
Streck Blood Collection Tube 	5.0 mL (13 x 75 mm)	16	16	–	–	–	36	40	–	–	–	–	
	10.0 mL (16 x 100 mm)	–	–	16	16	8	–	–	–	–	24	–	

\* Manufacturer adapter needed

## Centrifuge 5702 / 5702 R / 5702 RH

Centrifuge 5702/R/RH		Rotor A-4-38	Round bucket					Rectangular bucket				16 x 125 mm bucket
			5702 737 003	5702 719 005	5702 735 000	5702 724 009	5702 723 002	5702 763 004	5702 710 008	5702 711 004	5702 712 000	5702 764 000
Vessel type	BD Vacutainer® 	4.0 – 6.0 mL (ø 13 mm)	20	16	–	–	–	36	–	–	–	–
		8.0 – 11.0 mL (ø 16 mm)	–	–	16	16	–	–	–	–	24	–
		8.0 – 9.5 mL (conical, ø 16 mm)	–	–	–	–	8	–	–	–	24	–
	Greiner Vacuette® 	4.0 mL (ø 13 mm)	16	16	–	–	–	–	40	–	–	–
		6.0 – 10.0 mL (ø 16 mm)	–	–	16	16	–	–	–	–	24	–
	Sarstedt Monovette® 	3.2 – 6.0 mL (ø 13 mm)	16	16	–	–	–	–	40	–	–	–
		8.5 – 10.0 mL (ø 15 mm)	–	–	16	16	–	–	–	–	32	24
Test Tube 	15.0 mL	–	–	–	–	8	–	–	–	24	–	

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)