

Technical Specifications of MultiCon PC



Technical specifications

Processor	Intel Core® i5-6500TE (with Intel vPro)
Graphics	Intel® HD Graphics 530
Memory (RAM)	4GB 2133 MHz DDR4 SO-DIMM on 1 of 2 slots
Storage	128GB 2.5" SATA SSD 72 TBW
Chipset	Intel Q170 Express
Display	19.5" diagonal, Active Matrix TFT LCD (LED)
USB Ports	> 4 x USB 3.0 on I/O > 2 x USB 2.0 (USB Front Access)
Serial Ports	1 x Serial RS-232 (standard)
Ethernet Port	1 x LAN RJ45 (Gigabit)
Wireless	802.11 a/b/g/n/ac
Operating system	Microsoft® Windows® 10 IoT Enterprise LTSB 64-bit
Power Supply	> AC input voltage: 100-240 VAC > Input frequency: 50-60 Hz > Maximum output power: 150W
Power consumption	75W
Dimensions	492 mm x 357 mm x 252 mm
Weight	7.5 kg

Your local distributor: www.eppendorf.com/contact
 Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com