



Your Easy Choice

Ultracentrifuge Consumables

Consumables for Centrifuge CP-NX Series

Consumable		Fixed-angle rotors														
Type	Name	P21A2	P27A	P32AT	P42AT	P45AT	P50A3	P50AT2	P50AT4	P65A	P70AT	P70AT2	P80AT	P90AT	P100AT2	
Microtubes	1.5 mL Micro Tube									■ A						
	1.5 mL Micro Tube (C)						■									
Thick-walled tubes	0.9PC thick-walled tube									■ A		■ A	■ A	■ A		
	10PP thick-walled tube			■		■ A		■ A		■	■ A	■	■	■		
	10PC thick-walled tube			■		■ A		■ A		■	■ A	■	■	■		
	33PP thick-walled tube					■ A		■			■					
	33PC thick-walled tube					■ A		■			■					
	4.7PC thick-walled tube					■ A		■ A	■	■ A	■ A	■ A	■ A	■ A	■ A	■ C
Bottles	10PP bottle			■		■ A		■ A		■	■ A	■	■	■		
	10PC bottle			■		■ A		■ A		■	■ A	■	■	■		
	10PC bottle (B)			■		■ A		■ A		■	■ A	■	■	■		
	30PP bottle					■ A		■			■					
	30PP bottle (C)					■ A		■			■					
	30PC bottle (B3)					■ A		■			■					
	30PC bottle (C3)					■ A		■			■					
	70PP bottle (C)					■										
	70PC bottle (B)					■										
	70PC bottle (C)					■										
	160PC bottle			■ C												
160PP bottle			■ C													
230PP bottle	■ C															
Tubes	0.23PP tube				■											
	0.23PC tube				■											
	3PP tube					■ A		■ A		■ AC	■ A	■ AC	■ AC	■ AC	■ AC	
	3PC tube					■ A		■ A		■ A	■ A	■ A	■ A	■ A	■ A	
	4PP tube						■ AC		■ AC	■ AC	■ AC	■ AC	■ AC	■ AC	■ AC	
	4PET tube															
	5PP tube					■ AC		■ AC		■ AC	■ AC	■ AC	■ AC	■ AC	■ AC	
	5PET tube															
	6.5PP tube					■ AC		■ AC	■ C	■ AC	■ AC	■ AC	■ AC	■ AC	■ AC	■ C
	6.5PC tube					■ AC		■ AC	■ C	■ AC	■ AC	■ AC	■ AC	■ AC	■ AC	■ C
	6.5PET tube					■ AC		■ AC	■ C	■ AC	■ AC	■ AC	■ AC	■ AC	■ AC	■ C
	12PP tube			■ C		■ AC		■ AC		■ C	■ AC	■ C	■ C	■ C	■ C	
	12PC tube			■ C		■ AC		■ AC		■ C	■ AC	■ C	■ C	■ C	■ C	
	12PE tube			■ C		■ AC		■ AC		■ C	■ AC	■ C	■ C	■ C	■ C	
	12PET tube					■ AC		■ AC		■ C	■ AC					
	12SST tube					■ AC		■ AC		■ C	■ AC	■ C	■ C	■ C	■ C	
	12TI tube											■ C	■ C	■ C		
	13PP tube															
	13PET tube															
	16PP tube															
	16PET tube															
	40TI tube							■ C			■ C					
	40PP tube					■ AC		■ C			■ C					
	40PC tube															
	40PE tube					■ AC		■ C			■ C					
	40PET tube					■ AC		■ C			■ C					
	40SST tube (B)							■ C			■ C					
94PP tube					■ C											
94PC tube					■ C											
94SST tube					■ C											
Seal tubes	4PP Seal tube							■ AC							■ AC	
	5PP Seal tube							■ AC							■ AC	
	6PP Seal tube											■ AC	■ AC			
	6.5PP Seal tube							■ C							■ C	
	8PP Seal tube										■ AC	■ AC	■ AC			
	10PP Seal tube										■ C	■ C	■ C			
	12PP Seal tube								■ C		■ C	■ C	■ C			
	15PP Seal tube															
	20PP Seal tube										■ AC					
	30PP Seal tube							■ AC			■ AC					
	33PP Seal tube															
40PP Seal tube							■ C			■ C						
94PP Seal tube					■ C											

A: Requires adapter C: Requires cap. Technical specifications subject to change. Errors and omissions excepted.

Consumables for Centrifuge CP-NX Series

Consumable		Vertical rotors				Swing-bucket rotor						Neo-angle rotors			
Type	Name	P50VT2	P65VT2	P65VT3	P100VT	P28S	P32ST	P32ST2	P40ST	P55ST2	P56ST	P65ST	P65NT	P65NT2	P90NT
Microtubes	1.5 mL Micro Tube														
	1.5 mL Micro Tube (C)														
Thick-walled tubes	0.9PC thick-walled tube									■ A	■ A				
	10PP thick-walled tube														
	10PC thick-walled tube														
	33PP thick-walled tube														
	33PC thick-walled tube														
	4.7PC thick-walled tube										■ A				
Bottles	10PP bottle														
	10PC bottle														
	10PC bottle (B)														
	30PP bottle														
	30PP bottle (C)														
	30PC bottle (B3)														
	30PC bottle (C3)														
	70PP bottle (C)														
	70PC bottle (B)														
	70PC bottle (C)														
	160PC bottle														
160PP bottle															
230PP bottle															
Tubes	0.23PP tube														
	0.23PC tube														
	3PP tube										■ A	■ A			
	3PC tube										■ A	■ A	■		■
	4PP tube										■ A	■ A		■	■
	4PET tube														
	5PP tube										■ A		■		■
	5PET tube											■		■	
	6.5PP tube														
	6.5PC tube														
	6.5PET tube														
	12PP tube														
	12PC tube														
	12PE tube														
	12PET tube														
	12SST tube														
	12TI tube														
	13PP tube														
	13PET tube														
	16PP tube														
	16PET tube														
	40TI tube														
	40PP tube										■	■			
	40PC tube										■	■			
	40PE tube														
	40PET tube										■	■			
	40SST tube (B)										■	■			
94PP tube															
94PC tube															
94SST tube															
Seal tubes	4PP Seal tube														
	5PP Seal tube														
	6PP Seal tube										■ C			■ C	■ C
	6.5PP Seal tube														
	8PP Seal tube														
	10PP Seal tube														
	12PP Seal tube														
	15PP Seal tube														
	20PP Seal tube														
	30PP Seal tube														
	33PP Seal tube														
40PP Seal tube															
94PP Seal tube															

A: Requires adapter C: Requires cap. Technical specifications subject to change. Errors and omissions excepted.

Match Your Needs With Different Types of Rotors and Consumables

Fixed angle rotors:
Ideal for pelleting applications with 25-45° angles

Thick walled tubes:
Allow extra filling volume

Tubes:
From different materials for different applications

Swing-bucket rotors:
Perfect for density gradient applications

Special 1.5 mL micro tubes:
Convenient to handle and spin at up to 200,000 × *g* with Micro Tube (C)

Neo-angle rotors:
Short sedimentation distances for fast separations

Seal tubes:
Heat welded shut to enclose precious samples

Vertical rotors:
Rapid separation of e.g. nucleic acids

Bottles:
Easy to close and open multiple times

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.link/centrifuge-CP-NX

