



Your Easy Choice

Micro Ultracentrifuge Consumables

Consumables for Centrifuge CS-(F)NX Series

| Consumable | | Fixed-angle rotors | | | | | | | | | | | | | |
|--------------------|-------------------------|--------------------|-------|------|-------|--------|--------|---------|---------|--------|---------|---------|--------|--------|--|
| Type | Name | S50A | S55A2 | S58A | S70AT | S80AT2 | S80AT3 | S100AT3 | S100AT4 | S110AT | S120AT2 | S120AT3 | S140AT | S150AT | |
| Microtubes | Sampling Tube | | | | | | ■ | | ■ | ■ | | | | | |
| | 1.5 mL Micro Tube | | | | | | ■ | | ■ | ■ | | | | | |
| | 1.5 mL Micro Tube (C) | | ■ | | | | ■ | | ■ | ■ | | | | | |
| Thick-walled tubes | 0.9PC thick-walled tube | | | ■ A | | | | | | | | | | | |
| | 3PP thick-walled tube | | | | | | | | ■ | | | | | | |
| | 4PP thick-walled tube | | | | | | | | | ■ | | | | | |
| | 4.7PC thick-walled tube | | | ■ A | | | | | | | | | | | |
| | 6PP thick-walled tube | | | | | | ■ | | | | | | | | |
| | 6PC thick-walled tube | | | | | | ■ | | | | | | | | |
| | 10PP thick-walled tube | | | ■ | | | | | | | | | | | |
| | 10PC thick-walled tube | | | ■ | | | | | | | | | | | |
| | 25PP thick-walled tube | ■ | | | | | | | | | | | | | |
| | 25PC thick-walled tube | ■ | | | | | | | | | | | | | |
| | Bottles | 10PP bottle | | | ■ | | | | | | | | | | |
| | | 10PC bottle | | | ■ | | | | | | | | | | |
| 10PC bottle (B) | | | | ■ | | | | | | | | | | | |
| 20PP bottle (C) | | ■ | | | | | | | | | | | | | |
| 20PC bottle (C) | | ■ | | | | | | | | | | | | | |
| Tubes | 0.23PP tube | | | | | | | ■ | | | | | | | |
| | 0.23PC tube | | | | | | | ■ | | | | | | | |
| | 0.5PP tube | | | | | ■ | | | | | | ■ | | | |
| | 0.5PC tube | | | | | ■ | | | | | | ■ | | | |
| | 0.5PC tube (B) | | | | ■ | | | | | | | | | | |
| | 1PP tube | | | | | | | | | | ■ | | ■ | ■ | |
| | 1PC tube | | | | | | | | | | ■ | | ■ | ■ | |
| | 2.2PP tube | | | | | | | | | | | | | | |
| | 2.2PET tube | | | | | | | | | | | | | | |
| | 3PP tube | | | ■ AC | | | | | | | | | | | |
| | 3PC tube | | | ■ A | | | | | ■ | | | | | | |
| | 4PP tube | | | ■ AC | | | | | | | | | | | |
| | 4PC tube | | | | | | | | | ■ | | | | | |
| | 5PP tube | | | ■ AC | | | | | | | | | | | |
| | 5PET tube | | | | | | | | | | | | | | |
| | 5TI tube | | | | | | | | | | ■ | | | | |
| | 6.5PP tube | | | ■ AC | | | | | | | | | | | |
| | 6.5PC tube | | | ■ AC | | | | | | | | | | | |
| | 7PP tube (B) | | | | | | | | | | | | | | |
| | 7PC tube (B) | | | | | | | | | | | | | | |
| | 7PET tube | | | | | | | | | | | | | | |
| | 12PP tube | | | ■ | | | | | | | | | | | |
| | 12PC tube | | | ■ | | | | | | | | | | | |
| | 12PE tube | | | ■ | | | | | | | | | | | |
| | 12SST tube | | | ■ | | | | | | | | | | | |
| | 12TI tube | | | ■ | | | | | | | | | | | |
| | 30PP tube | ■ | | | | | | | | | | | | | |
| Seal tubes | 1.5PP Seal tube | | | | | | | | | | ■ AC | | ■ AC | ■ AC | |
| | 2 PP Seal tube | | | | | | | | | | ■ AC | | ■ AC | ■ AC | |
| | 2PP Seal tube | | | | | | | | ■ AC | | | | ■ AC | ■ AC | |
| | 4PP Seal tube | | | | | | | | | ■ AC | | | | | |
| | 5PP Seal tube | | | | | | | | | ■ AC | | | | | |
| | 6PP Seal tube | | | ■ AC | | | ■ AC | | | | | | | | |
| | 8PP Seal tube | | | ■ AC | | | ■ | | | | | | | | |
| | 12PP Seal tube | | | ■ | | | | | | | | | | | |
| | 20PP Seal tube | ■ AC | | | | | | | | | | | | | |
| | 30PP Seal tube | ■ | | | | | | | | | | | | | |

A: Requires adapter C: Requires cap. Technical specifications subject to change. Errors and omissions excepted.

Consumables for Centrifuge CS-(F)NX Series

| Consumable | | Swing-bucket rotors | | | Vertical rotors | Neo-angle rotors | | | |
|---------------------------|-------------------------|---------------------|-------|------|-----------------|------------------|--------|--------|--|
| Type | Name | S50ST | S52ST | S55S | S120VT | S100NT | S120NT | S140NT | |
| Microtubes | Sampling Tube | | | | | | | | |
| | 1.5 mL Micro Tube | | | | | | | | |
| | 1.5 mL Micro Tube (C) | | | | | | | | |
| Thick-walled tubes | 0.9PC thick-walled tube | | ■ | | | | | | |
| | 3PP thick-walled tube | | ■ | | | | | | |
| | 4PP thick-walled tube | | | | | | | | |
| | 4.7PC thick-walled tube | | | | | | | | |
| | 6PP thick-walled tube | | | | | | | | |
| | 6PC thick-walled tube | | | | | | | | |
| | 10PP thick-walled tube | | | | | | | | |
| | 10PC thick-walled tube | | | | | | | | |
| | 25PP thick-walled tube | | | | | | | | |
| | 25PC thick-walled tube | | | | | | | | |
| | Bottles | 10PP bottle | | | | | | | |
| | | 10PC bottle | | | | | | | |
| 10PC bottle (B) | | | | | | | | | |
| 20PP bottle (C) | | | | | | | | | |
| 20PC bottle (C) | | | | | | | | | |
| | | | | | | | | | |
| Tubes | 0.23PP tube | | | | | | | | |
| | 0.23PC tube | | | | | | | | |
| | 0.5PP tube | | | | | | | | |
| | 0.5PC tube | | | | | | | | |
| | 0.5PC tube (B) | | | | | | | | |
| | 1PP tube | | | | | | | | |
| | 1PC tube | | | ■ | | | | | |
| | 2.2PP tube | | | ■ | | | | | |
| | 2.2PET tube | | | ■ | | | | | |
| | 3PP tube | | | | | | | | |
| | 3PC tube | | ■ | | | | | | |
| | 4PP tube | | | | | | | | |
| | 4PC tube | | | | | | | | |
| | 5PP tube | | ■ | | | | | | |
| | 5PET tube | | ■ A | | | | | | |
| | 5TI tube | | | | | | | | |
| | 6.5PP tube | | | | | | | | |
| | 6.5PC tube | | | | | | | | |
| | 7PP tube (B) | ■ | | | | | | | |
| | 7PC tube (B) | ■ | | | | | | | |
| | 7PET tube | ■ | | | | | | | |
| | 12PP tube | | | | | | | | |
| | 12PC tube | | | | | | | | |
| | 12PE tube | | | | | | | | |
| | 12SST tube | | | | | | | | |
| | 12TI tube | | | | | | | | |
| | 30PP tube | | | | | | | | |
| | Seal tubes | 1.5PP Seal tube | | | | | | | |
| 2 PP Seal tube | | | | | | | | | |
| 2PP Seal tube | | | | | ■ | | | ■ | |
| 4PP Seal tube | | | | | | | ■ | | |
| 5PP Seal tube | | | | | | ■ | | | |
| 6PP Seal tube | | ■ | | | | | | | |
| 8PP Seal tube | | | | | | | | | |
| 12PP Seal tube | | | | | | | | | |
| 20PP Seal tube | | | | | | | | | |
| 30PP Seal tube | | | | | | | | | |

Match Your Needs With Different Types of Rotors and Consumables

Fixed angle rotors:
Ideal for pelleting applications with 25-45° angles

Thick walled tubes:
Allow extra filling volume

Tubes:
From different materials for different applications

Swing-bucket rotors:
Perfect for density gradient applications

Special 1.5 mL micro tubes:
Convenient to handle and spin at up to 200,000 × *g* with Micro Tube (C)

Neo-angle rotors:
Short sedimentation distances for fast separations

Seal tubes:
Heat welded shut to enclose precious samples

Vertical rotors:
Rapid separation of e.g. nucleic acids

Bottles:
Easy to close and open multiple times

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.link/centrifuge_CS-F-NX

