




# VisionNize<sup>®</sup> Lab Suite

Overview of compatible Eppendorf devices

# How Devices Can Be Connected



- 1 VisioNize box (hardware hub)
- 2 VisioNize sense gateway router
- 3 VisioNize sense environmental monitoring sensors
- 4 VisioNize sense external temperature sensor



# VisioNize Device Overview

## Overview of compatible Eppendorf devices for VisioNize Lab Suite

This list contains Eppendorf devices which can be connected to the platform. Visit [www.eppendorf.com/visionize](http://www.eppendorf.com/visionize) for the most current compatibility list.

### VisioNize touch enabled devices

> Eppendorf devices equipped with the VisioNize touch interface can be connected via ethernet cable:

#### Eppendorf ULT Freezers

CryoCube® F740i  
CryoCube® F740hi  
CryoCube® F740iw  
CryoCube® F740hiw

#### Eppendorf CO<sub>2</sub> Incubators

CellXpert® C170i

#### Eppendorf Shakers

Innova® S44i\*

\* only compatible after serial number S44IJH302546

#### Eppendorf Cyclers

Mastercycler® X50a  
Mastercycler® X50s  
Mastercycler® X50p  
Mastercycler® X50h  
Mastercycler® X50i  
Mastercycler® X50l  
Mastercycler® X50r  
Mastercycler® X50t

#### Bioprocess Systems

SciVario® twin

#### Eppendorf Centrifuge

Centrifuge 5910 Ri

### VisioNize box enabled

> Eppendorf devices which are connectable via VisioNize box to the VisioNize Lab Suite. The VisioNize box acts as a hardware hub and works as a translation system enabling data collection of the device data in VisioNize Lab Suite.

#### Eppendorf Centrifuges

Centrifuge 5910 R  
Centrifuge 5424 R  
Centrifuge 5425  
Centrifuge 5425 R  
Centrifuge 5427 R

#### Eppendorf Shakers

Innova® 40  
Innova® 42  
Innova® 43  
Innova® 44  
New Brunswick™ S41i

#### Eppendorf ULT Freezers

CryoCube® F440n  
CryoCube® F440h  
CryoCube® F570  
CryoCube® F570h  
CryoCube® F570hw  
CryoCube® FC660  
CryoCube® FC660h  
CryoCube® F740  
CryoCube® F740h  
CryoCube® F101h  
Innova® U101  
Innova® U535  
Innova® U570  
Innova® U725  
Innova® U725-G  
Innova® C760  
Innova® C585  
Premium U410  
Premium U570  
HEF U410

#### Eppendorf CO<sub>2</sub> Incubators

Galaxy® 48R  
CellXpert® C170

#### Eppendorf Thermomixer

Eppendorf ThermoMixer® C, F0.5, F1.5, F2.0, FP  
Eppendorf ThermoStat™ C

#### Eppendorf Cyclers

Mastercycler® Nexus  
Mastercycler® Nexus gradient  
Mastercycler® Nexus flat  
Mastercycler® Nexus X2  
Mastercycler® Nexus GX2  
Mastercycler® Nexus SX1  
Mastercycler® Nexus GSX1

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)

Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany  
eppendorf@eppendorf.com · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)