eppendorf




Get it Right!

Find the most suitable Pipetting System with the Eppendorf Liquid Handling Selection Chart

eppendorf


Liquid Handling Selection Chart 1.1: Small Volumes (µL up to 10 mL)



¹Liquids with viscosity, volatility, surface tension or density other than water as well as hot, cold or hazardous liquids. ²Complexity of task, e.g. pipetting of many different liquids or volumes.



Liquid Handling Selection Chart 1.2: Large Volumes (5 - 100 mL)



¹Complexity of task, e.g. pipetting of many different liquids or volumes.

²All vessels with single channel format like 50 mL conicals or cell culture vessels.

³Liquids with viscosity, volatility, surface tension or density other than water as well as hot, cold or hazardous liquids.

Liquid Handling Selection Chart 2: Associated Tips



Easypet®, Multipette®, Pipet Helper®, Varispenser®, Varitips®, Combitips®, epT.I.P.S.®, ep Dualfilt Liquid type Liquid handling device **Contamination risk** Recommended tip Available purity grade Varitips® S + Maxitip **Eppendorf Quality Aqueous** Varipette® 4720 No Volatile Varitips® S + Maxitip with valve **Eppendorf Quality** Other Challenging Varitips® P **Eppendorf Quality** Eppendorf®, the Eppendorf Brand Design, epMotion®, Move It®, Varipette®, Easypet®, Multipette®, Pipet Helper®, Varispenser®, Varitips®, Combiti ep Dualfilter T.I.P.S. SealMax® are registered trademarks of Eppendorf SE, Germany.
U.S. Design Patents are listed on www.eppendorf.com/ip. All rights reserved, including graphics and images. Copyright © 2023 by Eppendorf SE. Aqueous / Varitips® P Yes **Eppendorf Quality** Challenging Easypet® 3, Aqueous / Pipet Helper® **Eppendorf Serological Pipets** Sterile No / Yes (include a filter) Challenging¹ **Eppendorf Quality, PCR clean,** Combitips® advanced Multipette® **Avoids con-**Biopur®, Forensic DNA Grade tamination ViscoTip® **Eppendorf Quality** principle **Eppendorf Quality, PCR clean,** Combitips® advanced Biopur®, Forensic DNA Grade epT.I.P.S.®, **Eppendorf Quality,** Research plus, PCR clean, Biopur® epT.I.P.S.® Motion Reference 2, Xplorer, epT.I.P.S.®, **Eppendorf Quality, PCR clean** Challenging¹ epT.I.P.S.® Motion **Xplorer plus** Pipette Manager, epT.I.P.S.® Motion with filter, PCR clean, Move It® ep Dualfilter T.I.P.S.®, PCR clean + sterile. ep Dualfilter T.I.P.S.® SealMax® **Forensic DNA Grade** epMotion® 96, epMotion[®] epT.I.P.S.® Motion with filter, PCR clean, ep Dualfilter T.I.P.S.®, PCR clean + sterile ep Dualfilter T.I.P.S.® SealMax

¹Liquids with viscosity, volatility, surface tension or density other than water as well as hot, cold or hazardous liquids.