



S41i Incubator Shaker

S41i Incubator Shaker Performance Plans

The versatility of the New Brunswick S41i Incubator Shaker supports simultaneous culture of adherent and suspension cultures. The smallest deviations in the culturing environment can result in undesirable outcomes.

To ensure reproducible results and to protect your precious cultures from undesirable variables, we offer a range of maintenance services.

For your Quality Management requirements we also offer certification services that include Installation Qualification (IQ) and Operational Qualification (OQ) certification to provide optimal instrument performance and stability and accuracy of culturing parameters.

With our Performance Plans options we ensure that your New Brunswick S41i Incubator Shaker instruments are in good working condition, according to manufacturer specifications.

Optimized service—optimized yield—reliable results!

IQ/OQ Services

Installation Qualification & Operational Qualification services assure that instrumentation is installed and running according to manufacturer specifications. IQ/OQ documentation is supplied to help prepare for quality and regulatory audits.

■ ADVANCED MAINTENANCE

Preventive maintenance service to meet manufacturer specifications

■ ESSENTIAL CHECK

Check of all fundamental functions of the product

eppendorf

Performance tested on

Model: _____

Serial no.: _____

Serviced by: _____

Service no.: _____

Next service: _____

Date: _____

epServices
for premium performance

S41i Incubator Shaker Performance Plans

Performance Plans	ESSENTIAL CHECK	ADVANCED MAINTENANCE	PREMIUM SERVICE	Installation Qualification (IQ)	Operational Qualification (OQ)
Catalog No.	0082050003	0082050004	0082050005	0082050007	0082050008
External Checks and Maintenance					
Check installation environment	■	■	■	■	■
Check electrical/gas supply and fittings		■	■	■	■
Replace in-line CO ₂ filter and auto-zero filter		■	■		
Internal Equipment and Fittings					
Check security of shaker platform(s) and flask holders		■	■	■	■
Check security of incubator shelf/shelves and racking system		■	■	■	■
Replace platform support post boot seals		■	■		
Check condition of humidity trays		■	■		
Check condition of shaker drive belt		■	■		
Check condition of chamber CO ₂ and temperature sensors	■	■	■		■
Validation and Metrology					
Validation and adjustment of operating parameters to within Eppendorf specifications	■	■	■		■ ¹
Documentation					
Issue preventive maintenance Check List (signed document)	■	■	■		
Dated service sticker fixed to S41i Incubator Shaker to confirm Eppendorf service	■	■	■	■	■
Provide IQ Report and Certificate				■	
Provide OQ Report and Certificate					■
Supporting Information					
Travel	not included	not included	included	included	not included
Time/Labor	included	included	included	included	included
Additional parts	not included	not included	included	not included	not included

Additional catalog number for combined Installation and Operational Qualification (IQ/OQ) services is listed below.

Travel and labor are included. Additional parts are not included.
0082050009

¹ The detailed certification service contents are listed in the OQ report.

Field Service Only

Advantage OS	Catalog No.
S41i	ADSERV-S41I

Advantage OS

An easy-to-use value packaged service option to help plan in advance and protect the system. Advantage OS provides parts, labor, and travel for 2 visits; an Advanced Maintenance and a single repair visit within the year.

Eppendorf North America, Inc.
 Phone: 800-645-3050
 Email: info@eppendorf.com

www.eppendorf.com/epServices