

Eppendorf Certificate

Certificate of Quality Eppendorf Laboratory Consumables

STATEMENT ON BSE/TSE

As a leading manufacturer of biotech products, Eppendorf uses PP (polypropylene), PC (polycarbonate), and PE (polyethylene) granulates specifically suited for laboratory applications and the manufacturing process of laboratory consumables. PP, PC, and PE granulates may contain small amounts of materials derived from animals.

Eppendorf only works with granulate suppliers who guarantee that their animal components derive exclusively from countries without BSE (bovine spongiform encephalopathy) occurrences. Risk materials are not used. Thus, the requirements of the EU Regulation 1326/2001 and Commission Decision 2001/2/EC amending Decision 2000/418/EC regulating the use of materials presenting risks regarding transmissible spongiform encephalopathies (TSE) are fulfilled.

The granulate production includes – depending on the process – hydrolysis, esterification, or hydrogenation steps in different variations. The common features of these steps include processing conditions with temperatures above 235 °C and pressures above 3,000 kPa with retention times up to several hours. The final product is obtained through fractionation, neutralization, and purification. The subsequent extrusion (for the production of granulate) takes place at minimum 200 °C for several minutes.

Thus, the production chain of raw materials by far exceeds the stringent requirement of 200 °C for 20 minutes (Annex VI, chapter III of EU Regulation 1774/2002, in EU Directives 2000/6/EC and 1999/82/EC, referring to Document EMEA/410/01-Final, latest version: Rev. 3 – 05.03.2011, and in the Report WHO/CDS/VPH/95.145). Any virus, bacterium, or substance causing immunological diseases (TSE, BSE, CJD) is destroyed.

Eppendorf states that the materials and laboratory consumables are to be considered safe with respect to BSE and TSE transmission when used in consumer applications.

This certificate applies to the following Eppendorf Laboratory Consumables:

Pipette tips	epT.I.P.S. [®] epT.I.P.S. [®] 384 epT.I.P.S. [®] Long epT.I.P.S. [®] LoRetention epT.I.P.S. [®] Motion ep Dualfilter T.I.P.S. [®] ep Dualfilter T.I.P.S. [®] 384 ep Dualfilter T.I.P.S. [®] SealMax ep Dualfilter T.I.P.S. [®] LoRetention GELoader [®] Microloader™ Eppendorf Serological Pipets
--------------	--

Eppendorf Certificate

Positive displacement tips	Mastertip® Eppendorf Varitips® Combitips advanced® incl. adapters ViscoTip®
Eppendorf Tubes®	Eppendorf Tubes® 3810X/ Flex-Tube® Eppendorf Safe-Lock Tubes* Eppendorf DNA LoBind/ Protein LoBind Tubes, Eppendorf Tubes® 5 mL incl. adapters Eppendorf Conical Tubes 15 mL, 25 mL, 50 mL*
Eppendorf Plates®	Eppendorf Microplates* Eppendorf Deepwell Plates* Eppendorf DNA LoBind/ Protein LoBind Plates Eppendorf Assay/Reader Microplates
PCR Consumables	Eppendorf twin.tec® PCR Plates* Eppendorf twin.tec® PCR Plates LoBind Eppendorf twin.tec® <i>microbiology</i> PCR Plates Eppendorf twin.tec® <i>real-time</i> PCR Plates PCR Tube Strips Fast PCR Tube Strips <i>real-time</i> PCR Tube Strips PCR Cap Strips PCR Tubes PCR Films & Foils
Cell Culture Consumables	Eppendorf Cell Culture Dishes Eppendorf Cell Culture Plates Eppendorf Cell Culture Flasks Eppendorf Cell Imaging Dishes Eppendorf Cell Imaging Plates CCCadvanced® Cell Imaging Slides & Coverglasses
Cuvettes	UVette® Vis Cuvette
Sample Handling Consumables	Centrifuge Bottles

* applies also to the SafeCode variants

Hamburg, January 2020

Page 2 of 2

0015 033.510-01



Joana Tziolis
Expert Liquid Handling
LH-Consumables



Monika Schneider
Head of Global Quality Management &
Regulatory Affairs

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified

Eppendorf®, the Eppendorf Brand Design, epT.I.P.S.®, ep Dualfilter T.I.P.S.®, GELoader®, Mastertip®, Eppendorf Varitips®, CCCAdvanced®, Combitips advanced®, ViscoTip®, Eppendorf Tubes®, Flex-Tube®, Eppendorf twin.tec® and UVette® are registered trademarks of Eppendorf AG, Germany. Microloader™ is a trademark of Eppendorf AG, Germany. All rights reserved, incl. graphics and images. Copyright ©2020 by Eppendorf AG.