

# Eppendorf Certificate

## Certificate of Quality Combitips advanced®

### Quality Assurance

Combitips advanced are manufactured in a controlled area and are subject to regular quality inspections within statistical process control.

The Eppendorf Quality Management System is certified according to ISO 9001 and ISO 13485.

### Materials

Combitips advanced are made of virgin polypropylene and polyethylene of highest purity and quality. Our material suppliers do not use or intentionally incorporate the following agents into materials Eppendorf uses for the production of Combitips advanced:

- Slip agents (including oleamide, erucamide, stearamide)
- Biocides (including di(2-hydroxyethyl)methyldodecylammonium salts (DiHEMDA))
- Plasticizers (softeners)

The used raw materials for Combitips advanced comply with FDA regulations 21CFR § 177.1520 "Olefin Polymers" and 21CFR § 178.2010 "Antioxidants and Stabilizers for Polymers" for food containers.

During production of Combitips advanced no slip agents, biocides and plasticizers are used.

Combitips advanced and all product components do not contain or come in contact with:

- Metallic dyes
- Latex
- Glucan
- Cellulose
- Melamine

May 2015



Brigitte Klose  
Global Product Manager  
Liquid Handling Consumables/Instruments



Dr. Birgit Schreiber  
Vice President  
Quality Management & Regulatory Affairs

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
Eppendorf AG · 22331 Hamburg · Germany  
E-mail: [eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com)

ISO 9001  
Certified

ISO  
13485  
Certified

ISO  
14001  
Certified