

Eppendorf Certificate

Certificate of Quality

Eppendorf Tubes®

All Eppendorf Tubes are manufactured in an controlled ISO class 8 clean room environment according to ISO 14644-1 and are subject to regular quality inspections within a statistical process control. Tests being performed are e.g. vapor tightness, lid-opening force, wall strength, and centrifugal stability.

The Eppendorf Quality Management System is certified according to the international standards ISO 9001, ISO 13485, and ISO 14001.

Eppendorf Tubes are made of virgin polypropylene of highest purity and quality. Neither the material suppliers of Eppendorf nor Eppendorf use or intentionally incorporate the following agents into the materials Eppendorf is using for the production of Tubes:

- Slip agents (including oleamide, erucamide, stearamide)
- Biocides (including di(2-hydroxyethyl)methyldodecylammonium salts (DiHEMDA))
- Plasticizers

The used raw materials for Eppendorf Tubes comply with FDA regulations 21CFR 177.1520 “Olefin Polymers” and 21CFR 178.2010 “Antioxidants and Stabilizers for Polymers” for food containers.

July, 2019



Joana Tziolis
Expert Liquid Handling
LH-Consumables



Dr. Birgit Schreiber
Vice President Quality Management &
Regulatory Affairs

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified