

Eppendorf Certificate

Certificate of Quality Mastertip®

Mastertip® pipette tips are single-use plastic tips which, in combination with the Biomaster® pipette, form a perfectly functioning pipetting system.

Quality Assurance

Mastertip® pipette tips are manufactured in a controlled area and are subject to regular quality inspections within statistical process control. The finely extended tip guarantees accurate volume dispensing and assures reliable results and outstanding reproducibility.

The Eppendorf Quality Management System is certified according to ISO 9001 and ISO 13485.

Materials

Mastertip® pipette tips are made of virgin polypropylene and polyethylene of highest purity and quality. Our material suppliers do not use or intentionally incorporate the following agents into materials Eppendorf uses for the production of Mastertip® pipette tips:

- Slip agents (including oleamide, erucamide, stearamide)
- Biocides (including di(2-hydroxyethyl)methyldodecylammonium salts (DIHEMDA))
- Plasticizers (softeners)

The used raw materials for Mastertip® pipette tips comply with FDA regulations 21CFR § 177.1520 "Olefin Polymers" and 21CFR § 178.2010 "Antioxidants and Stabilizers for Polymers" for food containers.

During production of Mastertip® pipette tips no slip agents, biocides and plasticizers are used.

Mastertip® pipette tips and all product components do not contain metallic dyes or latex.

May 2015



Brigitte Klose
Global Product Manager
Liquid Handling Consumables/Instruments



Dr. Birgit Schreiber
Vice President
Quality Management & Regulatory Affairs

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · 22331 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO
13485
Certified

ISO
14001
Certified