

Eppendorf Certificate

Certificate of Quality – ep Dualfilter T.I.P.S.®

ep Dualfilter T.I.P.S.® filter tips do not differ in handling nor in specification for any of the properties listed below.

- > Orifice diameter, inner diameter, outer diameter
- > Length and geometry of the pipette tips
- > Weight / Attachment forces / Ejection forces
- > Compatibility with Eppendorf pipettes and non-Eppendorf pipettes
- > Quality of purity and sterility
- > Pipette tip material / Resistivity to chemicals
- > Position of the filter within the tip
- > Packing quantity, packing volume, packing weight

The optimization of filter properties (aerosol binding) refers exclusively to the pore sizes of the filter and the color marking of the upper filter layer with pigments.

ep Dualfilter T.I.P.S.® Racks	Order no. (international)	Order no. (North America)
0.1 - 10 µL S, 34 mm	0030 077.504	22491202
0.1 - 10 µL M, 40 mm	0030 077.512	22491211
0.5 - 20 µL L, 46 mm	0030 077.520	22491229
2 - 20 µL, 53 mm	0030 077.539	22491270
2 - 100 µL, 53 mm	0030 077.547	22491237
2 - 200 µL, 55 mm	0030 077.555	22491296
20 - 300 µL, 55 mm	0030 077.563	22491245
50 - 1000 µL, 76 mm	0030 077.571	22491253
50 - 1250 µL L, 103 mm	0030 077.750	22494002
100 - 5000 µL, 120 mm	0030 077.580	22491261
1 - 10 mL L, 243 mm	0030 077.598	22491288
ep Dualfilter T.I.P.S.® LoRetention Racks	Order no. (international)	Order no. (North America)
0.1 - 10 µL S, 34 mm	0030 077.610	22493000
0.5 - 20 µL L, 46 mm	0030 077.628	22493002
2 - 100 µL, 53 mm	0030 077.644	22493006
20 - 300 µL, 55 mm	0030 077.636	22493004
50 - 1000 µL, 76 mm	0030 077.652	22493008
ep Dualfilter T.I.P.S.® SealMax Racks	Order no. (international)	Order no. (North America)
0.1 - 10 µL S, 34 mm	0030 077.806	30077806
0.5 - 20 µL L, 46 mm	0030 077.814	30077814
2 - 100 µL, 53 mm	0030 077.822	30077822
2 - 200 µL, 55 mm	0030 077.830	30077830
20 - 300 µL, 55 mm	0030 077.849	30077849
50 - 1000 µL, 76 mm	0030 077.857	30077857

Date: Hamburg, September 2014

Brigitte Klose
Global Product Manager
Consumables

Dr. Birgit Schreiber
Vice President
Quality Management & Regulatory Affairs

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · 22331 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO
13485
Certified

ISO
14001
Certified