



Our Knowledge, Our Commitment

Pipette maintenance, calibration and adjustment services

epServices
for premium performance

Premium Service from the experts!

For more than, half a century Eppendorf has been an innovator and leading provider of liquid handling systems. We pursue a total approach philosophy extending beyond products with quality pipette calibration, adjustment, certification, maintenance and repair services.

Regular maintenance, calibration and adjustment by Eppendorf assure your pipettes and dispensers continue to generate reproducible results.

Pipette calibration facility at Eppendorf India in Chennai is ISO 17025 certified and accredited by NABL (National Accreditation Board for Testing and Calibration Laboratories)

Eppendorf offers pipette calibration at 9 different depots/centres in Chennai, Bengaluru, Hyderabad, Mumbai, Pune, Kolkata, New Delhi, Chandigarh and Lucknow. Calibrate or service any brand of pipettes at any of these locations near you.



Pipette Performance Plans

Eppendorf Service Levels

Piston operated pipette

	Standard Full Calibration in compliance with ISO 8655	Premium Full Calibration with certificate according to NABL
Initial inspection to determine pipette condition	✓	✓
Cleaning of pipette surface	✓	✓
Cleaning and lubrication of seals, O-rings and piston	✓	✓
Replacment of worn or defective parts, when necessary*1	✓	✓
Number of measurements per test volume (min, mid, max)	10	10
Adjustment according to ISO 8655	✓	✓
Calibration report provided	✓	✓
Calibration label pasted on pipettes	✓	✓
Calibration certificate according to NABL (ISO 17025) provided		✓

*1 Spare parts charges extra

Price details*2

Fixed-volume pipette	₹ 750	₹ 1000
Adjustable-volume pipette	₹ 1100	₹ 1450
8-channel pipette	₹ 1700	₹ 2350
12-channel pipette	₹ 2100	₹ 3100

*2 Prices mentioned above excluding taxes and shipping cost