



Pipette Calibration Service

Treat your pipettes well with premium services from the experts!

Precision and accuracy of your pipettes are important to ensure the quality of your results. Therefore, regular maintenance, calibration, and adjustment services by Eppendorf assure your pipettes and dispensers continue to generate reproducible results.

Our newly, clearly designed calibration service portfolio includes a thorough check of the pipette internally and externally, ensuring that wearable parts like O-rings are replaced in a timely manner.

Get what you want - Know what you get!

We offer a range of options from economic calibration to calibration in ISO 17025 accredited laboratories in our major markets.

Value of pipette service by Eppendorf

- > 75 years of manufacturer expertise for devices, consumables and services
- > Maintaining every pipette's conformity according to ISO 8655
- > ISO 17025 accredited certification in 19 countries
- > Transparent and modular service offers and no hidden costs

Watch video:
www.eppendorf.com/pipette-video



Pipette Calibration Service

Calibration Services – Ordering information									
Service	Calibration measurements "As found"	Calibration measurements "As Left"	Decontamination from outside ²⁾	Checking, cleaning, lubricating ⁴⁾ before calibration	Single channel pipette/dispenser	8-channel pipette	12-channel pipette	16-channel pipette	24-channel pipette
ISO 8655 including Preventive Maintenance		10 per test volume ¹⁾	included	included	0082 000 142 46,00 EUR	0082 000 165 122,00 EUR	0082 000 166 143,00 EUR	0082 000 185 175,00 EUR	0082 000 186 270,00 EUR

Calibration Service Packages – Ordering information									
Service	Calibration measurements "As found"	Calibration measurements "As Left"	Outside decontamination ²⁾ ; checking, cleaning, lubrication ⁴⁾	ISO 17025 Accreditation calibration certificate	Single channel pipette/dispenser	8-channel pipette	12-channel pipette	16-channel pipette	24-channel pipette
ISO 17025 including Preventive Maintenance		10 per test volume ¹⁾	included	included	0082 000 143 64,00 EUR	0082 000 167 149,00 EUR	0082 000 168 186,00 EUR	0082 000 187 215,00 EUR	0082 000 188 322,00 EUR
ISO 17025 Professional	10 per test volume ¹⁾	10 per test volume ¹⁾	included	included	0082 000 220 75,00 EUR	0082 000 221 160,00 EUR	0082 000 222 205,00 EUR	0082 000 223 240,00 EUR	0082 000 224 350,00 EUR

1) 10%, 50% and 100% of nominal volumes for variable volume instruments.
 2) Pipettes and dispensers and package equipment need to be decontaminated before sending in. Download decontamination form at www.eppendorf.com/epServices.
 3) ISO 17025 accreditation calibration certificate is included in the second calibration after the maintenance.
 4) Internal parts

We ensure our stringent and validated procedures are followed, for Eppendorf and other brand pipettes:



Qualification and Certification

ISO 8655

ISO 8655 lays out the standard criteria by which a pipette should be calibrated and what error limits should be used at each volume. ISO 8655 criteria includes measurements at 3 volumes, 10 measurements at each volume, and controls for room environment.

ISO 17025

The ISO 17025 accreditation certifies the high level of competence and compliance with strict quality rules. ISO 17025 criteria include a traceability chain for measurement equipment, measurement uncertainty values, and approved quality standards.

Your local distributor:

Eppendorf Ibérica S.L.U. · Avenida Tenerife 2 · Edificio 1
 28703 San Sebastián de los Reyes, Madrid · Spain
 Tel.: +34 91 651 76 94 · eppendorf@eppendorf.es · www.eppendorf.es

www.eppendorf.com/epServices