

Centrifuge 5702 / 5702 R / 5702 RH



General Benefits

- > **Small footprint:** Medium capacity (40 x blood tubes) at a very compact footprint
- > **High versatility:** Spin a large variety of different blood and urine tubes including pediatric tubes (Microtainer® or MiniCollect®)
- > **Optimal conditions:** Optimum pellet formation due to swing-out rotor for contamination-free removal of the supernatant
- > **Strong sample safety:** Timer function and program buttons (Centrifuge 5702 R) for reproducible results

		Centrifuge 5702	Rotor A-4-38	Round bucket					Rectangular bucket			16 x 125 mm bucket		
				5702 737 003	5702 719 005	5702 735 000	5702 724 009	5702 723 002	5702 763 004	5702 710 008	5702 711 004	5702 712 000	5702 764 000	
Blood sample	BD Vacutainer® 	Microtainer®	16*	16*	–	–	–	–	36*	40*	–	–	–	
		1.8 – 4.8 mL (13 x 75 mm)	16	16	–	–	–	–	–	36	–	–	–	–
		4.0 – 7.0 mL (13 x 100 mm)	16	16	–	–	–	–	–	–	40	32	–	–
		8.0 – 11.0 mL (16 x 100 mm)	–	–	16	16	–	–	–	–	–	–	24	–
		CPT™ 4.0 mL (13 x 100 mm)	16	16	–	–	–	–	–	–	40	–	–	–
		CPT™ 8.0 mL (16 x 125 mm)	–	–	–	–	–	–	–	–	–	–	–	4
	Greiner® Vacuette 	MiniCollect®	16	16	–	–	–	–	–	36*	40*	–	–	–
		0.9 – 4.0 mL (13 x 75 mm)	16	16	–	–	–	–	–	36	40	–	–	–
		3.5 – 6.0 mL (13 x 100 mm)	16	16	–	–	–	–	–	–	40	–	–	–
		6.0 – 10.0 mL (16 x 100 mm)	–	–	16	16	–	–	–	–	–	–	24	–
	Sarstedt Monovette® 	1.6 – 4.0 mL (ø 11/13 mm)	16	–	–	–	–	8	–	36	–	–	–	–
		4.0 – 5.0 mL (ø 11 mm)	16	–	–	–	–	–	–	36	–	–	–	–
		4.0 – 10.0 mL (ø 13 mm)	–	16	–	–	–	–	–	–	20	32	–	–
		4.0 – 10.0 mL (ø 15 mm)	–	–	16	16	–	–	–	–	–	–	24	–
Streck Blood Collection Tube 	5.0 mL (13 x 75 mm)	16	16	–	–	–	–	–	36	40	–	–	–	
	10.0 mL (16 x 100 mm)	–	–	16	16	8	–	–	–	–	–	24	–	

* Manufacturer adapter needed

The Centrifuge 5702/5702 R/5702 RH is a non-automatic centrifuge for the separation of liquid substance mixtures from the human body and is specifically intended to be used as an accessory together with an in vitro diagnostic device and to enable its use according to its intended purpose. For in vitro diagnostic use, Eppendorf centrifuges are exclusively intended for indoor use by trained specialists.

Centrifuge 5702 / 5702 R / 5702 RH

Centrifuge 5702		Rotor A-4-38	Round bucket					Rectangular bucket			16 x 125 mm bucket	
			5702 737 003	5702 719 005	5702 735 000	5702 724 009	5702 723 002	5702 763 004	5702 710 008	5702 711 004	5702 712 000	5702 764 000
Urine sample	BD Vacutainer® 	4.0 – 6.0 mL (ø 13 mm)	20	16	–	–	–	36	–	–	–	–
		8.0 – 11.0 mL (ø 16 mm)	–	–	16	16	–	–	–	–	24	–
		8.0 – 9.5 mL (conical, ø 16 mm)	–	–	–	–	8	–	–	–	24	–
	Greiner Vacuvette® 	4.0 mL (ø 13 mm)	16	16	–	–	–	–	40	–	–	–
		6.0 – 10.0 mL (ø 16 mm)	–	–	16	16	–	–	–	–	24	–
	Sarstedt Monovette® 	3.2 – 6.0 mL (ø 13 mm)	16	16	–	–	–	–	40	–	–	–
		8.5 – 10.0 mL (ø 15 mm)	–	–	16	16	–	–	–	–	32	24
	Test Tube 	15.0 mL	–	–	–	–	8	–	–	–	24	–

Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com