


# VisionNize<sup>®</sup> Lab Suite

Overview of compatible Eppendorf devices

# How Devices Can Be Connected



- 1 VisioNize box (hardware hub)
- 2 VisioNize sense gateway router
- 3 VisioNize sense environmental monitoring sensors
- 4 VisioNize sense external temperature sensor



# VisioNize

## Device Overview

### Overview of compatible Eppendorf devices for VisioNize Lab Suite

This list contains Eppendorf devices which can be connected to the platform. Visit [www.eppendorf.com/visionize](http://www.eppendorf.com/visionize) for the most current compatibility list.

#### VisioNize touch enabled devices

> Eppendorf devices equipped with the VisioNize touch interface can be connected via ethernet cable:

##### Eppendorf ULT Freezers

CryoCube® F740i\*  
 CryoCube® F740hi\*  
 CryoCube® F740iw\*  
 CryoCube® F740hiw\*

\*Out-of-the-box compatibility coming soon

##### Eppendorf CO<sub>2</sub> Incubators

CellXpert® C170i

##### Eppendorf Shakers

Innova® S44i\*

##### Eppendorf Centrifuge

Centrifuge 5910 Ri\*\*

\*\*launch planned for April 2021

##### Eppendorf Cyclers

Mastercycler® X50a\*  
 Mastercycler® X50s\*  
 Mastercycler® X50p\*  
 Mastercycler® X50h\*  
 Mastercycler® X50i\*  
 Mastercycler® X50l\*  
 Mastercycler® X50r\*  
 Mastercycler® X50t\*

#### VisioNize box enabled

> Eppendorf devices which are connectable via VisioNize box to the VisioNize Lab Suite. The VisioNize box acts as a hardware hub and works as a translation system enabling data collection of the device data in VisioNize Lab Suite.

##### Eppendorf ULT Freezers

CryoCube® F570  
 CryoCube® F570h  
 CryoCube® F570hw  
 CryoCube® FC660  
 CryoCube® FC660h  
 CryCube® F101h

Innova® U101  
 Innova® U535  
 Innova® U570  
 Innova® U725  
 Innova® U725-G  
 Innova® C760  
 Innova® C585  
 Premium U410  
 Premium U570  
 HEF U410

##### Eppendorf Shakers

Innova® 40  
 Innova® 42  
 Innova® 43  
 Innova® 44  
 New Brunswick™ S41i

##### Eppendorf Thermomixer

Eppendorf ThermoMixer® C, F0.5, F1.5, F2.0, FP  
 Eppendorf ThermoStat™ C

##### Eppendorf Cyclers

Mastercycler® Nexus  
 Mastercycler® Nexus gradient  
 Mastercycler® Nexus flat  
 Mastercycler® Nexus X2  
 Mastercycler® Nexus GX2  
 Mastercycler® Nexus SX1  
 Mastercycler® Nexus GSX1

##### Eppendorf Centrifuges

Centrifuge 5910 R  
 Centrifuge 5424 R  
 Centrifuge 5425  
 Centrifuge 5425 R

##### Eppendorf CO<sub>2</sub> Incubators

Galaxy® 48R  
 CellXpert® C170

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
 Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)