

Pipette Calibration Service

Calibration Services – Ordering information

Service	Calibration measurements »As Is«	Calibration measurements »As Left«	Outside decontamination ²	Checking, cleaning, lubricating ⁴ before calibration	Single channel pipette/dispenser	4-channel pipette	6-channel pipette	8-channel pipette	12-channel pipette	16-channel pipette	24-channel pipette
ISO 8655 including PM⁵		10 per test volume ¹	included	included	0082000142	0082000342	0082000343	0082000165	0082000166	0082000185	0082000186
Price[^]					\$115.00	\$200.00	\$220.00	\$244.00	\$288.00	\$360.00	\$432.00
3×4 including PM⁵		4 per test volume ¹	included	included	0082000141	0082000340	0082000341	0082000163	0082000164	0082000183	0082000184
Price[^]					\$94.00	\$130.00	\$150.00	\$173.00	\$216.00	\$288.00	\$324.00

Calibration Service Packages – Ordering information

Service	Calibration measurements »As Is«	Calibration measurements »As Left«	Outside decontamination ² checking, cleaning, lubricating ⁴	ISO 17025 Accreditation calibration certificate ³	Single channel pipette/dispenser	4-channel pipette	6-channel pipette	8-channel pipette	12-channel pipette	16-channel pipette	24-channel pipette
ISO 8655 Verified	4 per test volume ¹	10 per test volume ¹	included		0082000234	0082000324	0082000325	0082000235	0082000236	0082000237	0082000238
Price[^]					\$141.00	\$234.00	\$261.00	\$293.00	\$347.00	\$425.00	\$520.00
ISO 17025 Without PM⁵	10 per test volume ¹			included	0082000304	0082000305	0082000306	0082000300	0082000301	0082000302	0082000303
Price^{6,^}					\$95.00	\$162.00	\$172.00	\$185.00	\$201.00	\$216.00	\$286.00
ISO 17025 Advanced	4 per test volume ¹	10 per test volume ¹	included	included	0082000227	0082000322	0082000323	0082000228	0082000229	0082000230	0082000231
Price[^]					\$171.00	\$259.00	\$285.00	\$309.00	\$369.00	\$453.00	\$545.00
ISO 17025 Professional	10 per test volume ¹	10 per test volume ¹	included	included	0082000220	0082000320	0082000321	0082000221	0082000222	0082000223	0082000224
Price[^]					\$187.00	\$280.00	\$307.00	\$339.00	\$392.00	\$471.00	\$566.00

Additional Services – Ordering information

Advanced Documentation , Excel file for calibration information on all pipettes/dispensers per order	0082000212
Price[^]	\$50.00
On-site verification services in your own laboratory	Inquire

Additional Services Charges

Labor Charge of repair services	Inquire
--	---------

¹ 10%, 50% and 100% of nominal volumes for variable volume instruments. ² Pipettes and dispensers and package equipment need to be decontaminated before sending in. ³ ISO 17025 accreditation calibration certificate is included in the second calibration after the maintenance where applicable. ⁴ Internal parts. ⁵ Preventative Maintenance. ⁶ Only available for new pipettes accredited calibration. [^] All prices in AUD and exclude GST.

ISO 8655

ISO 8655 defines the standard methodology by which a pipette should be calibrated at 3 volumes (10%, 50% and 100% of the nominal volume) and 10 measurements at each volume, and controls for room environment.

Qualification and Certification

ISO 17025

The ISO 17025 accreditation certifies the high level of competence and compliance, providing a traceability chain for measurement equipment, measurement uncertainty values, and approved quality standards.

Download a decontamination form at www.eppendorf.com/epServices.

Your local distributor: www.eppendorf.com/contact

Eppendorf South Pacific Pty. Ltd. · Level 1, 97 Waterloo Road · Macquarie Park, NSW 2113 · Australia info@eppendorf.com.au · www.eppendorf.com

www.eppendorf.com/epservices