



eProcurement with Eppendorf

Bringing Quality and Value to Your Purchasing

Eppendorf has gained the trust of generations of researchers by providing innovative, high quality and reliable products. With an Eppendorf eProcurement catalog, users can access the products they need via a more efficient, simple and intelligent purchasing process. We share our customer's goal of streamlining the procure-to-pay process and support an user-friendly and customizable eProcurement environment.

Our capabilities include:

- > Hosted Catalogs
- > Punchout Catalogs
- > Level II and Live Price
- > cXML PO's and Invoices
- > Quote to Order
- > Stock Availability
- > Quick Order
- > Hotlists



Our dedicated digital business team will work with you to create customized data formats that best fits your eProcurement environment, your purchasing requirements and your business needs. We are connected with dozens of eProcurement systems, including Coupa®, Jaggaer® (SciQuest), Unimarket, Advanced, Prendio®, and many more. You can rely on Eppendorf's expertise to help you achieve a seamless connectivity experience.

An Eppendorf catalog can bring added value to your purchasing and procurement by reflecting custom account pricing, or you can opt in to discounts with one of Eppendorf's competitively bid Group Purchasing Contracts.

We value your time and your business and want to become a partner to help achieve your procurement goals. Contact us at eprocurement@eppendorf.com if you are interested in learning more.

Through Our Liquid, Sample, and Cell Handling Products, Eppendorf Helps Scientists to Research, Discover, and Publish

<div data-bbox="82 1473 513 1534" data-label="Section-Header"> <h2> LIQUID HANDLING</h2> </div> <ul style="list-style-type: none"> > Manual pipettes > Electronic pipettes > Pipette tips > Automated pipetting systems > Dispensers > Combitips > Bottletop dispensers 	<div data-bbox="571 1473 960 1534" data-label="Section-Header"> <h2> CELL HANDLING</h2> </div> <ul style="list-style-type: none"> > Fermentors and bioreactors > Single-use reactors > CO2 incubators > Biological shakers > Micromanipulators > Microinjectors > Electroporators 	<div data-bbox="1024 1473 1481 1534" data-label="Section-Header"> <h2> SAMPLE HANDLING</h2> </div> <ul style="list-style-type: none"> > Centrifuges and vacuum concentrators > Rotors and accessories > Thermomixer > PCR devices > Ultra-low temperature freezers > Multiwell plates > Tubes
---	---	---

Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eprocurement@eppendorf.com · www.eppendorf.com

www.eppendorf.com/eprocurement

Coupa® is a registered trademark of COUPA SOFTWARE INCORPORATED, USA. Jaggaer® is registered trademark of JAGGAER, LLC, USA. Prendio® is a registered trademark of Prendio, LLC, USA.
 ENA.C1.0292.B.US.21 Eppendorf® and the Eppendorf Brand Design are registered trademarks of Eppendorf SE, Germany. Product appearance, specifications, and/or prices are subject to change without notice. All rights reserved including graphics and images. Copyright © 2022 by Eppendorf SE.