



EA MLA Signatory
Český institut pro akreditaci, o.p.s.
Olšanská 54/3, 130 00 Praha 3

issues

according to section 16 of Act No. 22/1997 Coll., on technical requirements for products, as amended

CERTIFICATE OF ACCREDITATION

No. 406/2021

Eppendorf Czech & Slovakia s.r.o.
with registered office Voděradská 2552/16, 251 01 Říčany, Company Registration No. 27939031

to the Calibration Laboratory No. 2404
Calibration Laboratory

Scope of accreditation:

Calibration of piston volumetric devices to the extent as specified in the appendix to this Certificate.

This Certificate of Accreditation is a proof of Accreditation issued on the basis of assessment of fulfillment of the accreditation criteria in accordance with

ČSN EN ISO/IEC 17025:2018

In its activities performed within the scope and for the period of validity of this Certificate, the Body is entitled to refer to this Certificate, provided that the accreditation is not suspended and the Body meets the specified accreditation requirements in accordance with the relevant regulations applicable to the activity of an accredited Conformity Assessment Body.

This Certificate of Accreditation replaces, to the full extent, Certificate No.: 560/2020 of 11. 9. 2020, or any administrative acts building upon it.

The Certificate of Accreditation is valid until: 27. 7. 2026

Prague: 27. 7. 2021



Lukáš Burda
Director of the Department
of Testing and Calibration Laboratories
Czech Accreditation Institute
Public Service Company



Accredited entity according to ČSN EN ISO/IEC 17025:2018:

Eppendorf Czech & Slovakia s.r.o.
Calibration Laboratory
Voděradská 2552/16, 251 01 Říčany

CMC for the field of measured quantity: Volume

Ord. number	Calibrated quantity / Subject of calibration	Nominal range		Parameter(s) of the meas. quantity	Lowest expanded measurement uncertainty specified ^{2,4}	Calibration principle	Calibration procedure identification ³	Work- place
		min.	max.	unit				
1	Mechanical and electronic piston pipettes, pipettes, dispensers, burettes, dilutors	0.1	0.5	μl	0.03 μl	Gravimetric comparison with a standard (distilled water)	SOP.KL ECS 11 (ČSN EN ISO 8655-6, EURAMET cg-19)	
		0.5 10,000	to to	μl μl	0.15 % + 0.15 μl 0.06 %			

¹ Asterisk at the ordinal number identifies the calibrations, which the Laboratory is qualified to carry out outside the permanent laboratory premises.

² The expanded measurement uncertainty is in accordance with ILAC-P14 and EA-4/02 M, part of CMC, and it is the lowest value of the respective uncertainty. If not stated otherwise, its coverage probability is approx. 95%. If not stated otherwise, the uncertainty values stated without a unit are relative to the value measured. If the calibration is carried out outside the laboratory premises, the measurement uncertainty may be affected.

³ If the document identifying the calibration procedure is dated, only these specific procedures are used. If the document identifying the calibration procedure is not dated, the latest edition of the specified procedure is used (including any changes).

⁴ The lowest uncertainty includes the effect of the operator.

