

Press release

Eppendorf adds digital documentation to connected electronic pipette system

Hamburg, March 2023

The Pipette Manager serves as a wireless control panel for Eppendorf's connected electronic pipettes, adding new features and revolutionizing the speed and reproducibility of handheld pipetting.

With this innovative system, scientists can easily adjust connected pipette modes, settings, and volumes via a large touch display.

The system's Liquid Manager feature offers predefined pipette settings and technique hints for a variety of liquids, ensuring consistent and accurate results in several applications. Users can also define their own liquid types with respective pipette settings.

The Pipetting Records feature offers a solution for proper documentation of pipetting workflows, allowing users to automatically record activities and export recordings as a PDF file via USB.

Previously introduced as VisioNize® pipette manager in Europe, the US and Canada, the Pipette Manager has been renamed to distinguish it from other digital offerings from Eppendorf, such as VisioNize® Lab Suite.

Eppendorf is dedicated to continuous improvement and has added the new Pipetting Records feature as a free software update for existing customers in Europe, the US and Canada. The product is now also available for the first time in Asia and Latin America.

With the integration of Move It® pipettes in 2022, the Pipette Manager is now compatible with all Eppendorf Xplorer® (plus) pipettes and includes several other features designed to help scientists manage their pipette fleet and collaborate with ease. This essential, easy-to-use tool is a must-have for any laboratory looking to speed up and improve the efficiency and reproducibility of their pipetting workflows.

www.eppendorf.com/pipette-manager

Pressekontakt:

Bettina Grochulla
Barkhausenweg 1
22339 Hamburg

Telefon: +49 40 692106268
grochulla.b@eppendorf.de

www.eppendorf.com

About Eppendorf SE:

Eppendorf is a leading life science company that develops, produces and sells systems for use in laboratories around the world. The product portfolio is divided into four core business areas: Liquid Handling, Consumables, Separation & Instrumentation and Bioprocess and Workflow Solutions. It includes pipettes, pipette tips, centrifuges, mixers and ultra-low-temperature freezers. In addition, Eppendorf offers a wide variety of high-quality consumables.

Eppendorf products are used in academic and commercial research laboratories, e.g., in companies from the pharmaceutical and biotechnological as well as the chemical and food industries. They are also used in laboratories that perform clinical or environmental analysis, in forensic laboratories and in industrial laboratories where process analysis, production and quality assurance are performed.

Eppendorf has been in business since 1945, is headquartered in Hamburg, Germany, and employs about 5,000 people around the world. The company has subsidiaries in 33 countries and is represented on all continents and in all important markets by distributors.

Pressekontakt:

Bettina Grochulla
Barkhausenweg 1
22339 Hamburg

Telefon: +49 40 692106268
grochulla.b@eppendorf.de

www.eppendorf.com