

Eppendorf Certificate

Certificate of Quality Eppendorf Laboratory Consumables

STATEMENT ON PER- AND POLYFLUORINATED ALKYL SUBSTANCES (PFAS)

Eppendorf states that the materials and laboratory consumables listed below are safe with respect to per- and polyfluorinated alkyl substances (PFAS) transmission when used in laboratory applications.

PFAS are industrial chemicals that are used in numerous industrial processes, and consumer products due to their special technical properties.

Some PFAS are already considered substances of very high concern (SVHC) under REACH because they are very persistent in the environment, accumulate in organisms, and can be harmful to humans. Epidemiological studies have shown effects on the thyroid, immune system and reproduction at very low dose levels.

Eppendorf only works with granulate suppliers who guarantee according to their Absence of Substances List that their raw materials are free of PFAS.

This certificate applies to the following Eppendorf Laboratory Consumables:

Pipette tips	epT.I.P.S.® *1 epT.I.P.S.® 384 epT.I.P.S.® Long epT.I.P.S.® Motion ep Dualfilter T.I.P.S.® *1 ep Dualfilter T.I.P.S.® 384 ep Dualfilter T.I.P.S.® SealMax® *1 GELoader® Microloader Eppendorf Serological Pipets
Positive displacement tips	Mastertip® Varitips® Combitips® advanced incl. adapters ViscoTip®
Eppendorf Tubes®	Eppendorf Tubes® 3810X/ Flex-Tube® Eppendorf Safe-Lock Tubes* Eppendorf DNA LoBind®/ Protein LoBind® Tubes, Eppendorf Tubes® 5 mL*,*1 incl. adapters Eppendorf Conical Tubes 15 mL*,*1, 50 mL*,*1, SnapTec® 50 Eppendorf Conical Tubes 25 mL*1 incl. adapters

* applies also to the SafeCode variants

*1 applies also to the BioBased variants

Eppendorf Certificate

Eppendorf Plates®	Eppendorf Microplates* Eppendorf Deepwell Plates* Eppendorf DNA LoBind®/ Protein LoBind® Plates Eppendorf Assay/Reader Microplates
PCR Consumables	Eppendorf twin.tec® PCR Plates*.*1 Eppendorf twin.tec® PCR Plates LoBind®*1 Eppendorf twin.tec® <i>microbiology</i> PCR Plates Eppendorf twin.tec® <i>real-time</i> PCR Plates PCR Tube Strips Fast PCR Tube Strips <i>real-time</i> PCR Tube Strips PCR Cap Strips PCR Tubes PCR Films & Foils
Cuvettes	UVette® Vis Cuvette
Sample Handling Consumables	wide-neck bottles

* applies also to the SafeCode variants


*1 applies also to the BioBased variants

Hamburg, February 2024

Page 2 of 2



Joana Tziolis
Product Life Cycle Manager
Division Consumables



Thomas Hengstmann
Head of Global Quality Operations

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified

Eppendorf®, the Eppendorf Brand Design, epT.I.P.S.®, ep Dualfilter T.I.P.S.®, GELoader®, Mastertip®, Varitips®, CCCadvanced®, Combitips®, ViscoTip®, Eppendorf Tubes®, Flex-Tube®, Eppendorf twin.tec®, LoBind®, SealMax® and UVette® are registered trademarks of Eppendorf SE, Germany.
All rights reserved, incl. graphics and images. Copyright ©2024 by Eppendorf SE.