


Eppendorf Certificate

Certificate of Purity - Biopur®

This package contains a high-quality consumable manufactured under the Biopur® Purity Standard. The Biopur consumables are produced in a controlled environment according to ISO class 8 of ISO 14644-1 .

For this product Eppendorf certifies the following:

- > Sterile
- > Pyrogen/Endotoxin-free
- > Human DNA-free
- > Bacterial DNA-free
- > DNase-free
- > RNase-free
- > PCR inhibitor-free
- > ATP-free



These parameters are continuously monitored by an independent certified laboratory. Eppendorf guarantees the conformity within the following limits:

Sterility	in accordance with USP, Ph. Eur. 2.6.12
Pyrogens	< 0.001 EU/mL (tested according to Ph. Eur. 2.6.14 (LAL test))
Human DNA	< 2 pg
Bacterial DNA	< 50 fg
DNase	< 1.0 x 10 ⁻⁷ Kunitz units
RNase	< 1.0 x 10 ⁻⁹ Kunitz units
PCR inhibition	fewer than 10 targets amplifiable
ATP	< 5.5 x 10 ⁻¹² mg

Lot-specific certificates can be downloaded from the internet at www.eppendorf.com/certificates.

Hamburg, June 2024

Joana Tziolis
Product Life Cycle Manager
Division Consumables

Thomas Hengstmann
Head of Global Quality Operations

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified