

ACT[®] Ecolabel

Eppendorf

**epT.I.P.S.[®] Motion as Reload
System, with filter, PCR clean and
sterile, 10 µL, (24 trays x 96 tips)**

Pipette Tips

SKU: 0030014561

Oldenburg in Holstein, Germany

**Environmental
Performance Factor**

37.5

Certified April 2026



Environmental Performance

Product

Recycled/Renewable Content	0%	■
Chemicals of Concern	Yes	■
Electricity Consumed/Day	N/A	
Water Consumed/Day	N/A	
Refrigerant GWP	N/A	
Use Phase	Single Use	■
Recyclable Materials*	0%	■
Circularity Support	None	■

Packaging

Recycled/Renewable Content	40%	■
Shipping Method	Ambient	■
Recyclable Materials*	90%	■

Product Carbon Reporting

Reporting Framework	None	☑
Review	N/A	■
Product CO ₂ e*	Not Available	■

Manufacturing Facility

Best Practices	8.0/10	■
Renewable Electricity	65%	■
Renewable Energy	31%	■

Company GHG Reductions

Scope 1/2/3 Tracking	Yes/Yes/Yes	■
Near-Term Target	Yes	■
Net-Zero Target	No	■

Improvements

Near-Term Target Commitment

Audit Details



Third Party Reviewer: Verico

*See Definitions