

Eppendorf Certificate

General Certificate of Purity – Protein-free

This package contains a high-quality consumable manufactured under the protein-free Eppendorf Purity Standard. The protein-free consumables are produced in a controlled clean room environment according to class 8 of ISO 14644-1.

For this product Eppendorf certifies the following:
> Protein-free



These parameters are continuously monitored by an independent certified laboratory. Eppendorf guarantees the conformity within the following limits:

Protein < 1 µg/mL

The protein determination was performed according to the method of Bradford, M. M. (1976): A rapid and sensitive method for the quantitation of microgram quantities of protein utilizing the principle of protein-dye binding. Analytical Biochemistry, 72, 248–254.

Lot-specific certificates can be downloaded from the internet at www.eppendorf.com/certificates

Hamburg, July 2026


Marlene Jentzsch (Jul 29, 2026 14:53:24 GMT+2)

Dr. Marlene Jentzsch
Senior Vice President Core Business Portfolio



Katy Fischer
Director Quality Management Liquid Handling

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

Eppendorf® and the Eppendorf Brand Design are registered trademarks of Eppendorf SE, Germany.
All rights reserved, incl. graphics and images. Copyright ©2026 by Eppendorf SE.

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified