

# ACT<sup>®</sup> Ecolabel

## Eppendorf Conical Tubes 25 mL with SnapTec Cap, DNA LoBind, PCR clean

Centrifuge/Storage Tubes  
SKU: 0030108523  
Oldenburg in Holstein, Germany

Environmental  
Performance Factor

36.5

Certified December 2025

### Environmental Performance

#### Product

Recycled/Renewable Content	0%	■
Chemicals of Concern	Yes	■
Electricity Consumed/Day	N/A	
Water Consumed/Day	N/A	
Refrigerant GWP	N/A	
Use Phase	Single Use (Storage)	■
Recyclable Materials*	0%	■
Circularity Support	None	■

#### Packaging

Recycled/Renewable Content	60%	■
Shipping Method	Ambient	■
Recyclable Materials*	80%	■

#### Product Carbon Reporting

Reporting Framework	None	☑
Review	N/A	■
Product CO <sub>2</sub> e*	Not Available	■

#### Manufacturing Facility

Best Practices	8.0/10	■
Renewable Electricity	65%	■
Renewable Energy	31%	■

#### Company GHG Reductions

Scope 1/2/3 Tracking	Yes/Yes/Yes	■
Near-Term Target	Yes	■
Net-Zero Target	No	■

#### Improvements

Near-Term Target Commitment

#### Audit Details



Third Party Reviewer: Verico

\*See Definitions

