

ACT[®] Ecolabel

Eppendorf Conical Tubes 25 mL with SnapTec cap, amber (light protection), sterile, pyrogen- DNase- RNase- and DNA-free

Centrifuge/Storage Tubes
SKU: 0030108430
Oldenburg in Holstein, Germany

Environmental
Performance Factor

36.5

Certified December 2025



Environmental Performance

Product

Recycled/Renewable Content	0%	■
Chemicals of Concern	Yes	■
Electricity Consumed/Day	N/A	
Water Consumed/Day	N/A	
Refrigerant GWP	N/A	
Use Phase	Single Use (Storage)	■
Recyclable Materials*	0%	■
Circularity Support	None	■

Packaging

Recycled/Renewable Content	60%	■
Shipping Method	Ambient	■
Recyclable Materials*	80%	■

Product Carbon Reporting

Reporting Framework	None	☑
Review	N/A	■
Product CO ₂ e*	Not Available	■

Manufacturing Facility

Best Practices	8.0/10	■
Renewable Electricity	65%	■
Renewable Energy	31%	■

Company GHG Reductions

Scope 1/2/3 Tracking	Yes/Yes/Yes	■
Near-Term Target	Yes	■
Net-Zero Target	No	■

Improvements

Near-Term Target Commitment

Audit Details



Third Party Reviewer: Verico
*See Definitions