

ACT[®] Ecolabel

Eppendorf

Conical Tubes 25 mL with screw cap, Protein LoBind, PCR clean

Centrifuge/Storage Tubes

SKU: 0030122283

Oldenburg in Holstein, Germany

Environmental
Performance Factor

37.5

Certified December 2025

Environmental Performance

Product

Recycled/Renewable Content	0%	■
Chemicals of Concern	Yes	■
Electricity Consumed/Day	N/A	
Water Consumed/Day	N/A	
Refrigerant GWP	N/A	
Use Phase	Single Use (Storage)	■
Recyclable Materials*	0%	■
Circularity Support	None	■

Packaging

Recycled/Renewable Content	70%	■
Shipping Method	Ambient	■
Recyclable Materials*	80%	■

Product Carbon Reporting

Reporting Framework	None	<input type="checkbox"/>
Review	N/A	■
Product CO ₂ e*	Not Available	■

Manufacturing Facility

Best Practices	8.0/10	■
Renewable Electricity	65%	■
Renewable Energy	31%	■

Company GHG Reductions

Scope 1/2/3 Tracking	Yes/Yes/Yes	■
Near-Term Target	Yes	■
Net-Zero Target	No	■

Improvements

Near-Term Target Commitment

Audit Details



Third Party Reviewer: Verico

*See Definitions

